
4R10 Quickchange



DE	Gebrauchsanweisung	11
EN	Instructions for use	20
FR	Instructions d'utilisation	29
IT	Istruzioni per l'uso	39
ES	Instrucciones de uso	48
PT	Manual de utilização	57
NL	Gebruiksaanwijzing	66
SV	Bruksanvisning	76
DA	Brugsanvisning	85
NO	Bruksanvisning	94
FI	Käyttöohje	103
PL	Instrukcja użytkowania	112
HU	Használati utasítás	121
CS	Návod k použití	130
RO	Instrucțiuni de utilizare	139
HR	Upute za uporabu	148
SL	Navodila za uporabo	157
SK	Návod na používanie	166
BG	Инструкция за употреба	175
TR	Kullanma talimatı	184
EL	Οδηγίες χρήσης	193
RU	Руководство по применению	203
JA	取扱説明書	214
ZH	使用说明书	222
KO	사용 설명서	229

1 **5 Lieferumfang**
5 Scope of delivery

4R10=111

4R10=01


647G1248

2

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6 Technical data

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INFORMATION

Date of last update: 2020-04-29

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

INFORMATION

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings.

1.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

1.2 Meanings of pictograms in the illustration

	Numbering for the illustrations		Reference to the number of the corresponding section
	Numbering for a defined sequence		Numbering for the parts of an illustration

2 Product description

The product (4R10) consists of an upper section, the locking pin and a lower section. The 4R10=111 package contains two lower sections. This allows the patient to quickly change the distal prosthetic components without tools. For example, prosthetic feet for different heel heights and fields of application (everyday prosthesis, sports prosthesis – Required combinations: 4R10 + 1E95 CHALLENGER, bathing prosthesis) can be used.

3 Intended use

3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

3.2 Area of application

 CAUTION
Excessive strain on the product
Fall due to breakage of load-bearing components
► Only use the product according to its allowable field of application.

Allowable amputation types

Transtibial amputation, Knee disarticulation, Transfemoral amputation, Hip disarticulation

Allowable field of application

Everyday prosthesis, Bathing prosthesis, Sport prosthesis (Required combinations: 4R10 + 1E95 CHALLENGER)

Allowable body weight

4R10: 125 kg/275 lbs

4R10 + 1E95 CHALLENGER: 100 kg/220 lbs

Recommended mobility grade

2 – 4

Positioning in the prosthesis

Only according to the sections and illustrations in the instructions for use

3.3 Combination possibilities

CAUTION

Failure to observe the manufacturer requirements regarding combination possibilities

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- ▶ Contact the manufacturer with questions.

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

3.4 Environmental conditions

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Storage and transport

Temperature range -20 °C to $+60\text{ °C}$ (-4 °F to $+140\text{ °F}$), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature range: -10 °C to $+45\text{ °C}$ (14 °F to 113 °F)

Chemicals/liquids: fresh water, salt water, perspiration, urine, soapsuds, chlorine water

Allowable environmental conditions

Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions

Solids: Dust, occasional contact with sand

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see Page 27).

Unallowable environmental conditions

Solids: Highly hygroscopic particles (e.g. talcum), continuous contact with sand

Chemicals/liquids: Acids, continuous use in liquid media

3.5 Service life

CAUTION

Exceeding the service life

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the service life defined in this section is not exceeded.

This product has been tested by the manufacturer with 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

4 General safety instructions

CAUTION

Product contains a magnet

Influence on or damage to devices or items due to a strong magnetic field

- ▶ Keep the product away from devices and items that are sensitive to magnetic fields.

CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

5 Scope of delivery

The scope of delivery is shown in Figure 1 on page 2.

In the illustration, a reference number shown in bold stands for a package in which the product components shown in the frame can be reordered.

Only product components with reference numbers in the illustration can be reordered separately.

INFORMATION

After removing the 4X345 locking pin, the upper section and lower section can be separated from each other. This requires force because magnets are installed in the upper section.

6 Technical data

The technical data are shown in Figure 2 on page 3.

①	System height (SH): Without spacer plate: 12 mm With spacer plate: 15 mm	②	Allowable body weight: 125 kg/275 lbs
③	Mobility grade: 2 – 4	④	Weight: 220 g

⑤: Unallowable combinations

7 Preparing the product for use

7.1 Information on fabrication of a prosthesis

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

► Observe the alignment and assembly instructions.

With the 4R10, the patient can adapt the prosthesis to various fields of application by quickly changing the lower section (3 on page 4 – ①: **alignment examples**; ②: different heel heights; ③: various fields of application (Everyday prosthesis/Sport prosthesis); ④: with prosthetic knee joint). It is important to complete the bench alignment, optimisation of the static align-

ment and optimisation during the dynamic trial fitting for each field of application.

- ▶ Carefully plan the alignment of the prosthesis for the respective fields of application according to the sections that follow and the corresponding illustrations as well as the instructions for use of the prosthetic components that are used.

In doing so, determine the installation dimensions (system heights/structural heights) of the prosthetic foot and the other prosthetic components, and coordinate them with each other.

7.2 Bench alignment of the prosthesis

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions for thread lock.

Positioning the upper section and lower section (4 on page 5 – 1)

- ▶ For installation in the prosthesis, note the positioning of the upper section relative to the lower section (2: incorrect; 3: correct).

Prohibited positioning (5 on page 6 – 4)

- ▶ Do not position the product in the prosthesis as shown in Figure 5 on page 6 (5: prohibited; 6: proximal from the prosthetic knee joint).

Allowable positioning (6 on page 7 – 7)

- ▶ **INFORMATION:** For bench alignment, the 4R10 is positioned at the same height for each version. The distance to the prosthetic socket ($h_1 = h_2$) or to the alignment reference point of the prosthetic knee joint ($h_3 = h_4$) is always the same. The instructions for use of prosthetic sport feet provide information whether a different height has to be used for the bench alignment of a sports prosthesis due to the spring effect.

Position the product in the prosthesis according to the examples in Figure 6 on page 7 (9: allowable).

- 8: Recommended position (close to the TT prosthetic socket or prosthetic knee joint) – this positioning offers the most alignment possibilities.

- ⑩: Alternative position (close to the prosthetic foot) – this positioning shifts the centrifugal mass in the distal direction, but reduces the alignment possibilities and makes it more difficult to install a cosmetic cover.

Using the spacer plate (7 on page 8 – 11)

- **INFORMATION:** The system height can be changed with the spacer plate.

Only use or omit the spacer plate according to the illustrated, approved examples in reference to the socket adapter (12: flat contact surface; 13: allowable; 14: no flat contact surface; 15: prohibited).

Mounting the socket adapter (8 on page 9 – 16)

- **CAUTION!** Thread lock is only required when finishing the prosthesis. Install the socket adapter in the illustrated sequence.

7.3 Optimising the static alignment

- Optimise the static alignment of the prosthesis with its various distal prosthetic components according to the values in the instructions for use of the products that are used.

7.4 Optimising during dynamic trial fitting

- During the dynamic trial fitting, test and optimise the prosthesis with its various distal prosthetic components.

CAUTION! If the settings of the prosthesis with its various distal prosthetic components exhibit opposing effects, the safety-related aspects always have to take precedence over everything else in choosing a compromise setting.

7.5 Attaching the cosmetic cover

- Fabricate the cosmetic cover so that the quick-change function can be operated reliably.

7.6 Finishing the prosthesis

- To finish the prosthesis, all screw connections have to be tightened according to the prescribed specifications with regard to Nm and the use of thread lock, and the correct functioning of the prosthesis has to be tested with its various distal prosthetic components.

8 Use

CAUTION

Mechanical overload

Falling due to breakage of load-bearing components, impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Using the quick-change function

CAUTION

Using the product without tightening the locking pin

Risk of injury due to prosthetic components coming loose

- ▶ Before using the prosthesis, verify that the locking pin is tightened.
- ▶ Perform the quick change of the distal prosthetic components according to the sequence shown in Figure  on page 10.
 - ①: In doing so, tighten the locking pin by hand.

9 Cleaning

CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
- ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.

Cleaning light soiling

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

Cleaning heavier soiling

- > **Required materials:** cleaning cloth, 634A58 isopropyl alcohol, soft cloth
- ▶ **NOTICE!** Verify material compatibility.
Clean the product with a cleaning cloth and isopropyl alcohol.
- ▶ Dry the product with a cloth.

Cleaning with Disinfectants

- > **Required materials:** colourless, alcohol-free disinfectant (verify material compatibility!), soft cloth
- ▶ Disinfect the product with the disinfectant.
- ▶ Dry the product with a cloth.
- ▶ Allow to air dry in order to remove residual moisture.

Cleaning after contact with salt water or contaminated fresh water

- ▶ Rinse the product with clean, fresh water.
- ▶ Wipe the product with a damp, soft cloth.
- ▶ Dry the product with a soft cloth.
- ▶ Allow to air dry in order to remove residual moisture.

10 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.
- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

11 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

12 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

12.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

12.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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