**ottobock.**

***Skeo Sealing Liner.***

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DE

***Skeo Sealing liner.***

Notes for socket fabrication.

EN

***Manchon Skeo Sealing.***

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# General information.

## **For what residual limb shape is the *Skeo Sealing liner* suited?**

The *Skeo Sealing liner* for transfemoral amputees can be used for almost any residual limb shape, as long as there is sufficient soft tissue coverage. However, short and very conical residual limb shapes are contraindicated.

## **How do I choose the right ring height (10 or 17 cm)?**

For up to medium long residual limbs, please choose the 10 cm ring height (6Y110=\*-10/6Y111=\*-10); for longer residual limbs or knee disarticulation, use the 17 cm version (6Y110=\*-17/6Y111=\*-17).

## **What should I consider for long residual limbs and/or extremely tight sockets?**

Thanks to its silky-smooth outer coating, it is very easy to put on the *Skeo Sealing liner* and step into the socket. For this reason, a donning spray is generally not needed. However, for long residual limbs and/or extremely tight prosthetic sockets, we recommend using 640F18 donning spray when putting on the prosthesis. To do so, apply the donning spray to the sealing lip and the inside of the prosthetic socket to make the prosthesis easier to put on.

## **Can I use the *Skeo Sealing liner* in an existing socket?**

If the socket was originally designed for use with a different liner, the *Skeo Sealing liner* may not seal properly, since the material properties and wall thicknesses of other liners are not a 100% match. Therefore, the functionality of the *Skeo Sealing liner* should be checked before using it. A new socket may be required.

## Can I use the *Skeo Sealing liner* with waterproof walking aids?

Yes, the *Skeo Sealing liner* is suitable for use with waterproof walking aids.

# Preparing the plaster cast.

### Materials:

- Water basin (754W5)
- Scissors (719G1 or 719G2)
- Indelible pencil
- Cotton stockinette – width of 15 and 30 (623T1=15/30)
- Biplatrix plaster bandages – width of 15 and 20 (699G20=3X15/20)
- 642F18 Stretch film

- Trim the medial edge of the liner so it runs 1 cm below the perineum.
- The course of the lateral liner edge is at the same height as the edge of the prosthetic socket.
- Isolate the trimmed and donned liner on the outside with soap solution.
- Pull a cotton stockinette (1 – width of 15) over the residual limb. The distal half of the residual limb does not have a stockinette, so the distal part of the plaster cast will remain very smooth on the inner side. This is very important for modelling the plaster model.
- Cut an oval hole in the lower part of a cotton stockinette (width of 30)
- and put it on the patient as a protective layer (2). The cut hole serves as a reference line for the trim line in the posterior area below the gluteal muscles (3 – finger).
- Take measurements of the residual limb and document them.
- Secure the sealing lip with 642F18 stretch film.
- Important: Measure the first residual limb circumference 3 cm below the perineum. This measurement is crucial for modification of the plaster model.

Plaster the volume of the residual limb by wrapping it in circles up to the height of the trochanter.

To protect the distal end of the femur, mould the end of the residual limb slightly in the lateral and posterior direction by positioning your hands accordingly.

Hold the hand that is encompassing the residual limb from lateral to posterior above the end of the femur horizontal to facilitate better shaping of the soft tissue. Shape a slight support area.

The end result of the first plaster casting phase.

Place a four-layer plaster splint (longuette, width of 20) around the residual limb to take a cast of the socket entrance area, making sure that the ischial tuberosity and the lateral area of the subsequent prosthetic socket are included in the casting.


Grip the ischial tuberosity from the frontal direction (4) so your hand makes a 45° movement up towards the residual limb (5), and then execute two 45° tilting movements (6) to shape the tuberosity support.

The other hand lies flat against the femur, positioned directly below the greater trochanter.

## Modelling the plaster positive.


### View from medial:

Trim line for embedding the ramus and the adductor tendon are at the same height (7).


Reduce the circumference in the area below the medial trim line to the following values:

- Proximal area ⑩: residual limb circumference – max. 1 cm
- Distal area ⑨: do not reduce the residual limb circumference based on the residual limb measurements, only reduce at the height of the sealing ring

### View from posterior:

Sand the posterior course of the socket ⑩ so that it is nearly parallel to the anterior course of the


### View from anterior:

A flat shape of the medial trim line is important so the prosthetic socket is comfortable to put on (⑪ – height). Apply plaster accordingly as required ⑬. Flatten the anterior support area ⑭. To relieve the end of the femur, maintain the curve of the socket end in the lateral and posterior directions.


There should be a maximum of 2 cm between the medial socket brim and the tip of the greater trochanter **12**.

### View from lateral:

The anterior and posterior course of the socket are nearly parallel to each other.

Remove the bump caused by the lip **15**.

Reduce the residual limb length by 0.5–1 cm **16**:

- – 0.5 cm if there is virtually no soft tissue over the end of the femur
- – 1 cm if there is a long section of surplus soft tissue over the end of the femur

Otherwise, do not reduce the circumferences in the distal area.

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