Genium 3B1-2/3B1-2=ST

DE Gebrauchsanweisung (Benutzer)	3
EN Instructions for use (user)	51
FR Instructions d'utilisation (Utilisateur)	99
IT Istruzioni per l'uso (Utilizzatore)	149
ES Instrucciones de uso (Usuario)	197
PT Manual de utilização (Usuário)	245
NL Gebruiksaanwijzing (Gebruiker)	293
HR Upute za uporabu (Korisnik)	343
TR Kullanma talimatı (Kullanıcı)	389
JA 取扱説明書(ユーザー用)	435

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1	Important information.....	54
2	Foreword	55
3	Product description	55
3.1	Design	55
3.2	Function.....	55
4	Application.....	56
4.1	Indications for use	56
4.2	Area of application.....	56
4.3	Conditions of use	56
4.4	Indications.....	56
4.5	Qualification	57
5	Safety.....	57
5.1	Explanation of warning symbols.....	57
5.2	Structure of the safety instructions.....	57
5.3	General safety instructions.....	57
5.4	Information on the Power Supply/Battery Charging.....	60
5.5	Battery charger information.....	60
5.6	Information on Proximity to Certain Areas	61
5.7	Information on Use	62
5.8	Notes on the safety modes.....	63
5.9	Instructions for use with an osseointegrated implant system.....	64
5.10	Information on the use of a mobile device with the cockpit app.....	64
6	Scope of Delivery and Accessories.....	65
7	Charging the battery	65
7.1	Connecting the power supply and battery charger	66
7.2	Charging the prosthesis battery.....	66
7.3	Display of the current charge level	67
7.3.1	Display of battery charge level without additional devices	67
7.3.2	Display of the current charge level using the Cockpit app	67
7.3.3	Display of current charge level using remote control (optional accessory)	67
7.3.4	Display of the current charge level during the charging process	68
8	Cockpit app.....	68
8.1	System Requirements	68
8.2	Initial connection between cockpit app and prosthesis	69
8.2.1	Starting the cockpit app for the first time	69
8.3	Control elements for cockpit app.....	70
8.3.1	Cockpit app navigation menu.....	70
8.4	Managing components.....	70
8.4.1	Adding component	71
8.4.2	Deleting a component	71

8.4.3	Connecting component with multiple devices.....	71
9	Use	72
9.1	Movement patterns in basic mode (mode 1)	72
9.1.1	Standing	72
9.1.1.1	Stance function	72
9.1.2	Walking.....	73
9.1.3	Running short distances ("walk-to-run" function)	73
9.1.4	Sitting down	73
9.1.5	Sitting	73
9.1.5.1	Sitting function.....	74
9.1.6	Standing up	74
9.1.7	Walking up stairs step-over-step.....	74
9.1.8	Overcoming obstacles	75
9.1.9	Walking down stairs	76
9.1.10	Walking down a ramp	76
9.2	Changing prosthesis settings.....	76
9.2.1	Changing the prosthesis setting using the cockpit app.....	77
9.2.2	Overview of adjustment parameters in basic mode.....	77
9.2.3	Overview of adjustment parameters in MyModes	78
9.3	Turning Bluetooth on the prosthesis on/off.....	79
9.3.1	Switching Bluetooth off/on using the cockpit app	80
9.4	Querying the prosthesis status	80
9.4.1	Query status through cockpit app.....	80
9.4.2	Status display in the cockpit app	80
9.4.3	Status information on the display of the remote control (optional accessory)	80
9.5	Mute mode (silent mode)	81
9.5.1	Turning mute mode on/off using the Cockpit app	81
9.5.2	Turning mute mode on/off using the remote control	81
9.6	Deep sleep mode	81
9.6.1	Turning deep sleep mode on/off using the Cockpit app.....	82
9.6.2	Turning deep sleep mode on/off using the remote control.....	82
9.7	OPG function (Optimised Physiological Gait).....	82
10	MyModes	83
10.1	Switching MyModes with the cockpit app	83
10.2	Switching MyModes using motion patterns	84
10.3	Switching from a MyMode back to basic mode.....	85
11	Additional operating states (modes)	86
11.1	Empty battery mode	86
11.2	Mode for charging the prosthesis	86
11.3	Safety mode.....	86
11.4	Overheating mode	87
12	Storage and bleeding	87
13	Cleaning	87

- 14 Maintenance87**
- 15 Legal information88**
 - 15.1 Liability88
 - 15.2 Trademarks88
 - 15.3 CE Conformity88
 - 15.4 Local Legal Information88
- 16 Technical data89**
- 17 Appendices91**
 - 17.1 Symbols Used91
 - 17.2 Operating states/error signals92
 - 17.2.1 Signals for operating states92
 - 17.2.2 Warnings/error signals93
 - 17.2.3 Error messages while establishing a connection with the cockpit app95
 - 17.2.4 Status signals96

1 Important information

INFORMATION

Remote control function only possible with version 1.0.0.6 or above

Please note that you can only use this knee joint and the updated one with the **4X350 remote control in version 1.0.0.6 and higher.**

Therefore, check the version of the remote control as follows:

- ▶ With the remote control switched on, use the ▲ and ▼ keys to select the menu option **RC settings** and confirm with the ■ key.
- ▶ Use the ▲ and ▼ keys to select the menu option **Info** and confirm this with the ■ key.
- ▶ Use the ▲ and ▼ keys to select the **RC version: 1.0.0.4** menu option. If an older version appears than **1.0.0.6**, the remote control must be updated by an authorised Ottobock Service Centre.

The function of the remote control with this or the updated knee joint is only possible with version 1.0.0.6 or higher.

2 Foreword

INFORMATION

Date of last update: 2017-02-06

- ▶ Please read this document carefully before using the product.
- ▶ Have yourself instructed by qualified personnel in the proper and safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). The contact information for the manufacturer is found on the back page.
- ▶ Please keep this document in a safe place.


The product "Genium" is called the product/prosthesis/knee joint/component in the following. These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

3 Product description

3.1 Design

The product consists of the following components:

1. Proximal pyramid adapter
2. Optional flexion stops
3. Battery
4. Hydraulic unit
5. LED (blue) as indicator for the Bluetooth connection
6. Receiver of the inductive charging unit
7. Distal tube clamp screws
8. Connecting cable for tube adapter

3.2 Function

This product features microprocessor control of the stance and swing phase.

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the damping behaviour of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

Thanks to the microprocessor-controlled stance and swing phase, the system can be individually adapted to your needs.

For this purpose, the product is adjusted by qualified personnel using adjustment software.

The product features MyModes for special motion types (e.g. cross-country skiing, ...). These are pre-configured by the prosthetist using the adjustment software and can be activated with special motion patterns, the cockpit app and the remote control (optional accessory) (see Page 83).

In case of a system malfunction, safety mode makes restricted operation possible. Predefined resistance parameters are configured in the product for this purpose (see Page 86).

Empty battery mode permits safe walking when the battery is drained. Resistance parameters that are predefined by the product are configured for this purpose (see Page 86).

The microprocessor-controlled hydraulic unit offers the following advantages

- Approximation of the physiological gait pattern
- Stability while standing and walking
- Adaptation of product characteristics to various surfaces, inclines, gait situations and walking speeds

4 Application

4.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

4.2 Area of application

Area of application according to Ottobock's MOBIS classification:

Knee joint with attached 2R20 AXON tube adapter

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **150 kg (330 lbs)**.

Knee joint with attached 2R21 AXON tube adapter with torsion

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of up to **125 kg (275 lbs)**.

4.3 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see Page 89).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

4.4 Indications

- For users with knee disarticulation, transfemoral amputation and hip disarticulation (users with hip disarticulation or hemipelvectomy have to be fitted with the 7E10=* Helix^{3D}hip joint).
- For unilateral or bilateral amputation
- Dysmelia patients with residual limb characteristics corresponding to knee disarticulation, transfemoral amputation or hip disarticulation

- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

4.5 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

If the product is to be connected to an osseointegrated implant system, the qualified personnel must also be authorised for the connection to the osseointegrated implant system.

5 Safety

5.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

5.2 Structure of the safety instructions

 CAUTION The heading describes the source and/or the type of hazard The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible: > E.g.: Consequence 1 in case of failure to observe the hazard > E.g.: Consequence 2 in case of failure to observe the hazard ► This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

5.3 General safety instructions

 WARNING Non-observance of safety instructions Injury due to damage or unexpected product behaviour. ► Follow the safety instructions and the precautions specified in this accompanying document.

 WARNING Using the prosthesis while operating a vehicle Accidents due to unexpected behaviour of the prosthesis because of changed damping behaviour. ► All users are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre. ► Observe national legal regulations for retrofitting your vehicle in accordance with the type of fitting.

 WARNING Use of damaged power supply unit, adapter plug or battery charger Risk of electric shock due to contact with exposed, live components. ► Do not open the power supply unit, adapter plug or battery charger.
--

- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ CAUTION

Failure to observe warning/error signals

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see Page 93) and corresponding change in damping settings must be observed.

⚠ CAUTION

Failure to observe activated mute mode (silent mode)

Falling due to unexpected product behaviour because of changed damping behaviour.

The following feedback signals are deactivated when mute mode is activated:

- > Long vibration signal if the hydraulic unit overheats.
- > Beep and vibration signal to confirm that the movement pattern has been recognised (switching to a MyMode/basic mode with movement pattern).
- > Beep and vibration signal to indicate successful switching to a MyMode/basic mode.
- > Beep and vibration signal upon successfully switching to deep sleep mode.
- ▶ Before activating mute mode, note that these feedback signals will be deactivated. For more information about mute mode, see the section "Mute mode" (see Page 81).
- ▶ After switching to a MyMode/basic mode, verify the changed damping settings.
- ▶ Ensure that you stand securely during all switching processes.
- ▶ To deactivate mute mode, connect and then disconnect the battery charger.

⚠ CAUTION

Independent user manipulation of system components

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Mechanical stress on the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Use of the product when battery charge level is too low

Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the prosthesis if required.

- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

CAUTION

Risk of pinching in the joint flexion area

Injuries due to pinching of body parts.

- ▶ Ensure that fingers/body parts or soft tissue of the residual limb are not in this area when bending the joint.

CAUTION

Penetration of dirt and humidity into the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that neither solid particles, foreign objects nor liquids (such as body and/or wound fluids) penetrate into the product.
- ▶ The knee joint and AXON tube adapter are weatherproof but not resistant to corrosion. Therefore, the knee joint and AXON tube adapter should not come into contact with salt water, chlorinated water or other solutions (such as soap or shower gel, and body and/or wound fluids). Do not use the knee joint under extreme conditions like diving or jumping into water. The knee joint and the AXON tube adapter are protected against splashed water from all directions, but not designed for underwater use.
- ▶ After contact with water, remove the Protective Cover (if installed) and hold the prosthesis with the sole of the foot facing up until the water has drained from the knee joint and AXON tube adapter. Dry the knee joint and components with a lint-free cloth and allow the components to fully air dry.
- ▶ Should the knee joint or AXON tube adapter come into contact with **salt water, chlorinated water or other solutions** (such as soap or shower gel, and body and/or wound fluids), **promptly** remove the Protective Cover (if installed) and clean the **knee joint**. In order to do so, wipe knee joint and tube adapter with a damp cloth (fresh water) and let them dry.
- ▶ In case of a malfunction after drying, the knee joint and AXON tube adapter must be inspected by an authorised Ottobock Service Center. The O&P professional is your contact person.
- ▶ The knee joint and the AXON tube adapter are not protected against submersion and the penetration of jets of water and steam.

CAUTION

Signs of wear on system components

Falling due to damage or malfunction of the product.

- ▶ In the interest of the user's own safety and in order to maintain the operational reliability and the validity of the warranty, the specified service intervals must be observed.

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product with a damp cloth only (fresh water).

INFORMATION

Knee joint movement noise

When using exoprosthetic knee joints, servomotor, hydraulic, pneumatic or brake load dependent control functions can cause movement noise. This kind of noise is normal and unavoidable. It generally does not indicate any problems. If movement noise increases noticeably during the lifecycle of the knee joint, the knee joint should be inspected by an O&P professional immediately.

5.4 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the product without taking it off

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ For safety reasons, take the product off before charging the battery.

⚠ CAUTION

Charging the product with damaged power supply unit/charger/charger cable

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

5.5 Battery charger information

NOTICE

Penetration of dirt and humidity into the product

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Mechanical stress on the power supply/battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

NOTICE

Operating the power supply unit/charger outside of the permissible temperature range

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see Page 89).

NOTICE**Independent changes or modifications carried out to the battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

NOTICE**Contact of the battery charger with magnetic data storage devices**

Wiping of the data storage device.

- ▶ Do not place the battery charger on credit cards, diskettes, audio or video cassettes.

5.6 Information on Proximity to Certain Areas**⚠ CAUTION****Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.50 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.35 m
 - DECT cordless phones incl. base station: 0.18 m
 - WiFi (routers, access points,...): 0.11 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.11 m

⚠ CAUTION**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for people (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations...).
- If this cannot be avoided, ensure at least that you have a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected changes in the damping behaviour of the product.

⚠ CAUTION**Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)**

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of strong magnetic fields.
- ▶ Take off the product before entering a room or area with strong magnetic fields and store the product outside this room or area.
- ▶ Damage to the product caused by exposure to strong magnetic fields cannot be repaired.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see Page 89).

5.7 Information on Use

⚠ CAUTION

Walking up stairs

Falling due to foot placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking up stairs and place most of the area of the sole of your foot on the stair surface.
- ▶ Particular caution is required when carrying children up stairs.

⚠ CAUTION

Walking down stairs

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and roll over the edge of the step with the middle of the shoe.
- ▶ Observe the warning/error signals (see Page 93).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warning and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

⚠ CAUTION

Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- ▶ Be sure to pay attention when pulsating vibration signals start. They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you have to reduce the activity level so the hydraulic unit can cool down.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.
- ▶ If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. In this case, the product should be checked for damage by a prosthetist. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Overloading due to unusual activities

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, paragliding, etc.).

- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by a prosthetist. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Improper mode switching

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.
- ▶ Take the weight off the product and correct the switching, if required.

⚠ CAUTION

Improper use of the stance function

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using the stance function and check the lock of the knee joint before placing your full weight on the prosthesis.
- ▶ Make sure you have been instructed in the proper use of the stance function by a prosthetist or therapist. Information on the stance function see Page 72.

⚠ CAUTION

Quickly pushing the hip forward with the prosthesis extended (e.g. serve while playing tennis)

- > Falling due to unexpected activation of the swing phase.
- ▶ Note that the knee joint may flex unexpectedly when the hip is pushed forward quickly while the prosthesis is extended.
- ▶ Therefore, familiarise yourself with swing phase activation in such situations under secure conditions (e.g. while holding on to parallel bars) and with the instruction of trained, qualified personnel.
- ▶ For sports where this movement pattern may occur, use a corresponding pre-configured MyMode. For further information about the MyModes, see the section 'MyModes' (see Page 83).

5.8 Notes on the safety modes

⚠ CAUTION

Using the product in safety mode

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see Page 93) have to be observed.
- ▶ Particular caution is necessary when using a bicycle without a freewheel (with a fixed gear).

⚠ CAUTION

Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Do not continue using the defective product.

- ▶ Contact your prosthetist immediately.

⚠ CAUTION

Safety mode cannot be deactivated

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ If you cannot deactivate safety mode by recharging the battery, a permanent malfunction has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. The prosthetist is your contact person.

⚠ CAUTION

Safety signal occurs (ongoing vibration)

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warnings/error signals (see Page 93).
- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. In this case, the contact person is your prosthetist.

5.9 Instructions for use with an osseointegrated implant system

⚠ WARNING

High mechanical loads due to normal or unusual situations, such as falling

- > Overloading of the bone, which can lead to pain, loosening of the implant, necrosis or fracture among other things.
- > Damage or breakage of the implant system or its components (safety components...).
- ▶ Verify compliance with the fields of application, conditions of use and indications according to the information of the manufacturers, both for the knee joint and for the implant system.
- ▶ Note the instructions of the clinical personnel that indicated the use of the osseointegrated implant system.
- ▶ Note changes in your state of health that result in restrictions or doubt regarding the use of the osseointegrated connection.

5.10 Information on the use of a mobile device with the cockpit app

⚠ CAUTION

Improper use of the device

Falling due to altered damping behaviour as a result of unexpected switching into MyMode.

- ▶ Make sure you have been instructed in the proper use of the device with the cockpit app.

⚠ CAUTION

Independently applied changes or modifications made to the device

Falling due to altered damping behaviour as a result of unexpected switching into MyMode.

- ▶ Do not make any independent changes to the hardware of the device.
- ▶ Do not make any independent changes to the software/firmware of the device which are not included in the update function of the software/firmware.

⚠ CAUTION

Improper mode switching with the device

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter (beeper) and the display on the device.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.

NOTICE

Destruction of the device due to falling or penetration of water

Malfunction of the device.

- ▶ If necessary, let the device dry at room temperature (at least 1 day).
- ▶ If it should no longer be possible to switch back from a MyMode to basic mode, the component can only be switched back to basic mode by using a movement pattern (see Page 85) or by connecting/disconnecting the inductive battery charger.

NOTICE

Failure to observe the system requirements for the installation of the cockpit app

Malfunction of the device.

- ▶ Install the cockpit app only on the operating systems listed in the section "System requirements" (see Page 68). The tested devices are listed in this section as well.

6 Scope of Delivery and Accessories

Scope of Delivery

- 1 pc. 3B1-2=ST/3B1-2=ST-9.2 Genium (with threaded connector) or
- 1 pc. 3B1-2/3B1-2=9.2 Genium (with pyramid connector)
- 1 pc. 2R20 AXON tube adapter or
1 pc. 2R21 AXON tube adapter with torsion
- 1 pc. 647F542 prosthesis passport
- 1 pc. card holder for prosthesis passport
- 1 pc. 646D459, 646D459=1 instructions for use (user)
- 1 pc. 757L16* power supply
- 1 pc. 4E60* inductive charger
- 1 pc. cosmetic case for battery charger and power supply
- 1 pc. 646C107 Bluetooth PIN card
- Android app "Cockpit 4X441-Andr=V*" for download from the website: <http://www.ottobock.com/cockpitapp>

Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- 4X880=* Genium Protector
- 3S26 cosmetic foam cover
- 4X350 remote control
The user interface is available in the following languages: German, English, French, Italian, Spanish, Dutch, Swedish

7 Charging the battery

The following points must be observed when charging the battery:

- With average use, the capacity of the fully charged battery is sufficient for about 5 days.
- We recommend charging the product every day when used on a daily basis.

- When used daily, the complete charging unit (power supply – battery charger) may remain plugged into the wall socket.
- The battery should be charged for at least 3 hours prior to initial use.
- Note the permissible temperature range for charging the battery (see Page 89).
- Use the 757L16* power supply and 4E60* battery charger to charge the battery.
- The distance between the battery charger and the receiver on the product must not exceed 2 mm.

7.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
 - 2) Connect the round, **three-pin** plug of the power supply to the receptacle on the battery charger (see fig. 2) so that the plug locks into place.
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
 - 3) Plug the power supply unit into the outlet (see fig. 3).
 - The green LED on the back of the power supply lights up.
 - The LED ring (status indicator) on the rear of the charger lights up green to indicate the correct connection to the power supply.
- If the green LED on the power supply and the LED ring on the battery charger do not light up, there is an error (see Page 93).

7.2 Charging the prosthesis battery

- 1) Connect the inductive charger to the receiver of the charging unit on the rear of the product. The charger is held in place by a magnet.
 - The LED ring on the rear of the charger pulsates purple (4-second cycle).
 - If the LED ring lights up in a different colour, this indicates an error (see Page 93).
- 2) The charging process starts.
 - Once the product battery is fully charged, all LEDs on the side of the battery charger light up.
- 3) After the charging process is complete, hold the prosthesis still and remove the inductive charger from the receiver.
 - A self-test is performed. The joint is operational only after corresponding feedback (see Page 96).

7.3 Display of the current charge level

7.3.1 Display of battery charge level without additional devices

INFORMATION

The charge level cannot be displayed by turning the prosthesis over during the charging process.

- 1) Turn the prosthesis by 180° (the sole of the foot must face up).
- 2) Hold still for 2 seconds and wait for beep signals.

Beep signal	Battery charge level
5x short	more than 80%
4x short	60% to 80%
3x short	40% to 60%
2x short	20% to 40%
1x short	less than 20%

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app or using the remote control (see Page 76) or if mute mode (silent mode) is activated, there are no beep signals.

7.3.2 Display of the current charge level using the Cockpit app


Once the Cockpit app has been started, the current charge level is displayed in the bottom line of the screen:

1. 38% – Charge level of battery for currently connected component

7.3.3 Display of current charge level using remote control (optional accessory)


The current charge level is displayed in the status bar of the remote control when it is connected and switched on:

1. – Charge level of battery for currently connected prosthesis
 – Prosthesis is being charged
2. 68% – Charge level of battery for currently connected prosthesis, as a percentage

7.3.4 Display of the current charge level during the charging process

During the charging process, the current battery charge level is indicated by the LEDs lit on the side of the charger.

	Quantity	Battery charge level
	0	0%-10%
	1	10%-30%
	2	30%-50%
	3	50%-70%
	4	70%-90%
5	> 90%	

8 Cockpit app

The cockpit app enables switching from basic mode into the pre-configured MyModes. In addition, information about the product (step counter, battery charge level, etc.) can be called up.

The behaviour of the product can be changed to a certain extent on a day-to-day basis using the app (e.g. while becoming accustomed to the product). At the next appointment, the prosthetist can trace the changes using the adjustment software.

INFORMATION

The Cockpit app can be downloaded free of charge from the respective online store. For more information, please visit the website: <http://www.ottobock.com/cockpitapp>. To download the Cockpit app, the QR code on the supplied Bluetooth PIN card can be read with the mobile device (requirement: QR code reader and camera).

INFORMATION

The language of the Cockpit app user interface can only be changed by the prosthetist using the adjustment software.

8.1 System Requirements

The functioning of the cockpit app is assured on devices with the Android operating system version 4.0.3 or higher.

The functionality was tested on the following devices:

- Samsung Galaxy S5, Galaxy S4 mini, Galaxy SIII mini, Galaxy Note II, Galaxy Fame
- Sony Xperia Z, Xperia Z3, Xperia J
- HTC One, One mini
- LG Optimus L9, Optimus F5, Optimus G2
- Huawei Ascend P6
- Motorola Droid Razr Maxx, Moto X, Moto G
- Google Nexus 4, Nexus 5, Nexus 6

8.2 Initial connection between cockpit app and prosthesis

The following points should be observed before the initial connection:

- Bluetooth of the component must be switched on.
If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 79).
- Bluetooth of the device must be switched on.
- The device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The device must be connected to the internet.**
- The serial number and PIN code of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".

INFORMATION

If the Bluetooth PIN card with the PIN card and serial number of the component is lost, contact your prosthetist.


8.2.1 Starting the cockpit app for the first time

- 1) Tap the symbol of the Cockpit app ().
→ The end user license agreement (EULA) is displayed.
 - 2) Accept the end user license agreement (EULA) by tapping the **Accept** button. If the end user license agreement (EULA) is not accepted, the Cockpit app cannot be used.
→ The welcome screen appears.
 - 3) Tap the **Add component** button.
→ The "**Preparation**" screen appears. The serial number of the component should be entered here.
 - 4) Follow the subsequent instructions on the screen.
 - 5) After the PIN code is entered, a connection to the component is established.
→ While the connection is being established, 3 beep signals sound and the  symbol appears.
The  symbol is displayed when the connection has been established.
- Once the connection has been established, the data are read from the component. This process may take up to a minute.
Then the main menu appears with the name of the connected component.

INFORMATION


After the initial connection to the component has been established successfully, the app will connect automatically each time it is started. No further steps are required.

8.3 Control elements for cockpit app

1. ☰ Call up the navigation menu (see Page 70)
2. Product
The component name can only be changed with the adjustment software.
3. If connections to more than one component have been saved, you can switch between the saved components by tapping the entry **change**.
4. MyModes configured with the adjustment software.
Switching the mode by tapping the corresponding symbol and confirming by tapping "OK".
5. Currently selected mode
6. Charge level of the component.
 - Component battery fully charged
 - Component battery empty
 The current charge level is also displayed in %.
7. Display and name of the currently selected mode (e.g. **1. Basic Mode**)
8. Mute mode is activated
9. Connection to component has been established
 Connection to component has been interrupted. The app is attempting to re-establish the connection automatically.

8.3.1 Cockpit app navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

Product

Name of the connected component

MyModes

Return to the main menu to switch MyModes

Functions

Call up additional functions of the component (e.g. turn off Bluetooth) (see Page 79)

Settings

Change settings of the currently selected mode (see Page 76)

Status

Query status of the connected component (see Page 80)

Manage components

Add or delete components (see Page 70)

Imprint/Info

Display information/legal notices for the cockpit app

8.4 Managing components

Connections with up to four different components can be stored in the app. However, a component can only be connected to one device or remote control at a time.

The following points need to be observed before establishing the connection:

- Bluetooth of the component must be switched on.
If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 79).
- Bluetooth of the device must be switched on.
- The device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The device must be connected to the internet.**
- The serial number and PIN code of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".

INFORMATION

If the Bluetooth PIN card with the PIN card and serial number of the component is lost, contact your prosthetist.

8.4.1 Adding component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) In the following screen, tap the "**Add component**" button.
→ The "Preparation" screen appears. The serial number of the component should be entered here. The serial number begins with the letters "SN".
- 4) Follow the subsequent instructions on the screen.
- 5) After the PIN code is entered, a connection to the component is established.
→ While the connection is being established, 3 beep signals sound and the (📶) symbol appears.
The (📶) symbol is displayed when the connection has been established.
→ Once the connection has been established, the data are read from the component. This process may take up to a minute.
Then the main menu appears with the name of the connected component.

INFORMATION

If establishing a connection to a component is not possible, perform the following steps:

- ▶ Delete the component from the Cockpit app if applicable (see the section "Deleting a component")
- ▶ Add the component again in the Cockpit app (see the section "Adding a component")


8.4.2 Deleting a component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Tap the 🗑 symbol under the component you want to delete.
→ The component is deleted.

8.4.3 Connecting component with multiple devices

The connection for a component can be stored on more than one device. However, only one device or remote control can be connected to the component at one time.

If there is an existing connection between the component and a different device, the following information appears while the connection is being established with the current device:


- ▶ Tap the **OK** button.
- The connection to the last connected device is broken off and established with the current device.

9 Use

9.1 Movement patterns in basic mode (mode 1)

9.1.1 Standing

Knee control through high hydraulic resistance and static alignment. A stance function can be enabled by the prosthetist. Please see the following section for further information on the stance function.

9.1.1.1 Stance function

INFORMATION

To use this function, it must be enabled by the prosthetist. It also has to be activated using the cockpit app or remote control (see Page 77).

The stance function (standing mode) is a functional supplement to the basic mode (mode 1). This function makes it easier, for example, to stand on an inclined surface for a longer time. The joint is fixed in the flexion direction.

The standing function must be enabled by the prosthetist. In addition, the type of joint locking (conscious/intuitive) has to be established by the prosthetist. Changing the locking type using the Cockpit app or remote control is not possible.

Intuitive locking of the joint

The intuitive stance function recognises any situation that puts strain on the prosthesis in the flexion direction but where flexion is not permitted. Examples of this include standing on uneven or sloping surfaces. The knee joint is always locked in the flexion direction when the prosthetic leg is not fully extended, under some amount of load and at rest. Upon forward or backward rollover or extension, the level of resistance is immediately reduced to stance phase resistance again.

The knee joint is not locked when the above conditions are met and a sitting position is assumed (for example while driving).

Deliberate locking of the joint


- 1) Assume the desired knee angle.
- 2) Do not remove the entire load from the prosthesis.
- 3) Do not change the knee angle for a brief period (1/8 second). This time period prevents unintentional activation of the standing function.

→ The blocked joint can now be fully loaded in the flexion direction.

Deliberate unlocking of the joint

- ▶ By deliberately extending or unloading the knee joint, it is unlocked again.

9.1.2 Walking


Initial attempts at walking with the prosthesis always require the instruction of trained, qualified personnel.

The hydraulics stabilise the knee joint in the stance phase and release the knee joint in the swing phase so that the leg can swing forward freely.

Switching to the swing phase requires that the prosthesis roll over to the front out of the stride position.

9.1.3 Running short distances ("walk-to-run" function)

For covering short distances quickly, the knee joint detects a transition from walking to running in basic mode and automatically changes the following settings:

- The swing phase angle is increased
- Preflexion of 4° at heel strike (PreFlex) is reduced to 0°

The requirements to automatically switch to the running motion are fast forward movement of the prosthetic leg and high dynamic load on the knee joint. When stopping from the running motion, the changed settings are set back to the standard values.

9.1.4 Sitting down

The resistance in the prosthetic knee joint while sitting down ensures evenly bending into the sitting position.

- 1) Place both feet side by side at the same level.
- 2) While sitting down, distribute weight evenly between both legs and use armrests, if available.
- 3) Move the buttocks in the direction of the backrest and lean the upper body forward.

INFORMATION: Resistance while sitting down can be changed with the Cockpit app or the remote control via the parameter "Resistance" (see Page 77).

9.1.5 Sitting

INFORMATION

While sitting, the knee joint also switches to energy saving mode. This energy saving mode is activated regardless of whether the sitting function is activated or not.

If the patient is in a sitting position for more than two seconds (i.e. the thigh is close to horizontal and there is no load on the leg), the knee joint switches the resistance to a minimum in the extension direction.

A sitting function can be enabled by the prosthetist. For more information about the sitting function, see the following section.

9.1.5.1 Sitting function

INFORMATION

To use this function, it must be enabled by the prosthetist. It also has to be activated using the cockpit app or remote control (see Page 77).

In the sitting position, the resistance in the flexion direction is reduced in addition to the reduction of resistance in the extension direction. This makes it possible to swing the prosthetic leg freely.

9.1.6 Standing up

Flexion resistance is increased steadily while standing up.


- 1) Place the feet at the same level.
- 2) Lean the upper body forward.
- 3) Put the hands on armrests, if available.
- 4) Stand up with support from the hands while keeping weight evenly distributed on feet.

9.1.7 Walking up stairs step-over-step

INFORMATION

To use this function, it must be enabled by the prosthetist. It also has to be activated using the cockpit app or remote control (see Page 77).

Although the knee joint is passive, which means it cannot actively initiate movements, negotiating stairs step-over-step is possible.

This function must be practised and executed consciously.

- 1) Lift the extended prosthesis off the floor.
- 2) Immediately after lifting the extended leg off the floor, extend the hip briefly and then abruptly flex it. This requires a sufficiently secure hold in the socket and a certain level of residual limb strength.
 - This whip motion flexes the knee, because the knee joint automatically recognises the movement and sets the flexion resistance to minimum.

INFORMATION: Take note of people behind you before executing the whip motion.


- 3) When sufficient knee flexion has been achieved, the knee joint increases extension resistance so that there is enough time to position the foot on the next step before the knee joint is extended again.
- 4) Set the foot onto the next step.

The support area for the foot on the step must be large enough that the heel does not extend back too far over the edge. With too little support area, the lower leg would extend too early and position the leg too far backwards. In this phase, the knee joint has already set the flexion resistance to maximum (blocked). The knee joint can no longer be flexed but only extended. This ensures that the leg does not buckle if the hip strength is not sufficient for the extending motion.
- 5) Support yourself with your hand on the contralateral side. A smooth wall will also work. This lateral support is intended to prevent the residual limb from twisting in the socket. Twisting can lead to unpleasant surface tension between the skin and the socket. Lateral support also improves balance.
- 6) Bring the knee into extended position. When the knee joint is fully extended, the initial position has been reached.
- 7) You can climb the next step or continue walking normally.

9.1.8 Overcoming obstacles

INFORMATION

To use this function, it must be enabled by the prosthetist. It also has to be activated using the cockpit app or remote control (see Page 77).


The stair function can also be used to cross obstacles:

- 1) Lift the extended prosthesis off the floor.
- 2) Briefly extend the hip.
- 3) Quickly flex the hip. This causes the knee to flex.
- 4) With the knee flexed, step over the obstacle.

With sufficient knee flexion, the extension resistance is increased to allow enough time for crossing the obstacle.

9.1.9 Walking down stairs


This function must be practised and executed consciously. Only when the sole is properly positioned can the system react correctly and permit controlled flexion.

- 1) Hold the handrail with one hand.
- 2) Position the leg with the prosthesis on the step so that the foot projects halfway over the edge of the step.
→ This is the only way to assure a secure rollover.
- 3) Roll the foot over the edge of the step.
→ This flexes the prosthesis slowly and evenly in the knee joint.
- 4) Place the foot of the other leg onto the next step.
- 5) Place the foot of the prosthesis leg on the next step after that.

INFORMATION: The flexion speed of the knee joint can be changed using the Cockpit app or remote control via the parameter "Resistance" (see Page 77).

9.1.10 Walking down a ramp

Under increased flexion resistance, permit controlled flexion of the knee joint which lowers the body's centre of gravity.

INFORMATION: The flexion resistance at which the knee bends can be changed using the Cockpit app or remote control via the parameter "Resistance" (see Page 77).

9.2 Changing prosthesis settings

Once an active connection to a prosthesis has been established, the settings of the **mode currently active** can be changed using the cockpit app or the remote control (optional accessory).

For information on changing the prosthesis settings with the remote control (optional accessory), see the remote control instructions for use.

INFORMATION

Bluetooth on the prosthesis must be switched on to change the prosthesis settings.


If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. The connection must be established during this period.

Information for changing the prosthesis settings

- Before changing settings, always check the main menu of the cockpit app or the display of the remote control (optional accessory) to make sure the correct prosthesis has been selected. Otherwise, parameters could be changed for the wrong prosthesis.
- It is not possible to change prosthesis settings nor to switch to a different mode while the prosthesis battery is being charged. Only the status of the prosthesis can be called up. Instead of the  symbol, the  symbol appears in the bottom row of the screen in the cockpit app.
- The prosthetist's setting is located in the middle of the scale. After making adjustments, this setting can be restored by tapping the **"Standard"** button (Cockpit app) or by moving the slider control into the middle (remote control).

- Prosthesis settings should be optimised using the adjustment software. The cockpit app or remote control (optional accessory) is not intended for use by the prosthetist in order to set up the prosthesis. The behaviour of the prosthesis can be changed to a certain extent on a day-to-day basis using the remote control or app (e.g. while becoming accustomed to the prosthesis). The prosthetist can use the adjustment software to trace these changes at the next appointment.
- If the settings of a MyMode are to be modified, one must first switch to this MyMode.

9.2.1 Changing the prosthesis setting using the cockpit app

- 1) Once the component is connected and in the desired mode, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) Tap the menu entry "**Settings**".
→ A list appears with the parameters for the currently selected mode.
- 3) Change the setting of the desired parameter by tapping the "<", ">" symbols.

INFORMATION: The prosthetist's setting is marked and, after the setting has been changed, can be restored by tapping the "Standard" button.

9.2.2 Overview of adjustment parameters in basic mode

INFORMATION
When mute mode (silent mode) is activated, no beep and vibration signals are generated.

The parameters in basic mode describe the dynamic behaviour of the prosthesis in a normal gait cycle. These parameters act as basic settings for automatically adjusting the damping behaviour to the current motion situation (e.g. ramps, slow walking speed, etc.).

The stance function, sitting function and/or stairs and obstacles function can also be activated/deactivated. Further information on the stance function (see Page 72), sitting function (see Page 73), stairs and obstacles function (see Page 74).

The following parameters can be modified:

Parameter	Adjustment software range	App/remote control adjustment range	Meaning
Resistance	120 – 180	+/- 10	Resistance against flexion motion e.g. when walking down stairs or when sitting down

Parameter	Adjustment software range	App/remote control adjustment range	Meaning
Angle	55° – 70°	+/- 3°	Maximum flexion angle in the swing phase
Stance function	Deactivated Activated	0 - deactivated 1 - activated	Activation/deactivation of the stance function. For switching with the Cockpit app or remote control, this function has to be activated by the prosthetist. Further information (see Page 72).
Sitting function	Deactivated Activated	0 - deactivated 1 - activated	Activation/deactivation of the sitting function. For switching with the Cockpit app or remote control, this function has to be activated by the prosthetist. Further information (see Page 74).
Stair Function	Deactivated Activated	0 - deactivated 1 - activated	Activating/deactivating the stair and obstacle function. For switching with the Cockpit app or remote control, this function has to be activated by the prosthetist. Further information (see Page 74).
Pitch	1000 Hz – 4000 Hz	1000 Hz – 4000 Hz	Pitch of beep signal for confirmation tones
Volume	0 – 4	0 – 4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals in case of errors are still generated.

In addition to the parameters, the remote control (optional accessory) can also be used to activate deep sleep mode. Further information on deep sleep mode (see Page 81).

9.2.3 Overview of adjustment parameters in MyModes

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

The parameters in the MyModes describe the static behaviour of the prosthesis for a specific motion pattern such as cross-country skiing. Damping behaviour is not automatically controlled and adjusted in MyModes.

The following parameters can be modified in the MyModes (except MyMode "Brake knee"):

Parameter	Adjustment software range	App/remote control adjustment range	Meaning
Gain	0 – 100	+/- 10	Value for speed at which flexion resistance increases as the knee angle increases
Pitch	1000 Hz – 4000 Hz	1000 Hz – 4000 Hz	Pitch of beep signal for confirmation tones
Volume	0 – 4	0 – 4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals in case of errors are still generated.

The following parameters can be modified in MyMode "Brake knee":

Parameter	Adjustment software range	App/remote control adjustment range	Meaning
Stance function	0 – 180	+/- 20	Basic resistance for the knee joint without load
Pitch	1000 Hz – 4000 Hz	1000 Hz – 4000 Hz	Pitch of beep signal for confirmation tones
Volume	0 – 4	0 – 4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals in case of errors are still generated.

9.3 Turning Bluetooth on the prosthesis on/off

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the Bluetooth battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 79).

INFORMATION

Basic mode (mode 1) has to be active to turn off Bluetooth. If a MyMode is activated, one has to switch to basic mode to turn off Bluetooth.

9.3.1 Switching Bluetooth off/on using the cockpit app

Switching off Bluetooth

- 1) When the component is connected, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Deactivate Bluetooth**".
- 4) Follow the on-screen instructions.

Switching on Bluetooth

- 1) Turn the component over or connect/disconnect the battery charger.
→ Bluetooth is switched on for approx. 2 minutes. During this time, the app must be started in order to establish a connection to the component.
- 2) Follow the on-screen instructions.
→ If Bluetooth is switched on, the (↔) symbol appears on the screen.

9.4 Querying the prosthesis status

9.4.1 Query status through cockpit app

- 1) When the component is connected, tap the ☰ symbol in the main menu.
- 2) In the navigation menu, tap the entry "**Status**".

9.4.2 Status display in the cockpit app

Menu option	Description	Possible actions
Trip	Daily step counter (steps with the prosthesis side)	Reset the counter by tapping the " Reset " button.
Step	Total step counter (steps with the prosthesis side)	For informational purposes only
Service	Display of the next maintenance date	For informational purposes only
Batt.	Current prosthesis charge level, as a percentage	For informational purposes only
Stb/Act: 58/29	Estimated remaining operating time of prosthesis in hours. Rest mode (Stb) e.g. 58 hours, active use (Act) e.g. 29 hours	For informational purposes only

9.4.3 Status information on the display of the remote control (optional accessory)

Menu option	Description	Possible actions
Trip	Daily step counter (Steps with the prosthesis side)	Reset the counter by using the ■ key to confirm the menu option
Step	Total step counter (Steps with the prosthesis side)	For informational purposes only
Batt.	Current prosthesis charge level, as a percentage	For informational purposes only

Menu option	Description	Possible actions
Stb/Act: 58/29	Estimated remaining operating time of prosthesis in hours. Rest mode (Stb) e.g. 58 hours, active use (Act) e.g. 29 hours	For informational purposes only
Service	Display of the next maintenance date	For informational purposes only
Mute: On	Turn mute mode (silent mode) of the prosthesis on or off. The current status is displayed.	Mute mode can be turned on or off by confirming the menu point with the ■ key (see Page 81).
Bluetooth: On	Turning the Bluetooth function of the prosthesis on/off. The current status is displayed.	The Bluetooth function of the prosthesis can be turned on or off (see Page 79) by confirming the menu item with the ■ key.

9.5 Mute mode (silent mode)

Activating mute mode (silent mode) turns off the audible feedback signals and the vibration signals. However, warnings in case of component errors are still generated (see Page 93).

Mute mode can be activated/deactivated by using the Cockpit app or the remote control (optional accessory).

INFORMATION

Connecting the battery charger automatically deactivates mute mode again.

9.5.1 Turning mute mode on/off using the Cockpit app

- 1) When the component is connected, tap the **☰** symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Mute mode**".
- 4) Follow the on-screen instructions.

9.5.2 Turning mute mode on/off using the remote control

- 1) When a connection is established, use the **▲**, **▼** keys to select the "**Status**" menu option in the main menu and confirm this with the **■** key.
- 2) Use the **▲**, **▼** keys to select the **Mute** menu option. The current status is displayed.
- 3) By pressing the **■** key, mute mode is either turned on (On) or turned off (Off).
- 4) To leave the menu without saving changes, use the **▼** key to select the **⬅** symbol and confirm this with the **■** key.

9.6 Deep sleep mode

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app or using the remote control, there are no beep signals (see Page 76).

The cockpit app or remote control (optional accessory) can be used to place the knee joint in a deep sleep mode, in which power consumption is minimised. The knee joint offers no functionality in this mode. The safety mode damping values are activated.

It can be awakened from deep sleep mode with the cockpit app, remote control or by connecting the charger. Waking from deep sleep mode using the cockpit app or remote control can take up to 30 seconds.

After ending deep sleep mode, the knee joint is in basic mode again.

9.6.1 Turning deep sleep mode on/off using the Cockpit app

Activating deep sleep mode

- 1) When the component is connected, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
 - 2) In the navigation menu, tap the entry "**Functions**".
 - 3) Tap the entry "**Activate Deep Sleep**".
 - 4) Follow the on-screen instructions.
- The activation of deep sleep mode is indicated by a short beep signal and a short vibration signal, assuming that mute mode (silent mode) is not active.

Deactivating deep sleep mode

- 1) When deep sleep mode is active for the currently connected prosthesis, the **Exit deep sleep mode** button automatically appears when the Cockpit app is started.
- 2) Tapping this button establishes a connection to the prosthesis and deactivates deep sleep mode.

INFORMATION: Establishing a connection in deep sleep mode can take up to 30 seconds.

If a prosthesis is in deep sleep mode but not connected to the Cockpit app, a connection to the prosthesis has to be established (see Page 70).

9.6.2 Turning deep sleep mode on/off using the remote control

Activating deep sleep mode

- 1) Use the ▼, ▲ keys to select the menu option "**Device settings**" in the main menu and confirm this with the ■ key.
- 2) Use the ▼, ▲ keys to select the menu option **Deep sleep** and confirm this with the ■ key.
→ The activation of deep sleep mode is indicated by a short beep signal and a short vibration signal. Also the remote control is turned off.

Deactivating deep sleep mode

- 1) Use the ■ key to activate the remote control.
- 2) With the ▼ key, establish a connection to the prosthesis. Establishing the connection to the prosthesis deactivates deep sleep mode.

INFORMATION: Establishing a connection in deep sleep mode can take up to 30 seconds.

9.7 OPG function (Optimised Physiological Gait)

INFORMATION

To use this function, it must be enabled by the prosthetist.

The OPG function minimises the prosthesis wearer's prosthetic deviations from a harmonious gait pattern and promotes more biomechanically correct walking. Selection of this function enables the following features:

PreFlex

PreFlex ensures the knee is in 4° of flexion at the end of swing phase in preparation for loading response. This makes initiating stance phase flexion easier and forward movement is less restricted.

Adaptive yielding control

The knee joint has auto-adaptive stance and swing extension resistance. The stance flexion resistance experienced by the user is dependent on the slope or incline when walking downhill. When walking on a ramp, adaptive yielding control manages flexion depending on the angle of the ramp. The knee joint flexes slowly if the ramp is flat, and flexes quickly if the ramp is steep.

Dynamic stability control (DSC)

DSC ensures the knee will not release stance resistance during biomechanically unstable static and dynamic conditions. Constantly checking multiple parameters, DSC ensures the optimally timed decision for the knee to safely switch from stance to swing. Because DSC is always monitoring knee function, multi-directional movement and walking backward are also possible without risk of stance resistance releasing.

Adaptive swing phase control


Instantaneous adaptation to varied walking cadences and to changes of the pendular mass (e.g. varying footwear) ensures the knee always achieves the swing flexion target angle within (+/-) one degree. The swing phase extension and flexion resistance experienced by the user are auto-adaptive.

The flexed and partially loaded knee will also disable the stance phase on slopes and ramps to allow for greater knee flexion and more ground clearance in the swing phase.

10 MyModes

With the help of adjustment software, the prosthetist can activate and configure up to 5 MyModes in addition to the basic mode. These can be selected by using the Cockpit app or the remote control (optional accessory). Only the first 3 MyModes can be selected using movement patterns. Switching by using movement patterns has to be activated in the adjustment software by the prosthetist.

For information on switching MyModes with the remote control (optional accessory), see the remote control instructions for use.

The MyModes are intended for specific motion patterns or postures (e.g. inline skating, ...). Settings can be adjusted using the Cockpit app or remote control (see Page 78).

10.1 Switching MyModes with the cockpit app

INFORMATION


Bluetooth on the prosthesis must be turned on in order to use the Cockpit app.

If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the Bluetooth battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 79).

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app or using the remote control (see Page 76) or if mute mode (silent mode) is activated, there are no beep signals.

Once a connection to a prosthesis has been established, the cockpit app can be used to switch between the MyModes.

- 1) Tap the symbol of the MyMode (1) you want in the main menu of the app.
→ A security question for changing the MyMode appears.
- 2) If you want to change the mode, tap the "OK" button.
→ A beep signal sounds to confirm the switch.
- 3) After switching, a symbol (2) is displayed to identify the active mode.
→ The current mode is also indicated by the name on the lower edge of the screen (3).

10.2 Switching MyModes using motion patterns

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app or using the remote control, there are no beep signals (see Page 76).

Information on switching


- Switching and the number of motion patterns must be activated by the prosthetist in the adjustment software.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Switching using movement patterns has to be enabled by the prosthetist.
- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.
- Weight must be placed on the forefoot while bouncing.
- The load may not be taken off the foot fully when relieving the load.

Switching process

- 1) Position the prosthetic leg back slightly.
- 2) While maintaining constant contact with the floor, bounce on the forefoot with the leg extended a number of times in one second depending on the desired MyMode (MyMode 1 = 3 times, MyMode 2 = 4 times).
- 3) Keep the prosthetic leg still in this position with no load.
 - A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see Page 81).

- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.
 - A confirmation signal will sound to indicate that the prosthesis has successfully switched to the corresponding MyMode (2 times = MyMode 1, 3 times = MyMode 2).

INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process to correctly switch to the required mode. For more information about mute mode, see the section "Mute mode" (see Page 81).

10.3 Switching from a MyMode back to basic mode

Information on switching


- Regardless of the configuration of additional MyModes in the adjustment software, it is always possible to switch back to basic mode (mode 1) with a motion pattern.
- It is always possible to switch back to basic mode (mode 1) by connecting/disconnecting the battery charger.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Position the prosthetic leg back slightly and bounce on the forefoot with the leg extended while maintaining constant contact with the floor.
- Weight must be placed on the forefoot while bouncing.
- The load may not be taken off the foot fully when relieving the load.

Switching process

- 1) Position the prosthetic leg back slightly.

- 2) While maintaining constant contact with the floor and with the leg extended, bounce on the forefoot at 3 times or more.
- 3) Keep the prosthetic leg still in this position with no load.
 - A beep and vibration signal will sound to confirm that the movement pattern has been recognised.

INFORMATION: If this beep and vibration signal does not sound, the requirements were not observed when bouncing the foot or mute mode (silent mode) is activated. For more information about mute mode, see the section "Mute mode" (see Page 81).
- 4) After the beep and vibration signal sounds, keep the prosthetic leg extended and still for 1 second.
 - A confirmation signal will sound to indicate that the prosthesis has successfully switched over to basic mode.

INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly kept still or mute mode (silent mode) is activated. Repeat the process to correctly switch to the required mode. For more information about mute mode, see the section "Mute mode" (see Page 81).

11 Additional operating states (modes)

11.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level is 5% (see Page 93). During this time, damping settings are set to their safety mode values. The prosthesis is then switched off. You can switch back to basic mode (mode 1) from empty battery mode by charging the product.

11.2 Mode for charging the prosthesis

The product is non-functional during charging.

The product is set to the flexion resistance of safety mode. This may be low or high depending on the setting configured by the prosthetist.

11.3 Safety mode

The product automatically switches to safety mode if a critical system fault occurs (e.g. failure of a sensor signal). Safety mode remains in effect until the error has been rectified.

The switch to safety mode is indicated by beeps and vibration signals immediately prior to switching (see Page 93).

Safety mode can be disabled by connecting and disconnecting the battery charger. If the product switches into safety mode again, this means a permanent error exists. The product must be inspected by an authorised Ottobock Service Centre.

Depending on the type and severity of the error, different remaining functionality is offered in safety mode which enables limited walking by the user.

The following remaining functionality is available:

- Continuous stance phase flexion resistance is set, with the ability to initiate the swing phase.
- The swing phase control and stance phase extension resistance may be available or not depending on the type of error.
- Safety mode flexion resistance is configured. This may be low or high depending on the setting configured by the prosthetist.

The following functions are deactivated in safety mode:

- OPG function
- Stairs and obstacles function
- Stance function
- Sitting function

11.4 Overheating mode

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

When the hydraulic unit overheats due to uninterrupted, increased activity (e.g. extended walking downhill), damping is increased along with the rising temperature in order to counteract the overheating. When the hydraulic unit cools down, the product switches back to the damping settings that existed before the overheating mode.

Overheating mode is not activated in the MyModes.

Overheating mode is indicated by a long vibration every 5 seconds.

The following functions are deactivated in overheating mode:

- OPG function
- Stairs and obstacles function
- Stance function
- Sitting function
- Display of the battery charge level without additional equipment
- Switching to a MyMode

12 Storage and bleeding

Air may accumulate in the hydraulic unit if the product is stored for longer periods and not in an upright position. This is noticeable through sounds and irregular damping behaviour.

The automatic bleeding mechanism ensures that all functions of the product are again intact after approximately 10 - 20 steps.

Storage

- When storing the knee joint, press the top of the knee to the flexion stop.
- Avoid extended disuse of the product (use the product regularly).

13 Cleaning

- 1) Clean the product with a damp cloth (fresh water) when needed.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

14 Maintenance

INFORMATION

This component was tested for three million load cycles in accordance with ISO 10328.

Depending on the activity level, this corresponds to a service life of three to five years.

The duration of use can be individually extended depending on the intensity of use by performing regular service inspections.

Regular service inspections are recommended in the interest of the user's own safety, and in order to maintain operating reliability and protect the warranty. These service inspections include an inspection of the sensors and replacement of worn parts.

To have a service inspection carried out, please submit the product as well as the battery charger and power supply unit to the prosthetist.

INFORMATION

If a remote control was included in the scope of delivery as an optional product accessory, it must also be sent along with the product for service inspections.

15 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

15.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

15.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

15.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

This product meets the requirements of the European Directive 1999/5/EC for radio equipment and telecommunications terminal equipment. The conformity assessment was drawn up by the manufacturer in accordance with Annex III of the directive.

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

15.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

—Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA
Phone + 1-801-956-2400
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.
5470 Harvester Road
L7L 5N5 Burlington, Ontario
Canada
Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

16 Technical data

Environmental conditions	
Transportation in original packaging	-25°C/-13°F to +70°C/+158°F
Transportation without packaging	-25°C/-13°F to +70°C/+158°F Max. 93% relative humidity, non-condensing
Storage (≤3 months)	-20°C/-4°F to +40°C/+104°F Max. 93% relative humidity, non-condensing
Long-term storage (>3 months)	-20°C/-4°F to +20°C/+68°F Max. 93% relative humidity, non-condensing

Environmental conditions	
Operation	-10°C/+14°F to +60°C/+140°F Max. 93% relative humidity, non-condensing
Charging the battery	+10°C/+50°F to +45°C/+113°F

Product	
Reference number	3B1-2*/3B1-2=ST*
Mobility grade according to MOBIS	2 to 4
Maximum body weight with installed 2R20 tube adapter	150 kg
Maximum body weight with installed 2R21 tube adapter (torsion)	125 kg
Protection rating	IP54
Range of Bluetooth connection to remote control/mobile device	max. 10 m/32 ft
Weight of the prosthesis without protector and tube adapter	approx. 1400 g / 3.09 lbs

Data transfer	
Wireless technology	Bluetooth 2.0
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK
Data rate (over the air)	704 kbps
Maximum output power (EIRP)	-1.78 dBm (0.66 mW)

Tube adapter		
Reference number	2R20	2R21 (with torsion unit)
Weight	190-300 g/0.42-0.66 lbs	435-545 g/0.96-1.20 lbs
Material	Aluminium	
Max. body weight	150 kg	125 kg
Protection rating	IP54	

Prosthesis battery	
Battery type	Li-Ion
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500
Charge level after 1 hour charging time	30 %
Charge level after 2 hours charging time	50 %
Charge level after 4 hours charging time	80 %
Charge level after 8 hours charging time	Fully charged
Product behaviour during the charging process	The product is non-functional
Operating time of the prosthesis with new, fully charged battery at room temperature	Approx. 5 days with average use

Power supply	
Reference number	757L16*
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10 % to 93 % relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 90 % relative humidity, non-condensing
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ==

Battery charger	
Reference number	4E60*
Storage and transport in original packaging	-25 °C to 70 °C/-13 °F to 158 °F
Storage and transport without packaging	-25 °C to 70 °C/-13 °F to 158 °F Max. 93% relative humidity, non-condensing
Operation	5 °C to 40 °C/41 °F to 104 °F Max. 93% relative humidity, non-condensing
Protection rating	IP40
Input voltage	12 V ==

Cockpit app	
Article number	Cockpit 4X441-Andr=V*
Supported operating system	Android 4.0.3 or higher
Website for download	http://www.ottobock.com/cockpitapp

17 Appendices

17.1 Symbols Used

Legal manufacturer

Type BF applied part

Compliance with the requirements according to "FCC Part 15" (USA)

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Non-ionising radiation

IP40

Protection against penetration of solid foreign objects with a diameter greater than 1 mm, no protection against water

IP54

Protected against dust, protected against splashing water

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

SN YYY YW NNN

Serial number

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Lot number

Caution, hot surface

Protect from moisture

17.2 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

17.2.1 Signals for operating states

Battery charger connected/disconnected

Beep signal	Vibration signal	Event
	3x long	Charging mode started (3 sec. after connecting the battery charger)
1x short	1x before beep signal	Self-test completed successfully, product is operational

Mode switching

INFORMATION

When mute mode (silent mode) is activated, no beep and vibration signals are generated.

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app or using the remote control, there are no beep signals (see Page 76).

Beep signal	Vibration signal	Additional action performed	Event
1x short	1x short	Mode switching using the Cockpit app or remote control	Mode switching performed using the Cockpit app or remote control.
1x long	1x long	Bouncing on the forefoot followed by unloading the prosthetic leg	Bouncing pattern recognised.
1x short	1x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to basic mode (mode 1) carried out.
2x short	2x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 1 (mode 2) carried out.
3x short	3x short	Weight taken off prosthetic leg and leg kept still for 1 second	Switching to MyMode 2 (mode 3) carried out.



17.2.2 Warnings/error signals


Error during use

Beep signal	Vibration signal	Event	Required action
	1x long at intervals of approx. 5 seconds (if mute mode (silent mode) is activated, this signal is not generated)	Overheated hydraulic unit	Reduce activity.
	3x long	Battery charge level under 25%	Charge battery soon. Remaining operating time approx. 24 hours
	5x long	Battery charge level under 10%	Charge battery soon. Remaining operating time approx. 6 hours
5x long	5 x long repeated every 60 seconds	Moderately severe error/indication of safety mode activation e.g. a sensor is not operational Switch to safety mode (see Page 86).	Restricted walking possible. The change in flexion resistance has to be noted. The product must be inspected immediately by an O&P professional.

Beep signal	Vibration signal	Event	Required action
10x long	10x long	Battery charge level 5% After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.
30x long	1x long, 1x short repeated every 3 seconds	Severe error e.g. failure of valve drives Possibly no switching into safety mode (see Page 86).	Restricted walking possible. Please note the possible change in flexion/extension resistance. Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected immediately by an O&P professional.

Error while charging the product

LED on power supply	Status LED on battery charger	Error	Resolution
○		Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.
		Non-functional socket	Check socket with another electric appliance.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●		No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.

	Status LED	Charging status indicator (5 LEDs)	Error	Resolution
	The LED ring is lit in weak purple	No LED is lit	Distance between battery charger and receiver of the charging unit on the prosthesis too great. If the distance is more than 2 mm, the prosthesis cannot be charged.	Reduce distance between battery charger and receiver of the charging unit.
	The LED ring is lit up yellow	LED 2 and 4 are lit up	Charger excess temperature	Check whether the specified ambient conditions for charging the battery are met (see Page 89).
		LED 1, 3 and 5 are lit up	Excessively high or low temperature of the prosthesis	
		LED 3 is lit up	The prosthesis is not being charged Distance between battery charger and receiver of the charging unit too great.	Connection may be improved by reducing the distance between the charger and receiver of the charging unit.
	The LED ring is lit up green		The battery charger is functional but not connected to the receiver, or the distance between the battery charger and receiver of the charging unit is too great.	Connect the battery charger or reduce the distance between the charger and the receiver of the charging unit on the prosthesis.
	The LED ring is flashing red		The prosthesis is not being charged Defective battery charger.	Reset the error by disconnecting and connecting the power supply. If the error persists, the battery charger and power supply must be inspected by an authorised Ottobock Service Center.

17.2.3 Error messages while establishing a connection with the cockpit app

Error message	Cause	Correction
Device not found	The connection could not be established as no component with the entered serial number could be found.	Compare the serial number entered with that of the component, then try to establish the connection again.

Error message	Cause	Correction
Component was connected to another device. Establish connection?	The component was connected to another device/remote control.	To disconnect the original connection, tap the "OK" button. If the original connection is not to be disconnected, tap the "Cancel" button.
	A current connection to the prosthesis was interrupted.	Check the following points: <ul style="list-style-type: none"> • Distance from the prosthesis to the device • Battery charge level of the prosthesis • Bluetooth on the prosthesis switched on? (see Page 79) • If multiple prostheses were stored in the remote control, was the correct prosthesis selected?

17.2.4 Status signals

Battery charger connected

LED on power supply	Status LED on battery charger	Event
		Power supply and battery charger operational. Battery charger not yet connected to receiver.
		Battery charger connected to receiver, good connection. This display turns off automatically after approximately one minute to avoid bothersome light at night. The charging process is not affected.

Battery charger disconnected

Beep signal	Vibration signal	Event	Resolution
1x short	1x short	Self-test completed successfully. Product is operational.	

Beep signal	Vibration signal	Event	Resolution
3x short	3x short	Maintenance note: e.g. service interval has been exceeded, temporary disruption of a sensor signal	<ul style="list-style-type: none"> • Use the Cockpit app/remote control to check the next maintenance date for the prosthesis (see Page 80). Should the date lie within the next month, arrange a maintenance appointment with the prosthetist. For this appointment, the prosthetist has to be provided with the battery charger, power supply and remote control (if included) in addition to the prosthesis. • Conduct the self-test again by connecting/disconnecting the battery charger. • If the beep signal sounds again and the maintenance date has not yet been reached or exceeded, visit the prosthetist soon. If necessary, the prosthetist will forward the prosthesis to an authorised Ottobock Service Centre. • The product can be used without restrictions. However, vibration signals may not be generated.
5x long	5x long (every minute)	Self-test error	<ul style="list-style-type: none"> • Conduct the self-test again by connecting/disconnecting the battery charger. • If the beep/vibration signal occurs again, visit your prosthetist promptly. If necessary, the prosthetist will forward the prosthesis to an authorised Ottobock Service Centre.

Battery charge level

During the charging process, the current battery charge level is indicated by the LEDs lit on the side of the charger.

LEDs	0	1	2	3	4	5
Battery charge level	0%-10%	11%-30%	31%-50%	51%-70%	71%-90%	>90%

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Otto Bock Healthcare Products GmbH