






6A43

DE Gebrauchsanweisung	5
EN Instructions for Use	13
FR Instructions d'utilisation	21
IT Istruzioni per l'uso	30
ES Instrucciones de uso	38
PT Manual de utilização	47
NL Gebruiksaanwijzing	55
SV Bruksanvisning	63
DA Brugsanvisning	71
NO Bruksanvisning	79
FI Käyttöohje	86
PL Instrukcja użytkowania	94
HU Használati utasítás	102
CS Návod k použití	110
RO Instrucțiuni de utilizare	118
HR Upute za uporabu	126
SL Navodila za uporabo	134
SK Návod na používanie	142
BG Инструкция за употреба	150
TR Kullanma talimatı	158
EL Οδηγίες χρήσης	166
RU Руководство по применению	175
JA 取扱説明書	184
ZH 使用说明书	191
KO 사용 설명서	197

12

13

1 Produktbeschreibung

INFORMATION

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INFORMATION

Last update: 2014-02-24

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

1.1 Construction and Function

The 6A43 Lamination Adapter is laminated in the prosthetic socket. It serves to connect with the distal prosthetic components.

1.2 Combination Possibilities

6A43 Lamination Adapter	
Designation	Reference number
Shuttle Lock	6A40

For additional combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Application

2.1 Intended Use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of Application

Approved for a body weight of up to **125 kg**.

2.3 Environmental Conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of Warning Symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.

3.2 General Safety Instructions

 CAUTION
Using the product without following the instructions for use
Deteriorating health condition and damage to the product because of failure to follow the safety information
<ul style="list-style-type: none">▶ Observe the safety information in these instructions for use.▶ Make the patient aware of all safety instructions marked with "Inform the patient".

 CAUTION
Excessive strain on the product
Fall due to breakage of load-bearing components
<ul style="list-style-type: none">▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").▶ Inform the patient.

⚠ CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

⚠ CAUTION

Exceeding the service life and reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

⚠ CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).

- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).
- ▶ **Inform the patient.**

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of Delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering according to their respective symbols as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

6A43 MagnoFlex Lock Lamination Adapter				
Fig.	Item	Quantity	Designation	Reference number
–	–	1 (■)	Instructions for Use	647G1007
1	2	1 (■)	Lamination Adapter	6A43
1	5	1 (▲)	Screw	503S3
1	6	1 (■)	Vacuum Forming Dummy	–
1	7	4 (▲)	Snap Bushing	–
1	8	4 (▲)	Cap Screw	501T28=M6X8

5 Preparation for Use

⚠ CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- ▶ Observe the alignment and assembly instructions.

ℹ INFORMATION

Not all of the materials listed may be available in your country.

- ▶ If the materials are not available, please contact the local branch of the manufacturer to obtain information on alternative materials.

5.1 Fabricating the Prosthetic Socket

5.1.1 Preparing for Socket Fabrication

> Required tools and materials:

6A43 Lamination Adapter, 501T28=M6X8 Cap Screws, Vacuum forming dummy, 503S3 Screw, 633W8 Wax, 636K8=* Plastaband

- 1) Isolate the cap screws (thread and head from the outside) with wax.
- 2) Screw the cap screws into the lamination adapter (see Fig. 2).
- 3) Fill the heads of the cap screws with Plastaband (see Fig. 3).
- 4) Position the lamination adapter on the plaster positive along the longitudinal axis of the residual limb.
- 5) Place the vacuum forming dummy into the round opening of the lamination adapter and secure the two parts to the plaster positive with the flat head screw (see Fig. 4).
- 6) Isolate the head of the flat head screw with Plastaband.
- 7) **If there is an undercut between the lamination adapter and plaster positive:** Fill the undercut with wet plaster.

5.1.2 Fabricating a Check Socket

> Required tools and materials:

755T4=360 Outer Ring, 755X104=360 Vacuum Tube (with vacuum sealing ring), 755E9 Vacuum Pump, 756D2, 756B12=* or 756B20=* Oscillating Saw, 616T52=* ThermoLyn rigid or 616T83=* ThermoLyn clear, string

- 1) Prepare for fabrication of the prosthetic socket (Preparing for Socket Fabrication).
- 2) Position two pieces of string crosswise on the lamination adapter (see Fig. 5). The pieces of string form air channels that make it easier to mould the vacuum forming material to the contours.
- 3) Complete the vacuum forming process with ThermoLyn.
- 4) Finish the prosthetic socket (see Page 20).

5.1.3 Laminating a Transtibial Socket

> Required tools and materials:

99B81=70X19X5 and 99B81=100X19X5 PVA Bag, 623T3=8 or 623T3=10 Perlon Stockinette, 81A1=8 or 81A1=10 Tube Sock, 616B1=25x* Carbon Fibre Cloth Strap, 616G2 Carbon UD Hose, 616G15 Carbon Fibre Woven Hose, 617H119 Orthocryl Lamination Resin 80:20 PRO, string

- 1) Prepare for fabrication of the prosthetic socket (Preparing for Socket Fabrication).
- 2) Soak the shorter PVA bag and pull it over the plaster positive.

- 3) Pull a tube sock over the plaster positive.
- 4) Place a layer of carbon fibre cloth strap circularly around the MPT (middle patellar tendon) point (see Fig. 6).
- 5) Position a layer of carbon fibre cloth strap medially and laterally from the shuttle lock up to the condyles (see Fig. 7).
- 6) Place a layer of carbon fibre cloth strap anteriorly and posteriorly from the shuttle lock to the circular carbon fibre cloth strap (see Fig. 8).
- 7) Pull a tube sock over the plaster positive.
- 8) Loop a piece of string around the lamination adapter in the groove and tie tightly (see Fig. 9).
- 9) Cut a piece of carbon UD hose (1.5 times the length of the plaster positive).
- 10) Pull the carbon UD hose over the plaster positive up to the socket edge. Tie the upper section of the carbon UD hose and fold the remaining length over the plaster positive.
- 11) Loop a piece of string around the lamination adapter in the groove and tie tightly (see Fig. 10).
- 12) Pull a tube sock over the plaster positive.
- 13) Cut a piece of carbon fibre woven hose (1.3 times the length of the plaster positive).
- 14) Pull the carbon fibre woven hose over the plaster positive to the socket edge. Tie the upper section of the carbon fibre woven hose and fold the remaining length over the plaster positive (see Fig. 11).
- 15) Cut a piece of Perlon stockinette (twice the length of the plaster positive).
- 16) Pull the Perlon stockinette over the plaster positive up to the socket edge. Tie the upper section of the Perlon stockinette and fold the remaining length over the plaster positive.
- 17) Soak the longer PVA bag and pull it over the plaster positive.
- 18) Perform the casting process with Orthocryl.
- 19) Finish the prosthetic socket (see Page 20).

5.1.4 Laminating a Transfemoral Socket

> Required tools and materials:

99B81=100X26X5 and 99B81=130X26X5 PVA Bag, 623T3=12 or 623T3=15 Perlon Stockinette, 81A1=12 or 81A1=15 Tube Sock, 616B1=50x* Carbon Fibre Cloth Strap, 616G2 Carbon UD Hose, 616G15 Carbon Fibre Woven Hose, 617H119 Orthocryl Lamination Resin 80:20 PRO, string

- 1) Prepare for fabrication of the prosthetic socket (Preparing for Socket Fabrication).
- 2) Soak the shorter PVA bag and pull it over the plaster positive.
- 3) Pull a layer of Perlon stockinette over the plaster positive.
- 4) Place a layer of carbon fibre cloth strap circularly 3 cm below the perineum.
- 5) Place a layer of carbon fibre cloth strap medially and laterally as well as anteriorly and posteriorly from the shuttle lock up to the circular carbon fibre cloth strap.
- 6) Place a layer of carbon fibre cloth strap circularly 3 cm below the perineum.
- 7) Cut two pieces of carbon fibre cloth strap (length approx. 20 cm).
- 8) Form the carbon fibre cloth straps into Vs by folding them in the middle and position them in the area of the tuberosity support.
- 9) Pull a tube sock over the plaster positive.
- 10) Loop a piece of string around the lamination adapter in the groove and tie tightly.
- 11) Cut a piece of carbon UD hose (1.5 times the length of the plaster positive).
- 12) Pull the carbon UD hose over the plaster positive up to the socket edge. Tie the upper section of the carbon UD hose and fold the remaining length over the plaster positive.
- 13) Loop a piece of string around the lamination adapter in the groove and tie tightly.
- 14) Pull a tube sock over the plaster positive.
- 15) Cut a piece of carbon fibre woven hose (1.3 times the length of the plaster positive).
- 16) Pull the carbon fibre woven hose over the plaster positive to the socket edge. Tie the upper section of the carbon fibre woven tube and fold the remaining length over the plaster positive.
- 17) Cut a piece of Perlon stockinette (twice the length of the plaster positive).

- 18) Pull the Perlon stockinette over the plaster positive up to the socket edge. Tie the upper section of the Perlon stockinette and fold the remaining length over the plaster positive.
- 19) Soak the longer PVA bag and pull it over the plaster positive.
- 20) Perform the casting process with Orthocryl.
- 21) Finish the prosthetic socket (see Page 20).

5.1.5 Finishing the Prosthetic Socket

> Required materials:

Snap Bushings

- 1) Mark and trim the contour of the prosthetic socket.
- 2) Sand the distal end of the socket down to the vacuum forming dummy and the heads of the cap screws (see Fig. 12).
- 3) On a level surface, check whether the sanded surface is flat. Rework as needed.
- 4) Remove the vacuum forming dummy and the cap screws from the prosthetic socket.
- 5) Remove the prosthetic socket from the plaster positive.
- 6) Sand the contour of the prosthetic socket.
- 7) Insert the snap bushings into the openings from which the cap screws were removed (see Fig. 13).

5.2 Installing the Shuttle Lock and Socket Adapter

- ▶ Install the shuttle lock and the socket adapter according to the 647G931 Instructions for Use.

6 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

7 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this

document, particularly due to improper use or unauthorised modification of the product.

7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

8 Technical data

Reference number	6A43
Weight [g]	80
System height [mm]	7
Material	High-performance polyamide
Max. body weight [kg]	125

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Otto Bock HealthCare GmbH