
4R220 Dynamic Vacuum System

DE Gebrauchsanweisung	8	CS Návod k použití.....	102
EN Instructions for use	15	RO Instrucțiuni de utilizare.....	109
FR Instructions d'utilisation.....	22	HR Upute za uporabu	116
IT Istruzioni per l'uso	30	SL Navodila za uporabo	123
ES Instrucciones de uso	37	SK Návod na používanie	130
PT Manual de utilização.....	45	BG Инструкция за употреба	137
NL Gebruiksaanwijzing.....	52	TR Kullanma talimatı	145
SV Bruksanvisning	60	EL Οδηγίες χρήσης	152
DA Brugsanvisning	67	RU Руководство по применению	160
NO Bruksanvisning	73	JA 取扱説明書.....	169
FI Käyttöohje.....	80	ZH 使用说明书.....	176
PL Instrukcja użytkowania.....	87	KO 사용 설명서.....	181
HU Használati utasítás	95		




1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2021-01-22

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 4R220* Dynamic Vacuum System (DVS) enables an active vacuum prosthesis. It consists of a component with an integrated piston pump, a liner and a sealing sleeve.

The Dynamic Vacuum System uses the pistoning movements between the residual limb and prosthetic socket to generate the vacuum. The piston of the pump is equipped with magnets. A metallic counter-piece is installed on the liner. Air is suctioned out of the space between the liner and prosthetic socket into the cylinder during the swing phase, and this air is pressed out during the stance phase. The result is a permanent vacuum, with a level that adapts to the user's activity level.

1.2 Combination possibilities

INFORMATION

To ensure proper product function, only use the product with the components specified in this section.

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

Allowable liner: 6Y94

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthesis fittings.

2.2 Area of application

Product approved for max. body weight

- The maximum approved body weight is specified in the technical data (see page 22).

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

3 Safety

3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

3.2 General safety instructions

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- Use the product according to the specified area of application (see page 16).

Improper combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- Only combine the product with prosthetic components that are approved for that purpose.
- Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

Use under restricted environmental conditions

Risk of injury due to damage to the product

- Do not expose the product to restricted environmental conditions.
- If the product has been exposed to restricted environmental conditions, check it for damage.
- If damage is apparent or in case of doubt, do not continue using the product.
- Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Exceeding the lifetime and reuse

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the maximum lifetime is not exceeded.
- ▶ Only use the product for one person.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION

Product contains a magnet

Influence on or damage to devices or items due to a strong magnetic field

- ▶ Keep the product away from devices and items that are sensitive to magnetic fields.

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

4R220 Dynamic Vacuum System (see fig. 1)			
Item	Quantity	Designation	Reference number
	1	Instructions for use	–
1	1	Cylinder body	4X324=3

4R220 Dynamic Vacuum System (see fig. 1)			
Item	Quantity	Designation	Reference number
2	1	Piston	–
3	1	Fixation ring	4X339
4	1	Spacer plate	4X314
5	4	Snap bushing	5X163
6	4	Countersunk head screw	501S128=M6x22
	4		501S128=M6x25
7	1	Valve bushing	–
8	1	Duckbill valve	–
9	1	Plaster screw	503S3
10	1	Valve dummy	–
11	1	Silicone dummy	–
12	1	Piston dummy	–
13	1	Polylub GLY 801 special grease	633F30=2
14	1	Mounting wrench	4X338

Spare parts/accessories (not included in the scope of delivery)

Designation	Reference number
Single-component pack (1 piston, 1 fixation ring)	4X320=N
Single-component pack (1 valve bushing, 2 duckbill valves)	4X322
Dummy set (1 piston dummy, 1 silicone dummy, 1 valve dummy, 1 plaster screw)	4X326

5 Preparing the product for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

INFORMATION

- ▶ Not all of the materials may be available in your country. In this case, please contact the local branch of the manufacturer to obtain information on alternative materials.

5.1 Preparing for socket fabrication

INFORMATION

The Dynamic Vacuum System can be used with a total surface weight-bearing socket or a specific weight-bearing socket. Observe the following instructions:

- ▶ Mould the residual limb length with liner unchanged while taking the plaster cast.
- ▶ Do not shorten the plaster model during modelling.

> **Required materials:** piston dummy, silicone dummy, short or long valve dummy, 99B25 nylon stockinette, 99B81=70X19X5 PVA bag, string, 633W8 wax, 636K8* Plastaband

- 1) Distally flatten the plaster model so that the piston dummy can be positioned in the longitudinal residual limb axis (see fig. 2). Position the piston dummy on the plaster model and loosely secure it with a plaster screw.
- 2) Pull a nylon stockinette over the plaster model.
- 3) Tie off the nylon stockinette between the piston dummy and plaster model and trim off the excess nylon stockinette (see fig. 3).
- 4) Turn the plaster screw all the way in (see fig. 4).

- 5) **INFORMATION: No PVA bag is required for fabricating a check socket.**

Pull a soaked PVA bag over the plaster model.

- 6) Use string to tie off the PVA bag in the undercut of the piston dummy (see fig. 5).
- 7) Cut off the excess PVA bag with a scalpel.
- 8) Position the silicone dummy on the piston dummy. Check that no excess PVA bag is protruding during this process (see fig. 6).
- 9) **INFORMATION: Do not close the discharge opening in the cylinder body from the inside. Screw in the valve dummy after positioning the cylinder body.**
Position the cylinder body on the piston dummy. Perform the mediolateral alignment of the exhaust opening depending on the prosthetic alignment.
- 10) Isolate the heads of the cap screws with wax.
- 11) Fill the heads of the cap screws with Plastaband (see fig. 7).
- 12) Screw the valve dummy into the exhaust opening until the O-ring closes the opening on the valve dummy (see fig. 8).
- 13) Close the slot in the valve dummy with Plastaband.

5.2 Fabricating the prosthetic socket

5.2.1 Optional: fabricating the check socket

> **Recommended materials:** ThermoLyn (e.g. 616T52* ThermoLyn rigid, 616T83* ThermoLyn clear), 616S5* BetaSil

- 1) Prepare for fabrication of the prosthetic socket (see page 18).
- 2) Attach two strips of Plastaband lengthwise over the cylinder on both the anterior and the posterior side. This makes it easier to remove the cylinder after the trial fitting (see fig. 9).
- 3) Attach a strip of Plastaband circularly over the proximal edge of the cylinder to seal the prosthetic socket (see fig. 9).
- 4) Complete the vacuum forming process.
- 5) Finish the prosthetic socket (see page 20).
- 6) **If the prosthetic socket leaks:** seal the transition from the cylinder body to the inside of the prosthetic socket with BetaSil.

5.2.2 Laminating the transtibial socket

INFORMATION

The layout described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the layout.

- > **Required materials:** 623T3=8 or 623T3=10 Perlon stockinette, 81A1=8 or 81A1=10 nylon stockinette, 616B1=25x* carbon fibre cloth strap, 616G2 carbon UD stockinette, 616G15 carbon fibre stockinette, 99B81=100X19X5 PVA bag, 617H119 Orthocryl lamination resin 80:20 PRO, string, 627B40 polyethylene adhesive tape
- 1) Prepare for fabrication of the prosthetic socket (see page 18).
- 2) Cut a piece of Perlon stockinette (twice the length of the plaster model).
- 3) Pull the Perlon stockinette over the plaster model to the edge. Tie off the second half of the Perlon stockinette and fold it over the plaster model (see fig. 10).
- 4) Expose the valve dummy using a scalpel (see fig. 11).
- 5) Wrap string around the valve dummy and tie the Perlon stockinette off tightly in the distal groove of the cylinder body (see fig. 12).
- 6) **INFORMATION: Tie this off tightly so the prosthetic socket is airtight after the lamination process.**
Tie off the two layers of Perlon stockinette in the undercut of the cylinder body so they are up against the cylinder body (see fig. 13).
- 7) Position a layer of carbon fibre cloth strap from the medial condyle over the cylinder body to the lateral condyle (see fig. 14).
- 8) Position a layer of carbon fibre cloth strap from the MPT (mid-patella-tendon) point over the cylinder body to the opening for the back of the knee (see fig. 15).
- 9) Position a layer of carbon fibre cloth strap circularly at the height of the MPT (mid-patella-tendon) point (see fig. 16).

- 10) Pull a nylon stockinette over the plaster model.
- 11) Expose the valve dummy.
- 12) Wrap string around the valve dummy and tie the nylon stockinette off tightly in the distal groove of the cylinder body (see fig. 17).
- 13) Tie the nylon stockinette off tightly in the undercut of the cylinder body with Perlon string (see fig. 18).
- 14) Cut off a piece of carbon fibre stockinette (1.3 times the length of the plaster model).
- 15) Pull the carbon fibre stockinette over the plaster model to the edge (see fig. 19).
- 16) Tie off the excess carbon fibre stockinette at the distal end and fold it over the plaster model (see fig. 20).
- 17) Tie off the carbon fibre stockinette tightly in the distal groove of the cylinder body.
- 18) Tie off the carbon fibre stockinette tightly in the undercut of the cylinder body using Perlon string.
- 19) Expose the valve dummy.
- 20) Pull a nylon stockinette over the plaster model.

Prosthetic socket with carbon finish

- 1) Cut a piece of carbon UD stockinette (twice the length of the plaster model).
- 2) Pull the carbon UD stockinette over the plaster model up to the edge.
- 3) Tie off the excess carbon UD stockinette at the distal end and fold it over the plaster model.
- 4) Pull a nylon stockinette over the plaster model.
- 5) Soak the longer PVA bag and pull it over the plaster model.
- 6) Carry out the lamination process with Orthocryl.
- 7) **Once the lamination resin is sufficiently distributed:** shift excess lamination resin in the area of the cylinder body to the distal end by wrapping with polyethylene adhesive tape. Wrap tightly to make the prosthetic socket airtight.
- 8) Allow the lamination resin to cure.

- 9) Finish the prosthetic socket (see page 20).

Prosthetic socket without carbon finish

- 1) Cut a piece of carbon UD stockinette (1.5 times the length of the plaster model).
- 2) Pull the carbon UD stockinette over the distal half of the plaster model (see fig. 21).
- 3) Tie off the excess carbon UD stockinette at the distal end and fold it over the plaster model (see fig. 22).
- 4) Tie off the carbon UD stockinette in the undercut of the cylinder body (see fig. 23).
- 5) Cut a piece of Perlon stockinette (twice the length of the plaster model).
- 6) Pull the Perlon stockinette over the plaster model to the edge. Tie off the second half of the Perlon stockinette and fold it over the plaster model (see fig. 24).
- 7) Soak the longer PVA bag and pull it over the plaster model (see fig. 25).
- 8) Carry out the lamination process with Orthocryl.
- 9) **Once the lamination resin is sufficiently distributed:** shift excess lamination resin in the area of the cylinder body to the distal end by wrapping with polyethylene adhesive tape. Wrap tightly to make the prosthetic socket airtight.
- 10) Allow the lamination resin to cure.
- 11) Finish the prosthetic socket (see page 20).

5.3 Finishing the prosthetic socket

> **Required materials:** mounting wrench, snap bushings, 617H46 silicone bonding agent, 633F30=2 Polylub GLY 801

- 1) Mark and trim the contour of the prosthetic socket.
- 2) Sand free the valve dummy and remove the Plastaband.
- 3) Unscrew the valve dummy and retain it. After fabricating the check socket, the valve dummy is still required for the fabrication of the final prosthetic socket.
- 4) Remove the prosthetic socket from the plaster model.

- 5) Remove the piston dummy together with the silicone dummy from the plaster model and insert into the cylinder so it is closed (see fig. 26).
- 6) Close the opening in the piston dummy with Plastaband.
- 7) Screw in the valve dummy (see fig. 27).
- 8) Sand the contour of the prosthetic socket.
- 9) Sand the distal end of the socket down to the heads of the cap screws (see fig. 28).
- 10) Remove the cap screws.
- 11) Sand the distal end of the socket flat and check that its surface is even.
- 12) Remove any dirt from the 4 openings in the distal end of the socket. This ensures that the snap bushings will lie flat against the cylinder body.
- 13) **INFORMATION: The snap bushings must protrude approx. 1 mm from the laminate. This ensures forces are transferred to the Dynamic Vacuum System.**
Apply silicone bonding agent to the snap bushings and insert them into the openings with the groove end first (see fig. 29).
- 14) Unscrew the valve dummy.
- 15) Remove the piston dummy and silicone dummy from the prosthetic socket.
- 16) Wipe out the cylinder space with a lint-free cloth.
- 17) Grease the cylinder space and O-rings of the valve bushing with the supplied Polylub GLY 801.
- 18) **NOTICE! Do not allow the duckbill valve to come into contact with grease.**
Insert the duckbill valve into the valve bushing with the pointed end first (see fig. 30).
- 19) Screw the valve bushing into the prosthetic socket (torque value: **3 Nm**).
- 20) Insert the piston into the cylinder to the stop (see fig. 31).
- 21) Insert the fixation ring into the thread of the cylinder and tighten it with the mounting wrench (see fig. 32).

5.4 Fastening the socket adapter

The connection to the distal prosthetic components is established by a socket adapter. A spacer plate is positioned on the snap bushings in the distal end of the socket when mounting the socket adapter. To ensure forces are transferred to the Dynamic Vacuum System, the spacer plate must not lie against the laminate.

> Required tools and materials:

710D4 torque wrench, 636K13 Loctite®, spacer plate, 501S128=M6x22 or 501S128=M6x25 countersunk head screws

- 1) Set the spacer plate onto the prosthetic socket (see fig. 33).
- 2) **If the pyramid receiver or the pyramid adapter is rotatable:** Place the pressure plate on the adapter.
- 3) Apply the adapter to the prosthetic socket (see fig. 34).
- 4) **Optional:** Align the pyramid receiver or the pyramid adapter.
- 5) **NOTICE! Only use the screws specified here.**

Select suitable countersunk head screws.

Fixed pyramid adapter/pyramid receiver: 501S128=M6x22

Swivelling pyramid adapter/pyramid receiver:
501S128=M6x25

- 6) Secure the screws with Loctite®.
- 7) Screw in the 2 posterior countersunk head screws and tighten them (torque value: **12 Nm**).
- 8) Screw in the 2 anterior countersunk head screws and tighten them (torque value: **12 Nm**).

6 Use

⚠ CAUTION

Using the product without a sealing sleeve

Risk of injury due to unexpected loosening of the connection

- ▶ Always use the product with a sealing sleeve so that a vacuum can be generated.

6.1 Donning the prosthesis

- 1) Apply the liner to the residual limb so that the distal end is aligned with the longitudinal axis of the residual limb. Unroll the liner onto the residual limb so that it has no wrinkles, without shifting the soft tissue.
- 2) Slide into the prosthetic socket with the liner. The magnets in the piston of the Dynamic Vacuum System connect to the metallic counter-piece in the liner.
- 3) Unroll the sealing sleeve over the prosthetic socket onto the thigh.

6.2 Doffing the prosthesis

- 1) Unroll the sealing sleeve from the thigh onto the prosthetic socket.
- 2) Doff the residual limb and liner out of the prosthetic socket.

7 Cleaning

⚠ CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Dirt in the prosthetic socket or on the liner can impair the function of the Dynamic Vacuum System.

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

8 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

- To lubricate the cylinder or piston seal, use **Polylub GLY 801**.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Reference number	4R220
Weight [g]	210
System height [mm]	37
Material	Aluminium
Max. body weight [kg]	150

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