


## Kenevo Protector 4X840

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	4
<b>EN</b> Instructions for use (qualified personnel) .....	9
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	14
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	20
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	25
<b>PT</b> Manual de utilização (Pessoal técnico) .....	30
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	35
<b>SV</b> Bruksanvisning (Fackpersonal) .....	40
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	45
<b>NO</b> Bruksanvisning (Fagpersonell) .....	50
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	55
<b>PL</b> Instrukcja użytkowania (Personel fachowy) .....	60
<b>HU</b> Használati utasítás (szakszemélyzet) .....	66
<b>CS</b> Návod k použití (Odborný personál) .....	71
<b>SK</b> Návod na používanie (Odborný personál) .....	76
<b>TR</b> Kullanma talimatı (Uzman personel) .....	81
<b>RU</b> Руководство по применению (Квалифицированный персонал) .....	86

## INFORMATION

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### INFORMATION

Date of last update: 2021-03-01

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "4X840 Kenevo Protective Cover" is referred to simply as the Protective Cover/product below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The product consists of the following components:

1. Protective cover main unit
2. Upper closure (one side always stays on the main unit)
3. Lengthwise closure
4. Flexible cover cap

## 2.2 Function

The protective cover provides a means of enclosing the Kenevo electronic knee joint, including the tube adapter, to enhance its appearance and function. The housing of the knee joint is protected against external influences by the protective cover, e.g. scratching and surface contamination. The protective cover is not required for the proper and safe functioning of the joint.

## 2.3 Combination possibilities

Prosthetic components that are suited for combination with the product are listed below.

### Knee Joint Combination Possibilities

- Kenevo: 3C60, 3C60=ST

## 3 Intended use

### 3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

### 3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The fitting of patients with the product may only be carried out by qualified personnel (such as an O&P professional).

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

## 4.2 Structure of the safety instructions

### CAUTION

#### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

## 4.3 General safety instructions

### CAUTION

#### **Non-observance of safety instructions**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety instructions and the stated precautions in this accompanying document.

### CAUTION

#### **Improper disassembly, assembly and handling of the product**

- > Injury due to pinching of body parts.
- > Possible injury in case of parts with sharp edges.
- ▶ Use increased caution during disassembly/assembly and handling of the product.
- ▶ Instruct the patient in the proper disassembly, assembly and handling of the product.

### NOTICE

#### **Mechanical damage to the product**

Change in or loss of functionality due to damage.

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

## 5 Scope of delivery

- 1 pc. Kenevo Protective Cover main unit
- 1 pc. 4P840=R lengthwise closure
- 1 pc. 4P860=U upper closure
- 1 pc. 4P8 Distal cap
- 1 pc. Instructions for use (qualified personnel)
- 1 pc. Instructions for use (user)

## 6 Preparing the product for use

### 6.1 Required tools

- 718H5 Deburring knife
- 743B2 Tape measure
- 627B6=15 Crepe adhesive tape
- 636K13 adhesive (Loctite 241)
- Vibrating/band saw
- 710D20 Torque wrench

## 6.2 Shortening protective cover

- 1) Remove the prosthesis from the residual limb.
- 2) Check whether the closures on the protective cover are properly fastened and that the protective cover is not twisted.
- 3) Lay the protective cover next to the prosthesis.
- 4) Measure the length from the centre of the knee to the lower edge of the tube adapter on the prosthesis (see fig. 1).
- 5) Mark the measured length at at least three places all around the protective cover (see fig. 2).
- 6) **INFORMATION: For tall patients, shortening the protective cover may not be required.**
- 7) Apply crepe adhesive tape proximally all around at the mark.
- 8) Use a vibrating or band saw to shorten the protective cover along the bottom edge of the crepe adhesive tape (see fig. 3).  
**INFORMATION: Make sure that the closures are not shifted while sawing.**
- 9) Deburr the cut edge.
- 10) Using suitable tools, remove the profile rib on the interior of the protective cover (see fig. 4).

### NOTICE

#### Collision between protective cover and prosthetic foot

Damage to or functional impairment of the protective cover or prosthetic foot

- ▶ Carry out a function test with the user wearing the protective cover.
- ▶ If the protective cover collides with the prosthetic foot, shorten the protective cover again.

## 6.3 Protective Cover Installation

- 1) Remove the prosthetic foot from the tube adapter.  
In doing so, unscrew the two set screws that are screwed in the furthest but are not opposite each other. Mark the two set screws.
- 2) Slide the distal cap over the tube adapter (see fig. 5).
- 3) Install the prosthetic foot on the tube adapter (see fig. 6).
- 4) Apply Loctite 241 to the two set screws and screw them in.

### CAUTION

#### Improperly secured screws

Falling due to breakage of load-bearing components caused by loosened screw connections.

- ▶ The set screws in the tube adapter must be secured before they are tightened to the specified tightening torque.
- ▶ See the knee joint instructions for use for the specified tightening torque of the set screws.

- 5) Remove the crepe adhesive tape from the Protective Cover after final shortening.
- 6) Open or remove the closures.
- 7) Spread open the Protective Cover (see fig. 7).
- 8) Apply the Protective Cover to the knee.  
**INFORMATION: The retainers of the Protective Cover have to engage on the frame element (see fig. 8).**
- 9) Fasten the upper closure (see fig. 9) and lengthwise closure (see fig. 10) on the Protective Cover.
- 10) Fit the distal cap to the distal end of the Protective Cover main component. Align the lengthwise grooves in the extension of the lengthwise closure on the dorsal side of the Protective Cover (see fig. 11).

## NOTICE

### **Loss of stability due to inadequately installed retainers**

Protective Cover does not stay in place on the knee joint.

- ▶ Ensure that all retainers and closures are properly engaged and/or fastened.

## **6.4 Functional Test**

After assembly, the following points have to be checked:

- Has any collision between the protective cover and prosthetic foot been eliminated?
- Is the function of the rotation adapter impaired?

Material may be removed from the protective cover to improve the functionality.

**INFORMATION: If material is removed, the stability/flexibility of the protective cover on the knee joint may be impaired.**

## **7 Maintenance**

- ▶ During the standard inspection of the prosthesis components, also inspect the product.

## **8 Cleaning and Care**

### NOTICE

#### **Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

## **9 Legal information**

### **9.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### **9.2 CE conformity**

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

### **9.3 Trademarks**

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 10 Technical data

Ambient conditions	
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F max. 93% relative humidity, non-condensing
Operation	-10 °C/+14 °F to +40 °C/+104 °F max. 93% relative humidity, non-condensing

General	
Reference number	4X840
Product service life	Wear part subject to normal wear and tear
Protective cover weight (with closures)	380 g
Flexible cover cap weight	11 g

## 11 Appendices

### 11.1 Symbols Used

Manufacturer

Declaration of conformity according to the applicable European directives

Lot number (PPPP YYYY WW)  
PPPP – plant  
YYYY – year of manufacture  
WW – week of manufacture

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