
## 17B206, 17B206=7.5

<b>DE</b>	Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b>	Instructions for use (qualified personnel) .....	5
<b>FR</b>	Instructions d'utilisation (Personnel spécialisé) .....	7
<b>IT</b>	Istruzioni per l'uso (Personale tecnico specializzato) .....	10
<b>ES</b>	Instrucciones de uso (Personal técnico especializado) .....	12
<b>PT</b>	Manual de utilização (Pessoal técnico) .....	15
<b>NL</b>	Gebruiksaanwijzing (Vakmensen) .....	17
<b>SV</b>	Bruksanvisning (Fackpersonal) .....	20
<b>DA</b>	Brugsanvisning (Faguddannet personale) .....	22
<b>NO</b>	Bruksanvisning (Fagpersonell) .....	25
<b>FI</b>	Käyttöohje (Ammattihenkilöstö) .....	27
<b>CS</b>	Návod k použití (Odborný personál) .....	29
<b>SK</b>	Návod na používanie (Odborný personál) .....	32
<b>RU</b>	Руководство по применению (Квалифицированный персонал) ...	34
<b>ZH</b>	使用说明书 (专业人员) .....	37

1

## INFORMATION

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### INFORMATION

Last update: 2016-01-29

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the processing of the 17B206 and 17B206=7.5 medial support.

## 2 Product description

### 2.1 Components

Component	Article number
Medial support with dummy, 5° pre-flexion	17B206
Medial support with dummy, 7.5° pre-flexion	17B206=7.5

### 2.2 Design

The products 17B206 and 17B206=7.5 are medial supports. The medial supports are used in the medial position in the E-MAG system.

## 3 Application

### 3.1 Indications for use

The system knee joints are to be used **exclusively** for the orthotic fitting of the lower limbs on a patient.

The 17B206 medial support may only be used in an E-MAG Active orthosis for additional medial stabilisation. Using the medial support is mandatory for a patient weight over 85 kg up to a maximum of 100 kg. The 647G1165 instructions for use must be observed.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 General safety instructions

 <b>CAUTION</b>
<b>Incorrect alignment, assembly or adjustment</b> Risk of injury due to change in or loss of functionality <ul style="list-style-type: none"><li>▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.</li><li>▶ Observe the alignment, assembly and adjustment instructions.</li></ul>

## 5 Preparation for use

### 5.1 Mounting the accessories

- > **Additional mounting material and tools required:** 2x 17LS3=16(-T) lamination bars, Loctite® 241 adhesive.
- 1) Integrate the medial support in the orthosis being fabricated. Use the dummy for installation in the orthosis.
  - 2) For trial fitting, attach the medial support to the orthosis using the screws (torque value 6 Nm).
  - 3) Upon completion of the orthosis, attach the medial support to the orthosis using the screws and secure the screw connections with Loctite® 241 (torque value 6 Nm).

## 6 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer recommends a semi-annual inspection to verify functionality and check for wear of the system joints.

Follow the maintenance schedule included with the E-MAG instructions for use.

## 7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as class I accessories according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH