

5065N=* Omo Neurexa plus


DE Gebrauchsanweisung	5
EN Instructions for use	8
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	15
ES Instrucciones de uso	18
NL Gebruiksaanwijzing	22
SV Bruksanvisning	25
DA Brugsanvisning	28
NO Bruksanvisning	31
FI Käyttöohje	35
PL Instrukcja użytkowania	38
HU Használati utasítás	41
CS Návod k použití	45
SK Návod na používanie	48
TR Kullanma talimatı	51
RU Руководство по применению	54
JA 取扱説明書	58
ZH 使用说明书	61

	Größe/Size	Brustumfang/Chest circumference	
		cm	inch
	XXS	71 – 78	27 – 31
	XS	79 – 86	31 – 34
	S	87 – 94	34 – 37
	M	95 – 102	37 – 40
	L	103 – 110	40 – 43
	XL	111 – 118	43 – 47

Material	PE, PES, PU, Stahl / steel, Spiralfeder / coil spring
-----------------	---

INFORMATION

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INFORMATION

Date of last update: 2020-05-11

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on fitting and applying the 5065N Omo Neurexa plus orthosis.

2 Product description

2.1 Product overview

Item:	Description	Item:	Description
1	Shoulder cuff	2	Forearm cuff
3	Connecting straps, short	4	Strap guide loop
5	Connecting straps, long	6	Muscle stimulation pad
7	Axilla pad	8	Upper arm strap
9	Support straps		

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the shoulder and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Shoulder subluxation with neurological origin, or painful restriction of movement in the shoulder and arm musculature due to neurological causes

Indications must be determined by the physician.

3.3 Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distal to the support; sensory and circulatory disorders in the area of the arm and hand.

3.4 Effects

The orthosis relieves and supports the shoulder, and ensures slight outward rotation of the shoulder.

4 Safety

4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

4.2 General safety instructions

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Size Selection

- 1) Measure the chest circumference.
- 2) Determine the orthosis size (see size chart).

5.2 Fitting and Application

⚠ CAUTION

Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Open all hook-and-loop closures.
- 2) Guide the arm through the shoulder cuff (see fig. 2).
- 3) Apply the shoulder cuff to the shoulder (see fig. 3).
- 4) Fasten the hook-and-loop closure on the upper arm strap (see fig. 4). In doing so, make sure the red seam is below the shoulder joint gap.
- 5) Bring the chest strap around the chest (see fig. 5). In doing so, position the axilla pad under the axilla.
- 6) Fasten the chest strap buckle on the shoulder cuff.
- 7) Adjust the length of the chest strap with the hook-and-loop closure (see fig. 6). To do this, shorten the chest strap with scissors and fasten it with an Y-hook-and-loop.
- 8) Guide the forearm through the forearm cuff (see fig. 7). In doing so, the snap fasteners on the forearm cuff have to be on top.
- 9) Make sure that the snap fasteners of the forearm strap are close to the pivot point of the elbow joint without allowing the orthosis to cover the olecranon (see fig. 8). In doing so, the closure on the forearm cuff has to be at the front.

- 10) Fasten the hook-and-loop closure on the forearm strap (see fig. 9).
- 11) Choose the connecting straps and attach them to the snap fasteners on the shoulder cuff (see fig. 1 item: 3 and item: 5).
- 12) Attach the strap guide loops to the snap fasteners on the forearm cuff.
- 13) Guide the connecting straps through the strap guide loop and adjust the length with the hook-and-loop closure (see fig. 10). The strain on the shoulder joint should be relieved and the shoulder should be rotated slightly outwards.
- 14) Adjust the length of the support straps with the hook-and-loop closure. In doing so, make sure that the humeral head is raised.
- 15) Optional: Attach the muscle stimulation pad in the desired position with hook-and-loop (see fig. 11).
- 16) Check the fit of the orthosis. Adjust the length of the connecting and supporting straps as needed.

6 Cleaning

- 1) Close all hook-and-loop closures.
- 2) Hand wash the orthosis in warm water at 40°C with mild laundry detergent.
- 3) Rinse well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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