
## 28P30 Manu Neurexa plus, 28P31

<b>DE</b> Gebrauchsanweisung .....	10
<b>EN</b> Instructions for use .....	16
<b>FR</b> Instructions d'utilisation .....	22
<b>IT</b> Istruzioni per l'uso .....	28
<b>ES</b> Instrucciones de uso .....	35
<b>NL</b> Gebruiksaanwijzing .....	41
<b>SV</b> Bruksanvisning .....	48
<b>DA</b> Brugsanvisning .....	54
<b>NO</b> Bruksanvisning .....	59
<b>FI</b> Käyttöohje .....	65
<b>PL</b> Instrukcja użytkowania .....	71
<b>HU</b> Használati utasítás .....	77
<b>CS</b> Návod k použití .....	83
<b>SK</b> Návod na používanie .....	89
<b>TR</b> Kullanma talimatı .....	95
<b>RU</b> Руководство по применению .....	101
<b>JA</b> 取扱説明書 .....	108
<b>ZH</b> 使用说明书 .....	114

18

28P31

19

20


21

22


<b>Material</b>	PP, PA, PA/EL, PU, Aluminium
-----------------	------------------------------

**INFORMATION**

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### INFORMATION

Date of last update: 2022-03-10

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 28P30 Manu Neurexa plus wrist orthosis and the 28P31 wrist positioning splint.

## 2 Product description

### 2.1 Product overview

28P30 Manu Neurexa plus (see fig. 1)		
Item	Designation	Reference number
1	Torx wrench	
2	Self-adhesive hook-and-loop patches	
3	Hand support bar	
4	Orthosis	
5	Padded strap for hand support bar	29P30
6	Forearm straps	
7	Forearm pad	

<b>28P31 Hand positioning splint (see fig. 18)</b>		
<b>Item</b>	<b>Designation</b>	<b>Reference number</b>
1	Hand support pad	29P31
2	Y-hook-and-loop closure	
3	Finger closure strap	
4	Thumb-metacarpus closure strap	
5	Hand support	

### 3 Intended use

#### 3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the wrist joint and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

#### 3.2 Indications

Pareses of the forearm and hand musculature

Paralysed forearm and hand musculature with slight spasticity, for instance after:

- Stroke
- Brachial plexus injuries
- Cervical intervertebral disc prolapse

Indications must be determined by the physician.

#### 3.3 Contraindications

##### 3.3.1 Absolute Contraindications

- Moderate and severe spasticity of the forearm and hand musculature

##### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distal to the support; sensory and circulatory disorders in the area of the arm and hand.

#### 3.4 Lifetime

##### 28P30

The product is designed for a maximum lifetime of **2 years**.

##### 28P31

The product is designed for a maximum lifetime of **1 year**.

#### 3.5 Effects

The orthosis stabilises the wrist joint; the position of the wrist joint is continuously variable, and the movement of the hand is limited in the palmar or dorsal direction. The hand can be held with either palmar or dorsal support. The orthosis facilitates therapeutic treatment, hand function and physical activity. The 28P30 Manu Neurexa plus can be combined with the 28P31 hand positioning splint and also the 5065N Omo Neurexa plus shoulder support.

### 4 Safety

#### 4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

**NOTICE**

Warning regarding possible technical damage.

## 4.2 General safety instructions

**⚠ CAUTION****Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**⚠ CAUTION****Incorrect or excessively tight application**

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

**⚠ CAUTION****Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**NOTICE****Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

**NOTICE****Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

**NOTICE****Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

## 5 Handling

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

## 5.1 Size Selection

The product is supplied in a universal size and can be used on both sides. The orthosis length and forearm supports can be adapted.

## 5.2 Adaptation

### Adapting the orthosis length

The length of the orthosis is changed by sliding the forearm supports on the arm bar. The arm bar can be adapted to the requirements of the patient.

- 1) Loosen the screws for the bar (see fig. 1, item 4) with the wrench (see fig. 1, item 1).
- 2) Place the patient's elbow in a 90° flexion position.
- 3) Apply the orthosis laterally on the forearm.
- 4) Position the supports of the orthosis between the wrist joint (see fig. 2) and approx. 2 cm from the olecranon (see fig. 3). **INFORMATION: The distal support has to be applied on the ulnar head. In case of patient-specific pressure sensitivity at this point, the support can be applied directly before it. Ensure that the edge of the distal forearm support does not press on the back of the hand.**
- 5) Mark the length adjustment on the bar.
- 6) Take off the orthosis and tighten the screws with the wrench (torque value: **1 Nm**) (see fig. 4).

### Adapting the orthosis angle

The angle adjustment range is **30°** max. in palmar flexion or dorsiflexion, with one position respectively corresponding to **10°**. The positions are marked by holes in the orthosis joint.

- 1) Loosen the screw on the orthosis joint with the wrench (see fig. 7).
- 2) Place the patient's elbow in a 90° flexion position.
- 3) Apply the orthosis laterally on the forearm and align it.
- 4) Set the orthosis joint to the desired angle.
- 5) Mark the angle position.
- 6) Remove the orthosis and tighten the screw with the wrench (torque value: **3.6 Nm**).

### Adapting the hand support bar

- 1) **INFORMATION: Align the hand support bar with the curve facing up.**  
Insert the hand support bar (see fig. 1, item 3) into the orthosis (see fig. 8).
- 2) Place the patient's elbow in a 90° flexion position.
- 3) Apply the orthosis laterally on the forearm.
- 4) Slide the hand support bar and position it under the palm according to the individual conditions.
- 5) Take off the orthosis and tighten the screw with the wrench (torque value: **1 Nm**).

### Adapting the forearm supports

- 1) Place the patient's elbow in a 90° flexion position.
- 2) Apply the orthosis laterally on the forearm.
- 3) Mark the required lengths of the supports (see fig. 5).
- 4) Trim the supports with scissors (see fig. 6).
- 5) Sand the cut edges smooth.
- 6) Check the fit of the orthosis again and rework it as needed.

### Adapting the forearm pad

> **Prerequisite:** Take off the orthosis.

- 1) Apply the self-adhesive hook-and-loop patches to the inside of the supports (see fig. 9).
- 2) Position the pads inside the supports according to their shape and fasten with hook-and-loop.
- 3) Individually trim the pads with scissors, leaving a **5 mm** overhang along the edges of the supports (see fig. 10).

### Adapting the forearm straps

The forearm straps are attached in the integrated slots in the supports. With correct assembly, the pull loop of the strap has to be on top.

- 1) Pull the strap through the lower slot in the support and fasten it with the hook-and-loop closure.
- 2) Pull the strap through the upper slot in the support and fasten it at the hook-and-loop closure, individually adapting the strap circumference to the arm.

### **Adapting the pad of the padded strap for the hand support bar**

> **Prerequisite:** Take off the orthosis.

- 1) Open the padded strap (see fig. 1, item 5).
- 2) Cut the interior plastic tube to length with scissors (see fig. 11), taking the patient's hand width into account as the measurement.
- 3) **Optional:** Trim the shorter padded hook-and-loop closure with scissors (see fig. 12).
- 4) Slide the padded strap over the hand support bar and fasten the hook-and-loop closure (see fig. 13, see fig. 14).

### **Adapting the padded strap for the hand support bar**

> **Prerequisite:** The orthosis is attached to the forearm straps.

- 1) Lay the hand onto the padded hand support bar (see fig. 15).
- 2) Pull the strap over the hand between the thumb and fingers, and fasten it with the hook-and-loop closure under the padded strap (see fig. 16). **INFORMATION: Close the strap at the bottom or top, depending on the patient.**
- 3) **Optional:** Shorten the hook-and-loop with scissors.

### **Adapting the 28P31 hand positioning splint**

- 1) Position the patient's hand on the hand support (see fig. 18, item 5).
- 2) Mark the thumb side and the finger sizes on the hand support (see fig. 19).
- 3) Trim the hand support to the correct size along the slots using scissors (see fig. 20).
- 4) Sand the cut edges smooth.
- 5) Starting on the bottom, thread the thumb-metacarpus closure strap (see fig. 21), positioning the strap guide loop on the thumb side. Slide the strap opening over the safety closure for the hand support.
- 6) Attach the self-adhesive hook-and-loop patch to the bar in the area of the centre of the thumb and use it to secure the strap (see fig. 22).
- 7) Shorten the closure strap so that there is sufficient material for attaching the Y-hook-and-loop closure. Attach the Y hook-and-loop closure. **Optional:** Sew the Y-hook-and-loop closure to the strap.
- 8) Starting on the bottom, thread the finger closure strap into the bar (see fig. 23, see fig. 24, see fig. 25).
- 9) Connect the hand support to the hand support bar (see fig. 26).
- 10) Fasten the safety closure (see fig. 27).
- 11) Position the patient's hand on the hand support.
- 12) Fasten the split closure strap on the thumb (see fig. 28).
- 13) Fasten the split closure strap on the metacarpus.
- 14) Fasten the finger closure strap (see fig. 29).

### **Applying the hand support pad**

- 1) Attach the self-adhesive hook-and-loop patches evenly distributed to the top of the hand support.
- 2) Put the pad onto the hand support.
- 3) Use scissors to cut the pad to size.

## **5.3 Application**

### **Application of the 28P30 Manu Neurexa plus**

> **Prerequisite:** All hook-and-loop closures are open.

> **Prerequisite:** Place the orthosis on a table in front of the patient, opened to the top.

- 1) Insert the forearm into the orthosis until the hand is lying fully on the hand support.

- 2) Fasten the straps.
- 3) Make sure that the straps fit closely without constriction.

#### **Application of the 28P30 Manu Neurexa plus with 28P31 hand positioning splint**

- > **Prerequisite:** All hook-and-loop closures are open.
  - > **Prerequisite:** Place the orthosis on a table in front of the patient, opened to the top.
- 1) Insert the forearm into the orthosis until the hand is lying fully on the hand positioning splint.
  - 2) Fasten the straps.
  - 3) Make sure that the straps fit closely without constriction.

#### **Application of the 28P30 Manu Neurexa plus with 5065N Omo Neurexa plus**

- > **Prerequisite:** The covers of the snap fasteners on the 28P30 Manu Neurexa plus have been removed.
  - > **Prerequisite:** The 28P30 Manu Neurexa plus has been put on.
  - > **Prerequisite:** The 5065N Omo Neurexa plus has been put on.
  - > **Prerequisite:** The suitable connecting straps have been attached to the snap fasteners of the shoulder cuff.
  - > **Prerequisite:** The forearm cuff of the 5065N Omo Neurexa plus has been removed.
- 1) Open the hook-and-loop closures of the connecting straps.
  - 2) Attach the strap guide loops to the snap fasteners on the wrist orthosis (see fig. 39, see fig. 40), connecting the snap fasteners of the same colour to each other (black-black/red-red).
  - 3) Guide the connecting straps through the strap guide loops from the inside to the outside and adjust the length with the hook-and-loop closure (see fig. 41, see fig. 42). The strain on the shoulder joint should be relieved and the shoulder should be rotated slightly outwards.
  - 4) Check the relief function on the shoulder. Adjust the length of the straps if necessary.

### **5.4 Retrofitting**

#### **Changing the hand support of the 28P31 hand positioning splint to the hand support bar of the 28P30 Manu Neurexa plus**

- 1) Open all straps (see fig. 30, see fig. 31).
- 2) Open the safety closure (see fig. 32).
- 3) Remove the hand support from the hand support bar (see fig. 33, see fig. 34).
- 4) Attach the hand support pad to the hand support bar with the fixation strap (see fig. 35, see fig. 36) and fasten the strap (see fig. 37).

### **5.5 Removal**

- 1) Open all straps.
- 2) Remove the orthosis from the forearm.

## **6 Cleaning**

### **NOTICE**

#### **Use of improper cleaning agents**

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

#### **Textile pads**

- 1) Remove all textile pads from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag or net.
- 4) Wash the textile pads in warm water at **30 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 6) Reattach the textile pads to the orthosis.

## Orthosis frame

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 7 Disposal

Dispose of the product in accordance with national regulations.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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