

743A7

DE Gebrauchsanweisung	3
EN Instructions for Use	4
FR Instructions d'utilisation	5
IT Istruzioni per l'uso	6
ES Instrucciones de uso	7
SV Bruksanvisning	8
NL Gebruiksaanwijzing	9
PT Manual de utilização	10
HU Használati utasítás	11
PL Instrukcja użytkowania	12
TR Kullanma talimatı	13

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Achtung!

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Attention!

Date of the last update: 2013-07-12

- Please read this document carefully.
- Follow the safety instructions.

743A7 Pivot point adjustment aid for alignment and positioning of joint axes for knees and ankles on plaster negatives or plaster positives.

2 Components (fig. 1)

- (a) Square bar, stainless steel
- (b) Holding fork for ankle and knee joint, plastic
- (c) Centering arbor, stainless steel
- (d) Thumb screw
- (e) Joint plates, plastic
- (f) Foot plate, aluminum
- (g) Thumb screw for ankle rotation
- (h) Fixation arbor

3 Function and Description

To mount the foot plate, insert the foot plate's stud into the square bar, and make sure that the bore holes are positioned correctly. Secure with the thumb screw (fig. 2).

The pivot point adjustment aid can be placed in the vise vertically (fig. 3).

The foot plate allows the plaster negative or plaster positive to be placed on it vertically (fig. 4).

The red thumb screws allows to rotate the ankle (fig. 4).

Align the plaster model in the holding forks using the fixation arbors and mark the pivot points with the centering arbors (fig. 5). For fixing the centering arbors, the thumb screws can be used.

The joint plates provided can be used to create parallel joint surfaces on the medial and lateral side of the plaster positive. Roughen the plaster in the joint areas. Slightly lubricate the joint plates with 640Z5 Plaster Parting Agent Cream. Apply either 636K17 Light Putty or plaster to the joint plates. Next, put the joint plates on the centering arbors and press against the plaster model.

4 Cleaning

If necessary, all components can be cleaned with soap solution and water and then wiped dry. From time to time, lubricate the fixation arbors and centering arbors with 633G6 Ottobock Special Lubricant.

5 Security

Attention!

Risk of injury, while working with the fixation arbors and the centering arbors!!

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HealthCare Inc.
South East Asia Co., Ltd.
Otto Bock HealthCare GmbH