

50A3=7 Epi Forsa Plus

DE Gebrauchsanweisung.....	3
EN Instructions for use.....	7
FR Instructions d'utilisation	12
IT Istruzioni per l'uso.....	17
ES Instrucciones de uso	22
PT Manual de utilização	27
NL Gebruiksaanwijzing	31
SV Bruksanvisning	36
DA Brugsanvisning.....	40
FI Käyttöohje	45
PL Instrukcja użytkowania	49
HU Használati utasítás	54
CS Návod k použití	59
HR Upute za uporabu	64
SK Návod na používanie.....	68
BG Инструкция за употреба	73
TR Kullanma talimatı.....	78
EL Οδηγίες χρήσης	82
RU Руководство по применению	87
JA 取扱説明書	93
ZH 使用说明书.....	97

Material	PE, PA, PU, EL, PES
-----------------	---------------------

INFORMATION

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INFORMATION

Date of last update: 2022-03-31

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50A3=7 Epi Forsa Plus elbow support.

2 Intended use

2.1 Indications for use

The support is intended **exclusively** for treating the upper limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

2.2 Indications

- Epicondylitis humeri radialis (tennis elbow)
- Epicondylitis humeri ulnaris (golfer elbow)

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent, swollen scars, reddening and hyperthermia in the relevant area; pronounced varicose veins, especially with impaired blood return, lymph drainage disorders and unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the forearm and hand, e.g. associated with diabetic neuropathy.

2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

2.5 Mechanism of Action

The support relieves the muscles of the lower arm and reduces pain.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

 CAUTION
Contact with heat, embers or fire
Risk of injury (such as burns) and risk of product damage
▶ Keep the product away from open flames, embers and other sources of heat.

 CAUTION
Reuse on other persons and improper cleaning
Skin irritation, formation of eczema or infections due to contamination with germs
▶ The product may be used by one person only.
▶ Clean the product regularly.

NOTICE
Contact with oils, salves, lotions or other products that contain oils or acids
Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

The product is available in a universal size.

4.2 Adaptation

The contour of the pad is marked on the support.

- 1) **Treating the left arm:** Affix the pad to the support. The wide side of the pad is on the **left**.
- 2) **Treating the right arm:** Affix the pad to the support. The wide side of the pad is on the **right**.

4.3 Application

Determine the pad position with epicondylitis humeri radialis.

- 1) Make the fist into a ball.
- 2) Raise the hand. This tenses the muscle on the outside of the forearm on which the pad will be located.

Determine the pad position with epicondylitis humeri ulnaris.

- ▶ Flex the hand downwards This tenses the muscle on the inside of the forearm on which the pad will be located.
- > The hook-and-loop closure of the support is open.
- 1) Pull the support over the forearm.
- 2) Position the pad on the muscle concerned. In doing so, place the lower edge of the support about **4 to 5 cm** below the elbow joint.
- 3) Pull the hook-and-loop strap through the strap guide loop and fasten it. Make sure that the support sits tightly but does not pinch.

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.
- 1) Fasten the hook-and-loop closure.
 - 2) **Recommendation:** Use a laundry bag or net.
 - 3) Wash the support in warm water at **30 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
 - 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage

caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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