
## 3R66

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	6
<b>FR</b> Instructions d'utilisation .....	8
<b>IT</b> Istruzioni per l'uso .....	11
<b>ES</b> Instrucciones de uso .....	14
<b>PT</b> Manual de utilização .....	17
<b>NL</b> Gebruiksaanwijzing .....	20
<b>SV</b> Bruksanvisning .....	22
<b>PL</b> Instrukcja użytkowania .....	25
<b>HU</b> Használati utasítás.....	28
<b>EL</b> Οδηγίες χρήσης.....	30
<b>RU</b> Руководство по применению.....	33

①


②

③

④

⑤

## INFORMATION

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## INFORMATION

Date of last update: 2020-12-07

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

## 1 Description and function

### 1.1 Application

The 3R66 Modular Knee Joint is to be used for the prosthetic fitting of the lower limb, **exclusively** with children from the ages of 3 to 6 years. It is limited to children weighing less than 35 kg and having a foot size not exceeding 21 cm.

The 3R66 Modular Knee Joint is not suitable for hip disarticulation.

#### Attention – Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.

### 1.2 Construction and function

The 3R66 Modular Polycentric Knee Joint is made of aluminum. Its main components are –

- upper joint section with pyramid
- mid-joint section with adjustable stop
- lower joint section with extension assist mechanism and clamp shell for tube adapter and integrated rotation unit.

The sections are connected to each other by four axes and two polycentric linkage bars\*.

This polycentric joint provides a rotational and translational movement and a shifting of its center of rotation during flexion and extension. In the extended position the center of rotation (instantaneous center of rotation) is placed above the pyramid and behind the weight bearing line, thus providing greater knee stability during stance phase.

The extension assist can be adjusted for control during swing phase. Reversal of the extension assist action occurs when flexion reaches approximately 90°.

A new technical feature of this knee is its **integrated rotation** capability which allows for the actual rotation of the prosthetic foot (fig. 1, 2), if desired.

Combined with a generous flexion angle of approx. 165° (fig. 2) this enables the child to easily move into a squatting, kneeling or sitting position with the foot rotated to the outside.

The prosthetic foot will automatically turn back to the neutral position when unloaded.

### 1.3 Lifetime

This product has been load-tested by the manufacturer. The maximum lifetime is 2 years.

\* the enclosed decorative decals can be attached here.

## 2 Alignment and assembly (fig. 5)

The starting point for achieving correct static alignment is to determine the individual weight bearing line and to transfer the measurements to the prosthetic components.

The position of the axes influences the joint's function. The joint is factory adjusted to the basic alignment, where the reference line runs in front of the anterior upper axis (see diagram).

The use of the 1K10 Dynamic Foot for Children or the 1K30/ 1S30 SACH Foot makes for an excellent combination with this knee joint. The tube adapter connects the knee joint to the prosthetic foot. Torque moment for tightening the clamp screw: 5 Nm (3.7 ft.lbs.) using the 710D1 Torque Wrench. Torque moment for tightening the adjustment screws of the adapter: 9 Nm (6.6 ft. lbs.) using the 710D1 Torque Wrench.

Secure the screws using 636K13 Loctite.

## 3 Adjusting stance and swing phase control

Before making any changes to the factory-settings, the following must be kept in mind:

### 3.1 Stance phase stability

This polycentric knee joint with its special joint geometry and arrangement of axes has been designed for the child patient to comply with his/her requirements regarding stability. **The joint is factory-set to maximum stability.** It can be adjusted by turning the stop screw using a 3 mm wrench (fig. 3).

Turning to the **right = less stability** = (flexion can be easily initiated).

Turning to the **left = more stability** = (flexion is more difficult).

### 3.2 Adjusting tension of the extension assist

The extension assist facilitates extension movement during swing phase. **The knee joint has been factory-set to maximum tension**, so that, when performing first trial walking, the joint will safely come into contact with the extension stop. By loosening the adjustment screw, tension can be reduced (fig. 4).

## 4 Maintenance instructions

After the knee joint has been worn for several weeks as well as **every 3 months** at least, Otto Bock recommends to inspect the knee joint and make any necessary adjustments. Particular attention must be paid to the resistances to movement and abnormal noises. It must be ensured that full extension is reached.

**During each maintenance of the prosthesis, the 4Z75=8x10 Extension Stop Bumper should be replaced** (available as single component). Use 617H46 silicone bonding agent to glue the extension stop bumper in the upper section.

**Do not disassemble the knee joint! In case of any problems, please return the entire knee joint to Otto Bock.**

### **Attention – Please inform your patients!**

The knee joint function may be affected by both environmental and service conditions. To reduce risk for the patient, the knee joint should not be used when noticeable functional changes occur. Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

### Action to be taken

Have the prosthesis examined at a specialist workshop.

**When using the knee joint do not reach into the joint mechanism since there is the risk of injury.**

### Attention

To eliminate noise in the cosmetic foam cover, use 519L5 Silicone Spray. **Do not use talcum powder!** Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of falling. Using this medical product after application of talcum powder will render all claims against Otto Bock HealthCare null and void.

Apply 519L5 Silicone Spray directly on the friction surfaces of the cosmetic foam cover!

## 5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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