












ottobock.

9X52, 9X53

| | | |
|-----------|----------------------------------|----|
| DE | Gebrauchsanweisung | 4 |
| EN | Instructions for use | 12 |
| FR | Instructions d'utilisation | 19 |
| IT | Istruzioni per l'uso | 27 |
| ES | Instrucciones de uso | 34 |
| PT | Manual de utilização | 42 |
| NL | Gebruiksaanwijzing | 50 |
| SV | Bruksanvisning | 58 |

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MD

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Date of last update: 2021-03-22

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

1 Application

The 9X52 Linear Transducer and the 9X53 4-Step Transducer are to be used **solely** in exoprosthetic fittings of the upper extremities for myoelectric control of arm prostheses.

Only start up the product in accordance with the information contained in the accompanying documents supplied.

Instruct the patient in the proper use and care of the product. Handover to the patient without prior instruction is not permitted.

2 Contraindications

All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3 Safety information

- When servicing the Transducers check the pull-cable as well as the housing for damage.
- When mounting a cable with cable housing, make sure there are no tight angles or kinks. Make sure the cable runs smoothly.
- Carefully read and observe the lamination instructions.

- When mounted into a provisional socket, the Transducers may only be worn by an amputee within an orthopedic workshop or as an in-patient in a rehabilitation facility.
- When installing the electrode cable, make sure there are no tight angles or kinks.
- Please note that the cable travel of the Transducers is limited. To prevent excessive strain on the Transducer, make sure the harness is correctly fit and accomodates the patient's movements.
- Instruct the patient in the correct use of the Transducer (see chapter 3 "Patient Information").
- The product and any damaged components may only be opened and/or repaired by certified Ottobock Myo-Service technicians.

4 Patient information

- Do not increase your tension on the cable once you have felt the cable come to a stop within the Transducer.
- Handle the Transducer with care and check the product regularly for damages (e.g. cracks or breakage). Should the Transducer be damaged or not function properly, please contact your Ottobock Service facility.
- Driving a vehicle or operating machinery is strictly prohibited while wearing a provisional socket fitted with a Transducer.
- Make sure that neither dirt nor fluids penetrate the Transducer. Do not expose the Transducer to intensive smoke or dust or to mechanic vibrations or shocks or high temperatures.

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

5 Range of application

5.1 9X52 Linear Transducer

The 9X52 Linear Transducer is suitable for patients who cannot generate enough muscle signal to control a myoelectric prosthesis but who want to take advantage of the high rehabilitation value of myoelectric prostheses.

The 9X52 Linear Transducer can be used to control proportional 1-channel and 2-channel electrode controls of prosthetic components such as the Ottobock System Electric Hand, System Electric Greifer or the Ottobock Eletrohand 2000 (Electronic Control 9E369, 9E370 or 9E420).

A combination of Linear Transducer and MyoBock Electrodes is also possible. By installing the Linear Transducer between the provisional socket and the shoulder harness the patient can begin learning how to control motion sequences during the first stages of rehabilitation and thus change more smoothly to the permanent prosthesis. The Linear Transducer is suitable for both right and left-sided amputees as well as double-sided amputees.

5.2 9X53 4-Step Transducer

Along with a properly configured prosthesis, the 4-Step Transducer allows a patient to switch between different components.

6 Description

6.1 9X52 Linear Transducer

By pulling on the cable of the 9X52 Linear Transducer a proportional control signal is created, comparable to the signal in an Ottobock MyoBock Electrode. The Linear Transducer control signal is at 0 Volts when released; when the Linear Transducer is pulled to the maximum, the control signal is equivalent to that of the maximum electrode signal from an Ottobock MyoBock Electrode. The advantages of the 9X52 Linear Transducer are:

- Control of the myoelectric prosthesis without a muscle signal or electrodes.

- DMC-control thorough proportional outgoing signal
- Minimal construction size
- Can be laminated into socket of upper arm or forearm

6.2 9X53 4-Step Transducer

Pulling the cable of the 9X53 4-Step Transducer creates a control signal. Depending on the step selected, signals are output at varying strengths. If the prosthesis has been properly configured, the 4-Step Transducer will switch between the various components. The 4-Step Transducer control signal is at 0 Volts when released.

Example: Switching the DynamicArm 12K100

| | |
|---------|----------------------------------|
| Step 1: | The Transducer releases |
| Step 2: | System Electric Hand |
| Step 3: | Electric wrist rotator mechanism |
| Step 4: | DynamicArm |

7 Contents of delivery (Fig. 1)

| | | | |
|-------------|-------|---|-------|
| Item 1 | 1 pc. | 9X52 Linear Transducer or 9X53 4-Step Transducer | |
| Item 2 | 1 pc. | Stirrup with Wedge Lock | 21Y79 |
| Item 3 | 1 pc. | Housing Assembly | 9X54 |
| Item 4 | 1 pc. | Lamination Pattern | 9X55 |
| no illustr. | 1 pc. | Instruction Manual | |

8 Installing the electrode cable (Fig. 2 to Fig. 5)

An electrode cable **is not** included in the Transducer delivery. It must be ordered separately and can be shortened to meet individual requirements. Electrode cables such as 13E129=G300, 13E129=G600 and 13E129=G1000 are compatible.

- 1) Liberally apply 633E11 Silicone Grease to the inner side of the Transducer housing around the contact pins. The grease insulates the connector against humidity (Fig. 2).
- 2) Guide the electrode cable (with the gray side facing upward) as far as it will go into the connector (Fig. 3). Remove any excessive grease.

- 3) With the help of a screwdriver or similar tool, press the connector with the inserted electrode cable as deeply as possible into the Transducer (Fig. 4). The connection is now active.

The linkage can be detached by pulling the electrode cable upwards at a right angle to the Transducer (Fig. 5). If the cable is to be pressed back into place again, cut off approx. 5 mm from the loose end of the electrode cable and repeat the procedure beginning with step 2.

9 Assembly of pull-cable and stirrup (Fig. 6, Fig. 7)

The Transducer is delivered with a pre-assembled pull-cable. When mounting the Transducer to a provisional socket or permanent prosthesis you must connect the free end of the pull-cable to the stirrup (Fig. 1, item 2).

- 1) Determine the necessary cable length, but do not yet shorten the cable.
- 2) Form a loop and feed the cable into the stirrup (Fig. 6).
- 3) Place the wedge into the stirrup and tighten the pull-cable (Fig. 7).

Then check the length of the cable and adjust if necessary. Remove wedge for corrections. When you have the correct length of cable, you may cut off the free end.

10 Lamination and assembly information (Fig. 8 to Fig. 15)

- 1) Transfer the position of the Transducer from the provisional socket to the permanent prosthesis (Fig. 8).
- 2) Form the foam core for the permanent outer socket with the help of the 743A23 Alignment Tool for ErgoArm following the 647H447 Instruction Manual (Fig. 9).
- 3) Affix the 9X55 Lamination Pattern on the predetermined position of the permanent prosthesis with a double-sided adhesive tape (Fig. 10).
- 4) Cover the socket with a 99B25 Nylon Stockinette first (Fig. 11-13), followed by a 99B81 PVA Bag.
- 5) Cover the socket with six layers of 623T3 Perlon Stockinette white (Fig. 14).
- 6) Cover the socket with the 99B81 PVA Bag.
- 7) Tint the 617H19 Orthocryl Lamination Resin with 617Z2 Pigment Paste and laminate the socket (Fig. 15).

Carefully follow the outline of the Lamination Pattern in order to prevent the Transducer from slipping out from between the outer socket and the Thermoplastic Inner Socket.

11 Examples

11.1 Transducer with provisional socket (Fig. 16, Fig. 16a)

The Transducer can be attached to the provisional socket during rehabilitation. Slide the Transducer into the 9X54 Housing Assembly (Fig. 1, item 3) and mount it in the provisional socket using a 503F3 Bolt. The Transducer is actuated by the patient's shoulder harness. Make sure that the pull-cable of the Transducer is aligned with the fitted harness and that there are no kinks or sharp angles present.

11.2 Transducer with permanent prosthesis (Fig. 17, Fig. 18)

If the patient is being fitted with a permanent prosthesis you can integrate the Transducer into the socket. The integration into the socket permits a clean cosmetic design. Drill a hole into the permanent outer socket to lead the pull-cable of the Transducer to the outside.

Attention!

Position and drill the hole so that the pull-cable, when actuated, does not come into contact with the permanent outer socket. The pull-cable of the Transducer can be fed through a cable housing to the shoulder harness. When positioning the cable housing, make sure there are no tight angles or edges and that the cable runs smoothly.

12 Technical data

| | |
|--------------------------------------|---|
| Cable Travel | 8 mm |
| Max. Operating Force | 10 N |
| Output Voltage | 0 – 2.5 V |
| Power Supply | 757B20/757B21 Ottobock EnergyPack MyoEnergy Integral 757B25=* 757B15 X-ChangePack |
| Service life | 5 years |
| Storage (with and without packaging) | +5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing |

| | |
|--|--|
| Transport (with and without packaging) | -20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing |
| Operation | -5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing |

13 Cleaning and care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Ensure that no liquid penetrates into the system component. Dry the product with a soft cloth afterwards.

14 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

15 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices. The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>.

16 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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