















8R1=L / R 7 1/4, 8R1=L / R 7 3/4, 8R1=L / R 8 1/4



DE Gebrauchsanweisung	7
EN Instructions for use	13
FR Instructions d'utilisation	19
IT Istruzioni per l'uso	25
ES Instrucciones de uso	31
NL Gebruiksaanwijzing	37
SV Bruksanvisning	43
PT Manual de utilização	48

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Date of the last update: 2020-11-30

- Please read this document carefully before using the product and observe the safety notices.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

1 Application

The 8R1=* Endoskeletal Adapter is to be used **solely** for the prosthetic fitting of the upper extremities in conjunction with the 8E44 Transcarpal-Hand.

Only start up the product in accordance with the information contained in the accompanying documents supplied.

Instruct the patient in the proper use and care of the product. Handover to the patient without prior instruction is not permitted.

2 Safety information

- To ensure a solid bond, please use acetone to remove all grease from any areas to which adhesive will be applied. Once you have cleaned the area, do not touch it, and immediately apply the adhesive!
- Follow all fabrication and safety considerations of the adhesive manufacturer.
- Please follow all recommend torque values.
- The 8R1=* Endoskeletal Adapters may only be removed and repaired by an authorized Ottobock technician.
- Use the Patient Information section of this pamphlet (Section 3) to inform your patient about the function and use of the Endoskeletal Adapter.

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3 Patient information

- The Endoskeletal Adapter has been developed for everyday life activities and must not be used for unusual activities such as extreme sports (free climbing, parachute jumping, paragliding, etc.).
- Careful handling of the prosthesis and its components not only increases their service life but, above all, serves the patient's safety!
- If the prosthesis has been subjected to unusual stresses (e.g. caused by a fall or similar), immediately contact a prosthetist and have the prosthesis examined for any damage. If necessary, the prosthetist will pass the prosthesis on to the Ottobock Myo-Service.
- Repair of any damaged components may only be carried out by the certified Ottobock Myo-Service.

Operating a motor vehicle

Whether or not the user of an arm prosthesis can safely operate a vehicle is a complex question. It depends on the type of amputation (amputation height, unilateral or bilateral, residual limb condition, type of prosthesis) as well as the individual capabilities of the prosthesis user.

All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center.

In general, Ottobock recommends that the vehicle to be used is outfitted with the necessary requirements (such as a steering fork) by a specialty outfitter. The user must be able to operate the vehicle with the electro-hand turned off. Operating the vehicle with the hand turned on can lead to steering problems, and can endanger not only the driver but also others nearby.

Caution!

Before operating a motor vehicle, the Quick Disconnect Unit of the 8R1=* must be so positioned that a slight turn of it, which can occur during steering, cannot disconnect the hand from the prosthesis!

4 Components (Fig. 1)

Pos. 1	1 Pc.	9S263=*	Plate
Pos. 2	1 Pc.	9S264	Threaded Connector with Tube Adapter
Pos. 3	1 Pc.	9S265	Threaded Connector
Pos. 4	1 Pc.	9S267=*	Retaining Ring
Pos. 5	2 Pc.	9E167	Cable Connector (1 Pc. Reserve)
Pos. 6	1 Pc.	9E388	Cable Seal
Pos. 7	1 Pc.	9S266	Quick Disconnect Unit
Pos. 8	1 Pc.	9E397=8	Coaxial Bushing, White
Pos. 9	1 Pc.		Coding Ring, Brown
wo. Fig.	1 Pc.	636W23	UHU plus endfest 300 Adhesive
wo. Fig.	1 Pc.	647H501	Instructions for Use

5 Description

The 8R1=* Endoskeletal Adapter is designed to work with the 8E44 Transcarpal-Hand, and can be used to fit both long and short residual limbs.

Because the Transcarpal-Hand is so light, and more of the weight is distributed proximally, patients feel noticeable relief on their residual limb.

The Ottobock 10S17 Electric Wrist Rotator can also be installed in the prosthesis, and the integrated quick disconnect allows for the pronation and supination of the prosthesis and a quick change to a Greifer (with Quick Disconnect). The 8R1=* Endoskeletal Adapter can be set up for use with either right or left sided fittings.

6 Endoskeletal adapter for fitting short residual limbs

Before finishing a short socket, please be aware of the limited amount of space inside the socket. Position the battery receptacle such that, if a 10S17 Electric Wrist Rotator is used, they cannot come in contact with each other.

6.1 Determining the tube length (Fig. 2 – Fig. 3)

1. To determine the tube length, separate the Transcarpal-Hand from the lamination plate (as shown in Instructions for Use 647H398). Attach the Endoskeletal Adapter to the 8E44 Transcarpal-Hand, without securing the threads and without bonding the Tube in place. For a temporary attachment of the Threaded Connector (Fig. 1, Pos. 2), to the Threaded Connector with Tube (Fig. 1, Pos. 3), secure the parts with adhesive tape. Pull the System Inner Hand over the Transcarpal-Hand and insert the assembly into the patient's socket.
2. Measure the patient's hand with a laser measuring tool (such as the Ottobock 743L20=230 LaserLine) to determine (Fig. 2) the required tube length (Fig. 1, Pos. 2)

Position the Lasermark on the tip of the thumb of the opposing hand and copy the Laser- mark onto the prosthesis, measuring the distance between the Lasermark on the System Inner Hand to the thumb tip of the 8E44 Transcarpal-Hand (Fig. 3).

The calculated measurement is the length to which the tube should be cut. Finally, remove the assembly from the patient's socket. To unscrew the Quick Disconnect Unit, use a pair of pliers with plastic jaws.

6.2 Assembly (Fig. 4 – Fig. 20)

Use the 719R2 Tube Cutter to cut the determined length from the Threaded Connector with Tube Adapter (Fig. 1, Pos. 2). Next, use the 718R1 Pipe Reamer to file the inside and outside of the tube. To guarantee a solid bond, clean the areas on both the Threaded Connector with Tube Adapter (Fig. 1, Pos. 2) and the Threaded Connector (Fig. 1, Pos. 3) with a cotton swab dipped in acetone to remove any grease. Do not touch the areas cleaned with acetone, and immediately coat with adhesive (Fig. 4–Fig 6).

Coat the degreased inside of the tube and the unanodized, degreased outside of the cylinder of the Threaded Connector with an approximately 10 mm wide band of adhesive (such as UHU-plus, endfest 300 (636W23) (Fig. 6).

Connect the two parts and promote the adhesive area with a maximum temperature of 70 °C. The components should harden in a vertical position to guarantee that the adhesive creates a consistent bond. Please follow the adhesive manufacturer's directions concerning varying cure times and how they may change depending on temperature.

The Wrist Connector Unit (Fig. 1, Pos. 7) is delivered with a coaxial bushing (Fig. 1, Pos. 8) and a white coding ring installed. The white coding ring specifies that Digital Twin control is in use.

If you are using DMC plus control, you must install a brown coding ring (Fig. 1, Pos. 9). To do this, remove the coaxial bushing and remove the white coding ring from the coaxial housing with a screwdriver or similar tool (Fig. 7).

Place the side of the brown coding ring with the raised letters face down on a hard surface and press the coaxial housing straight down on it with your thumb until you hear a 'click' (Fig. 8). The coding ring must completely lock-in.

Finally, clean all threads (except the two bolts in the plate) with acetone and then coat with Loctite 245 (636K12): Plate (Fig. 1, Pos. 1), Threaded Connector with Tube (Fig. 1, Pos. 2), Threaded Connector (Fig. 1, Pos. 3), and Wrist Connector Unit (Fig. 1, Pos. 7) (Fig. 9).

Feed the four-pole cable through the tube adapter until it comes out the other end (Fig. 10), then use pliers with plastic jaws to screw together the Quick Disconnect Unit (Fig. 1, Pos. 3) and the Threaded Connector (Fig. 1, Pos. 3) until they are hand tight (about 8 Nm).

Feed the Four-Pole Cable through the Plate (Fig. 1, Pos. 1), then screw the Threaded Connector with Tube and the Plate together until they are handtight (8 Nm). Make sure that the cable is not damaged during the assembly (Fig. 11).

The assembled and glued Endoskeletal Adapter (Fig. 12).

Slide the Cable Seal (Fig. 1, Pos.6) onto the Four-Pole Cable. Make sure that the cable in the tube is **not stretched tight**. Leave a length of cable in the tube for later repairs. Press the Cable Seal into the Plate. Besides its function as a seal, the Cable Seal acts as a strain-relief for the cable (Fig.13).

Mount the 8E44 Transcarpal-Hand to the Plate (Fig. 1, Pos. 1), but do not secure the nuts with Loctite yet (Fig. 14). Make sure that the cable can not be pinched.

To determine the appropriate length for the cable, temporarily connect it to the electronic circuit on the 8E44 Transcarpal-Hand (Fig. 15). Mark the necessary lengths.

Disconnect the 8E44 Transcarpal-Hand from the Plate (but do not disconnect Plate from the Threaded Connector) and cut the cable to length (Fig.16).

Remove approximately 2 mm of the cable housing and prepare the exposed wires with solder (Fig. 17).

Solder the four-pole cable to the Cable Connector (Fig. 1, Pos.5) (Fig. 18).

INFORMATION

Polarity of the cables: While three of the housing contacts have been spaced equally apart, another contact has been set apart from them with a greater distance. Solder the end gray wire of the Four-Pole Cable to this contact (Arrows, Fig. 18). Pay attention to the order of the wires and make sure they are connected to the correct contacts!

Mount the 8E44 Transcarpal-Hand to the Plate as shown in the Instructions for Use 647H398 (Fig. 1, Pos. 1). Secure the nuts with Loctite 245 (636K12). Make sure that the cable is connected to the correct contacts (Fig. 19).

Plug the connection piece into the electronic circuit (Fig. 19).

The Ottobock 757T10 MyoSimulator can be used to test this assembly. To do this read the Instructions for Use 647H452.

Pull the Inner Hand Shell over the Transcarpal-Hand and fix in place with the Retaining Ring (Fig. 1, Pos. 4).

Check the hand and the prosthesis of the patient with a Laser measuring device (such as the Ottobock 743L20=230 LaserLine) (Fig. 20).

For the following adjustment of the 8E44 Transcarpal-Hand and for handling and use by the patient please read the Instruction Manual 647H398.

7 Endoskeletal adapter for long residual limbs (Fig. 7 – 26)

7.1 Determining socket length (Fig. 21)

To determine the length of the socket, measure the patient's contralateral limb with a laser measuring tool (such as the Ottobock 743L20=230 LaserLine).

Then determine the distance from the elbow joint of the opposing arm to the tip of the thumb in millimeters (Fig. 21, Pos.A).

From this length, subtract 115 mm (the distance from the thumb tip of the System Inner Hand of the Transcarpal-Hand to the distal end of the socket). The resulting length is the necessary socket length (i.e., the distance of the elbow to the thumb tip (370 mm) minus 115 mm = 225 mm).

7.2 Assembly

The Wrist Connector Unit (Fig. 1, Pos. 7) is delivered with a coaxial bushing (Fig. 1, Pos.8) and a white coding ring, which specifies that Digital Twin control is being used.

If DMC plus control is being used, you must install a brown coding ring (Fig. 1, Pos.9).

To do this, disassemble the coaxial bushing and remove the white coding ring from the coaxial housing with a screwdriver or similar tool (Fig.7).

Place the side of the brown coding ring with the raised letters face down on a hard surface and press the coaxial bushing straight down on it with your thumb until you hear a ,click' (Fig.8). The coding ring must completely lock-in.

Finally, clean all threads (except the two bolts in the plate) with acetone and then coat with Loctite 245 (636K12): Quick Disconnect Unit (Fig. 1, Pos. 7) and Plate (Fig. 1, Pos. 1; Fig. 7) (Fig.9). Feed the Four-Pole Cable through the Plate and use pliers with plastic jaws to screw the two components together handtight (8 Nm) (Fig. 22). Make sure that cable is not damaged in this process.

Slide the Cable Seal (Fig. 1, Pos.6) over the Four-Pole Cable and press the Cable Seal into the Plate. Besides its function as a seal, the Cable Seal acts as a strain-relief for the cable. (Fig. 23).

Attach the Plate (Fig. 1, Pos. 1) to the 8E44 Transcarpal-Hand , but do not secure the nuts with Loctite.

Determine the appropriate cable length for connecting it to the electronic circuit on the Transcarpal-Hand (Fig. 24).

Disconnect the 8E44 Transcarpal-Hand from the Plate and cut the cable to length (Fig.16).

Remove approximately 2 mm of the cable housing and prepare the wire with solder (Fig. 17).

Solder the four-pole cable to the electronic circuit (Fig. 1, Pos.5) (Fig. 18).

INFORMATION

Polarity of the cables: While three of the housing contacts have been spaced equally apart, another contact has been set apart from them with a greater distance. Solder the end gray wire of the Four-Pole Cable to this contact (Arrows, Fig. 18). Pay attention to the order of the wires and make sure they are connected to the correct contacts!

Mount the 8E44 Transcarpal-Hand to the Plate (Fig. 1, Pos. 1) as shown in Instruction Manual 647H398. Secure the nuts with Loctite 245 (636K12). Make sure that the cable is connected to the correct contacts (Fig. 25).

Plug the connection piece into the electronic circuit. Fig. 25 shows the assembled Endosketetal Adapter with the 8E44 Transcarpal-Hand.

The Ottobock 757T10 MyoSimulator can be used to test the functionality of this assembly. To do this read the Instruction Manual 647H452.

Pull the System Inner Hand over the Transcarpal-Hand and fix in place with the Retaining Ring (Fig. 1, Pos. 4).

Check the prosthesis of the patient with a Laser measuring device (such as the Ottobock 743L20=230 LaserLine) (Fig. 26).

For the following adjustment of the 8E44 Transcarpal-Hand and for handling and use by the patient please read the Instruction Manual 647H398

8 Cleaning and care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Ensure that no liquid penetrates into the system component.

Then, dry the prosthesis component with a soft cloth.

9 Technical data

All data refer to the 8R1=* Endoskeletal Adapter, without attached 8E44 Transcarpal-Hand.

Operating Temperature: 0–70°C

Endoskeletal adapter for short residual limb fittings

Max. Length (w/out Transcarpal-Hand): 180 mm

Max. Weight: 140 g

Endoskeletal-adapter for long residual limb fittings

Length (w/out Transcarpal-Hand): 17 mm

Weight: 87 g

Service life of the component: 6 years

10 Symbols used

Medical device

11 Legal information

11.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

11.2 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address:
<http://www.ottobock.com/conformity>

11.3 Trademarks

All product names mentioned in this accompanying document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners. Should trademarks used in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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