



5R9, 4R60

DE Gebrauchsanweisung	3	CS Návod k použití.....	65
EN Instructions for use	7	RO Instrucțiuni de utilizare.....	69
FR Instructions d'utilisation.....	12	HR Upute za uporabu	74
IT Istruzioni per l'uso	17	SL Navodila za uporabo	79
ES Instrucciones de uso	22	SK Návod na používanie	83
PT Manual de utilização.....	27	BG Инструкция за употреба	88
NL Gebruiksaanwijzing.....	32	TR Kullanma talimatı	93
SV Bruksanvisning	37	EL Οδηγίες χρήσης	98
DA Brugsanvisning	41	RU Руководство по применению	103
NO Bruksanvisning	46	JA 取扱説明書.....	108
FI Käyttöohje.....	50	ZH 使用说明书.....	113
PL Instrukcja użytkownika.....	55	KO 사용 설명서.....	117
HU Használati utasítás	60		

1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2020-07-06

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 5R9 socket attachment block is bonded and laminated to the foamed distal end of the prosthetic socket. The 4R60 socket adapter serves to connect with the distal prosthetic components and is fastened to the 5R9 socket attachment block.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors for children. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Environmental conditions

Storage and transport

Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature range: -10 °C to +45 °C

Moisture: relative humidity: 20% to 90%, non-condensing

Unacceptable environmental conditions

Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water

Solids: dust, sand, highly hygroscopic particles (e. g. talcum)

2.2 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.3 Area of application

Approved for a body weight up to **max. 45 kg.**

Approved for a body height up to **max. 145 cm.**

Approved for a foot size up to **max. 21 cm.**

2.4 Service life

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum service life of 5 years.

3 Safety

3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

3.2 General safety instructions

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 8).

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION**Exceeding the service life**

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

⚠ CAUTION**Use under unallowable environmental conditions**

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION**Reuse on another patient**

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

⚠ CAUTION**Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.

- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

5R9			
see fig. 1, item	Quantity	Designation	Reference number
–	1	Instructions for use	–
1	1	Socket attachment block	5R9
2	1	Lamination dummy	4X8
–	4	Countersunk head screw	501S41=M5x1-6

4R60			
see fig. 1, item	Quantity	Designation	Reference number
–	1	Instructions for use	–
3	1	Socket adapter	4R60
–	4	Set screw	506G3=M6X10

5 Preparing the product for use**⚠ CAUTION****Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

► Observe the alignment and assembly instructions.

INFORMATION

The layout described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the layout.

5.1 Fabricating the prosthetic socket

First lamination process

- > **Required materials:** 99B81 PVA bag, 616G6 Dacron® felt, 623T3 Perlon stockinette, 616G12 carbon fibre cloth, 617H119 Orthocryl lamination resin 80:20 PRO
- 1) Pull a soaked PVA bag over the model.
 - 2) Place a layer of Dacron® felt over the entire model.
 - 3) Place a layer of Dacron® felt around the proximal edge of the model.
 - 4) Pull 2 layers of Perlon stockinette over the model.
 - 5) Place 3 layers of carbon fibre cloth (e.g. 15 cm x 15 cm) on the distal end of the model, offsetting the fibre alignment.
 - 6) Pull 2 layers of Perlon stockinette over the model.
 - 7) Pull a soaked PVA bag over the model.
 - 8) Complete the lamination process with lamination resin.
 - 9) Once the lamination resin has cured, remove the PVA bag.

Foaming and fastening the socket attachment block

- > **Materials required:** 617H32 Pedilen rigid foam 300, 617H21 Orthocryl sealing resin
- 1) Sand down the surplus lamination resin on the distal end of the prosthetic socket.
 - 2) Foam the distal end of the prosthetic socket using rigid foam.
 - 3) Trim the rigid foam as much as possible and sand it smooth in the desired position.
 - 4) Fasten the socket attachment block to the rigid foam with sealing resin.

- 5) Mill a transition between the socket attachment block and the rigid foam. Do not sand the shape or size of the adapter connection surface, which is determined by the socket attachment block.
- 6) Insert the lamination dummy in the threaded holes of the socket attachment block (see fig. 2).
- 7) Apply coloured sealing resin to the outside of the prosthetic socket and to the lamination dummy (see fig. 3).

Second lamination process

- > **Required materials:** 99B81 PVA bag, 623T3 Perlon stockinette, 616B1=50 carbon fibre cloth strap, 616G12 carbon fibre cloth, 616G15 woven carbon fibre stockinette, 617H119 Orthocryl lamination resin 80:20 PRO
- 1) Pull a Perlon stockinette over the model.
 - 2) Apply 2 strips of carbon fibre cloth strap in a cross over the socket attachment block up to the middle of the model.
 - 3) Place 2 layers of carbon fibre cloth (e.g. 15 cm x 15 cm) on the distal end of the model, offsetting the fibre alignment.
 - 4) Pull a Perlon stockinette over the model.
 - 5) Pull 2 layers of woven carbon fibre stockinette over the model. Do not allow the stockinette to protrude beyond the socket edge.
 - 6) Pull 2 layers of Perlon stockinette over the model.
 - 7) Pull a soaked PVA bag over the model.
 - 8) Complete the lamination process with lamination resin.
 - 9) Once the lamination resin has cured, remove the PVA bag.

5.2 Final assembly

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- Clean the threads before every installation.
- Apply the specified torque values.

- Follow the instructions regarding the length of the screws and about how to secure the screws.

Fastening the socket adapter

> **Materials required:** sanding tool, 710D20 torque wrench

- 1) **CAUTION! Do not sand down the laminate around the edge of the product.**

Sand down the laminate on the distal end of the prosthetic socket to the surface of the lamination dummy and remove the lamination dummy (see fig. 4).

- 2) **CAUTION! Do not use thread lock (e.g. Loctite®) to fasten the socket adapter. Chemical thread locks damage the laminate.**

Fasten the socket adapter to the socket attachment block using the countersunk head screws:

Screw in 2 countersunk head screws on the posterior side (**8 Nm**).

Screw in 2 countersunk head screws on the anterior side (**8 Nm**).

Connecting the pyramid adapter and pyramid receiver

The pyramid adapter is fixed with the set screws of the pyramid receiver.

> **Materials required:** 710D4 torque wrench, 636K13 Loctite®

- 1) **Trial fitting:** Insert the pyramid adapter into the pyramid receiver and tighten the set screws (**7 Nm**).
- 2) **Definitive mounting:** Insert the pyramid adapter into the pyramid receiver, secure the set screws with Loctite® and tighten them (**7 Nm**).
- 3) Replace any set screws that are protruding or are recessed too much with suitable set screws (see selection chart).

Selection chart for set screws

6 mm	12 mm	16 mm
506G3=M6X6	506G3=M6X12	506G3=M6X16

Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

► Observe the following maintenance instructions.

- A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- Have the prosthetic components inspected every 3 months.
- Inspect the entire prosthesis for wear during normal consultations.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Reference number	5R9	4R60
Material	Plastic	Aluminium
System height [mm]	30	33
Weight [g]	125	45
Max. body weight [kg]	45	
Max. body height [cm]	145	
Max. foot size [cm]	21	

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