## 4R63, 4R63=5.1, 4R68, 4R100

<b>DE</b> Gebrauchsanweisung .....	3	<b>CS</b> Návod k použití.....	56
<b>EN</b> Instructions for use .....	7	<b>RO</b> Instrucțiuni de utilizare.....	60
<b>FR</b> Instructions d'utilisation.....	11	<b>HR</b> Upute za uporabu .....	64
<b>IT</b> Istruzioni per l'uso .....	15	<b>SL</b> Navodila za uporabo .....	68
<b>ES</b> Instrucciones de uso .....	19	<b>SK</b> Návod na používanie .....	71
<b>PT</b> Manual de utilização.....	24	<b>BG</b> Инструкция за употреба .....	76
<b>NL</b> Gebruiksaanwijzing.....	28	<b>TR</b> Kullanma talimatı .....	80
<b>SV</b> Bruksanvisning .....	32	<b>EL</b> Οδηγίες χρήσης .....	84
<b>DA</b> Brugsanvisning .....	36	<b>RU</b> Руководство по применению .....	88
<b>NO</b> Bruksanvisning .....	40	<b>JA</b> 取扱説明書.....	93
<b>FI</b> Käyttöohje.....	44	<b>ZH</b> 使用说明书.....	97
<b>PL</b> Instrukcja użytkowania.....	47	<b>KO</b> 사용 설명서.....	100
<b>HU</b> Használati utasítás .....	52		

1

# 1 Produktbeschreibung

Deutsch

## INFORMATION

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### INFORMATION

Date of last update: 2021-02-16

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The 4R63, 4R63=5.1, 4R68 and 4R100 lamination anchors are laminated into a TT prosthetic socket. They provide the connection with the distal prosthetic components and are equipped with a pyramid. The 4R63, 4R68 and 4R100 lamination anchors have a bore hole in the pyramid for inserting the pin of a shuttle lock. The pyramid of the 4R63=5.1 does not have a bore hole.

### 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 2.2 Area of application

The product must be used in TT prostheses only.

- **4R68, 4R100:** Approved for a body weight of up to **max. 100 kg**.
- **4R63, 4R63=5.1:** Approved for a body weight of up to **max. 136 kg**.

### 2.3 Environmental conditions

#### Storage and transport

Temperature range  $-20\text{ °C}$  to  $+60\text{ °C}$  ( $-4\text{ °F}$  to  $+140\text{ °F}$ ), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

#### 4R63, 4R63=5.1, 4R68

#### Allowable environmental conditions

**Temperature range:**  $-10\text{ °C}$  to  $+45\text{ °C}$

**Moisture:** relative humidity: 20% to 90%, non-condensing

**Chemicals/liquids:** fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)

**Solids:** dust

#### Prohibited environmental conditions

**Chemicals/moisture:** salt water, perspiration, urine, acids, soapsuds, chlorine water

**Solids:** dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

## 4R100

### Allowable environmental conditions

**Temperature range:** -10 °C to +45 °C (14 °F to 113 °F)

**Chemicals/liquids:** fresh water, salt water, perspiration, urine, soapsuds, chlorine water

**Moisture:** Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions

**Solids:** Dust, occasional contact with sand

**Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage** (see page 10).

### Unallowable environmental conditions

**Solids:** Highly hygroscopic particles (e.g. talcum), continuous contact with sand

**Chemicals/liquids:** Acids, continuous use in liquid media

## 2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION!**

### Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 7).
- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Observe the maximum lifetime of the product.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

### Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

## 4 Scope of delivery

4R63, 4R63=5.1, 4R68, 4R100 Lamination anchor				
Fig.	Item	Quantity	Designation	Reference number
-	-	1	Instructions for use	-
1	-	1	Lamination anchor	-
1	2	-	Lamination dummy	4X3
The following parts are only for the 4R63, 4R63=5.1 and 4R100:				
1	1	1	Lamination dummy	4X52

## 5 Preparing the product for use

### ⚠ CAUTION

#### **Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### ⚠ CAUTION

#### **Improper assembly of the screw connections**

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

### 5.1 Laminating the transtibial socket

> **Required materials:** 99B81\* PVA bag, 623T3\* Perlon stockinette, 616B1\* carbon fibre cloth strap, 699B1\* fibreglass roving, 616G12\* carbon fibre cloth, 616G15\* carbon fibre woven stockinette, 617H119\* Orthocryl lamination resin 80:20 PRO

- 1) Pull a soaked PVA bag over the model. The PVA bag should be sealed by heat welding at the distal end.
- 2) Pull two layers of Perlon stockinette over the model.
- 3) Place three layers of woven carbon fibre (e.g. **15 cm x 15 cm**) with offset fibre alignment on the distal end of the model.
- 4) Pull two layers of Perlon stockinette over the model.
- 5) **CAUTION! Do not bend the anchor arms of the lamination anchor.**  
Move the lamination anchor arms to the correct **a-p** and **m-l** position.
- 6) Set the 4X3 lamination dummy onto the pyramid of the lamination anchor.

- 7) **For 4R63\*, 4R100:** Place the 4X52 lamination dummy on the proximal side of the lamination anchor. The lamination dummy remains in the prosthetic socket.
- 8) Place the lamination anchor at the socket end as required for proper alignment and bring it into position. Use the marking on the 4X3 lamination dummy to do so.
- 9) Stabilise and fix the lamination anchor with a sealing resin/talcum powder mixture.
- 10) Pull the fibreglass roving through the holes in the anchor arms and allow it to hang in a loop. If there are no holes, allow the fibreglass roving to hang around the anchor arms in a loop.
- 11) Place three layers of carbon fibre cloth (e.g. **15 cm x 15 cm**) in an offset arrangement over the lamination anchor arms.
- 12) Pull one layer of Perlon stockinette over the model.
- 13) Pull two layers of carbon woven stockinette over the model.
- 14) Pull two layers of Perlon stockinette over the model and tie off in the groove of the lamination anchor.
- 15) Pull a soaked PVA bag over the model.
- 16) Complete the lamination process with Orthocryl lamination resin.
- 17) After the lamination resin has cured, remove the 4X3 lamination dummy.

### 5.2 Assembling the distal prosthetic components

#### **Connecting the pyramid adapter and pyramid receiver**

The pyramid adapter is fixed with the set screws of the pyramid receiver.

- > **Required materials:** torque wrench (e.g. 710D20), 636K13 Loc-tite 241
- 1) **Fitting:**  
Screw in the set screws.  
Use the torque wrench to tighten the set screws (**10 Nm**).

- 2) **Definitive assembly:**  
Secure the set screws with Loctite.  
Screw in the set screws.  
Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten them (**15 Nm**).
- 3) Replace any set screws that are protruding or are recessed too much with appropriate set screws (see selection chart).

Selection table for set screws	
Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

### Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

### Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

### 6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

### 7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

### 8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

### 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

#### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

#### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

#### 9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

### 10 Technical data

Reference number	4R63	4R63=5.1	4R68	4R100
Weight [g]	85	95	70	40
System height [mm]	-7			
Build height [mm]	11			

Reference number	4R63	4R63=5.1	4R68	4R100
Material	Stainless steel		Aluminium	Titanium
Max. body weight [kg]	136		100	

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