6A30=20N

DE Gebrauchsanweisung	5	CS Návod k použití.....	60
EN Instructions for use	9	RO Instrucțiuni de utilizare.....	64
FR Instructions d'utilisation.....	13	HR Upute za uporabu	68
IT Istruzioni per l'uso	18	SL Navodila za uporabo	72
ES Instrucciones de uso	22	SK Návod na používanie	76
PT Manual de utilização.....	26	BG Инструкция за употреба	81
NL Gebruiksaanwijzing.....	31	TR Kullanma talimatı	85
SV Bruksanvisning	35	EL Οδηγίες χρήσης	89
DA Brugsanvisning	39	RU Руководство по применению	94
NO Bruksanvisning	43	JA 取扱説明書.....	99
FI Käyttöohje.....	47	ZH 使用说明书.....	103
PL Instrukcja użytkowania.....	51	KO 사용 설명서.....	106
HU Használati utasítás	55		

1


2





3

4

1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2020-12-22

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 6A30=20N shuttle lock serves to secure a suitable liner in the prosthetic socket. The shuttle lock is laminated distally into the prosthetic socket.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The product is embedded in the prosthetic socket and is no structural component. The maximum body weight is determined by the prosthetic socket's laminate and the prosthetic components used.

For use with 5R2

Approved for a body weight of up to **125 kg**.

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60 °C

Storage and transport conditions: -10° C to +60° C, relative humidity: no restrictions

Allowable environmental conditions

Humidity: fresh water, salt water – cleaning and drying is required after contact and use.

Allowable immersion depth: 1 m

Unallowable environmental conditions

Sauna, extreme watersports

Mechanical vibrations or impacts

Perspiration, urine, acids or bases

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Lifetime

Depending on the patient's level of activity, the service life of the product is 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 Warning regarding possible risks of accident or injury.

3.2 General safety instructions

Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 9).
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.

- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

6A30=20N Shuttle Lock		
Quantity	Designation	Reference number
1	Instructions for use	–
1	Shuttle lock housing	5X120
1	Ratchet unit	6A52=K
1	Pin	6Y13=1
1	Set screw	506G21=M4x-10
1	Dummy set Consisting of: 1 plate dummy, 1 ratchet unit dummy, 1 countersunk head screw, 4 nails	5X108

6A30=20N Shuttle Lock		
Quantity	Designation	Reference number
1	Dummy set Consisting of: 1 cap screw, 1 sleeve dummy, 2 O-rings	6A94

Spare parts/accessories (not included in the scope of delivery)		
Quantity	Designation	Reference number
–	Pin, short	6Y13=2
–	Pin, long	6A13=L1

5 Preparing the product for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

5.1 Preparing for socket fabrication

> **Required materials:** Plate dummy, ratchet unit dummy, shuttle lock housing, silicone grease, Plastaband 636K8*

- 1) **NOTICE! Align the shuttle lock with the longitudinal axis of the residual limb to prevent unnecessary load and wear (see fig. 4).**

Position the plate dummy on the distal surface of the plaster model along the longitudinal axis of the residual limb and secure it with 4 nails (see fig. 5).

- 2) Grease the thread of the ratchet unit dummy and screw the dummy into the shuttle lock housing.
- 3) Fill the slot in the head of the ratchet unit dummy with Plastaband.

5.2 Fabricating the prosthetic socket

Use without lamination disc

> **Required materials:** Sleeve dummy, cap screw from 6A94 dummy set

- 1) Insert the sleeve dummy into the shuttle lock housing.
- 2) Secure the shuttle lock housing and sleeve dummy to the plate dummy with the cap screw.
- 3) Fabricate the prosthetic socket as described in the instructions for use of the socket attachment block or the lamination anchor.
- 4) Perform final assembly (see page 12).

INFORMATION

The layup described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the layup.

Use with lamination disc

> **Required materials:** 5R2 lamination disc, 636K8* Plastaband, 501S74=3.5x9.5 fillister head self-tapping screws, countersunk head screw from 5X108 dummy set, silicone grease, 99B81 PVA bag, 623T9 Nylglass stockinette, 616G4 fibreglass mat, 616G12 carbon fibre cloth, 616G14* woven carbon fibreglass stockinette, 616G3 fibreglass woven hose, 617H119 Orthocryl lamination resin 80:20 PRO

- 1) Fill the 4 screw holes on the underside of the lamination disc with Plastaband (see fig. 6).
- 2) Position the lamination disc on the shuttle lock so that the marking hole is aligned with the ratchet unit (see fig. 7).
- 3) Secure the lamination disc on the shuttle lock housing with 6 fillister head self-tapping screws (see fig. 8).
- 4) Secure the shuttle lock housing and lamination disc on the plate dummy with the countersunk head screw.
- 5) Apply grease to the contact surface of the lamination dummy and put on the lamination dummy (see fig. 9, see fig. 10).

- 6) Pull a soaked PVA bag over the plaster model (see fig. 11).
- 7) Cut a piece of Nylglass stockinette (2 times the length of the plaster model). Pull the Nylglass stockinette halfway over the plaster model.
- 8) Apply a strip of carbon fibre cloth (approx. **8 cm** wide) in a circle on the lamination disc.
- 9) Fold the second half of the Nylglass stockinette over the plaster model and tie it off with perlon stockinette.
- 10) Lay one layer of carbon fibre cloth mediolaterally over the plaster model, up to approximately **2 cm** distally from the socket brim.
- 11) Attach a strip of carbon fibre cloth (approx. **8 cm** wide) in a circle distally from the MPT line.
- 12) Trim a layer of woven carbon fibreglass stockinette (0.3 times the length of the plaster model). Pull the woven carbon fibreglass stockinette over the plaster model in the area of the shuttle lock.
- 13) Cut a piece of Nylglass stockinette (2 times the length of the plaster model). Pull the Nylglass stockinette halfway over the plaster model.
- 14) Tie off the Nylglass stockinette in the groove of the lamination disc. Fold the excess Nylglass stockinette over the plaster model.
- 15) Cut a piece of fibreglass stockinette (1.5 times the length of the plaster model). Pull the fibreglass stockinette over the distal half of the plaster model.
- 16) Fold the second half of the fibreglass stockinette over the plaster model.
- 17) Reinforce the distal area of the socket around the lamination disc and shuttle lock with carbon fibre cloth.
- 18) Pull 2 layers of Nylglass stockinette over the plaster model.
- 19) Pull a soaked PVA bag over the plaster model.
- 20) Complete the lamination process with Orthocryl.
- 21) Perform final assembly (see page 12)

5.3 Final assembly

- > **Prerequisite:** The prosthetic socket has been finished.
 - > **Required materials:** 710D4 torque wrench, torque wrench (adjustable to 0.5 Nm), 506G21=M4x10 set screw
- 1) Grind away the ratchet unit dummy and remove the dummy (see fig. 12).
 - 2) Screw the ratchet unit into the shuttle lock (torque value: **10 Nm**).
 - 3) Screw the set screw into the shuttle lock on the inside of the prosthetic socket (torque value: **0.5 Nm**). This secures the ratchet unit (see fig. 13).

6 Use

CAUTION

Pin not fully inserted into the lock

Risk of injury due to loss of connection to the prosthesis

- ▶ Guide the pin into the lock in accordance with the instructions.
- ▶ Always verify that the pin is engaged in the lock before using the prosthesis.

The liner and shuttle lock are connected to each other with the pin. The pin engages in the shuttle lock and holds the liner in place.

Putting on the Liner

- 1) Check the pin and liner for damage.
- 2) **NOTICE! Align the pin with the longitudinal axis of the residual limb to prevent damage.**
Roll up the liner and place it on the end of the residual limb.
- 3) Unroll the liner onto the residual limb so that it has no wrinkles, without shifting the soft tissue.
- 4) Check the fit and the alignment of the liner.

Donning the Prosthetic Socket

- 1) Slide into the prosthetic socket with the liner until the pin slides into the opening of the shuttle lock.
- 2) Insert the pin fully into the shuttle lock.

- 3) Before using the prosthesis, verify that the pin is fully engaged in the shuttle lock.

Doffing the Prosthetic Socket

- ▶ Press and hold the side button and pull the residual limb with the liner out of the prosthetic socket.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9 Technical data

Reference number	6A30=20N
Weight [g]	125
Overall height [mm]	42
Material	Plastic

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