
## 17Y32, 17Y35

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	7
<b>FR</b> Instructions d'utilisation .....	11
<b>IT</b> Istruzioni per l'uso .....	16
<b>ES</b> Instrucciones de uso .....	20
<b>PT</b> Manual de utilização .....	25
<b>NL</b> Gebruiksaanwijzing .....	29
<b>SV</b> Bruksanvisning .....	34
<b>DA</b> Brugsanvisning .....	38
<b>NO</b> Bruksanvisning .....	42
<b>PL</b> Instrukcja użytkowania .....	46
<b>HU</b> Használati utasítás .....	51
<b>CS</b> Návod k použití .....	55
<b>TR</b> Kullanma talimatı .....	59
<b>EL</b> Οδηγίες χρήσης .....	64
<b>RU</b> Руководство по применению .....	69

1

## INFORMATION

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### INFORMATION

Date of last update: 2021-05-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the processing of the 17Y32 and 17Y35 system extension assist.

## 2 Product description

### 2.1 Available sizes

System extension assist	
Reference number	For orthotic knee joint
17Y32	17B26*, 17B3
17Y35	17B26*, 17B3

### 2.2 Components

Scope of delivery (see fig. 1)			
Item	Quantity [pcs.]	Designation	Reference number
1	2	Flat head screw	501S43=M3x4 ( <b>17Y32, 17Y35</b> )
2	2	Oval countersunk screw	501S27=M3x4 ( <b>17Y32, 17Y35</b> )
3	2	Rubber ring	18Z1=20.8x14.5 ( <b>17Y32</b> ) 18Z1=24.6x18.3 ( <b>17Y32</b> ) 18Z2=21.0x12.9 ( <b>17Y35</b> ) 18Z2=26.8x17.7 ( <b>17Y35</b> )
4	2	Special screw	501S79=M5x4.1x16.1 ( <b>17Y32, 17Y35</b> )
5	1	Serrated lock washer	507S11=5.3 ( <b>17Y32, 17Y35</b> )

Not included in scope of delivery			
Item	Quantity [pcs.]	Designation	Reference number
6	1	Orthotic knee joint, lower section	17B26*, 17B3*
7	1	Orthotic knee joint, upper section	17B26*, 17B3*

### 2.3 Function

The 17Y32 and 17Y35 system extension assists are suitable for leg prostheses with an extension assist function, and are unilaterally or bilaterally mounted on the 17B26\* or 17B3\* orthotic knee joints. The stronger rubber rings of the 17Y35 result in a more powerful extension assist function.

## 3 Intended use

### 3.1 Indications for use

The system extension assists are to be used **exclusively** for the unilateral or bilateral orthotic treatment of the lower limbs on one patient.

### 3.2 Indications

- Partial or total paralysis of the leg muscles
- Indications must be determined by the physician.

## 4 Safety

### 4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 Safety instructions

 **CAUTION**

#### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

#### **Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

 **CAUTION**

#### **Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

### CAUTION

#### **Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

### CAUTION

#### **Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### NOTICE

#### **Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

## **5 Preparing the product for use**

- 1) Flex the orthotic knee joint.

**INFORMATION: When the orthotic knee joint is flexed, the rubber rings of the system extension assist are tension free.**

- 2) Slide the serrated lock washer onto the special screw.
- 3) Install the flat head screw into the outer threaded bores of the orthotic knee joint upper section and lower section.
- 4) Secure the screw connection with Loctite 241.

## **6 Cleaning**

Promptly clean the product after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the product with clear fresh water.

- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e. g. from an oven or radiator).

## 7 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.

## 8 Disposal

Dispose of the product in accordance with national regulations.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medicinal devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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