C-Leg compact Prosthetic System

DE	Gebrauchsanweisung (Fachpersonal)	2
EN	Instructions for use (qualified personnel)	23
FR	Instructions d'utilisation (Personnel spécialisé)	44
ES	Instrucciones de uso (Personal técnico especializado)	66

Adapter	 <p>4R104=60 4R104=75</p> <p>4R57 4R57=ST</p> <p>4R41 4R43 4R89</p> <p>4R111=N 4R111 4R116</p> <p>4R40</p> <p>4R118</p>	 <p>757L16-*</p> <p>4E50-*</p>	Lade- / Netzgerät
C-Leg compact Kniegelenk	<p>mit Justierkern</p>  <p>3C96-1</p>  <p>4X250</p> <p>mit Gewindeanschluss</p>  <p>3C86-1</p>  <p>4X250</p>	 <p>4X160=1.2</p>  <p>4X160=5.6</p> <p>3S26 (o. Abb)</p> <p>3R59 (o. Abb)</p>	Kosmetik / Protector
Rohradapter	 <p>2R80</p>  <p>2R81</p>		
Prothesenfuß	 <p>1C40</p> <p>1E56</p> <p>1E57</p> <p>1C30</p> <p>1D35</p> <p>1A30</p> <p>1D10</p> <p>1M10</p>	 <p>4X180=*</p>  <p>60X3</p> <p>60X5</p>	C-Soft / Bionclink
	 <p>4X78/ 4X79</p> <p>4X74</p>  <p>4X77=D</p>	<p>4X83=430-USB</p>	Zubehör

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Deutschland: (+49) 55 27-848 13 33

Österreich: (+43) 1-5 23 37 86 430

Adapter	 <p>4R104=60 4R104=75</p> <p>4R57 4R57=ST</p> <p>4R41 4R43 4R89</p> <p>4R111=N 4R111 4R116</p> <p>4R40</p> <p>4R118</p>	 <p>757L16-*</p> <p>4E50-*</p>	Charger / AC adapter
C-Leg compact Knee Joint	<p>with pyramid adapter</p>  <p>3C96-1</p> <p>4X250</p> <p>with threaded connector</p>  <p>3C86-1</p> <p>4X250</p>	 <p>4X160=1.2</p>  <p>4X160=5.6</p> <p>3S26 (not illustrated)</p> <p>3R59 (not illustrated)</p>	Cosmetic Cover / Protector
Tube Adapter	 <p>2R80</p>  <p>2R81</p>		
Prosthetic Foot	 <p>1C40</p> <p>1E56</p> <p>1E57</p> <p>1C30</p> <p>1D35</p> <p>1A30</p> <p>1D10</p> <p>1M10</p>	 <p>4X180=*</p> <p>60X3</p> <p>60X5</p>	C-Soft / BionicLink
	 <p>4X78/ 4X79</p> <p>4X74</p>  <p>4X77=GB</p>	<p>4X83=430-USB</p>	Accessories

Table of Contents

1 Foreword	25
Explanation of symbols	25
Recommended accessories for prosthetists	25
2 Important C-Leg compact information – Intended use	25
2.1 Medical purpose	25
2.2 Field of application	25
2.3 Conditions of use	26
2.4 Functionality	26
2.5 Qualification of the prosthetist	26
2.6 Triton combination possibilities	26
2.7 Safety instructions	27
3 Alignment and Adjustment	31
3.1 Connecting the Tube Adapter	31
3.2 Bench Alignment Using an Alignment Tool such as L.A.S.A.R. Assembly	32
3.2.1 Checking the Socket after Bench Alignment	33
3.3 Torque Values of the Screw Connections	34
3.4 Charging the C-Leg compact	34
3.5 Static Alignment Optimization (with 743L100 L.A.S.A.R. Posture)	35
3.6 BionicLink / optional: Serial Data Cable	36
3.7 C-Soft	36
3.8 Optional: Parameter Settings with Foam Cover and C-Soft	36
3.9 Remote Control	37
3.9.1 Modes and Use of the Remote Control	37
3.9.2 Attachment of the Remote Control	38
3.9.3 Pairing (Establishing the Communication between the C-Leg compact and Remote Control)	39
3.9.4 Changing the Battery of the Remote Control	39
3.9.5 Switching between Compact Mode and Therapy Mode	39
3.10 Finishing the Prosthesis	39
3.11 Important User Instructions	40
4 Additional Information	41
4.1 Service Intervals	41
4.2 Damage Event	41
4.3 Technical Information	41
4.3.1 Symbols on the Knee Joint	41
4.3.2 Symbols on remote control	42
4.3.3 Symbols on the battery charger	42
4.4 Transportation	42
4.5 Guaranty	42
4.6 Liability	42
4.7 Trademarks	42
4.8 CE Conformity	43

INFORMATION

Last update: 23.03.2015

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

1 Foreword

Only put the product into use in accordance with the information contained in the accompanying documents supplied. Instruct the user in the proper and safe use of the product.

Explanation of symbols

▲ DANGER Warnings regarding directly impending risks of severe accident or injury.

▲ WARNING Warnings regarding possible risks of severe accident or injury.

▲ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

Recommended accessories for prosthetists

710D4	Torque Wrench	
710Y2=5	Hex Bit	
743L100	L.A.S.A.R Posture	(647H189 Instructions for Use)
743L200 / 743L300	L.A.S.A.R Assembly	(647H193 Instructions for Use)
60X3	BionicLink	(647G192 Instructions for Use)
4X180=*	C-Soft	(647G268 Instructions for Use)

2 Important C-Leg compact information – Intended use

INFORMATION

Before using the C-Leg compact, please read these Instructions for Use thoroughly! Please pay special attention to the safety instructions!

The patient must be taught how to handle, care for and operate his/her prosthesis properly. Please see the following Sections: 1.3 Conditions of Use; 1.6 Safety Instructions; 2.4 Charging the C-Leg compact; 2.11 Important User Instructions; 3.1 Service Intervals; 3.3 Technical Information; 3.6 Liability.

2.1 Medical purpose

The **3C96-1** and **3C86-1** C-Leg compact knee joints are to be used **exclusively** for the exoprosthetic fitting of amputations of the lower limb.

2.2 Field of application

The C-Leg compact was designed for a broad field of application and can be used by transfemoral amputees as well as by hip disarticulation amputees. For knee disarticulation amputees, a version for long residual limbs is available (3C86-1).

Field of application according to the **Ottobock MOBIS® Mobility System:**


Recommended for **Mobility Grades 2 and 3**

(restricted outdoor walkers and unrestricted outdoor walkers).

A therapy mode enlarges the field of application and allows targeted use of the C-Leg compact in rehabilitation.

Approved for a patient weight of up to **125 kg / 275 lbs.**

Exception: 2R80=110 Tube Adapter is approved for a patient weight of up to 100 kg / 220 lbs.

2.3 Conditions of use

INFORMATION

Advise your patients of the information in this section.

The C-Leg compact was developed for everyday activities and not for extreme sports such as free climbing, parachuting or paragliding. For the necessary environmental conditions, please see Section 3.3 “Technical Information”.

The C-Leg compact knee joint is designed **exclusively** for use on one patient. Use of the product by other persons is not approved by the manufacturer.

2.4 Functionality

The 3C96-1/3C86-1 C-Leg compact is a monocentric knee joint with microprocessor-controlled stance phase and hydraulically controlled swing phase.

Strain gauges in the tube adapter measure the anterior and posterior flexion moment, while a knee angle sensor measures the flexion angle and angular speed of the knee joint. These measured variables are transmitted to the microprocessor, which then calculates the necessary movement resistances for the stance phase. Servomotors open and close hydraulic valves to provide the required damping during the stance and swing phases. This allows the C-Leg compact to adapt to the individual requirements and the activities of the prosthesis wearer. Different modes also allow, in addition to everyday use, for use in connection with early mobilisation following an amputation.

2.5 Qualification of the prosthetist

The fitting of a patient with a C-Leg compact may only be carried out by prosthetists who have been authorized with the corresponding Ottobock training.

2.6 Triton combination possibilities

The combination of the prosthetic foot Triton 1C60 / Triton Vertical Shock 1C61/ Triton Harmony 1C62/ Triton Low Profile 1C63/Triton Heavy Duty 1C64 with C-Leg compact 3C96-1/ 3C86-1 is only allowed in the following ranges (see chart below).

For the combination of the Triton prosthetic feet with or C-Leg compact 3C96-1/ 3C86-1 outside of the approved ranges please contact Ottobock Customer Service.

C-Leg compact 3C96-1/ 3C86-1																								
Tube Adapter		Body Weight		Foot Size																				
				21 cm		22 cm		23 cm		24 cm		25 cm		26 cm		27 cm		28 cm		29 cm		30 cm		
kg	lbs	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	
2R82=110	up to 100	up to 220																						
2R82=120	2R82=/2R81=160	up to 100	up to 220																					
2R82=/2R81=200	2R82=/2R81=240	101 – 125	221 – 275																					

2.7 Safety instructions

INFORMATION

Advise your patients of the information in this section.

⚠ CAUTION

Non-observance of safety instructions. Failure to follow the below-mentioned safety instructions can lead to a faulty control or malfunction of the C-Leg compact and result in risk of injuries for the patient as well as damage to the C-Leg compact.

General safety instructions

⚠ CAUTION

Alignment and adjustment error. During alignment and adjustment of the prosthesis, errors can occur resulting in malfunction of the joint up to loss of function due to structural failure. This can cause the patient to fall.

- Participation in an Ottobock product training course for the C-Leg compact is obligatory before the first fitting. Additional product training courses may become necessary to qualify for fitting product updates.
- The tubes must not be shortened as that may damage the integrated strain gauges.
- During the data transfer (PC to C-Leg compact), the patient must remain standing or sitting still and the Bionic-Link or the communication cable must not be removed.
- If the prosthesis wearer uses crutches or walking canes during adjustment, readjustment is required as soon as he stops using these walking aids.

⚠ CAUTION

Use of inappropriate prosthetic components. If inappropriate prosthetic components are installed in the prosthesis, malfunction of the joint can occur up to loss of function due to structural failure. This can cause the patient to fall.

Combine the C-Leg compact only with adapters and feet approved by Ottobock (see component overview on the folded front page).

⚠ CAUTION

Manipulations on system components. Any changes or modifications you make to system components on your own initiative can lead to malfunction of the joint up to loss of function due to structural failure. This can cause the patient to fall.

- Any changes or modifications to the device may limit its use.
- The opening and repairing of the joint may only be performed by authorized Ottobock technicians, and the handling of the battery may only be carried out by Ottobock Service Centres (exchanges are not permissible).

⚠ CAUTION

Incorrect battery charging. Charging the battery with battery chargers that have not been approved by Ottobock can lead to defects and result in malfunction of the joint. This can cause the patient to fall.

Advise the patient about the charging procedure and refer to the Patient Information. Advise the patient of the following patient instructions.

⚠ CAUTION

Improper use of the joint. Any kind of overloading or excessive strain as well as improper use can lead to defects and result in malfunction of the joint up to loss of function due to structural failure. This can cause the patient to fall.

Advise the patient of the proper use of the C-Leg compact as well as of the following patient instructions.

⚠ CAUTION

Improper use of the remote control. The remote control can get damaged by improper use. This can lead to malfunction of the remote control and result in unexpected actions of the joint. This can cause the patient to fall.

Advise the patient about the proper use of the remote control and refer to the Patient Information. Advise the patient of the following patient instructions.

⚠ CAUTION

Transport damage. Mechanical impact or stress during transportation of the joint such as shocks and vibrations can lead to

- defects and result in malfunction of the joint;
- defects on the battery and hydraulic damper resulting in leakage of liquid; or
- loss of function due to structural failure.

This can cause the patient to fall as well as result in skin irritation.

Always use the transport packaging for transport.

⚠ CAUTION

Results of product deterioration. Wear and tear on system components can lead to malfunction of the joint. This can cause the patient to fall.

In the interest of the prosthesis wearer (maintenance of operational safety and guaranty), the specified service intervals must be complied with.

⚠ CAUTION

Using the tube adapter without calibration with 4X180=* C-Soft. If a tube adapter is used without subsequent calibration of the joint with 4X180* C-Soft or Slider software, the damping behaviour of the joint cannot be properly controlled. This can cause the patient to fall in certain cases.

When using a new tube adapter, calibration (zero setting) must be carried out with 4X180=* C-Soft or Slider software.

Patient instructions

⚠ WARNING

Risk of accident when driving a motor vehicle. The ability of leg prosthesis users to drive a motor vehicle is determined on a case-by-case basis. Criteria include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's individual abilities. All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test centre. Ottobock recommend that the motor vehicle be professionally retrofitted to the user's individual needs (e.g. automatic shift). Risk-free driving must be ensured even when the prosthesis is not functioning.

⚠ CAUTION

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WLAN devices)

Unexpected behaviour of the prosthesis due to interference with internal data communication. This can cause a fall in certain situations.

Therefore, maintain the following minimum distance from these HF communication devices:

- Mobile phone GSM 850/GSM 900: 0.99m
- Mobile phone GSM 1800/GSM 1900/UMTS: 0.7m
- DECT cordless phones incl. base station: 0.35m
- WLAN (routers, access points,...): 0.22m
- Bluetooth devices (third-party products not approved by Ottobock): 0.22m

⚠ CAUTION

Magnetic interference fields. Wearing the joint in the vicinity of high-voltage lines, transmitters, transformers, computer tomographs or other sources of strong electromagnetic radiation (e.g. retail security systems in department stores or body scanners at airports) can result in malfunctions of the joint. This can cause the patient to fall.

Avoid remaining in the vicinity of sources of strong magnetic and electrical interference (e.g. transformer stations, transmitters). If remaining in these areas cannot be avoided, make sure to at least walk/stand in a safe manner (e.g. with a handrail or support from another person).

⚠ CAUTION

Thermal overloading. Extended exposure to high temperatures can lead to defects and result in malfunction of the joint up to loss of function due to structural failure. This can cause the patient to fall.

Avoid areas with extreme temperatures (see Section 3.3 “Technical Information”).

⚠ CAUTION

Mechanical overloading. Exterior mechanical impact or stress such as shocks and vibrations can lead to

- short circuits in the electronics and battery and result in malfunction of the joint;
- defects on the battery and hydraulic damper resulting in leakage of liquid; or
- loss of function due to structural failure.

This can cause the patient to fall as well as result in skin irritation.

Do not expose system components to mechanical vibrations or shocks.

⚠ CAUTION

Penetration of dirt and humidity. Penetration of dirt and humidity into the system components can lead to

- short circuits in the electronics and battery and result in malfunction of the joint;
- defects on the hydraulic damper resulting in leakage of liquid; or
- loss of function due to structural failure.

This can cause the patient to fall as well as result in skin irritation.

- Do not let foreign particles or liquids enter the system components. Should the joint come into contact with liquid, remove the cosmetic cover and let the components dry. The joint must then be sent to an authorized Ottobock Service for inspection. Your prosthetist is the contact person.
- Always close the plastic cover of the protective elastomer cap after use.
- If the C-Leg compact comes into contact with salt water, immediately clean it with a cloth moistened with fresh-water and let it dry. The joint must then be sent to an authorized Ottobock Service for inspection. Your prosthetist is the contact person.

⚠ CAUTION

Improper use of the joint. Any kind of overloading or excessive strain as well as improper use can lead to

- defects and result in malfunction of the joint;
- loss of function due to structural failure; or
- defects on the battery and hydraulic damper resulting in leakage of liquid.

This can cause the patient to fall as well as result in skin irritation.

- The C-Leg compact was designed for everyday activities and must not be used for unusual activities such as extreme sports (i.e. free climbing, paragliding, etc.).
- Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures your personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service Centre.

⚠ CAUTION

Overheating of the hydraulic unit. Extended, continuous use (e.g. lengthy downhill walks) can lead to

- overheating of the hydraulic unit and result in malfunction of the joint; or
- defects on the hydraulic damper resulting in leakage of liquid.

This can cause the patient to fall as well as result in skin irritation. Touching overheated components can cause burns.

- Pay attention to the vibration signals that will occur in such cases to alert you that there is a risk of overheating. As soon as these vibrations begin, all activities must be stopped and the hydraulic unit be allowed to cool down. You may resume your activities once the vibration signals stop.
- If activities are continued despite the vibration signals, the hydraulic element may overheat and, in extreme cases, lead to a damage to the C-Leg compact. The joint should then be sent to an authorized Ottobock Service for inspection.

⚠ CAUTION

Risk of falling when walking backwards. When putting down the toe first when walking backwards, the C-Leg compact can switch from the high stance phase resistance to swing phase resistance. When the patient actively flexes the hip joint at this point in time, this can cause the patient to fall. When walking backwards, secure the joint actively with the residual limb muscles.

⚠ CAUTION

Risk of falling when going downstairs. When walking downstairs, the banister or handrail should always be used and the prosthetic foot should be placed on the step so that the heel (max. centre of the foot) is close to the edge of the step to facilitate rollover.

- Stop walking downstairs immediately whenever the sound warning beeps. Make careful tests to verify if the stance phase stabilization is active (see Section 2.11 “Important User Instructions”).
- Pay attention to the vibration and sound warnings (beeps) of the C-Leg compact.
- Special caution is required when walking downstairs while carrying children.

⚠ CAUTION

Non-active safety mode. If the safety mode can no longer be activated, there is the risk that the patient will fall. If the C-Leg compact does not switch to the safety mode (e.g. due to short-circuit due to water penetration), the amputee must actively stabilize the C-Leg compact at heel strike with his/her residual limb muscles until a prosthetist can be reached or a prosthesis replacement be accomplished.

⚠ CAUTION

Self-discharge of the battery. When the joint is not in use for an extended period of time, self-discharge of the battery will result. This can lead to insufficient power supply to the electronics of the joint and result in undefined conditions. This can cause the patient to fall.

- We therefore recommend recharging prior to every use.

⚠ CAUTION

Results of product deterioration. Wear and tear on system components can lead to malfunction of the joint. This can cause the patient to fall.

For your own safety, comply with the specified service and inspection intervals (maintenance of operational safety and guaranty).

⚠ CAUTION

Malfunction of the joint. Malfunctions of the joint can cause the prosthesis wearer to fall.

Pay attention to the vibration and sound warnings (beeps) of the C-Leg compact.

INFORMATION

When using exoprosthetic knee joints, servomotoric, hydraulic, pneumatic or brake load dependent control functions can cause movement noise. This kind of noise is normal and unavoidable. Usually, it does not cause any problems.

If the movement noise noticeably increases in the knee joint's life cycle, the joint should be inspected immediately by a prosthetist.

⚠ CAUTION


Risk of pinching where the joint bends. Ensure that fingers and other body parts are not in this area when bending the joint.

3 Alignment and Adjustment

3.1 Connecting the Tube Adapter

- 1 Mount the foot to the tube adapter (lightly tighten the screws).
- 2 Plug the tube adapter plug into the knee joint (only one polarity possible).
- 3 Push the protruding cable into the tube and push the tube adapter into the C-Leg compact. In doing so, take care to achieve the minimum insertion depth (see following table).

⚠ CAUTION Error during alignment of the prosthesis. Errors during the alignment of the prosthesis can lead to malfunction of the joint up to loss of function due to structural failure. This can cause the patient to fall. With the 3C96-1 and 3C86-1 models, the scaling must be positioned in front of the tube. The 2R80=110 Tube Adapter is approved for a maximum weight of 100 kg / 220 lbs.

- 4 Turn the foot slightly outward (adjustment range $\pm 30^\circ$ to suit the needs of the patient) and pre-fasten the screws lightly at the tube clamp.


3.2 Bench Alignment Using an Alignment Tool such as L.A.S.A.R. Assembly

INFORMATION

To align the prosthesis please proceed in two steps:

First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.

In a second step, the static alignment is optimized with the 743L100 L.A.S.A.R. Posture (see Section 2.5 "Static Alignment Optimization").

A correct bench alignment (e.g. using the 743L200 L.A.S.A.R. Assembly) ensures that the user can benefit from all the advantages of the C-Leg compact. The optimal residual limb position must be anticipated when positioning the socket connector. Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

- Position the middle of the foot (MF) approx. 30 mm anterior to the alignment reference line (A). This applies to all feet that are recommended for use with the C-Leg compact, independently of the specifications in the Instruction for Use of those feet!
Add 5 mm to the required heel height. Set correct outward rotation of the foot (Fig. 1).
- Clamp the knee joint with the mounted tube adapter into the alignment tool. Place the alignment reference point (=knee axis) approx. 0 – 5 mm anterior to the alignment reference line (Fig. 2). Consider the knee-floor distance and outward rotation of the knee (approx. 5° are provided for by the adapter insert in the L.A.S.A.R. Assembly). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- Connect the foot to the modular knee joint using a tube adapter. To do so, tilt the joint in the correct position and set the required tube length (Fig. 2).
- Mark the centre of the socket proximally (M) and distally on the lateral side. Draw a line through both marks from the socket brim to the distal end of the socket (Fig. 3).
Now position the socket such that the alignment reference line (A) passes through the proximal centre mark (M). Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity-to-floor distance (Fig. 4).


CAUTION Error during alignment of the prosthesis. Errors during the alignment of the prosthesis can lead to malfunction and insufficient fitting results of the joint up to loss of function due to structural failure. This can cause the patient to fall. If the residual limb flexion is not taken into account, the joint will be positioned too far to the front.

- Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41 Socket Adapter, etc.). For alignment correction, use the 4R112 Sliding Adapter (only for temporary use). When using socket adapters with a 4-hole connection (e.g. 5R1=*, 5R6=*), the 4R118 Sliding Adapter can be used for a permanent, additional posterior placement of 10 – 25 mm.


CAUTION Error during alignment of the prosthesis. If incorrectly positioned, the support for the lamination anchor is not provided. This puts too much strain on the screwed connection adapter and may cause damage and subsequently cause the patient to fall. The rotation adjustment **must not** be used for length adjustment. Screw in the thread as far as possible and clamp with the cap screw. In case of the **4R43** Lamination Anchor with threaded connector, position the **clamp laterally or medially**. In case of the **4R111=N** Lamination Anchor with threaded connector, the **anchor arms** must be positioned **in the anterior/posterior and medial/lateral directions**. Here, the clamp must be positioned to show in the anterior/medial or anterior/lateral direction.

- Control of the bench alignment and all measurements.

3.2.1 Checking the Socket after Bench Alignment

CAUTION Error during alignment of the prosthesis. Errors during the alignment of the prosthesis can damage the hydraulic unit of the joint. A defective hydraulic unit can lead to malfunction of the joint up to loss of function due to structural failure. This can cause the patient to fall. The hydraulic unit may become damaged by contact with the socket! At maximum flexion and insofar as contact with the frame of the C-Leg compact cannot be avoided (in case of voluminous residual limbs), the socket must lie flat on the frame. Soft cushioning at the socket will assist in maintaining the socket flat. At maximum flexion, it is essential that the minimum distance of 3 mm (1/8") is maintained between the hydraulic unit and the socket.

CAUTION Error during alignment of the prosthesis. Errors during the alignment of the prosthesis can damage the electronics of the joint. Defective electronics can lead to malfunction of the joint. This can cause the patient to fall. The electronic unit may become damaged by contact with the socket! At maximum extension (is reached under full capacity load!) the minimum distance of 5 mm (1/5") must be maintained between the electronic unit and the socket at the beginning of the dampening phase. For fittings with the 4X160=* C-Leg Protector, the minimum distance between the electronic unit and the socket must be at least 10 mm (3/8").

CAUTION

Error during alignment of the prosthesis. If the socket drain or air expulsion hole is incorrectly positioned, the draining sweat may drop into the electric connector and lead to short-circuits in the electronics and malfunction of the knee joint which may subsequently cause the patient to fall.

Make sure the socket drain or air expulsion hole is not positioned frontally above the electric connector. Position the hole on the side or back of the socket.

3.3 Torque Values of the Screw Connections

INFORMATION

The torque values for additional as well as optional system components (in item 3 below) are indicated in the Instructions for Use of those products.

Using the 710D4 Torque Wrench with the 710Y2=5 Hex Bit, turn the screws alternatively to the prescribed torque so that the torque is gradually increased:

- 1 Tube adapter: 15 Nm
- 2 Tube clamp: 8 Nm

Fitting for short residual limb

- 3 Rotation adapter or sliding adapter: 15 Nm

Fitting for long residual limb

- 4 Lamination anchor with threaded connector: 10 Nm

3.4 Charging the C-Leg compact

INFORMATION

Advise your patients of the information in this section.

⚠ CAUTION

Incorrect behaviour while charging the joint. If the patient walks while the battery charger is connected to the prosthesis, he can get caught on the cable and fall.

Take off the prosthesis prior to charging.

⚠ CAUTION

Incorrect handling while charging the joint. If the joint is not flexed completely during charging, the charging plug can get damaged and the joint can no longer be charged. This can lead to insufficient power supply to the electronics of the joint and result in undefined conditions. This can cause the patient to fall.

During charging, keep the C-Leg compact fully flexed.

INFORMATION

Prior to charging with the 4E50-* Charger, read the corresponding 647G262 Instructions for Use.

The battery is empty upon delivery. The battery of the C-Leg compact must therefore be charged prior to the first fitting of the patient.

- 1 Bend the C-Leg compact until it hits the stop.
- 2 Open the plastic cover of the protective elastomer cap.
- 3 Connect the 4E50-* Charger and verify the charging process (see following figures).


			
Both LEDs are off.	Both LEDs are lit.	Yellow LED flashes; green LED is lit.	Yellow LED is off; green LED is lit.
No power supply (or defective charger)	Battery is being charged. Battery capacity is lower than 50%.	Battery is being charged. Battery capacity is above 50%.	Battery is fully charged (or connection with C-Leg compact is interrupted).

- 4 After charging, remove the charging plug (C-Leg compact performs self-test) and close the plastic cover. Only afterwards, bring the prosthesis into extension.

INFORMATION

- The capacity of a fully charged battery is sufficient for one full day. We recommend charging the battery over night when using the prosthesis on a daily basis.
- Charging is only possible at temperatures above 0 °C.

3.5 Static Alignment Optimization (with 743L100 L.A.S.A.R. Posture)

The static alignment can be substantially improved using L.A.S.A.R. Posture 743L100. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- 1 After the self-calibration of the L.A.S.A.R. device, the load line can be measured. For that, the transfemoral amputee steps onto the force sensing plate with the C-Leg compact and onto the height compensation plate with the other leg. The prosthetic side should be sufficiently loaded (> 35% of the body weight).
- 2 The alignment is adapted exclusively by modifying the plantar flexion so that the load line/laser line runs approximately 30 mm in front of the alignment reference point (=knee axis) of the C-Leg compact.
- 3 After adjusting the C-Leg compact with C-Soft (see Section 2.7 "C-Soft") perform dynamic optimizations during trial walking.

3.6 BionicLink / optional: Serial Data Cable

INFORMATION

Prior to connecting the 60X3 BionicLink, read the corresponding 647G192 Instructions for Use.

- 1 Remove the plastic cover of the protective elastomer cap.
- 2 Connect the wireless adapter cable (BionicLink) (optional: connect serial data cable: To do so, the protective elastomer cap must be removed beforehand, Fig. 4.).
- 3 Snap the BionicLink into place and attach it medially (for left-side fittings) or laterally (for right-side fittings) to the C-Leg compact.

3.7 C-Soft

INFORMATION

Prior to the parameter adjustment with 4X180=* C-Soft (compatible from version 2.0), read the corresponding 647G268 Instructions for Use.

INFORMATION

Ottobock Data Station is the platform for Ottobock applications such as C-Soft and others.

Ottobock C-Soft cannot be used without the Ottobock Data Station. The Ottobock Data Station is a standard component of C-Soft and is automatically installed when installing C-Soft.

3.8 Optional: Parameter Settings with Foam Cover and C-Soft

The 3S26 and 3R59 Foam Covers can be used for the cosmetic enhancement of the leg prosthesis.

For that purpose, a 4X78 Charger Extension Cable can be integrated into the C-Leg compact fitting. If the charger extension cable is not used, a connection cap made of 2 mm Pedilin® makes it possible to pull down the foam cover in order to reach the charging receptacle at the C-Leg compact.

Allow a 60 mm "compression allowance" to minimize the effect of the protective foam cover on the knee function. When determining the proper cover length, increase the thigh dimension by 30 mm and the shin dimension by the same amount. The mounted foam cover changes the damping behaviour. With the 60X3 BionicLink (wireless programming interface), which is attached below the pre-shaped foam block, the damping behaviour can be completely taken into account by changing the parameter settings.

That process furthermore requires the 60X5 (60X4) BionicLink PC for the data transfer and the 4X180=* C-Soft (from version 2.0) for the parameter adjustment (see Section 2.7 "C-Soft").

3.9 Remote Control

INFORMATION

Advise your patients of the information in this section.

The C-Leg compact has **two modes**: The **compact mode** is intended for everyday use; in this mode, an additional standing function can optionally be activated and deactivated with the remote control.

The **therapy mode** enlarges the field of application of the C-Leg compact and offers additional functionality for early mobilization after amputation. This mode allows to permanently lock the knee joint in full extension. By pressing a button on the remote control, the joint can be unlocked, e.g. for sitting down.

INFORMATION

The mode can only be changed with C-Soft (from version 2.0); with previous versions this is not possible.

⚠ CAUTION

Improper switching between modes in C-Soft and use of the remote control. When switching between modes in C-Soft as well as when using the remote control, the functionality of the joint changes. In certain situations, this can cause the patient to fall.

When switching between modes in C-Soft as well as when using the remote control, the patient must always stand securely. Then he should try carefully, whether the desired function has been activated.

⚠ CAUTION

Penetration of water in the remote control. The remote control is not waterproof. If water penetrates the remote control, the device may become damaged (guaranty will become void). This can lead to malfunction of the remote control, result in unexpected actions of the joint, and cause the patient to fall.

Should water have penetrated the remote control, dry the device for at least one day at room temperature. Before starting to reuse the remote control, return it to an authorized Ottobock Service for inspection.

⚠ CAUTION

Manipulation on the remote control. Any changes or modifications you make to the remote control on your own initiative can lead to malfunction and result in unexpected actions of the joint. This can cause the patient to fall. Any changes or modifications to the device may limit its use.

INFORMATION

Important user instructions on the use of the remote control. Keep the C-Leg compact and the residual limb still (no bending or extending!).

For safety purposes, the working range of the remote control is limited to approximately 70 cm. However, the remote control will only function properly if held at least 30 cm away from the joint. If the activation/deactivation per remote control was not performed, ensure that the remote control is placed between 30 and 70 cm of the C-Leg compact and repeat the command. If you are still not successful, perform a pairing process (see Section 2.9.3).


3.9.1 Modes and Use of the Remote Control

Compact mode:

Upon delivery the knee joint is in the compact mode. If required, the joint can be locked or unlocked in the flexion direction at a flexion angle between 0° and 30° (standing function) with the remote control. For safety reasons, movement in the extension direction is always possible; in full extension position the joint remains locked (locking function).

To activate the standing or locking function (Fig. 1):

Press the remote control button for 2 seconds (the joint emits 2 beep sounds at intervals of one second). Two long beep sounds and two vibration signals follow to confirm the activation.

INFORMATION Operating the remote control at a flexion angle larger than 30° exclusively activates the locking function: **The joint automatically locks once full extension is reached.**

To deactivate the standing or locking function:

Again press the remote control button for 2 seconds. One long beep sound and one vibration signal follow to confirm unlocking.

INFORMATION The standing function can only be activated and deactivated, when the knee joint is kept still during operation of the remote control.

INFORMATION If the standing or locking function has been activated and the remote control is defective or the battery is empty, the knee joint can be unlocked by connecting the battery charger to the C-Leg compact. This works only in the compact mode.


Therapy mode:

The therapy mode can be activated in the C-Soft adjustment software (from version 2.0). Once activated, the joint is locked in full extension, for example for therapeutic applications for early mobilization after amputation. If required, the joint can be unlocked with the remote control (e.g. for sitting down). The hydraulic damping of the system supports the prosthesis wearer when sitting down or standing up. The prosthesis can be loaded.

To unlock the knee joint:

Press the remote control button for 2 seconds (the joint emits 2 beep sounds at intervals of one second). One long beep sound and one vibration signal follow to confirm unlocking. Now the knee joint is unlocked and can be flexed, e.g. for sitting down (Fig. 1). The joint remains in stance phase damping until it is brought into extension again.

As soon as the knee angle during extension is less than 7°, the joint automatically activates the lock in the flexion direction. This is confirmed by 2 long beep sounds and 2 vibration signals (Fig. 2).


INFORMATION If the knee joint is not flexed more than 7° within 7 seconds after unlocking, it automatically locks again for safety reasons. During this period, the remote control cannot be used.

3.9.2 Attachment of the Remote Control

The remote control can be:

- 1 attached to a key lanyard (Fig. 1);
- 2 attached to the belt of the prosthesis wearer with a clip (Fig. 2); or
- 3 attached to a holder provided on the socket (integrated in the laminate, attached with buckles or glued in place; Fig. 3).

3.9.3 Pairing (Establishing the Communication between the C-Leg compact and Remote Control)

⚠ CAUTION

Fault while connecting the C-Leg compact and remote control (pairing). If several joints are present in close proximity, this can result in unwanted connection of the remote control with another joint (pairing). This can cause the patient to fall.

Given that only one joint may be paired with a remote control, it must be ensured that no other joint is in the circumference of 3 m during the pairing.

ℹ INFORMATION

If the C-Leg compact still is in the delivery state, pairing is not required.

Pairing must be done

- as soon as a new remote control (replacement) is used with the C-Leg compact.
- if a switch cannot be effected with fully charged batteries (of both the joint and the remote control) within the working range of the remote control (30 – 70 cm).
- upon receipt of the remote control with a new battery.

Execution:

Briefly press the start button by poking a thin object (e.g. a paper clip) through the small hole of the remote control. The joint will confirm the successful pairing with 5 short beep signals.

3.9.4 Changing the Battery of the Remote Control

⚠ CAUTION

Impermissible replacement of the battery of the remote control. Replacing the battery on your own initiative can result in defects of the remote control. Malfunction and subsequent unexpected actions of the joint are then possible. This can cause the patient to fall.

The battery of the remote control may only be changed by an authorized Ottobock Service. The battery is automatically replaced at the biennial service intervals.

To change a battery outside of the biennial service intervals, send the remote control to an authorized Ottobock Service. Upon receipt of the remote control with a new battery, the remote control must be paired with the C-Leg compact (see Section 2.9.3).

3.9.5 Switching between Compact Mode and Therapy Mode

Switching from the compact mode to the therapy mode and vice versa is made with the C-Soft adjustment software (from version 2.0) and may only be carried out by an authorized prosthetist.

3.10 Finishing the Prosthesis

Upon finalizing all settings, all screws must be fastened and tightened to the proper torque.

- 1 Optional: Remove the cosmetic foam cover (and also the BionicLink, if present).
- 2 Secure all screws required for the prosthesis alignment (see Section 2.3) with 636K13 Loctite.
- 3 Tighten all screws to the proper torque (see Section 2.3) with the 710D4 Torque Wrench.
- 4 Optional: Apply the cosmetic foam cover or the 4X160 C-Leg Protector.

3.11 Important User Instructions

INFORMATION

Advise your patients of the information in this section.

Safety mode

⚠ CAUTION

Non-active safety mode. If the safety mode can no longer be activated, there is the risk that the patient will fall. If the C-Leg compact does not switch to the safety mode (e.g. due to short-circuit due to water penetration), the amputee must actively stabilize the C-Leg compact at heel strike with his/her residual limb muscles until a prosthetist can be reached or a prosthesis replacement be accomplished.

⚠ CAUTION

Danger when activating the safety mode. Upon activation of the safety mode, the damping behaviour of the joint changes. In certain situations, this can cause the patient to fall. As soon as the sound and vibration signals go off simultaneously, the patient must stop all activities with the leg prosthesis. After approximately 10 seconds, and from a secure standing position, check to see if the safety mode with the high flexion damping has been activated by slightly bending the C-Leg compact repeatedly under controlled weight bearing.

⚠ CAUTION

Safety mode cannot be deactivated. If the joint is exclusively in the safety mode, the joint has a defect. When in this condition, malfunctions can occur. This can cause the patient to fall. If the safety mode could not be deactivated by recharging the battery, an error has occurred. Contact your prosthetist for troubleshooting.

Aside from the normal operation mode, the C-Leg compact has a safety mode. As soon as a critical error occurs in the system or when the battery is almost empty, the C-Leg compact automatically switches to the safety mode. In that mode, the prosthesis activates a high flexion damping (high safety/reduced comfort), which enables the patient to walk safely despite a non-active system.

The switch to the safety mode is announced immediately prior to the switch with sound and vibration signals. If the safety mode was prompted by an empty battery, recharge the battery. The joint will then switch back automatically to the operation mode.

Battery capacity

During normal operation, the battery capacity is indicated with different vibration signals:

Battery capacity	Vibration impulse
Approx. 1 hour of operation remains	3 ×
Approx. 30 min. of operation remains	5 ×
When turning off	10 ×, then turns off

Safety signal for incorrect connection to the tube adapter

⚠ CAUTION

Safety signal occurs. The safety signal indicates an error of the joint. Malfunctions may occur in this situation. This can cause the patient to fall.

If the safety signal occurs, discontinue operating the C-Leg compact and resolve the fault by connecting the tube adapter and re-calibrating (zero setting) the joint (by the prosthetist or service technician). If a contact fault (between the tube adapter plug and C-Leg compact) has occurred (e.g. following fitting by the prosthetist), send the knee joint (incl. tube adapter) to an authorised Ottobock Service Centre.

If the tube adapter is somehow disconnected from the C-Leg compact joint while in operation, the C-Leg compact will emit short beeps and, for approximately five minutes, vibrate with slow, pulsing signals. The C-Leg compact will activate the safety mode and maintain it until the tube adapter cable is correctly attached (provided the battery capacity remains sufficient).

Ventilation

Air may accumulate in the hydraulic unit if the C-Leg compact is not stored in an upright position. This is noticeable through sounds and irregular damping behaviour.

The automatic ventilation mechanism ensures that all functions of the C-Leg compact are again intact after approximately 10 – 20 steps.

4 Additional Information

4.1 Service Intervals

INFORMATION

A maintenance service should be obligatorily carried out every two years from the date of delivery.

To ensure proper function and the safety of the patients, a maintenance service is required every two years. This service includes the inspection of the sensors and the replacement of worn parts.

To obtain the service, send in the knee joint with mounted C-Leg tube adapter and remote control as well as the complete charger unit including the AC adapter.

4.2 Damage Event

To simplify problem resolution, describe damage events in detail and document patient statements.

Unless the cause of a damage is very apparent, send in all accessories.

4.3 Technical Information

Operating and storage temperature range	-10 °C to +60 °C
Weight of the joint	1210 g
Weight of the 2R80 Tube Adapters	178 g to 256 g
Weight of the 2R81 Tube Adapters	438 g to 482 g
Max. flexion angle	125 °
Operating voltage of the battery charger:	100 V - 260 V AC
Operating frequency of the battery charger:	50 Hz - 60 Hz
Battery charging temperature:	> 0 °C
C-Leg compact service life:	8 years if service intervals are maintained
Battery service life:	500 cycles
Battery type:	lithium-ion battery
Charging time:	max. 6 hours

4.3.1 Symbols on the Knee Joint

Legal manufacturer

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Serial number

Caution, hot surface

4.3.2 Symbols on remote control

Legal manufacturer

[LOT]PPPP YYYY WW

Lot number

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Non-ionising radiation

Compliance with the requirements according to "FCC Part 15" (USA)

Compliance with the requirements according to the "Radiocommunications Act" (AUS)

4.3.3 Symbols on the battery charger

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

[LOT]PPPP YYYY WW

Chargen-Nummer

4.4 Transportation

Always use the **X-3C100 Shipping Case** when transporting your electronic joint system and the 3C96-1=S/ 3C86-1=S Service Joint.

4.5 Guaranty

Ottobock offers an extensive guaranty service on the basis of the current C-Leg compact guaranty concept (see 646D201=* C-Leg compact Service Card). Please comply with the prescribed service intervals (see Section 3.1 "Service Intervals"). For your own interest, we also recommend obtaining a confirmation from your patients for the provision of instructions and the handing over of the device.

4.6 Liability

The device is only to be used under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

4.7 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

4.8 CE Conformity

The electronic 3C96-1/3C86-1 C-Leg compact Prosthesis System meets the requirements of the 93/42/EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

The electronic 3C96-1/3C86-1 C-Leg compact Prosthesis System also meets the requirements of the 1999/5/EC guidelines for radio equipment and telecommunications terminal equipment. Conformity assessment was carried out by Ottobock according to appendix IV of the guidelines. A copy of the declaration of conformity can be requested from the manufacturer (see back side).

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España: (+34) 91 806 30 00
 (+34) 609 41 09 80

FCC-Statement: This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care

Two Carlson Parkway North, Suite 100

55447 Minneapolis, Minnesota, USA

Phone + 1-763-553-9464

IC-Statement:

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation. The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website <http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.

5470 Harvester Road

L7L 5N5 Burlington, Ontario

Canada

Patents:

European Patent No. 0 549 855 (B, D, GB, F, I, L, NL, A, P, S, CH, E),

patented in Canada 1991 No. 2 057 108, patented in Russia No. 2 089 138, R.O.C.

Invention Patent No. 076 088,

patented in Japan No. 3 131 933, patented in Korea No. 176 977

Otto Bock Healthcare Products GmbH

Brehmstraße 16 · 1110 Wien · Austria

T +43-1 523 37 86 · F +43-1 523 22 64

info.austria@ottobock.com · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.