



2R57/2R58, 2R76/2R77 4R91, 4R82, 4R82=P

DE	Gebrauchsanweisung	4
EN	Instructions for Use	6
FR	Instructions d'utilisation	8
IT	Istruzioni per l'uso	11
ES	Instrucciones de uso	13
PT	Manual de utilização	16
NL	Gebruiksaanwijzing	18
SV	Bruksanvisning	20
DA	Brugsanvisning	23
PL	Instrukcja użytkowania	25
HU	Használati utasítás	28
CS	Návod k použití	30
HR	Upute za uporabu	32
BS	Uputstvo za upotrebu	35
TR	Kullanma talimatı	37
EL	Οδηγίες χρήσης	39
RU	Руководство по применению	42
JA	取扱説明書	44
ZH	使用说明书	46

● 4D4	▲	■
(D) Einzelteile-Pack	Mindestmenge	Einzelteile
(GB) Single-Component Pack	Minimum order quantity	Single components
(F) Kit de pièces de rechange	Quantité minimum	Pièces à l'unité
(I) Imballo dei singoli componenti	Quantità minima	Singoli componenti
(E) Kit componentes	Cantidad mínima	Componentes
(P) Conjunto de peças de reposto	Quantidade mínima	Componentes
(NL) Onderdelenpakket	Minimikvantitet	Onderdelen
(S) Service-Set	Minimum aantal	Separata delar
(DK) Komponentpakke	Mindste mængde	Komponenter
(PL) Zestaw naprawczy	Ilość minimalna	Podzespoły
(H) Alkatrész-csomag	Minimum mennyiség	Alkatrészek
(CZ) Sada jednotlivých dílů	Minimální množství pro objednání	Jednotlivé díly
(RO) Set componente individuale	Cantitate minimă	Componente individuale
(HR) Pakiranje pojedinačnih dijelova	Minimalna količina	Pojedinačni dijelovi
(TR) Müfritler parçaları ambalajı	Minimum miktar	Müfritler parçaları
(GR) Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα	Μεμονωμένα εξαρτήματα
(RUS) Комплект деталей	Минимальное количество	Отдельные детали
(J) 交換パーツセット	最低注文個数あり	個別販売
(中国) 维修组件	最低起订量	单个部件

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English

Date of the last update: 2012-12-03

- Please read this document carefully.
- Follow the safety instructions.

1 Scope of Delivery

1.1 4D4 Single-Component Pack (for Tube Clamp Adapters) (ill. 5+6, ●)

Includes 2 straight pins (2) und 1 cap screw (3).

1.2 Minimum Order Quantity (for Tube Clamp Adapters) (ill. 5+6, ▲)

- 506G3=M8×14 (5) (also for Tube Adapters 2R57, 2R58, 2R76, 2R77)
- 4X28=3 (4)

1.3 Single Components (ill. 5+6, ■)

None

2 Description

2.1 Application

The Tube Adapters 2R76/2R77, 2R57/2R58 and the Tube Clamp Adapters 4R91, 4R82 and 4R82=P are to be **exclusively** used for the prosthetic fitting of the lower limb.

2.2 Field of Application

Field of application according to the **Ottobock Mobility System MOBIS®**

4R91, 4R82, 4R82=P, 2R57, 2R76, 2R58, 2R77:

Approved for a patient with a weight of up to 150 kg/330 lbs.

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water and acids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

2.3 Function

The Tube Adapter is to be used in conjunction with the Tube Clamp Adapter as an adjustable and length-adaptable connection element for Ottobock Modular Lower Extremity Prostheses.

For maximum structural strength the Tube Clamp Adapter must be placed proximally (as shown in the cover illustration). The clamping slot must always be positioned anteriorly.

3 Technical Data

	2R57	2R76	2R58	2R77	4R82	4R82=P	4R91
Weight	220 g	260 g	330 g	370 g	95 g	90 g	140 g
Diameter	34 mm	34 mm	34 mm	34 mm	34 mm	34 mm	34 mm
System height	-	-	-	-	33 mm	-12 mm	33 mm
Min.system height	77 mm	77 mm	77 mm	77 mm	-	-	-

	2R57	2R76	2R58	2R77	4R82	4R82=P	4R91
Max.system height	282 mm	282 mm	472 mm	472 mm	–	–	–
Material	Titanium	INOX*	Titanium	INOX*	Titanium	Titanium	INOX*
Maximum patient weight	150 kg	150 kg	150 kg	150 kg	150 kg	150 kg	150 kg

* stainless steel

4 Handling

4.1 Connecting the Internal Adapter with Tube to the Tube Clamp Adapter

The Tube Adapter must be shortened as necessary based on the amputee's measurements. This must be done correctly, as noted below, to avoid distorting the tube. Due to the close tolerances of quality Ottobock adapters, a distorted tube will not fit easily into the Clamp Adapter.

Attention!

Do not clamp the tube in a vise, or it may become permanently distorted!

The cut surface must be perpendicular to the long axis of the tube (ill. 1). Use of the 719R2 Tube Cutter or equivalent is recommended.

All inner and outer burrs must be removed from the cut edge (ill. 2). Use of the 718R1 Reamer or equivalent is recommended.

Loosen the Cap Screw (ill. 3/5/6) of the Tube Clamp Adapter by two turns, check, whether the plastic ring has been inserted into the Tube Clamp Adapter completely and with the chamfer on the ring's outer diameter first, and then insert the tube 5 mm deep into the Tube Clamp Adapter while slightly rotating the tube. Next, place the Tube Clamp Adapter with the tube inserted vertically on the worktop and fully insert the Tube Adapter into the Adapter until the tube is in contact with the base rim inside the clamp.

Attention!

The interior surface of the Tube Clamp Adapter and the exterior surface of the Tube Adapter must be cleaned with 634A3 (Aceton) in the insertion area prior to final assembly (insertion of the Tube Adapter into the Tube Clamp Adapter).

For definitive use, the cap screw must be **torqued to 5 Nm/3.7 ft. lbs/44 In. lbs** with 710D1 Torque Wrench.

Attention!

For reasons of structural strength, the residual pieces of cut tube adapters must not be used as components in prostheses.

4.2 Adjustment Screws (ill. 3/4/5)

Adjustments may be made through adjusting the four 506G3 Adjustment Set Screws during alignment, fitting and at any time after the prosthesis has been completed. The optimal position of the components may be preserved when exchanging parts or when servicing the prosthesis by removing only the two most deeply set adjacent screws.

Once positioning has been determined, replace any set screws that are too long (to prevent foam cover damage) or too short (lacking strength).

Screws are available in the following sizes:

Article no.	Length	Article no.	Length
506G3=M8x10	10 mm	506G3=M8x14	14 mm
506G3=M8x12-V	12 mm	506G3=M8x16	16 mm

Use a 710D4 Torque Wrench to tighten the set screws.

Attention!

(Applies to the application of the Tube Adapters 2R57; 2R58; 2R76; 2R77 and Tube Clamp Adapters 4R82; 4R91)

For testing: with the torque wrench 710D4, screw in the set screws 506G3 of the Adapter clockwise at **10 Nm (7.4 lbf. ft / 88.5 lbf. in)**.

For finalizing the prosthesis: secure the set screws 506G3 with Loctite 636K13, tighten the set screws clockwise at **10 Nm (7.4 lbf. ft / 88.5 lbf. in)** first and then tighten completely at **15 Nm (11.1 lbf. ft / 133 lbf. in)**.

4.3 Maintenance Instructions

Note:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

5 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

6 Warranty

Ottobock warrants this product/these products from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local Ottobock distribution company for more information.

7 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Ottobock HealthCare GmbH