




ottobock.

4R91, 4R82, 2R57, 2R76, 4R82=P, 2R58, 2R77

DE	Gebrauchsanweisung	4
EN	Instructions for use	9
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	23
PT	Manual de utilização.....	28
NL	Gebruiksaanwijzing.....	32
SV	Bruksanvisning.....	37
DA	Brugsanvisning.....	42
PL	Instrukcja użytkowania.....	46
HU	Használati utasítás	51
CS	Návod k použití.....	55
HR	Upute za uporabu	60
TR	Kullanma talimatı.....	64
EL	Οδηγίες χρήσης	68
RU	Руководство по применению	73
JA	取扱説明書.....	78
ZH	使用说明书.....	82

	● 4D4	▲	■
DE	Einzelteile-Pack	Mindestmenge	Einzelteile
EN	Single-Component Pack	Minimum order quantity	Single components
FR	Kit de pièces de rechange	Quantité minimum	Pièces à l'unité
IT	Imballo dei singoli componenti	Quantità minima	Singoli componenti
ES	Kit componentes	Cantidad mínima	Componentes
PT	Conjunto de peças de reposto	Quantidade mínima	Componentes
NL	Onderdelenpakket	Minimikvantitet	Onderdelen
SV	Service-Set	Minimum aantal	Separata delar
DA	Komponentpakke	Mindest mængde	Komponenter
PL	Zestaw naprawczy	Ilość minimalna	Podzespoły
HU	Alkatrész-csomag	Minimum mennyiség	Alkatrészek
CS	Sada jednotlivých dílů	Minimální množství pro objednání	Jednotlivé díly
RO	Set componente individuale	Cantitate minimă	Componente individuale
HR	Pakiranje pojedinačnih dijelova	Minimalna količina	Pojedinačni dijelovi
TR	Münferit parça ambalajı	Minimum miktar	Münferit parçalar
EL	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα	Μεμονωμένα εξαρτήματα
RU	Комплект деталей	Минимальное количество	Отдельные детали
JA	交換パーツセット	最低注文個数あり	個別販売
ZH	维修组件	最低起订量	单个部件

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Date of the last update: 2014-12-10

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1 Scope of delivery

1.1 4D4 Single-component pack (for Tube Clamp Adapters) (ill. 5+6, ●)

Includes 2 straight pins (2) und 1 cap screw (3).

1.2 Minimum order quantity (for Tube Clamp Adapters) (ill. 5+6, ▲)

- 506G3=M8×14 (5) (also for Tube Adapters 2R57, 2R58, 2R76, 2R77)
- 4X28=3 (4)

1.3 Single components (ill. 5+6, ■)

None

2 Description

2.1 Indications for use

The Tube Adapters 2R76/2R77, 2R57/2R58 and the Tube Clamp Adapters 4R91, 4R82 and 4R82=P are to be **exclusively** used for the prosthetic fitting of the lower limb.

2.2 Field of application

Field of application according to the **Ottobock Mobility System MOBIS**

4R91, 4R82, 4R82=P, 2R57, 2R76, 2R58, 2R77:

Approved for a patient with a weight of up to 150 kg/330 lbs.

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water and acids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

2.3 Function

The Tube Adapter is to be used in conjunction with the Tube Clamp Adapter as an adjustable and length-adaptable connection element for Ottobock Modular Lower Extremity Prostheses.

For maximum structural strength the Tube Clamp Adapter must be placed proximally (as shown in the cover illustration). The clamping slot must always be positioned anteriorly.

3 Technical data

	2R57	2R76	2R58	2R77	4R82	4R82=P	4R91
Weight	220 g	260 g	330 g	370 g	95 g	90 g	140 g
Diameter	34 mm	34 mm	34 mm	34 mm	34 mm	34 mm	34 mm
System height	–	–	–	–	33 mm	-12 mm	33 mm
Min.system height	77 mm	77 mm	77 mm	77 mm	–	–	–
Max.system height	282 mm	282 mm	472 mm	472 mm	–	–	–
Material	Titanium	INOX*	Titanium	INOX*	Titanium	Titanium	INOX*
Maximum patient weight	150 kg	150 kg	150 kg	150 kg	150 kg	150 kg	150 kg

* stainless steel

4 Handling

4.1 Connecting the internal adapter with tube to the tube clamp adapter

The Tube Adapter must be shortened as necessary based on the amputee's measurements. This must be done correctly, as noted below, to avoid distorting the tube. Due to the close tolerances of quality Ottobock adapters, a distorted tube will not fit easily into the Clamp Adapter.

Attention!

Do not clamp the tube in a vise, or it may become permanently distorted!

The cut surface must be perpendicular to the long axis of the tube (ill. 1). Use of the 719R2 Tube Cutter or equivalent is recommended.

All inner and outer burrs must be removed from the cut edge (ill. 2). Use of the 718R1 Reamer or equivalent is recommended.

Loosen the Cap Screw (ill. 3/5/6) of the Tube Clamp Adapter by two turns, check, whether the plastic ring has been inserted into the Tube Clamp Adapter completely and with the chamfer on the ring's outer diameter first, and then insert the tube 5 mm deep into the Tube Clamp Adapter while slightly rotating the tube. Next, place the Tube Clamp Adapter with the tube inserted vertically on the worktop and fully insert the Tube Adapter into the Adapter until the tube is in contact with the base rim inside the clamp.

Attention!

The interior surface of the Tube Clamp Adapter and the exterior surface of the Tube Adapter must be cleaned with 634A3 (Aceton) in the insertion area prior to final assembly (insertion of the Tube Adapter into the Tube Clamp Adapter).

For definitive use, the cap screw must be **torqued to 5 Nm/3.7 ft. lbs/44 In. lbs** with 710D4 Torque Wrench.

Attention!

For reasons of structural strength, the residual pieces of cut tube adapters must not be used as components in prostheses.

4.2 Adjustment screws (ill. 3/4/5)

Adjustments may be made through adjusting the four 506G3 Adjustment Set Screws during alignment, fitting and at any time after the prosthesis has been completed. The optimal position of the components may be preserved when exchanging parts or when servicing the prosthesis by removing only the two most deeply set adjacent screws.

Once positioning has been determined, replace any set screws that are too long (to prevent foam cover damage) or too short (lacking strength).

Screws are available in the following sizes:

<i>Article no.</i>	<i>Length</i>	<i>Article no.</i>	<i>Length</i>
506G3=M8x10	10 mm	506G3=M8x14	14 mm
506G3=M8x12-V	12 mm	506G3=M8x16	16 mm

Use a 710D4 Torque Wrench to tighten the set screws.

Attention!

(Applies to the application of the Tube Adapters 2R57; 2R58; 2R76; 2R77 and Tube Clamp Adapters 4R82; 4R91)

For testing: with the torque wrench 710D4, screw in the set screws **506G3** of the Adapter clockwise at **10 Nm** (7.4 lbf. ft / 88.5 lbf. in).

For finalizing the prosthesis: secure the set screws 506G3 with Loctite 636K13, tighten the set screws clockwise at **10 Nm** (7.4 lbf. ft / 88.5 lbf. in) first and then tighten completely at **15 Nm** (11.1 lbf. ft / 133 lbf. in).

4.3 Maintenance instructions

Note:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

5 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

5.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH