(D)	Wichtiger Hinweis	B2
(GB)	Important	B3
(F)	Consigne importante	B4
(I)	Nota	B5
(E)	Indicación importante	B6
(P)	Recomendações importantes	B7
(NL)	Belangrijk	B8
(S)	Viktiga tips	B9
(PL)	Ważna wskazówka	B10
(CZ)	Důležité upozornění	B11
(TR)	Önemli	B12
(GR)	Σημαντική υπόδειξη	B13

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For optimal use of the Otto Bock 757T13 MyoSelect, we recommend a careful and thorough reading of this supplementary sheet and of the instructions for use. Understanding the device and its functions will let you to utilize all the features of this device!

Important

- With *MyoSelect*, only one electronic MYOBOCK component can be identified and adjusted at a time.
- Should a prosthesis contain more than one MYOBOCK component (e.g. *13E195 Four-channel Processor II* and *12K50 ErgoArm® Electronic plus*), only the component that is directly connected to the *MyoSelect* can be identified and adjusted at a time (here: *Four-channel Processor II*).

Some components may have to be dismantled to allow for the identification and adjustment of the another ones. (For example, dismantle the *Four-channel Processor II* and then connect the *ErgoArm Electronic plus* directly to the *MyoSelect*.)

- To adjust a 13E195 Four-channel Processor II which is used in combination with an ErgoArm® Electronic plus (version 1.0 or 1.2), the power supply must be interrupted first. For this purpose, dismantle the Four-channel Processor II and pull the cable off the battery connection of the Four-channel Processor II.

Alternative method:

It is also possible to adjust a 13E195 Four-channel Processor II which is used in combination with an ErgoArm® Electronic plus (version 1.0 or 1.2) without dismantling the Four-channel Processor II. In this case, the following steps must be followed **exactly**:

- 1) Remove the battery from the prosthesis.
- 2) Connect the deactivated MyoSelect (The MyoSelect automatically turns off after approximately 15 seconds if no entry or connection to a Myo component takes place. The MyoSelect **MUST BE AND REMAIN** turned off!)
- 3) Insert the battery of the prosthesis: The MyoSelect will automatically turn on and establish a connection to the Four-channel Processor II.

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Otto Bock Healthcare Products GmbH