







## 3R41

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	12
<b>FR</b> Instructions d'utilisation .....	19
<b>IT</b> Istruzioni per l'uso .....	27
<b>ES</b> Instrucciones de uso .....	35
<b>PT</b> Manual de utilização .....	43
<b>NL</b> Gebruiksaanwijzing .....	51
<b>SV</b> Bruksanvisning .....	59
<b>DA</b> Brugsanvisning .....	66
<b>NO</b> Bruksanvisning .....	73
<b>FI</b> Käyttöohje .....	80
<b>HU</b> Használati utasítás .....	88
<b>CS</b> Návod k použití .....	95
<b>RO</b> Instrucțiuni de utilizare .....	103
<b>HR</b> Upute za uporabu .....	111
<b>TR</b> Kullanma talimatı .....	118
<b>EL</b> Οδηγίες χρήσης .....	125
<b>RU</b> Руководство по применению .....	134
<b>JA</b> 取扱説明書 .....	142
<b>ZH</b> 使用说明书 .....	149

**INFORMATION**

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## INFORMATION

Last update: 2021-02-18


- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

## 1 Indications for use

### 1.1 Medical purpose

The 3R41 is to be used **exclusively** for the exoprosthetic fitting of amputations of the lower limb. The 3R41 can be used for unilateral or bilateral amputees.

### 1.2 Application

Field of application according to the Ottobock MOBIS Mobility System.

The 3R41 is intended exclusively for applications within mobility grade 1 (indoor walkers).

Approved for a **patient weight of up to 125 kg/275 lbs.**

### 1.3 Application conditions

The 3R41 has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

For the required environmental conditions, please refer to the Technical Data section.

### 1.4 Environmental conditions

#### ⚠ CAUTION

##### Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table „Unallowable environmental conditions“ in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

##### Allowable environmental conditions

###### Temperature range:

Use: -10 °C – +45 °C

Storage, Transportation: -20 °C – +60 °C

Relative humidity: 20 % – 90 % (Nicht kondensierend)

Commercially available, solvent-free cleaning agents

## Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

### 1.5 Reuse and lifetime

#### Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

### 1.6 Prosthetist qualifications

The fitting of patients with the 3R41 may only be carried out by a prosthetist.

## 2 Safety instructions

#### Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- Observe the safety information in these instructions for use.
- Make the patient aware of all safety instructions marked with „Inform the patient“.

#### Exceeding the service life and reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Ensure that the approved service life is not exceeded (see section „Service Life“).
- Only use the product for a single patient.
- Inform the patient.

#### Excessive strain on the product

Fall due to breakage of load-bearing components

- Use the prosthesis components according to the MOBIS classification (see section „Area of Application“).
- Inform the patient.

#### Mechanical damage to the product

Injuries due to changes in or loss of functionality

- Use caution when working with the product.
- Check a damaged product for proper functionality and readiness for use.
- In case of changes in or loss of functionality, do not continue using the product (see „Signs of changes in or loss of functionality during use“ in this section).
- Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- Inform the patient.

## Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthesis components relative to each other and by the development of noises.

### 3 Patient information

- The knee joint should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage.
- The knee joint function may be affected by both environmental and service conditions. To reduce the risk for the patient, the knee joint should not be used when noticeable functional changes occur. Noticeable functional changes may include different locking and unlocking behavior, poor response, incomplete extension, play in the bearings, and abnormal noises. If you notice any damage or functional changes, please have the prosthesis inspected by a specialist workshop.
- The knee joint may only be repaired by an authorized Ottobock service center.
- Regularly clean the product with a cloth slightly dampened with soap solution. Cleaning with chemicals or solvents is not permissible, since there is a risk of material damage and loss of functionality.
- The knee joint must not be lubricated or greased, since there is a risk of material damage and loss of functionality.
- Always make use of the handrail when negotiating stairs.
- Do not unlock the knee joint while walking or under flexion moment loading.
- Put weight on the knee joint only when it is fully extended and has been locked. Be sure you hear the lock latch's acoustic feedback when locking the knee joint.
- If the pull cable becomes entangled, this can lead to unlocking of the knee joint and cause the standing or walking patient to fall. In such a case, apply a slight flexion moment to test whether the knee joint has locked again.
- When handling the knee joint without the cosmetic foam cover, be careful not to get your fingers jammed.
- When extending and locking the knee joint, ensure that the cosmetic foam cover does not get caught since this could prevent the joint from locking correctly.

### 4 Maintenance instructions

Ottobock recommends readjusting the knee joint's settings once the patient has gotten used to the prosthesis. This period of time varies depending on individual patient characteristics.

It is recommended to have the 3R41 examined by a prosthetist for visible signs of wear and tear as well as damage within the scope of a routine inspection (at least once a year). As a precaution, retighten the Allen screw (wrench size 4 mm) of the tube clamp to 10 Nm. If necessary, the O-ring and the elastomer elements can be replaced at that time (see section 10). Note the correct mounting orientation during installation.

Special attention should be paid to the movement, unlocking, and locking resistance, the bearing clearance, and to unusual noise. Full flexion and extension must be guaranteed.

## 5 Scope of delivery and spare parts

Scope of delivery	Spare parts (Fig. 1)
1 pc. of 3R41 Modular Knee Joint with Lock	2 pcs of Elastomer Elements (Item 2) Order no. 4Z105
1 pc. of Instructions for Use (647G340=1)	1 pc. of O-Ring (Item 3) Order no. 4Z104
1 pc. of 4F18=N Lock Slide	1 pc. of Lock Slide, complete (no illustration) Order no. 4F18=N

## 6 Description and function

The 3R41 is a monocentric knee joint with lock. By using fibre-reinforced polyamide we managed to produce a lightweight knee joint which is easy to maintain.

### 6.1 Joint lock

The lock mechanism, which is located in the upper joint section, automatically secures the joint in the extended position; you can hear and feel the latch lock in place.

### 6.2 Flexion release/unlocking

Flexion release is activated via the pull cable, which is attached to the perlon cord loop. The lock unlocks with a stroke length of 15 mm and a force of up to 35 N. Another way to achieve flexion release is to press the lock latch (Fig. 2, item 1, arrow). For safety reasons, unlocking of the knee joint becomes more difficult with increasing flexion moment loading.

The patient may unlock the joint before or after sitting down. To unlock the 3R41, the patient either relieves most of the load on the prosthesis or applies load on the toes (extension moment) while unlocking the joint. While holding the unlocked state (pressed lock latch or pulled "lock"), the patient initiates the flexion movement.

### 6.3 Locking

After standing up, the patient automatically locks the joint with an extension moment with full extension up to the extension stop.

### 6.4 Plastic flap

A movably attached plastic flap (Fig. 2, item 3) covers the anterior joint space over the entire flexion range and thus protects the cosmetic foam cover and the fingers from getting jammed.

### 6.5 Connection

The pyramid adapter connects the knee joint proximally (to the socket). Distal connection is provided by the clamp of the lower joint section into which the 30 mm tube of the tube adapter is inserted.

## 7 Alignment recommendations and assembly

### Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- Observe the alignment and assembly instructions.

The cap screw (see Fig. 1, item 1) secures the pyramid adapter.

Do not loosen the cap screw.

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis. Function of the joint is influenced by the position

of the joint axis. In the extended position, the joint axis is located posterior to the alignment reference line (Fig. 9), which ensures sufficient stability and requires little effort for unlocking (Fig. 9).

**The optimal residual limb position must be anticipated when positioning the socket connector** (Fig. 8). Plumb lines in the frontal and sagittal planes (drawn from the hip joint's center of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

**To align the prosthesis please proceed in two steps:**

- 1) First, bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.
- 2) Then static alignment optimization using 743L100 L.A.S.A.R. Posture.

**7.1 Bench Alignment with Alignment Tool** (the following steps refer to Fig. 8)

- ① Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ② Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- ③ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, **place the joint axis (alignment reference point) 15 mm posterior to the alignment reference line.** At that point the pyramid base should be horizontal. Pay attention to the knee ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point is on the level of the medial tibial plateau.
- ④ Connect the foot to the modular knee joint using a tube adapter. Also see section 7.5.
- ⑤ Mark the center of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- ⑥ Now position the socket such that the alignment reference line passes through the proximal center mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip flexion contracture) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.
- ⑦ Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R41, 4R51, 4R55, 4R95, 4R111, 5R1=\* Socket Adapter). Also see section 7.6.

**7.2 Static alignment optimization using 743L100 L.A.S.A.R. posture**  
(the following steps refer to Fig. 9)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability under load and when unlocking, please proceed as follows:

- ① To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded. For the user's safety, this should take place between parallel bars or with the aid of forearm crutches.
- ② Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The load line (laser line) should be approx. **45 mm anterior to the joint axis (alignment reference point).**

**7.3 Combination possibilities**

**Failure to observe the manufacturer requirements regarding combination possibilities**

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- ▶ Contact the manufacturer with questions.

## INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

### Unallowable combinations

Prosthetic hip joints

7E10\*

#### 7.4 Mounting the tube adapter to the knee joint

- 1) Determine the tube length and cut off the tube. Use a tube cutter (e.g. 719R2 Ottobock Tube Cutter) to shorten the tube. Do not clamp the tube into a vise to avoid deformation and surface damage of the tube! The cutting plane must be level and at a right angle to the longitudinal axis of the tube.
- 2) Debur the inside and outside of the tube using an interior and exterior pipe reamer (e.g. 718R1 Ottobock Interior and Exterior Pipe Reamer).
- 3) Check the outside of the tube for damage. The surface must be smooth and free of burs. If a damaged tube is inserted into the knee joint, this can impair the strength of the knee joint.
- 4) Lift the plastic flap slightly and **insert the tube up to the stop**.
- 5) Using a 710D4 torque wrench, tighten the Allen screw (wrench size 4 mm) of the tube clamp (Fig. 2, item 2) to 10 Nm.

#### 7.5 Mounting the knee joint to the socket

The contact surfaces of the pyramid adapter and pyramid receiver must be undamaged, clean and smooth. Static corrections are possible at any time during alignment and trial walking as well as after completion of the prosthesis using the four adjustment screws.

Make sure that the set screws of the adapters used have the appropriate length. Set screws which are too short may not provide adequate stability for the prosthesis. Screws which are too long can damage the cosmetic foam cover.

#### Note!

Apply 636K13 Loctite to the set screws, screw in and slightly tighten them. Then adjust the adapter and tighten the set screws to the prescribed torque. For the tightening torques please refer to the Instructions for Use of the corresponding adapters. Use a 710D4 torque wrench.

To test the functionality, unlock the knee joint and flex it fully. The socket must contact the tube clamp cover (Fig. 3).

#### Attention!

The upper joint section must not contact the lower joint section, since this could put too much strain on the knee joint.

## 8 Finishing the prosthesis

Upon delivery, the perlon cord loop and pull cable are connected to the 3R41. Cut off the pull cable such that it is slightly taut but **does not yet cause the knee joint to unlock**. The pull cable should only be unlocked after a pull stroke of approx. 10 mm. The pull cable should be on the lateral side of the prosthesis.

### Note!

If the perlon cord loop has been cut by mistake, do not pull the cord through the knee joint. The cord may only be threaded by an Ottobock service technician after disassembly of the knee joint.

## 9 Cosmetic foam cover

Use the 3S107 Foam Cover for the 3R41.

## 10 Replacing the o-ring (Fig. 7)

- 1) Remove the screw (item 1).
- 2) Bend up the tube clamp cover slightly (item 2) and remove the O-ring 4Z104 (item 5).

### Note!

Make sure that the pin with the bore (item 3) and the pin with the internal thread (item 4) do not fall out!

If the pins fall out:

- Completely remove the tube clamp cover.
  - Position the pins on the left and right side of the thread. Be sure to have the correct mounting orientation of the pin with the bore (item 3).
  - Put the tube clamp cover over them.
- 3) Put the new O-ring into the notch provided.
  - 4) Lower the tube clamp cover and tighten the screw.

## 11 Technical data

Proximal connection	Pyramid adapter
Distal connection	Tube clamp Ø 30 mm
Knee flexion angle	150 °
Weight	385 g
System height	24 mm
Proximal system height from alignment reference point	12 mm
Distal system height up to alignment reference point	12 mm
Material	Fibre-reinforced polyamide
Admissible temperature range	-10 °C / +40 °C
Max. user weight	125 kg / 275 lbs
Mobility Grade	1

## 12 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

## 12.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

## 12.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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