


MyoWrist Transcarpal 10V38=*

DE	Gebrauchsanweisung	5
EN	Instructions for use	12
FR	Instructions d'utilisation	18
IT	Istruzioni per l'uso	25
ES	Instrucciones de uso	31
PT	Manual de utilização.....	38
NL	Gebruiksaanwijzing.....	44
SV	Bruksanvisning.....	51
TR	Kullanma talimatı.....	57

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1 Indications for use

Last update: 2020-12-10

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The MyoWrist Transcarpal 10V38=* **is to be used exclusively** for exoprosthetic fitting of the upper limb in combination with a 8E44=* Transcarpal-Hand.

It provides individual flexion and extension of the wrist joint as well as fixation in 5 different positions, in 20° increments. This way it supports the physiologically correct posture when working with the prosthesis and it particularly reduces the need for unnatural compensating movements when performing longer monotonous works.

For **fitting long residual limbs** the MyoWrist Transcarpal 10V38=* is connected to the **9S266 Chassis with Quick-Disconnect Wrist**.

For the **fitting of short and medium residual limbs** combination with the **8R1 Endoskeletal-Adapter** is recommended.

In both cases **passive pro- and supination** is possible **thanks to the Quick-Disconnect Wrist**. In many cases **electromotoric pro- and supination** recommends in combination with the optional **10S17 Electric Wrist Rotator**.

2 Safety instructions

Attention!

- To guarantee a solid bond, clean the bonding areas with acetone before applying the adhesive. Once you have degreased the areas with acetone, do not touch them, but immediately coat them with adhesive as soon as the acetone has dried.
- The Processing and Safety Instructions of the adhesive producer must be observed.

- The indicated turning moments must be observed.
- Disassembly and repair of the MyoWrist Transcarpal 10V38=* must be carried out solely by authorized Ottobock Service technicians.
- The patient must be instructed in the proper use according to section 3 of the "Patient Informati".

3 Patient information

The MyoWrist Transcarpal 10V38=* has been developed for everyday activities and must not be used for unusual activities such as extreme sports (free climbing, skydiving, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the prosthetist will pass the prosthesis on to Ottobock Myo-Service. Any damaged components may only be repaired by the certified Ottobock Myo-Service.

Attention!

- Kindly note that when carrying heavy objects in case of full flexion or extension this might cause a torque moment that unintentionally makes the hand rotate.
- Make sure the hand can flex unintentionally, if the flexion switch has been pressed erroneously.

3.1 Operating a motor vehicle

An upper limb amputee's ability to drive a motor vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork). The user must be able to operate the vehicle without any risk with the electric hand turned off. Driving with the electric hand turned ON may present a risk to other persons if the hand inadvertently opens due to unintentional muscle contraction or other causes.

**Caution:**

Before operating a motor vehicle, the Quick-Disconnect Unit and flexion angle of the MyoWrist Transcarpal 10V38=* must be so positioned that a slight turn of it, which can occur during steering, cannot disconnect the hand from the prosthesis!

4 Scope of delivery (display of title)

MyoWrist Transcarpal 10V38=*

9S267=* Retaining Ring

9E388 Cable Seal

5 Description and function

The MyoWrist Transcarpal 10V38=* allows for the fitting of both, long and short residual limbs, in combination with a Ottobock 8E44=* Transcarpal Hand and in addition optionally with a 8R1 Endoskeletal-Adapter.

The Ottobock MyoWrist Transcarpal 10V38=* allows for flexion and extension settings in different positions.

The wrist unit is unlocked by pressing the black release knob (Fig. 2). The integrated ratchet system holds the hand in the desired position until it is locked and thus facilitates exact positioning. Flexion and extension take place passively. Stop pressing the release knob and the joint will relock.

The wrist units are obtainable in the following sizes:

10V38=L7 1/4 10V38=R7 1/4

10V38=L7 3/4 10V38=R7 3/4

10V38=L8 1/4 10V38=R8 1/4

6 MyoWrist Transcarpal 10V38=* for the fitting of short residual limbs

Combination with a 8R1=* Endoskeletal-Adapter is necessary for this fitting (Fig. 1).

Before finishing a short socket, please be aware of the limited amount of space inside the socket.

6.1 Determining socket length

Attach the 8E44=* Transcarpal Hand to the MyoWrist* Transcarpal 10V38=*, but do not secure the nuts with Loctite. (Fig. 3).

Next, pre-assemble the tube of the Endoskeletal-Adapter to the MyoWrist Transcarpal 10V38=* , without bonding the tube in place. For further proceeding, please refer to the Instructions Manual of the Endoskeletal-Adapter 647H501.

6.2 Assembly instructions

After having cut the tube to the desired length and having assembled the Endoskeletal-Adapter, please proceed as follows (Fig. 4):

Clean all threads (except the two bolts in the plate) with acetone. Do not touch the cleaned areas! After the acetone has evaporated, coat the following threads with 636K12 Loctite 245: Threaded Connector with Tube, Threaded Connector. Coat the Chassis with Quick-Disconnect Wrist with 636K13 Loctite 241.

Feed the four-pole cable through the tube adapter until it comes out the other end. Screw the Chassis with Quick-Disconnect Wrist together with the threaded connector (hand tight, approx. 8 Nm). Use pliers with plastic jaws for this purpose. Feed the four-pole cable through the MyoWrist Transcarpal 10V38=* and tighten (approx. 8 Nm) to the Endoskeletal-Adapter by hand.

Please be aware, that the cable can become twisted in this process (Fig. 5). The cable must straightly pass through the wrist unit in order not to be damaged (Fig. 6).

Slide the Cable Seal over the Four-Pole Cable (Fig. 7) and press the Cable Seal into the plate of the MyoWrist Transcarpal (Fig. 8).

Make sure that the cable **is not** stretched tight, since the flexion movement of the wrist will require a certain amount of cable length. Besides functioning as a seal, the Cable Seal acts as a strain-relief for the cable. It is also important that enough cable is left in the tube for later repairs.

Attach the Transcarpal-Hand 8E44=* to the plate as shown in the Instructions Manual of the Endoskeletal-Adapter 647H501.

The cable is connected as described in the Instructions Manual of the Endoskeletal-Adapter 647H501.

7 MyoWrist Transcarpal 10V38=* for the fitting of long residual limbs (Fig. 2)

For this fitting the unit must be combined with the Chassis of the Quick-Disconnect Wrist 9S266 (Fig 9). The scope of delivery does NOT include the coaxial cable for attachment of a Transcarpal Hand. The latter must be ordered separately. The DMC plus control requires 9E397=7; the Digital Twin control requires 9E397=8.

7.1 Assembly instructions

Clean all threads (except the two bolts in the plate) with acetone.

Do not touch the areas cleaned with acetone and after the acetone has evaporated, immediately coat with Loctite 241 (636K12).

Feed the four-pole cable through the flexed wrist until it comes out the other end and screw the components together hand tight (about 8 Nm).

Please be aware, that the cable can become twisted in this process (Fig. 5). The cable must straightly pass through the wrist unit in order to prevent it from being damaged (Fig. 6).

Slide the Cable Seal over the Four-Pole Cable (Fig. 7) and press it into the plate of the MyoWrist Transcarpal (Fig. 8).

Make sure that the cable **is not stretched** tight, since the flexion movement of the wrist will require a certain amount of cable length. Besides functioning as a seal, the Cable Seal acts as a strain-relief for the cable. It is also important that enough cable is left in the tube for later repairs

Attach the 8E44=* Transcarpal Hand to the MyoWrist* Transcarpal 10V38=*, but do not secure the nuts with Loctite. (Fig. 3).

Determine the appropriate cable length for connecting it to the electronic circuit on the Transcarpal-Hand (Fig. 11). Mark the necessary length. Disconnect the Transcarpal-Hand 8E44=* from the wrist joint and cut the cable to length (Fig. 10).

Remove approximately 2 mm of the cable housing and prepare the wire with solder (Fig. 12). Solder the four-pole cable to the Cable Connector (Fig. 13).

Polarity of the cables:

While three of the housing contacts have been spaced equally apart, another contact has been set apart from them at a greater distance. Solder the outermost gray wire of the Four-Pole Cable to this contact (Arrows, Fig. 13). Pay attention to the order of the wires and make sure they are connected to the correct contacts!

Attach the Transcarpal-Hand 8E44=* to the plate of the MyoWrist Transcarpal 10V38=* as described in the Instruction Manual 647H398. Secure the nut with 636K12 Loctite 245. Make sure that the cable is properly connected to the electronic circuit (Fig. 14)!

Pull the System Inner Hand over the Transcarpal-Hand and fix in place with the Retaining Ring.

8 Technical data

Setting positions: 5 setting positions in 20° increments

Neutral or zero setting 0°

Extension setting: -20° / -40°

Flexion setting: 20° / 40°

Length: 24,7 mm

Outside Ø: 36 mm

Weight: 7 1/4 85 g
7 3/4 87 g
8 1/4 89 g

Maximum load capacity: 32,2 Nm ~27 kg

9 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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