


642V15=1, 87A5

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English

Date of the last update: 2014-01-10

- Please read this document carefully.
- Follow the safety instructions.

1 Instructions for Use

Impresil is a silicone material designed for the use in orthotics. For this reason, it must be used only by qualified professionals who have been instructed with regard to its application, for example in an Ottobock seminar. Prior to first-time use of the silicone material, these instructions for use must be read thoroughly. The manufacturer recommends the product be used only under the specified conditions and for the intended purposes. The product must be maintained and processed in accordance with these instructions for use. The manufacturer will not assume any liability whatsoever for damages that occur as a result of non-compliance with the instructions for use.

2 Safety Advice

The silicone material Impresil is suitable for direct application to the intact surface of the skin (no open wounds). It is used for taking impressions of parts of the body, face and residual limbs, and is designed for external use only. Please make sure that the silicone material does not get into body orifices, such as eyes, mouth, nostrils, ears or genitals. For safety reasons, you should perform a test prior to actual application, for example on the inner side of the wrist. Do not use if erythema or skin rash occur!

3 Properties of the Material

- Impresil is a two-component silicone material which crosslinks at room temperature (polyaddition crosslinking).
- It is thixotropic so the material should not flow downwards even on vertical surfaces when applied with a thickness of less than 5 mm.
- The material is provided in cartridges, each containing 75 ml of component A and 75 ml of component B.
- The Shore A hardness of the vulcanized silicone amounts to approx. 30°.

- Impresil does not contain any irritating or allergenic components. It is suitable for application to the intact surface of the skin (no open wounds).
- Impresil always provides the right mixture and thus cures reliably.
- It maintains a high degree of detail in the impressed area (skin details, scars etc.).
- Impresil is easy to transport (using the supplied kit) and can be used immediately wherever necessary.
- Impresil can be processed cleanly under normal environmental conditions.
- In contrast to plaster bandages, Impresil leaves no residues on the skin.
- It bonds with other silicones (Pastasil 85H11 or Chlorosil 85P31) and sticks to other plastics, provided these have a rough surface. When using the material on smooth surfaces, such as plastics or glass, a self-adhesive strap can be attached prior to silicone application in order to prevent silicone detachment.
- Impressions can be made more cost-effectively because they consume less time and material.
- Less stress and strain for the patient.
- Impresil does not stick to the skin and is thus easy to release.

4 Storage (Fig. 1)

The material may be stored at room temperature until the indicated use-by date if the cartridges are fitted with plugs. If stored at too high a temperatures (above 25 °C), the material may become stiff, which prevents appropriate mixing.

5 Processing Time

At normal room temperature (23 °C), the processing time of Impresil is about 6 minutes. This period will be prolonged if the silicone is stored in a cool place. If stored in a warm place or heated after application (using a hot air gun), this period is shortened accordingly.

At the end of this time, vulcanization is not yet complete, but the silicone will have become so hard that it is not possible to further modify. In order to prevent deformation, the silicone should be left on the patient's skin for a few more minutes before the impression is removed.

6 Properties of Impressions

- Impresil impressions do not shrink after release. They can be stored and shipped for an unlimited period of time.
- They are not altered by external impacts (temperature, atmospheric pressure, air humidity).
- Impresil impressions can be filled with any mold-making material.
- When released carefully, they can be filled several times.

Even in completely vulcanized state, Impresil will crosslink with itself again. If an impression is to be extended, just apply a new portion of Impresil. However, the new material should overlap the old by approx. 1 cm. The shorter the time between the first and second application, the tighter the two parts will bond together.

7 Handling

Use the dispenser (mixing gun) (Fig. 2) and the static mixing tubes (Fig. 3) to mix the components homogeneously and in the correct proportion. They can be applied immediately thereafter. If work interruptions are shorter than the processing time, you may keep on working on a continuous basis. If you should require more than one Impresil cartridge for a large-scale impression, the static mixer can be reused.

Unlock and pull off protective cap by turning it counter-clockwise 90 degrees (Figs. 4–6).

Put on static mixer in correct position and lock by turning it clockwise 90 degrees (Fig. 7).

Pull out dispenser plunger completely and put cartridge into dispenser with notch at bottom. Lock dispenser and push plunger into cartridge with the operating handle several times until static mixer is filled with silicone.

Avoid air pockets when applying the material. If required, use silicone spatula 89G10 to distribute and smooth the material. When doing so, make sure that you apply the silicone with an even thickness of at least 2 mm.

When removing the static mixer and putting on the cover cap again, the rest of the silicone in the cartridge will remain unmixed and can be used for further applications. Only the silicone in the static mixer will vulcanize and thus be lost. (It can be easily removed after hardening) (Fig. 8).

8 Skin Isolation

Application of a thin layer of vaseline will facilitate release, in particular if the impression is not to be cut open. In case of very hairy skin, it is important to apply a sufficient layer of vaseline.

9 Suitability and Consumption (Fig. 9)

Impresil can be used on the patient in many areas of application to take impressions and to make corrections during trial fitting:

Finger prostheses	¼ cartridge per finger
Partial hand prostheses, scar compression gloves	1– 2 cartridges per hand
Orthoses	1 – 1 ½ cartridges for one wrist orthosis for children
Liners for scar contractures of a transfemoral residual limb	½ – 1 cartridges
Scar compression masks	2 – 3 cartridges for one face
Socket volume adjustment for transtibial and transfemoral prostheses (depending on volume)	½ – 1 cartridge

These figures are approximate only since the consumption of the material is subject to considerable fluctuations depending on the size, thickness and length of the impression to be made.

When making impressions, a minimum thickness of 2 mm must be ensured in order to provide sufficient stability of the impression. Large-scale impressions (masks, pads etc.) must be stabilized by applying outside reinforcements made of silicone, plaster bandages or other materials.

If one part of the impression is to be made of silicone while the remainder will consist of plaster bandages (e.g. a liner impression with scar contractures), the silicone section must be molded first. While still unvulcanized, the silicone back must be fitted with a piece of stockinette, which will then attach to the plaster bandages and thus fix the silicone part in the plaster negative (Fig. 9).

10 Release

When releasing – even extensive – impressions, it is often not necessary to cut them open since the silicone is elastic compared to a plaster model and can thus expand. Should it become necessary to cut or tear it open, please ensure that the impression is glued together again immediately after release. For this purpose, ensure correct alignment of cutting edges and use Impresil to glue from outside.

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt/Germany
T +49 5527 848-0 · F +49 5527 72330
healthcare@ottobock.de · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.