
# ottobock.

## 6A30=10N

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for use .....	7
<b>FR</b>	Instructions d'utilisation .....	10
<b>IT</b>	Istruzioni per l'uso .....	14
<b>ES</b>	Instrucciones de uso .....	19
<b>PT</b>	Manual de utilização.....	23
<b>NL</b>	Gebruiksaanwijzing.....	27
<b>SV</b>	Bruksanvisning .....	31
<b>DA</b>	Brugsanvisning.....	34
<b>PL</b>	Instrukcja użytkowania.....	38
<b>CS</b>	Návod k použití.....	42
<b>RO</b>	Instrucțiuni de utilizare.....	46
<b>HR</b>	Upute za uporabu .....	50
<b>EL</b>	Οδηγίες χρήσης .....	54
<b>TR</b>	Kullanma talimatı.....	58
<b>RU</b>	Руководство по применению .....	62
<b>JA</b>	取扱説明書.....	66
<b>ZH</b>	使用说明书.....	70

**Achtung!**

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**Caution!**

Date of last update: 2021-02-22

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

## 1 Indications for use

The 6A30=10N Shuttle Lock is to be used **solely** for the prosthetic fitting of the lower limb. It is not suitable for being used in bathing prostheses.

**Caution!**

If using unsuitable liners you can cause skin irritations. For further information, please see Liner Selection Tool 646D252=GB.

## 2 Description

The Shuttle-Lock is used to connect the prosthetic socket of a transtibial or transfemoral prosthesis to the socket liner and is integrated into the laminated socket during the lamination process.

A release button on the ratchet unit can be used to unlock the Pin from the Shuttle-Lock.

**Caution!**

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water and acids.

**Please inform your patients.**

### 3 Component parts (fig. 1)

The Shuttle Lock Set consists of:

- (1) 6Y13=1 Pin
- (2) Shuttle Lock Housing
- (3) 6A52 Ratchet Unit  
and 5X125 Dummy Set with screw
- (4) Dummy Sleeve
- (5) Dummy Plate with four screws
- (6) Dummy for Ratchet Unit, long
- (7) Dummy for Ratchet Unit, short

### 4 Assembling and laminating


#### Caution!

Installation of the pin in a liner that is not approved

Falling due to loosening of the screw connection

- Only install the pin in liners with a metal thread and observe the specified tightening torque.
- Inform the patient.

For more details about the fabrication of a transtibial prosthesis with Liner and 6A30=10N Shuttle Lock refer to the Technical Information Sheet 2.1.6.

#### Caution!

Sometimes the alignment line (solid line) differs from the long axis of the stump (broken line). The shuttle-lock must follow the broken line.

- Use a rasp to flatten the distal end of the plaster positive to obtain an even surface for the dummy plate (5).
- Secure the dummy plate with the four screws provided with the kit (fig. 2).
- Pull over a 99B81 PVA Bag. Tie off distally in the groove of the dummy plate and cut off.
- Apply silicone grease to the thread of the short dummy for the ratchet unit (7) and insert it in the shuttle lock housing.
- Place the shuttle lock housing (2) with the dummy sleeve (4) and secure with the screw with the dummy plate (5) (fig. 3).

**Note: For trial fitting, slightly grease the thread of the ratchet unit (3) with 633F11 Silicone Grease.**

- **Degrease the ratchet unit (3) and the shuttle lock housing, secure it using 636K13 Loctite and screw it into the shuttle lock housing. Tighten using 710D21 Torque Wrench. Torque: 20 Nm.**

**Caution!**

Improper tools may damage the ratchet unit or impair its function!

Never use an open-end wrench for unscrewing and retightening.  
Use only the tool indicated in the Instructions for Use.

Use the following tools:

- 710D21 Torque Wrench
- the suitable 709S26=1/4"×1/2" adapter piece and
- the S6A20-2 Socket Wrench Kit.

**Caution!**

Improper tools may damage the push button!

To remove the push button from the ratchet unit, please use a screw driver with 2 mm tip or a coin of comparable thickness.

Whenever grinding inside the socket, protect the shuttle lock from dirt by sealing the lock opening with tape or a piece of clean cloth.

Please find further processing instructions in the Technical Information 2.1.6.

Secure the Pin with Loctite 241 (636K13), screw it into the liner's thread and tighten to 3 Nm (2.2 lbf.ft., 27 lbf.in.) using the 710D21 Torque Wrench or equivalent.

### **Caution!**

The pretension of the spring influences the rotation of the ratchet unit during donning the prosthesis. If the spring is removed (no audible feedback) or the pretension is reduced, the rotation resistance of the ratchet unit is also reduced.

## **5 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **5.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **5.2 CE conformity**

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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