


## 1M10 Adjust

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	9
<b>FR</b> Instructions d'utilisation .....	15
<b>IT</b> Istruzioni per l'uso .....	21
<b>ES</b> Instrucciones de uso .....	26
<b>PT</b> Manual de utilização .....	32
<b>NL</b> Gebruiksaanwijzing .....	38
<b>SV</b> Bruksanvisning .....	44
<b>DA</b> Brugsanvisning .....	49
<b>NO</b> Bruksanvisning .....	54
<b>FI</b> Käyttöohje .....	60
<b>PL</b> Instrukcja użytkowania .....	65
<b>HU</b> Használati utasítás .....	71
<b>CS</b> Návod k použití .....	76
<b>RO</b> Instrucțiuni de utilizare .....	82
<b>HR</b> Upute za uporabu .....	88
<b>SL</b> Navodila za uporabo .....	93
<b>SK</b> Návod na použitie .....	98
<b>BG</b> Инструкция за употреба .....	104
<b>TR</b> Kullanma talimatı .....	110
<b>EL</b> Οδηγίες χρήσης .....	115
<b>RU</b> Руководство по применению .....	121
<b>JA</b> 取扱説明書 .....	127
<b>ZH</b> 使用说明书 .....	133
<b>KO</b> 사용 설명서 .....	137

## INFORMATION

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### INFORMATION

Date of last update: 2020-05-18

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The 1M10 Adjust prosthetic foot is suitable for use in modular prostheses. The prosthetic foot distinguishes itself with a high level of safety while standing and walking. This is ensured by a multi-

axial joint in the load line. The heel characteristics can be individually adapted with a function module.

### 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 1 (indoor walker) and mobility grade 2 (restricted outdoor walker).

Maximum body weight				
Size [cm]	22 to 23	24 to 25	26 to 27	28 to 30
<b>Stiffness 1</b>	52 kg	58 kg	72 kg	77 kg
<b>Stiffness 2</b>	68 kg	76 kg	95 kg	100 kg
<b>Stiffness 3</b>	80 kg	100 kg	125 kg	125 kg

### 2.3 Environmental conditions

Storage and transport
Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions
<b>Temperature range:</b> -10 °C to +45 °C
<b>Moisture:</b> relative humidity: 20% to 90%, non-condensing
<b>Chemicals/liquids:</b> fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)
<b>Solids:</b> dust

Prohibited environmental conditions
<b>Chemicals/moisture:</b> salt water, perspiration, urine, acids, soapsuds, chlorine water
<b>Solids:</b> dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

Other environmental conditions	
<b>UV light</b>	Not resistant
<b>Hydrolysis</b>	Not resistant

### 2.4 Service life

This product has been tested by the manufacturer with 2 million load cycles according to ISO 22675. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

#### **CAUTION**

##### **Exceeding the service life and reuse on another patient**

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

#### **CAUTION**

##### **Use under unallowable environmental conditions**

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

#### **CAUTION**

##### **Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification.

#### **CAUTION**

##### **Unallowable combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

#### **NOTICE**

##### **Mechanical damage to the product**

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**NOTICE****Mechanical overload**

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**Signs of changes in or loss of functionality during use**

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

**4 Scope of delivery**

Quantity	Designation	Reference number
1	Instructions for use	-
1	Prosthetic foot	-
1	Protective sock	SL=Spectra-Sock-7

**Spare parts/accessories (not included in the scope of delivery)**

Fig.	Reference number	Reference number
1	Footshell	2C1
1	Insert for function module	2D11

**5 Preparing the product for use****⚠ CAUTION****Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

**INFORMATION**

- ▶ Not all of the materials may be available in your country. In this case, please contact the local branch of the manufacturer to obtain information on alternative materials.

**5.1 Alignment****5.1.1 Applying/removing the footshell****Putting on the footshell**

- 1) Slide the prosthetic foot into the footshell (see fig. 3, see fig. 4).
- 2) Press the heel of the prosthetic foot into the footshell until it engages (see fig. 5).

**Removing the footshell**

- 1) Push the lock of the footshell to the rear and pull up the heel of the prosthetic foot.
- 2) Remove the prosthetic foot from the footshell.

**5.1.2 Bench Alignment****TT bench alignment****Bench alignment process**

**Required materials:** 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)

<b>Bench alignment process</b>	
Perform the assembly and alignment of the prosthetic components in the alignment apparatus according to the following specifications:	
<b>Sagittal plane</b>	
①	Heel height: <b>Effective heel height</b> (shoe heel height - sole thickness of forefoot) + <b>5 mm</b>
②	Exterior foot rotation: <b>approx. 5°</b>
③	Anterior placement, middle of the prosthetic foot to the alignment reference line: <b>30 mm</b>
④	Connect the prosthetic foot and prosthetic socket using the chosen adapters. The instructions for use of the adapters must be observed.
⑤	Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: <b>Individual residual limb flexion + 5°</b>
<b>Frontal plane</b>	
⑥	Alignment reference line of prosthetic foot: <b>between big toe and second toe</b> Alignment reference line of prosthetic socket: <b>along the lateral patella edge</b>
⑦	Observe the abduction or adduction position.

## TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

### 5.1.3 Static Alignment

- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219\***, TT modular leg prostheses: **646F336\***) may be requested from Ottobock.
- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.

### 5.1.4 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.

#### 5.1.4.1 Optimising the heel characteristics

#### INFORMATION

Changing the configuration of the function module has no effect on the load line while standing.

The heel characteristics are adjusted with the function module. The following settings are possible:

Heel strike	Selected configuration (see fig. 2)	Position of insert
Soft	1	Insert removed, opening posterior
Moderate	2 (delivery condition)	Insert posterior
Firm	3	Insert anterior

- Turn the function module to the desired configuration (see fig. 6).

## 5.2 Optional: Installing the foam cover

The foam cover sits between the prosthetic socket and prosthetic foot. It is cut longer in order to compensate for the movements of the prosthetic foot and prosthetic knee joint. During flexion of the prosthetic knee joint, the foam cover undergoes posterior compression and anterior elongation. The foam cover should be stretched as little as possible in order to increase its service life.

There is a connecting element (such as a connection plate, connection cap or connection cover) on the prosthetic foot.

- > **Required materials:** degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive
- 1) Measure the length of the foam cover on the prosthesis and add the length allowance.  
**TT prostheses:** Distal allowance for movement of the prosthetic foot.  
**TF prostheses:** Allowance proximal of the knee rotation point for flexion of the prosthetic knee joint and distal allowance for movement of the prosthetic foot.
- 2) Cut the pre-shaped foam cover to length and fit it in the proximal area on the prosthetic socket.
- 3) Pull the foam cover over the prosthesis.
- 4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or rests on the foot adapter.
- 5) Install the prosthetic foot on the prosthesis.
- 6) Mark the outer contour of the connecting element on the distal face of the foam cover.
- 7) Disassemble the prosthetic foot and remove the connecting element.
- 8) Clean the connecting element using a degreasing cleaner.
- 9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.
- 10) Let the glue dry (approx. **10 minutes**).
- 11) Install the prosthetic foot and adapt the exterior cosmetic shape. Take into account compression caused by cosmetic stockings or SuperSkin.

## 6 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 7 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## 9 Technical data

Sizes [cm]	22	23	24	25	26	27	28	29	30
Heel height with 2C1=*N [mm]	10 ±5								

Sizes [cm]	22	23	24	25	26	27	28	29	30	
Heel height with 2C1=*S [mm]	20 ± 5						–			
System height with 2C1=*N [mm]	46	48	53		57	59	64			
System height with 2C1=*S [mm]	38		44		49	–				
Product weight without footshell [g]	255		320		385		465			
Max. body weight [kg]	80		100		125					
Mobility grade	1 + 2									
Colours	beige (4), light brown (15)									

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