15K10

DE Gebrauchsanweisung	3
EN Instructions for use	5
FR Instructions d'utilisation.....	8
IT Istruzioni per l'uso	10
ES Instrucciones de uso	13
PT Manual de utilização.....	16
NL Gebruiksaanwijzing.....	18
SV Bruksanvisning	21
DA Brugsanvisning	23
NO Bruksanvisning.....	26
FI Käyttöohje.....	28
PL Instrukcja użytkowania.....	31
HU Használati utasítás	33

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INFORMATION

Date of last update: 2020-04-03

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 15K10 modular foam kit serves to create an attractive exterior contour for modular prostheses from the forearm to the shoulder girdle level. It contains the 13R8=64 mounting flange.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is to be used solely for upper limb prosthetic fittings.

2.2 Area of application

The product serves to create passive arm prostheses.

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +45 °C

Relative humidity: max. 90%, non-condensing

Moisture: fresh water and salt water as splashing water, sweat

Unallowable environmental conditions

Mechanical vibrations or impacts

Acids

Dust, sand, highly hygroscopic particles (e.g. talcum), smoke

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

3.2 General safety instructions

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

CAUTION

Reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

Signs of changes in or loss of functionality during use

Instability or lack of joint function indicate a loss of function.

4 Scope of delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

15K10 modular foam kit				
Fig.	Item	Quantity	Designation	Reference number
-	-	1(■)	Instructions for use	647G456
1	-	1(■)	Forearm moulded foam part	-

15K10 modular foam kit				
Fig.	Item	Quantity	Designation	Reference number
1	–	1()	Upper arm moulded foam part	–
1	–	1()	Perlon frizz stockinette	–
1	–	1()	Mounting flange	13R8=64

5 Preparing the product for use

5.1 Preparation

- > The modular prosthesis has been prepared.
- INFORMATION: In both areas, add an allowance of 2 to 3 cm for compression.**
Forearm area: determine the length of the prosthesis from the terminal device to the modular elbow joint.
Upper arm area: determine the length of the prosthesis from the modular elbow joint to the prosthesis socket/prosthetic shoulder joint.
 - Determine the circumferences on the contralateral side. For a bilateral fitting, pay attention to the overall aesthetic impression.
 - Cut the foam moulded parts to length.

5.2 Assembling

- > **Materials needed:** 636W17 plastic adhesive, 13R9* connection disc with groove (for connection with a wrist), 13R8=64 connection flange (for connection to the lamination ring), 623T8=9 Perlon frizz stockinette
- Bond the required connection components to the moulded foam parts with plastic adhesive.
 - Push the moulded foam parts onto the modular prosthesis.
 - Optional:** Bond the moulded foam parts to each other.
 - Grind the contour of the cosmetic foam cover.
 - Pull a double layer of Perlon frizz stockinette over the cosmetic foam cover and secure it to the connection disc.

6 Cleaning

- Clean the product with a damp, soft cloth.
- Dry the product with a soft cloth.
- Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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