




743A18

DE	Gebrauchsanweisung (Fachpersonal)	5
EN	Instructions for Use (Qualified Personnel)	9
FR	Instructions d'utilisation (Personnel spécialisé)	14
IT	Istruzioni per l'uso (Personale tecnico specializzato)	19
ES	Instrucciones de uso (Personal técnico especializado)	24
PT	Manual de utilização (Pessoal técnico)	29
NL	Gebruiksaanwijzing (Vakmensen)	34
SV	Bruksanvisning (Fackpersonal)	38
DA	Brugsanvisning (Faguddannet personale)	43
NO	Bruksanvisning (Fagpersonell)	47
PL	Instrukcja użytkowania (Personel fachowy)	51
HU	Használati utasítás (szakszemélyzet)	56
CS	Návod k použití (Odborný personál)	60
RO	Instrucțiuni de utilizare (Personal de specialitate)	65
TR	Kullanma talimatı (Uzman personel)	70
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό)	74
RU	Руководство по применению (Квалифицированный персонал)	79
JA	取扱説明書 (有資格担当者)	84
ZH	使用说明书 (专业人员)	88

INFORMATION

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English

INFORMATION

These instructions for use were last updated on 02.11.2020.

INFORMATION

Please read the instructions for use carefully.

1 Single Components

1.1 Single Components

For available single components refer to the catalogue.

1.2 Accessories

- 743Y42=24 Lamination Dummy for 8K5=* Physolino Babyhand
- 11S12 Attachment Piece for laminated inner sockets
- 99B83=* Tube Dummies for the fabrication of vacuum formed inner sockets
- 507S15 Steel Washer with hole

The accessories are not included in the delivery. They must be ordered separately.

2 Description

2.1 Intended Use

The 743A18 Alignment Tool is to be used **solely** for the fabrication of prosthetic sockets for the exoprosthetic fitting of the upper extremities.

2.2 Field of Application

The 743A18 Alignment Tool is used as an alignment aid for forearm prostheses.

2.3 Function

The 743A18 Alignment Tool consists of the following single components:

- 743Y41 Alignment Rod
- 743Y42=40 Lamination Dummy
- 743Y42=45 Lamination Dummy
- 743Y42=50 Lamination Dummy
- 743Y42=54 Lamination Dummy

For a myoelectrically controlled prosthesis, the following lamination dummies are to be used:

- 743Y42=45 Lamination Dummy for hand size 7 and 7 ¼
- 743Y42=50 Lamination Dummy for hand size 7 ¾ and 8E33=* System Electric Greifer
- 743Y42=54 Lamination Dummy for hand size 8 ¼

For all other fittings, the following applies:

- 743Y42=40 Wrist diameter of 40 mm
- 743Y42=45 Wrist diameter of 45 mm
- 743Y42=50 Wrist diameter of 50 mm

For vacuum formed inner sockets, the tube dummies are offered in two different sizes matching the 99B13=* PVC Connection Tube:

- 99B83=16 for children
- 99B83=21 for adults

3 Scope of delivery

1 pc. 743A18 Alignment Tool

1 pc. 647G474 Instructions for Use

4 Handling

4.1 Fabrication of the Inner Socket

The later position of the 743Y41 Alignment Rod and the position of the 99B13=* PVC Connection Tube is determined with the 11S12 Attachment Piece and 99B83=* Tube Dummy respectively. The residual limb and fitting conditions of the patient (unilateral or bilateral) must be taken into account. In order to achieve the greatest possible degree of freedom for the patient, it is recommended for positioning the 11S12 Attachment Piece and 99B83=* Tube Dummy to attach them on the plaster positive – starting from the neutral O-position – in a slightly adducted and flexed direction.

4.2 Laminated Inner Socket

In case of a laminated inner socket, the 11S12 Attachment Piece is to be attached to the distal end of the plaster positive with light putty. After laminating the inner socket, the 743Y41 Alignment Rod is screwed into the Attachment Piece.

4.3 Vacuum Formed Inner Socket

- 1) Vacuum form the inner socket with 616T53=* or 616T69=* ThermoLyn Soft over the plaster positive and Tube Dummy (Fig. 1).
- 2) Attach matching 743Y42=* Lamination Dummy onto the 743Y41 Alignment Rod. Make sure that the neck of the Lamination Dummy is mounted in the direction of the inner socket.
- 3) Screw the 743Y41 Alignment Rod into the 99B83=* Tube Dummy.

4.4 Fabrication of the Shaped Part for the Outer Socket

- 1) Determine the position of the Lamination Dummy on the 743Y41 Alignment Rod. For this purpose, move the 743Y42=* Lamination Dummy on

the 743Y41 Alignment Rod until the upper edge reaches approx. 5 mm beyond the required forearm length to be achieved (Fig. 2). To make it easier to find the stud bolt of the 743Y42=* Lamination Dummy after the foaming process, position the stud bolt in the same direction as the angled portion of the 743Y41 Alignment Rod. Fix the 743Y42=* Lamination Dummy with the stud bolt.

- 2) Pull a 99B71=* PVC Bag over the inner socket and 99B83=* Tube Dummy (Fig. 2). Tie off the 99B71=* Bag to the vacuum tube above the 99B83=* Tube Dummy (Fig. 3). Switch on vacuum.
- 3) Glue on 627B5=* Neoprene Cell Rubber Tape (Fig. 4). Wrap 616T3=1 Trolene sheeting conically around the inner socket. Mix 617H12=* Pedilen Rigid Foam and 617P21=* Pedilen Hardener thoroughly as described in the enclosed processing instructions, and pour into the gap between 616T3=1 Trolene wrapper and model (Fig. 5). Remove 616T3=1 Trolene wrapper once the reaction time has elapsed.
- 4) Cut in the cured foam part and lay open the stud bolt of the 743Y42=* Lamination Dummy (Fig. 6). Unscrew the bolt and remove the 743Y41 Alignment Rod. Take off the foam part and grind it so that you get a smooth transition from the 743Y42=* Lamination Dummy to the inner socket (Fig. 7).

4.5 Fabrication of the Outer Socket

- 1) Cover the bore hole and stud bolt of the 743Y42=* Lamination Dummy using 636K6 Plastilin Clay and mask with adhesive tape (Fig. 8). Apply talcum powder to the shaped foam part. Pull a soaked 99B81 PVA Bag over the positive model and tie off at both ends without creating wrinkles. Switch on vacuum (Fig. 9).
- 2) Reinforce with 4 layers of Perlon stockinette. In order to achieve secure force transmission between the outer socket and inner socket during high patient activity, use 616G12 Carbon Fibre Cloth in the transition area (Fig. 10). In the area of the fixing holes, use the 507S15 Steel Washer with hole. Pull a soaked 99B81 PVA Bag over the positive model (Fig. 11). Mix 617H19=* Orthocryl Lamination Resin and 617P37=* Hardening Powder in accordance with the manufacturer's specifications. Add 617Z2=* Resin Colour Paste as required and pour the lamination resin into bag.
- 3) Tie off the 99B81 PVA bag above resin level and work the lamination resin into the reinforcement layer. Spread the lamination resin, trying to keep

the wall thickness as thin as possible. In order to create a high-quality socket, there must be no air bubbles inside the lamination resin (Fig. 12).

- 4) Mark the later contours of the socket and fixing holes on the outer socket and cut free proximally (Fig. 13). Create fixing holes for the inner and outer socket. Remove the plaster from the inner socket and pull the inner socket out of the outer socket. Remove the 99B83=* Tube Dummy out of the inner socket.
- 5) Lay open the socket at the distal end until the 743Y42=* Lamination Dummy becomes visible. Remove the 743Y42=* Lamination Dummy (Fig. 14). Grind the outer socket until you have achieved the intended contours (Fig. 15). Remove the PVA bag.
- 6) Drill holes into the inner socket for the 29C5=* Setting Nuts. Connect the inner and outer socket with 503F3=* Saucer Head Screw. Determine the dimensions of the outlet bore for the 99B13=* PVC Connection Tube. The 99B13=* PVC Connection Tube is attached to the inner socket in place of the 99B83=* Tube Dummy. Drill a hole into the outer socket using a 726W9=* Conical Drill and mount the 99B13=* PVC Connection Tube (Fig. 16).

5 Liability

Otto Bock Healthcare GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise.

The device may only be opened and repaired by authorised Ottobock technicians.

6 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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Otto Bock Healthcare Products GmbH