50S9 Malleo Direxa Stirrup

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	23
PT	Manual de utilização	28
NL	Gebruiksaanwijzing	33
SV	Bruksanvisning	39
DA	Brugsanvisning	44
NO	Bruksanvisning	49
FI	Käyttöohje	54
PL	Instrukcja użytkowania	59
HU	Használati utasítás	64
CS	Návod k použití	69
HR	Upute za uporabu	74
SK	Návod na používanie	79
RU	Руководство по применению	84
JA	取扱説明書	90
ZH	使用说明书	95

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INFORMATION

Date of last update: 2020-05-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the application of the 50S9 Malleo Direxa Stirrup ankle brace.

2 Intended use

2.1 Indications for use

The product is intended **exclusively** for orthotic fittings of the lower limbs. The brace must be used in accordance with the indications.

2.2 Indications

The indication must be determined by the physician.

- Functional treatment of ligament instability
- Post-operative treatment of ligament ruptures
- Ankle sprains
- Functional, conservative therapy of lateral ligament ruptures

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the foot area.

2.4 Further usage restrictions

The product is not suitable for use in sports.

2.5 Mechanism of Action

This brace stabilizes the ankle with lateral elements and straps. Eversion and inversion of the ankle joint are restricted without negatively affecting plantar and dorsal flexion.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

 **CAUTION**

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorized, qualified personnel.

⚠ CAUTION**Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE**Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

4.1 Selecting the Size

The brace is available in a universal size and no size selection is required.

4.2 Adaptation

The brace is delivered ready for trial fitting. Adjustment is not required.

4.3 Application

CAUTION

Improper molding or application

Injuries or damage to the brace due to overloading of the material and improper fit of the brace due to breakage of load-bearing components

- ▶ Do not make any improper changes to the brace.
- ▶ Always apply the brace according to the information in the instructions.

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

CAUTION

Incorrect use of the cooling elements

Local freezing due to direct contact with the skin

- ▶ Do not use the cooling elements directly on the skin.

- 1) For application, completely open all hook-and-loop closures and position the foot in the middle of the brace (see fig. 1).
- 2) **Optional:** Apply the textile pads (see fig. 4) or the cooling elements directly to the side elements.
- 3) Position the side elements medially and laterally on the ankle joints and close the lower hook-and-loop closure (see fig. 2).
- 4) Fasten the upper hook-and-loop closure (see fig. 3).
- 5) If necessary, adjust the lower hook-and-loop closure so the brace fits well.
- 6) Ensure that the brace fits tightly without constriction.

Cooling elements

During the acute phase, patients can use the cooling elements for cooling as long as this is comfortable for them. Otherwise use the textile pads.

- ▶ Charge the cooling elements in the freezing compartment of a refrigerator for at least **6 hours** (overnight if possible, no colder than **-18 °C**).

4.4 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Side elements

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

Textile pads

- 1) Remove the textile pads from the brace.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the textile pads in warm water at **30 °C** with standard mild detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

Cooling elements

- 1) Rinse by hand with warm water at **30 °C** as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

6.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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