

## Michelangelo Hand 8E500

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for use (qualified personnel) .....	14
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	25
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	38
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	49
<b>PT</b> Manual de utilização (Pessoal técnico) .....	61
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	72
<b>SV</b> Bruksanvisning (Fackpersonal) .....	84
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	95
<b>NO</b> Bruksanvisning (Fagpersonell) .....	106
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	117
<b>CS</b> Návod k použití (Odborný personál) .....	128
<b>JA</b> 取扱説明書 (有資格担当者) .....	140

## INFORMATION

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### INFORMATION

Last update: 2017-03-01

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The Axon-Bus prosthetic system is referred to below simply as the product.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

## 2 Product description

### 2.1 Function

The Michelangelo Hand is a myoelectrically controlled gripping component in the modular Ottobock Axon-Bus prosthetic system. The complex gripping kinematics combined with the anatomical appearance and low weight support the patient in his or her daily activities, thereby providing the highest rehabilitation value.

The Michelangelo Hand is equipped with two drives to achieve a natural motion pattern. The main drive is responsible for the gripping movements and gripping force. The thumb drive facilitates the two gripping modes 'Opposition Mode' and 'Lateral Mode'. Actively driven components are thus the thumb, index finger and middle finger. The ring finger and little finger passively follow the movements.

#### 2.1.1 Wrist unit

Starting at the neutral position, the joint can be flexed by approximately 75° in 4 locking positions; extension is 45° with 3 locking positions. Flexion and extension are performed passively.

A flexible and a rigid mode are available according to the application situation.

#### Flexible mode

Flexible mode simulates the natural movement characteristics of a relaxed wrist joint. This flexible condition closely approximates the physical movement characteristics of the hand and wrist joint.

To select the flexible mode, push the unlock lever to the stop (see fig. 2) until it engages. The joint can be moved without engaging at the ratchet positions. Pushing the lever again terminates the flexible mode and the wrist unit engages at the next available position in the rigid mode.

#### Rigid mode

Various everyday conditions require an individually adjustable wrist unit of the gripping component in the rigid mode. When the unlock lever is only pushed lightly and not to the stop (see fig. 2), the wrist unit can be moved to the desired position. When the unlock lever is released, the wrist unit engages at the next available position.

### 2.2 Combination possibilities

This product is to be combined exclusively with components in the Axon-Bus prosthetic system. Components of the Ottobock MyoBock system or components from other manufacturers cannot be used with this product.

## 3 Application

### 3.1 Indications for Use

The Michelangelo Hand 8E500 is to be used **exclusively** for exoprosthetic fittings of the upper limbs.

### 3.2 Use/Field of Application

The Michelangelo Hand 8E500 is suitable for unilateral or bilateral amputations starting with the transradial/transhumeral amputation level or, in case of dysmelia, for forearm or upper arm fittings.

### 3.3 Conditions of use

The Axon-Bus prosthetic system was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist unit (pushups, downhill, mountain biking) or extreme sports (free climbing, paragliding, etc.). Furthermore, the Axon-Bus prosthetic system should not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The prosthesis is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see Page 24).

### 3.4 Qualification

Fitting a patient with the 8E500 Michelangelo Hand may only be carried out by CPOs who have been authorised by Ottobock after completion of a corresponding training course.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 <b>CAUTION</b>
<p><b>The heading describes the source and/or the type of hazard</b></p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none"><li>&gt; E.g.: Consequence 1 in case of failure to observe the hazard</li><li>&gt; E.g.: Consequence 2 in case of failure to observe the hazard</li><li>▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.</li></ul>

### 4.3 General safety instructions

 <b>WARNING</b>
<p><b>Use of damaged power supply unit, adapter plug or battery charger</b></p> <p>Risk of electric shock due to contact with exposed, live components.</p> <ul style="list-style-type: none"><li>▶ Do not open the power supply unit, adapter plug or battery charger.</li><li>▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.</li><li>▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.</li></ul>

 <b>WARNING</b>
<p><b>Operating the prosthetic system near active implanted systems</b></p> <p>Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the prosthetic system.</p> <ul style="list-style-type: none"><li>▶ When operating the prosthetic system in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.</li><li>▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.</li></ul>

 <b>CAUTION</b>
<p><b>Independent manipulation of the product</b></p> <p>Injury due to malfunction and resulting unexpected prosthesis actions.</p>

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

#### **⚠ CAUTION**

##### **Signs of wear on the product**

Injury due to faulty control or malfunction of the product.

- ▶ In the interest of the patient's safety and in order to maintain operating reliability, the product should be serviced at regular intervals.

#### **4.4 Information on Proximity to Certain Areas**

#### **⚠ CAUTION**

##### **Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the prosthetic system.

#### **⚠ CAUTION**

##### **Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.50 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.35 m
  - DECT cordless phones incl. base station: 0.18 m
  - WiFi (routers, access points, ...): 0.11 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.11 m

#### **⚠ CAUTION**

##### **Wearing in extreme ambient temperatures**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see Page 24).

#### **4.5 Information on Alignment/Adjustment**

#### **⚠ CAUTION**

##### **Operator errors during the adjustment process with the adjustment software**

Injury due to unexpected product behaviour.

- ▶ Participation in an Ottobock product training course is mandatory prior to using the product. During this product training course you will receive a password giving you access to the adjustment software. Additional product training courses may become necessary to qualify for software updates.
- ▶ Do not share your unlock PIN.
- ▶ Use the online help which is integrated into the software.

### **⚠ CAUTION**

#### **Incorrect electrode settings**

Injury due to unexpected product behaviour.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.
- ▶ Set the electrode gain as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

### **INFORMATION**

Do not use silicone spray when putting on the prosthetic glove. This impairs the secure fit of the glove.

Observe the instructions for use provided for the prosthetic glove and the process for putting it on and taking it off described in the same.

## **4.6 Information on Use**

### **⚠ CAUTION**

#### **Improper use**

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Instruct the patient in the proper use of the product.

### **⚠ CAUTION**

#### **Mechanical stress on the product**

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

### **⚠ CAUTION**

#### **Operation outside of the permissible temperature range**

Injury due to faulty control or malfunction of the product.

- ▶ Avoid operation in areas outside the permissible temperature range (see Page 24).

### **⚠ CAUTION**

#### **Penetration of dirt and humidity into the product**

Injury due to unexpected behaviour of the prosthesis system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

### **⚠ CAUTION**

#### **Changing Axon-Bus gripping components when switched on**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Power down the Axon-Bus prosthetic system by pressing the button in the charging receptacle before changing any Axon-Bus components (e.g. Axon-Bus gripping component).

### **⚠ CAUTION**

#### **Unintentional unlocking of the Axon-Bus gripping component**

Injury due to releasing the Axon-Bus gripping component from the forearm (e.g. while carrying objects).

- ▶ The two release buttons should only be activated in order to change the Axon-Bus gripping component, deliberately and with due consideration of the respective situation.

### **⚠ CAUTION**

#### **Storing the 8E500 Michelangelo Hand in the closed position**

Injury due to faulty control or malfunction because of damage to the sensors and mechanism.

- ▶ Store the 8E500 Michelangelo Hand in the open state (neutral position) only.

### **NOTICE**

#### **Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Use only the following products for cleaning/disinfecting the inner socket:  
**Cleaning:** 453H10=1 Ottobock DermaClean  
**Disinfection:** 453H16 Ottobock DermaDesinfect.

### **INFORMATION**

Switching off the Axon-Bus prosthetic system during extended passive pauses (e.g. air or rail travel, visit to a theatre or cinema, etc.) will make the battery last longer before requiring a recharge. Only the entire Axon-Bus prosthetic system with all connected components can be switched off. Individual Axon-Bus components cannot be switched off separately.

## **4.7 Information on the Power Supply/Battery Charging**

### **⚠ CAUTION**

#### **Charging the product with soiled or damaged contacts**

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

### **NOTICE**

#### **Use of incorrect power supply unit/battery charger**

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

#### NOTICE

##### **Contact of the charging plug with magnetic data carriers**

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

#### INFORMATION

The Axon-Bus prosthetic system cannot be used as long as the charging plug is connected to the charging receptacle. The Axon-Bus prosthetic system is deactivated for the duration of charging.

## 5 Delivery/Storage

### 5.1 Scope of Delivery

- 1 Michelangelo Hand 8E500
- 1 Instructions for Use '647G587'

### 5.2 Storage

Observe the specifications in the chapter "Technical Data".

## 6 Preparation for use

### 6.1 Donning the Prosthetic Glove

For daily use, the Michelangelo Hand must be worn with the AxonSkin prosthesis glove. The glove protects the mechanism against environmental influences such as moisture, dirt and dust. For information about the use (putting on and taking off) and care of the prosthesis glove, please refer to the instructions for use included with the prosthesis glove.

#### **The following prosthesis gloves are available:**

- 8S501=\* AxonSkin Natural for men (skin colour)
- 8S502=\* AxonSkin Natural for women (skin colour)
- 8S500=\* AxonSkin Visual (sheer, translucent)
- 8S511=\* AxonSkin silicone for men (skin colour)
- 8S512=\* AxonSkin silicone for women (skin colour)

### 6.2 Settings


The Michelangelo Hand settings can be adjusted by means of Bluetooth data transfer and the AxonSoft 560X500=\* software. For this purpose, a wireless connection must be established between the AxonMaster 13E500 and the PC with the help of the BionicLink PC 60X5.

For further information, see the instructions for use provided with the AxonMaster 13E500 and the AxonSoft 560X500=\* software.

## 7 Application


The Michelangelo Hand is equipped with different gripping modes.

### 7.1 Neutral Position


Natural, physiological appearance in the rest position (Fig. 1).

### 7.2 Lateral Power Grip


The thumb moves laterally to the index finger; this laterally fixates medium-sized objects when the thumb position is semi-open (Fig. 2).

### 7.3 Lateral Pinch


The thumb moves laterally to the index finger; this laterally fixates flat objects when the thumb position is closed (Fig. 3).

### 7.4 Open Palm


In the Open Palm position the thumb is positioned far palmar; this achieves a flat hand position with a completely opened thumb position (Fig. 4).

### 7.5 Opposition Power Grip


The opening width makes it possible to hold cylindrical objects with large diameters and with an open thumb position (Fig. 5).

### 7.6 Tripod Pinch


Together with the index finger and the middle finger, the thumb forms a three-point support; small objects can then be fixated with the thumb position closed (Fig. 6).

### 7.7 Finger Abduction/Adduction


By spreading the fingers, a flat, thin object (< 3 mm/< 0.12 inch) can be held between the fingertips when closing the hand (Fig. 7).

### 7.8 Turning the Michelangelo Hand On and Off

The Michelangelo Hand is turned on and off through pressing of the button in the charging receptacle on the outside of the socket (Fig. 8).

## 7.9 Open the Axon-Bus Gripping Component in an Emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open. A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

## 8 Maintenance

### 8.1 Service Information

Since all moving mechanical parts are subject to wear and tear, regular service is required within the guarantee period. In the course of service, Ottobock Service (Ottobock Myo-Service) checks the entire product. Wear parts are replaced as required. Omitting a scheduled service inspection voids the guarantee.

## 9 Legal information

### 9.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

### 9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 9.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

## 10 Appendices

### 10.1 Symbols Used

Legal manufacturer

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

SN YYYY WW NNN Serial number

### 10.2 Technical data

<b>Environmental conditions</b>	
Storage and transport in original packaging	-20 °C/-4 °F to +40 °C/+104 °F
Storage and transport without packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 80% relative humidity, non-condensing
Operation	-10 °C/+14 °F to +60 °C/+140 °F max. 80% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F

<b>General</b>	
Reference number	8E500
Opening width	120 mm/4.72 inch
Weight of the Michelangelo Hand alone without AxonRotation adapter and prosthesis glove	approx. 510 g / 18 oz
Flexion of the wrist unit	75°/4 ratchet positions
Extension of the wrist unit	45°/3 ratchet positions
Product service life	5 years

The following gripping forces and load limits apply only when the battery of the Axon-Bus prosthetic system is fully charged and at room temperature.

<b>Maximum gripping forces</b>	
Gripping force "Opposition Mode"	Approx. 70 N
Gripping force "Lateral Mode"	Approx. 60 N
Gripping force "Neutral Mode"	Approx. 15 N

<b>Load limits</b>	
Maximum vertical load on the palm with locked wrist joint (e.g. when holding a sphere)	10 kg / 22.1 lbs

<b>Load limits</b>	
Maximum load on the actively driven fingers (index finger, middle finger) in the fully open hand position (e.g. holding a plate)	10 kg / 22.1 lbs
Maximum load on the actively driven fingers (index finger, middle finger) in the closed hand position (e.g. carrying bags)	20 kg / 44.1 lbs
Maximum vertical load on the knuckles (e.g. bracing oneself on the fist)	60 kg / 132 lbs
Weight of objects (typical diameter 19 mm/0.75 inch) before they slide out of the hand ("Power Grip" grip type)	18 kg / 39.6 lbs

<b>Prosthesis battery</b>	
Battery type	Li-Ion
Battery service life	2 years
Charging time until battery is fully charged	4 hours
Behaviour of the product during the charging process	The product is non-functional
Operating time of the product with fully charged battery	1 day with average use

### 10.3 Accessories

The following additional components are necessary for operating the 8E500 Michelangelo Hand:

- 1 pc. 757B500 AxonEnergy Integral or 757B501 AxonEnergy Integral
- 1 pc. 757L500 AxonCharge Integral
- 1 pc. 13E500 AxonMaster
- 1 pc. 560X500=\* AxonSoft
- 1 pc. Adapter for hand rotation, e.g. 9S501 AxonRotation Adapter
- 1 pc. 8S501=\*/8S502=\* AxonSkin Natural or 8S500=\* AxonSkin Visual
- 1 pc. 8S511=\*/8S512=\* AxonSkin silicone

MyoBock components such as suction socket electrodes, electrode cables, pull-in tube and suction valve.

### 10.4 Glossary

#### 10.4.1 Axon

The term "Axon" stands for **A**daptive **eX**change **O**f **N**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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**The 8E500 Michelangelo is covered by the following patents:**

Australia:	AU 2006 332 253; AU 2006 332 292; AU 2006 332 315; AU 2006 332 317; AU 2006 332 318; AU 2006 332 316;
Canada:	CA 2 631 970; CA 2 631 982; CA 2 632 241; CA 2 632 551; CA 2 632 240; CA 2 676 196; CA 2 678 987; CA 2 631 966
China:	CN 101 340 864; CN 101 340 865; CN 101 340 866; CN 101 340 867; CN 101 346 107; CN 101 346 109; CN 101 681 704
Japan:	JP 4 843 055; JP 5 242 409; JP 5 123 206; JP 5 074 414; JP 5 155 183; JP 5 389 667; JP 5 074 415
Mexico	MX 285 670; MX 286 932; MX 291 871; MX 291 872; MX 292 784
Russia:	RU 2 387 412; RU 2 414 871; RU 2 416 379; RU 2 423 952; RU 2 427 348; RU 2 429 803; RU 2 469 429
South Korea	KR 101 131 692; KR 101 169 834; KR 101 178 679; KR 101 178 646; KR 101 509 265; KR 101 509 264; KR 101 353 867; KR 101 265 934
Taiwan:	R.O.C. Invention Patent No. I421884
USA:	US 7 867 287; US 8 016 969; US 8 257 446; US 8 188 835; US 8 579 991; US 8 663 339; US 8 690 963
European Patent	EP 1962732 in AT, DE, FR, GB, IT, NL, SE, TR; EP 1962734 in AT, DE, FR, GB, IT, NL, SE, TR; EP 1962737 in AT, DE, FR, GB, IT, NL, SE, TR; EP 1962738 in DE, FR, GB, TR; EP 1971297 in AT, DE, FR, GB, IT, NL, SE, TR; EP 2129340 in DE, FR, GB, IS, IT, NL, SE, TR; EP 2115752 in AT, CH, DE, ES, FR, GB, IT, NL, SE, TR; EP 2528549 in DE, FR, GB, IS, PL, TR

Patents pending in Brazil, EPA, Germany, India, Mexico and USA.

**The 8E500 Michelangelo or parts thereof are covered by the following registered designs and design patents:**

Australia:	317212; 317213; 317214; 317789; 317790; 317791; 317792; 317866; 317867; 317868; 317869; 317870; 317871
Canada:	©Ottobock No. 122162; 122163
China:	ZL 200730154423.X; ZL 200730154429.7; ZL 201130050582.1; ZL 201130050654.2; ZL 200730154424.4; ZL 200730154425.9
European Design:	No.000786421; No.000786694; No.001824004
Germany:	40701345.8; 40701357.1
India:	212295; 212296; 212297; 212298; 212299; 212300; 212301; 212302; 212303
Japan:	Registered Design No. 1 365 277; 1 365 278
Russia:	69 071; 70 542
Taiwan:	R.O.C. Design Patent D 128 170; D 128 171
USA:	Patent US D 595,854; D 597,672; D 694,189; D697 030

Michelangelo is a tradename of Otto Bock Healthcare GmbH.

Michelangelo is a registered trademark in many countries of the world, beside others registered at the US Patent and Trade-mark Office, Reg.–No. 4008171.

Otto Bock Healthcare Products GmbH  
Brehmstraße 16 · 1110 Wien · Austria  
T +43-1 523 37 86 · F +43-1 523 22 64  
info.austria@ottobock.com · www.ottobock.com