


## 452A1=\*

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	5
<b>FR</b> Instructions d'utilisation .....	8
<b>IT</b> Istruzioni per l'uso .....	11
<b>ES</b> Instrucciones de uso .....	14
<b>PT</b> Manual de utilização .....	17
<b>NL</b> Gebruiksaanwijzing .....	20
<b>SV</b> Bruksanvisning .....	23
<b>DA</b> Brugsanvisning .....	25
<b>NO</b> Bruksanvisning .....	28
<b>FI</b> Käyttöohje .....	31
<b>PL</b> Instrukcja użytkowania .....	33
<b>HU</b> Használati utasítás .....	36
<b>CS</b> Návod k použití .....	39
<b>RO</b> Instrucțiuni de utilizare .....	42
<b>HR</b> Upute za uporabu .....	44
<b>SL</b> Navodila za uporabo .....	47
<b>SK</b> Návod na používanie .....	50
<b>BG</b> Инструкция за употреба .....	52
<b>EL</b> Οδηγίες χρήσης .....	55
<b>RU</b> Руководство по применению .....	58
<b>JA</b> 取扱説明書 .....	61
<b>ZH</b> 使用说明书 .....	64
<b>KO</b> 사용 설명서 .....	66

## INFORMATION

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### INFORMATION

Date of last update: 2021-12-16

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

## 1.1 Construction and Function

The 452A1 ProSeal ring is intended for the proximal sealing of a prosthetic socket for the thigh with the use of active and passive vacuum systems. The prosthetic socket has to be custom fabricated for the use of a ProSeal ring. The ring is glued into the prosthetic socket.

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.  
The product must be used in TF prostheses only.

### 2.2 Lifetime

The product is a wear part, which means it is subject to normal wear and tear.

## 3 Safety

### 3.1 Explanation of warning symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION!**

#### Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Make sure the product is applied properly and fits correctly. Improper or excessively tight application can cause pressure points and constriction of blood vessels and nerves.
- ▶ Do not continue to use the product if it does not fit properly.

#### Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may only be used by one person; reuse on another person is prohibited.
- ▶ Observe the cleaning instructions.

 **NOTICE!**

#### Risk of product damage and limited functionality

- ▶ Do not use talcum on the product or other prosthetic components. Talcum can plug the valves of vacuum systems and deplete lubricants in prosthetic joints.

## 4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	ProSeal ring	–
1	Lamination dummy	–
1	ProSeal mounting ring	–

## 5 Preparing the product for use

### 5.1 Preparing a plaster model

> **Requirement:** The socket encompasses the ischium.

- 1) Measure the circumference of the residual limb after putting on the liner, prepare the plaster cast and fabricate the plaster model.
- 2) **Based on the circumference measurements taken on the patient:** Reduce the circumference of the plaster model **5 cm** distal to the perineum by **1 cm** (see fig. 2).
  - Create a soft, round shape in the course of this reduction.
  - Do not reduce the remaining plaster model.

### 5.2 Size selection

- 1) Measure the circumference of the plaster model **3 cm** to **5 cm** distal to the perineum.
  - The ProSeal ring can be inserted into the prosthetic socket at a slight angle to the longitudinal axis.
  - Medially (in the area of the perineum), the ring should be at least **3 cm** to **5 cm** below the socket brim (see fig. 3).
- 2) Select the appropriate ProSeal ring. If the measurement is between two sizes, select the larger ring.

### 5.3 Fabricating the prosthetic socket

#### INFORMATION

**The lamination dummy may be used for a complete fitting process. Observe the information for use and dispose of the dummy after fabricating the definitive prosthetic socket.**

- ▶ **Prosthetic socket made of thermoplastic material:** Place a textile material (e.g. 99B25) under the lamination dummy during fabrication. This ensures the optimal fit of the ring in the recess of the prosthetic socket.
- ▶ **Soft-walled inner socket:** Note that the thermoplastic material of the inner socket must exhibit virtually no shrinkage.
- ▶ **Prosthetic socket with lamination technique:** Pull a moistened PVA bag over the lamination dummy.

### Positioning the lamination dummy

- 1) Align the lamination dummy. The triangular mark has to face the distal socket end.
- 2) Position the lamination dummy on the plaster model and check that it fits perfectly against the plaster model.
- 3) Fabricate the prosthetic socket.
- 4) Remove the dummy from the groove and put it aside.

**NOTICE! Dispose of the dummy after fabricating the definitive prosthetic socket.**

### Installing the ProSeal ring

- 1) Clean the ProSeal ring and the groove in the prosthetic socket with isopropyl alcohol.
- 2) Distribute silicone bonding agent in the groove in the prosthetic socket.

- 3) Align the ProSeal ring with the size marking facing in the proximal direction (see fig. 4) and insert it into the groove in the prosthetic socket.  
**TIP:** Pull the ring into the groove from distal to proximal.
- 4) Remove excess bonding agent and allow the bonding agent to cure.

## 6 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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