



17PK1=*

DE	Gebrauchsanweisung	4
EN	Instructions for use	10
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	23
ES	Instrucciones de uso	29
PT	Manual de utilização	36
NL	Gebruiksaanwijzing	42
SV	Bruksanvisning	48
DA	Brugsanvisning	54
NO	Bruksanvisning	60
PL	Instrukcja użytkowania	66
HU	Használati utasítás	72
CS	Návod k použití	79
HR	Upute za uporabu	85
TR	Kullanma talimatı	91
RU	Руководство по применению	97


1

2

743L100

INFORMATION

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INFORMATION

Date of last update: 2020-04-02

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the installation and processing of the 17PK1=* plastic orthotic knee joints.

2 Product description

2.1 Components/design

Item	Component	Reference number			
		17PK1=14		17PK1=20	
1	Countersunk head screw	29PK1=14		29PK1	
2	Lock cover				
3	Unlock cable (perlon)				

Item	Component	Reference number			
		17PK1=14		17PK1=20	
4	Compression spring	29PK1=14		29PK1	
5	Threaded sleeve		21A12		21A12
6	Lock				
7	Upper joint section				
8	Receiver for set screw				
9	M4x12 set screw	506G21=M4x12		506G21=M4x12	
10	Bearing bushing	29PK2=14		29PK2	
11	Joint bolt				
12	Countersunk head screw (17PK1=14) Cap screw (17PK1=20)				
13	Joint nut with M6 internal thread				
14	Lower joint section				
15	System bar receiver				
16	Countersunk head screw (17PK1=14) Cap screw (17PK1=20)	29PK4=14		29PK4	
17	Parallel truing tube adapter (17PK1=14 only)				
18	Parallel truing tube				
19	System bar (not included in the scope of delivery!)				

3 Intended use

3.1 Indications for use

The orthotic knee joints are to be used **exclusively** in pairs to provide orthotic treatment for a patient's lower limbs.

Reference number	Max. body weight
17PK1=14	45 kg
17PK1=20	100 kg

3.2 Indications

Partial or total paralysis of the leg muscles.
Indications must be determined by the physician.

3.3 Further Usage Restrictions

The product is not suitable for use in sports.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

 **CAUTION**

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

 **CAUTION**

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **CAUTION**

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

⚠ CAUTION**Unallowable use of lubricants**

Risk of injury and loss of functionality as well as damage to the product

- ▶ The product is designed to be free of lubricants. Avoid contact with lubricants (e.g. grease, oil).
- ▶ After contact with lubricants, promptly clean the product or parts of the product with a degreasing cleaning agent (e.g. acetone or isopropyl alcohol).

NOTICE**Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)**

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

NOTICE**Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparation for use

INFORMATION

- ▶ Check the exact parallel alignment of the system joints before delivering the orthosis.

INFORMATION

In order to facilitate the parallel alignment of the system joints, each pair of system joints is supplied with an adjustment aid. The adjustment aid consists of 2 screws, 1 parallel truing tube and, depending on the size of the system joints, 2 adapter pieces.

Due to varying alignment parameters, such as the generated forefoot leverage or the flexion position of the orthosis, the system knee joints are subjected to increased strain.

Assembling the lock

The lock is pre-assembled in the orthosis joint upon delivery. The 29PK1* set is required to replace the lock. Insertion aids are attached to the locks so that trying them out is easier.

- 1) Insert the bearing bushings (10) into the upper joint section (7).
- 2) Slide the lower joint section (14) into the upper joint section (7) and secure it with the joint bolt (11), the joint nut with internal thread (13) and the cap screw (12).
- 3) Insert the lock (6) in size 3 into the lock channel of the joint up to the stop.
- 4) Check the play of the joint.

INFORMATION: If the joint has play or the lock cannot be inserted up to the stop, exchange the lock for a different size (1, 2, 4 or 5).

- 5) Break the chosen lock off the insertion aid.
- 6) Screw the threaded sleeve (5) onto the unlock cable (3).
- 7) Insert the unlock cable (3) with threaded sleeve (5) into the chosen lock (6).
- 8) Guide the unlock cable (3) through the compression spring (4) and lock cover (2).
- 9) Insert the lock (6) and the compression spring (4) with the unlock cable (3) into the lock channel.
- 10) Attach the lock cover (2) in the upper joint section (7) with the counter-sunk head screw (1).

Installing the system bar in the system knee joint

- 1) Insert the system bar (19) into the system bar receiver (15) and mark the bore on the system bar (19).
- 2) Remove the system bar (19) from the system bar adapter (15) and drill a **Ø 4 mm** bore in the system bar (19).
- 3) Secure the system bar (19) and the system knee joint with the set screw (9).

Gluing the system bar

For stability reasons, the system bar (19) must be glued into the system bar receiver (15) of the system knee joints using 636W28 Ottobock special adhesive.

- 1) Use a degreasing cleaner (such as 634A1) to degrease the bonding surfaces of the system bar (19) and the system bar receiver (15).

- 2) Mix the special adhesive and hardener thoroughly and apply it to the components being glued.
- 3) Insert the system bar (19) and secure it with the set screw (9). The final bonding strength of the adhesive is reached after **16 h**.
- 4) Prior to delivering the orthosis:
Secure the joint bolt (11) and the joint nut (13) with 636K13 Loctite 241.
Verify the proper functionality of the product.

Removing the system bar

- 1) Unscrew the set screw (9).
- 2) Tap out and remove the system bar (19).

6 Cleaning

After contact with water containing salt, chlorine or soap, or when they get dirty, the system joints must be cleaned promptly:

- 1) Rinse the system joints with clean fresh water.
- 2) Dry the system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

The system joints must be cleaned promptly after coming into contact with lubricants (e.g. grease, oil):

- 1) Moisten a cloth with a degreasing cleaning agent (e.g acetone or isopropyl alcohol).
- 2) Wipe down the system joints or components with the cloth.

7 Maintenance

⚠ CAUTION

Use of an incorrect lock

Injuries from falling due to an incorrect lock being installed

- ▶ Only use locks of the same size when reinstalling or replacing.

ℹ INFORMATION

The lock is chosen according to the joint and may not be used for installation on other joints.

- ▶ Observe the information in the service set when replacing parts.

The manufacturer recommends a semi-annual inspection to verify functionality and check for wear of the system joints.

Spare parts are listed in the section "Components/design".

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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