


1C61 Triton Vertical Shock Prosthetic Foot


| | | |
|-----------|---------------------------------|-----|
| DE | Gebrauchsanweisung | 7 |
| EN | Instructions for use | 17 |
| FR | Notice d'utilisation..... | 27 |
| IT | Istruzioni d'uso | 37 |
| ES | Instrucciones de uso | 47 |
| PT | Manual de instruções..... | 57 |
| NL | Gebruiksaanwijzing | 67 |
| SV | Bruksanvisning | 77 |
| DA | Brugsanvisning | 86 |
| NO | Bruksanvisning | 95 |
| FI | Käyttöohje | 104 |
| PL | Instrukcja obsługi | 114 |
| HU | Használati utasítás..... | 124 |
| CS | Návod k použití | 133 |
| RO | Instrucțiuni de utilizare | 142 |
| HR | Upute za korištenje..... | 152 |
| SL | Navodila za uporabo | 162 |
| SK | Návod na použitie..... | 171 |
| BG | Инструкция за употреба | 181 |
| TR | Kullanım kılavuzu..... | 191 |
| EL | Οδηγίες χρήσης..... | 200 |
| RU | Руководство по применению..... | 211 |
| JA | 取扱説明書 | 222 |
| ZH | 使用说明 | 231 |
| KO | 사용설명서 | 239 |

| | ■ | ● | ▲ |
|----|------------------------|---------------------------------|----------------------------------|
| DE | Einzelteile | Einzelteile-Pack | Mindestmenge |
| EN | Single components | Single-Component Pack | Minimum order quantity |
| FR | Pièces à l'unité | Kit de pièces de rechange | Quantité minimum |
| IT | Singoli componenti | Imballo dei singoli componenti | Quantità minima |
| ES | Componentes | Kit componentes | Cantidad mínima |
| PT | Componentes | Conjunto de peças de reposto | Quantidade mínima |
| NL | Onderdelen | Onderdelenpakket | Minimikvantitet |
| SV | Separata delar | Service-Set | Minimum aantal |
| DA | Komponenter | Komponentpakke | Mindste mængde |
| PL | Podzespoły | Zestaw naprawczy | Ilość minimalna |
| HU | Alkatrészek | Alkatrész-csomag | Minimum mennyiség |
| CS | Jednotlivé díly | Sada jednotlivých dílů | Minimální množství pro objednání |
| RO | Componente individuale | Set componente individuale | Cantitate minimă |
| HR | Pojedinačni dijelovi | Pakiranje pojedinačnih dijelova | Minimalna količina |
| TR | Münferit parçalar | Münferit parça ambalajı | Minimum miktar |
| EL | Μεμονωμένα εξαρτήματα | Σετ μεμονωμένων εξαρτημάτων | Ελάχιστη ποσότητα |
| RU | Отдельные детали | Комплект деталей | Минимальное количество |
| JA | 個別販売 | 交換パーツセット | 最低注文個数あり |
| ZH | 单个部件 | 维修组件 | 最低起订量 |

INFORMATION

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INFORMATION

Last update: 2020-06-30

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

Explanation of symbols

CAUTION Warning messages regarding possible risks of accident or injury.

NOTICE Warning messages regarding possible technical damages.

INFORMATION Additional information on the fitting/application.

1 Single components

1.1 Scope of Delivery

The 1C61 Triton Vertical Shock is delivered with the following components:

1. 1C61 Triton Vertical Shock foot module (incl. 4X260 Functional Ring*)
2. 4X260* Functional ring (replacement)
3. 2F60 Heel Wedge Kit
4. SL = Spectra-Sock Protective Sock
5. 647G674 Instructions for Use
6. Pre-compression kit

* Additional versions of the 4X260 Functional Ring can be ordered separately (for selections refer to Section 4.3.1).

1.2 Accessories ● (Fig. 1)

- (1) Footshell (includes connection cap)

2C6=*S Slim Footshell

2C6=*N Normal Footshell

1.3 Single Components ■ (Fig. 1; for reordering)

- (2) Connection Cap for Footshell (can be ordered individually)

2C20= Slim Connection cap

2C19= Normal Connection cap

- (3) SL = Spectra Sock Protective Sock

1.4 Single component pack ● (Fig. 1; for reordering)

- (4) 2F60 = Triton Heel Wedge Kit

2 Product description

2.1 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2.2 Indications for use

The 1C61 Triton Vertical Shock is a Carbon Fiber Foot with integrated Shock Absorber and Torsion unit.

The product is to be used solely for lower limb prosthetic fittings.

2.3 Service life

INFORMATION

This component has been tested according to ISO 22675 standard for two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks.

2.4 Field of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

1C61 = 21 – 24 cm

Recommended for amputees with **Mobility Grade 3 and 4** (unrestricted outdoor walkers and unrestricted outdoor walkers with particularly high demands).

Approved for patient weight of up to 100 kg (220 lbs).

1C61 = 25 – 30 cm

Recommended for amputees with **Mobility Grade 4** (unrestricted outdoor walkers and unrestricted outdoor walkers with particularly high demands).

Approved for patient weight of up to 125 kg (275 lbs).

1C61 = 25 – 30 cm

Recommended for amputees with **Mobility Grade 3** (unrestricted outdoor walkers).

Approved for patient weight of up to 150 kg (330 lbs).

Please refer to the selection chart below to determine the appropriate stiffness required according to Ottobock recommendations.

Stiffness chart (spring – functional ring)

| Patient Weight | | Foot Size | | | | | | | | | | | | | | | | | | | | | |
|----------------|--------|---|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|---|---|
| | | 21cm | | 22cm | | 23cm | | 24cm | | 25cm | | 26cm | | 27cm | | 28cm | | 29cm | | 30cm | | | |
| kg | lbs | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | | |
| 40-47 | 88-103 | 1-0 Special Order – please contact Customer Service | | | | | | | | | | | | - | - | - | - | - | - | - | - | - | - |

Stiffness chart (spring – functional ring)

| Patient Weight | | Foot Size | | | | | | | | | | | | | | | | | | | |
|----------------|---------|-----------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|
| | | 21cm | | 22cm | | 23cm | | 24cm | | 25cm | | 26cm | | 27cm | | 28cm | | 29cm | | 30cm | |
| kg | lbs | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal | slim | normal |
| 48-55 | 104-121 | 1-1 | - | 1-1 | - | 1-1 | - | 1-1 | 1-1 | 1-1 | 1-1 | 1-1 | 1-1 | - | - | - | - | - | - | - | - |
| 56-65 | 122-143 | 2-2 | - | 2-2 | - | 2-2 | - | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 | 2-2 |
| 66-75 | 144-165 | 2-3 | - | 2-3 | - | 2-3 | - | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 | 2-3 |
| 76-87 | 166-192 | 3-4 | - | 3-4 | - | 3-4 | - | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 | 3-4 |
| 88-100 | 193-220 | 3-5 | - | 3-5 | - | 3-5 | - | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 | 3-5 |
| 101-112 | 221-247 | - | - | - | - | - | - | - | - | 4-6 | 4-6 | 4-6 | 4-6 | 4-6 | 4-6 | 4-6 | 4-6 | - | 4-6 | - | 4-6** |
| 113-125 | 248-275 | - | - | - | - | - | - | - | - | 4-7 | 4-7 | 4-7 | 4-7 | 4-7 | 4-7 | 4-7 | 4-7 | - | 4-7 | - | 4-7** |
| 126-137 | 276-302 | - | - | - | - | - | - | - | - | 5-8 | 5-8 | 5-8 | 5-8 | 5-8* | 5-8* | - | 5-8* | - | 5-8* | - | 5-8* |
| 138-150 | 303-330 | - | - | - | - | - | - | - | - | 5-9 | 5-9 | 5-9 | 5-9 | 5-9 | 5-9 | 5-9 | 5-9 | - | 5-9 | - | 5-9 |

* When combining this configuration with the C-Leg, please contact Ottobock Customer Service.

** When combining this configuration with the C-Leg or C-Leg compact, please contact Ottobock Customer Service.

2.5 Safety instructions

Please inform your patients about these safety issues!

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- Use caution when working with the product.
- If the product is damaged, check it for proper function and readiness for use.
- In case of changes in or loss of functionality, do not continue using the product (see “Signs of changes in or loss of functionality during use” in this section).
- Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer’s customer service, etc.).
- **Inform the patient.**

NOTICE

Damage due to incorrect environmental conditions.

Prosthetic components must not be exposed to the following elements:

- corrosive elements (such as fresh water, salt water, and acids)
- abrasive substances (e. g. sand)

The footshell material is not UV- or water-resistant.

Do not grind the surface of this prosthetic foot, as this can lead to premature wear.

Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.

INFORMATION

Risk of injury as a result of excessive strain.

Excessive impact loading typical with competitive sports beyond running or jumping activities, may weaken the spring elements.

To avoid the risk of accident and possible injury to the patient, the foot must no longer be used after a noticeable loss of function has occurred.

This noticeable function loss may be the result of reduced spring effect or forefoot resistance, in some cases accompanied with abnormal noises.

Measure: Have the prosthetic foot examined and, if necessary, replaced at a specialist workshop.

The Functional Ring is to be protected against mechanical damage. The Functional Ring must be replaced by qualified personnel if it shows any signs of damage (see Sections 4.3 and 4.6).

2.6 Function

2.6.1 Description of the Components (Fig. 2)

- (A) Carbon Fiber Forefoot Spring
- (B) Carbon Fiber Heel Spring
- (C) High Performance Polymer Base Spring
- (D) Carbon Fiber Attachment Spring
- (E) Aluminum Housing
- (F) Optional Heel Wedge
- (G) Functional Ring
- (H) Shaft Assembly

2.6.2 Carbon Fiber Foot with integrated Shock Absorber and Torsion unit

The main features of the 1C61 Triton Vertical Shock carbon fiber foot are the lightweight construction and the integrated shock absorber with rotation functionality. This system permits amputees with MG 3 or 4 the ability to walk effortlessly in everyday life situations including participation in high energy activities.

The integrated Shock Absorber and Torsion unit helps to reduce vertical and torsional forces in the extremities during walking. The result is an improvement for both the function and loading situation of the prosthesis wearer.

The 1C61 Triton Vertical Shock is designed to be used with the footshell and the detachable Foam Connection Cap. The Foam Connection Cap provides effective connection between the foot and cosmetic skin foam cover.

The combination of the carbon fiber foot together with the integrated Shock absorber and Torsion unit offers the following functional advantages:

- Secure heel strike with noticeable plantar flexion
- Progressive movement of the ankle in the mid-stance phase for the natural rollover of the foot
- Dynamic transition from the stance to swing phase
- Smooth dynamic forefoot response and energy return
- Compensation for uneven terrain
- Axial shock absorption
- Compensation of torsion moments

3 Technical data

| | | | | | | | | | | |
|---|----------------|-----|-----|-----|----------------|-----|------|------|------|------|
| Sizes (cm) | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 |
| Heel height with slim footshell (mm) | 15 | 15 | 15 | 15 | 15 | 15 | 15 | – | – | – |
| Heel height with normal footshell (mm) | – | – | – | 10 | 10 | 10 | 10 | 10 | 10 | 10 |
| System height with slim footshell (mm) | 163 | 164 | 166 | 167 | 175 | 177 | 177 | – | – | – |
| System height with normal footshell (mm) | – | – | – | 173 | 175 | 177 | 181 | 183 | 189 | 191 |
| Weight with slim footshell (g) | 790 | 800 | 860 | 890 | 940 | 960 | 1020 | – | – | – |
| Weight with normal footshell (g) | – | – | – | 900 | 950 | 980 | 1030 | 1050 | 1110 | 1140 |
| Body weight MG3 up to | 100 kg/220 lbs | | | | 150 kg/330 lbs | | | | | |
| Body weight MG4 up to | 100 kg/220 lbs | | | | 125 kg/275 lbs | | | | | |
| Mobility Grade | 3+4 | | | | | | | | | |

4 Handling

⚠ CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- Observe the alignment and assembly instructions.

4.1 Setup

4.1.1 Footshell application/removal

NOTICE

Risk of damage as a result of improper use.

The foot is designed to be used with the footshell. Please avoid frequent removal from and insertion into the footshell as this can damage the carbon fibre elements.

NOTICE

Risk of damage as a result of improper footshell assembly.

Do not use a screwdriver for assembly and disassembly as this can cause damage to the carbon elements and/or footshell. Ottobock recommends use of 2C100 Footshell Insertion/Removal tool (accessories).

NOTICE

Risk of damage as a result of improper use.

The footshell must not be dusted with talcum powder as this would increase abrasion. Use silicone spray for lubrication if needed. Replace the footshell in case it shows any damage or wear.

Inserting the foot into the footshell:

- Before inserting into the footshell, attach the tube adapter to the prosthetic foot. This facilitates the procedure.
- Use the 2C100 Footshell Insertion/Removal tool provided by Ottobock to insert the carbon foot into the footshell.
- Make sure that the heel of the polymer base spring locks in under the footshell heel tab.

Removing the footshell:

- Insert the 2C100 Footshell Insertion/Removal tool under the heel of the polymer base spring which is located under the heel tabs in the posterior part of the footshell and pry the heel out of the shell. It is then possible to pull the foot out of the footshell.

4.1.2 Bench alignment

The prosthetic foot with cover is clamped in the alignment apparatus.

TT-fitting

| Foot size [cm] | Mid-foot anterior to the alignment reference line [mm] | Heel height in alignment apparatus |
|----------------|--|---|
| 21 – 25 | 25 | Effective heel height (shoe heel height – sole thickness of forefoot) + 5mm |
| 26 – 28 | 30 | |
| 29 – 30 | 35 | |

TF-fitting

Observe the alignment instructions for the chosen prosthetic knee joint.

| Effective shoe heel height [mm] | Suitable footshell |
|---------------------------------|--------------------|
| 5 – 15 | Normal |
| 10 – 20 | Slim |

4.1.3 Static alignment optimization

If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

Use the L.A.S.A.R. Posture for static alignment optimisation of the prosthesis, if available.

4.1.4 Dynamic alignment optimization

Adapt the prosthetic alignment in the frontal plane (ML) and in the sagittal plane (AP) by angle adjustments and/or sliding adjustments, in order to ensure correct heel contact, easy roll-over, and optimal shifting of weight onto the contralateral side. For transtibial amputees please make sure that there is physiological knee flexion during stance phase.

4.2 Using the optional heel wedges

Risk of falling due to change of the knee dynamics when using heel wedges.

The use of heel wedges can facilitate stance phase flexion when bearing weight on the prosthesis. Therefore, particular caution is advised during the dynamic alignment optimization – especially when a knee joint without stance phase control is used.

4.2.1 Heel wedge application

Two heel wedges are included to each 1C61 Triton Vertical Shock prosthetic foot:

- „Soft“ Heel Wedge: translucent
- „Stiff“ Heel Wedge: graphite color

Should the patient feel the heel strike too soft or the knee in hyperextension during the mid-stance phase, the heel can be stiffened by using one of the two heel wedges.

4.2.2 Inserting a heel wedge

INFORMATION

Ottobock recommends to test the foot first without using a heel wedge.

- Orient the heel wedge (Fig. 3, Pos. F) for insertion so the locking recess (Fig. 3, Pos. L) in the wedge is positioned in the proximal and posterior direction.
- Slide the heel wedge into position from the posterior side of the prosthetic foot between the attachment spring (Fig. 3, Pos. D) and heel spring (Fig. 3, Pos. B).

INFORMATION

- The heel wedge is properly installed when the locking recess (Fig. 3, Pos. L) in the heel wedge is fully seated around the mating stop (Fig. 3, Pos. K) on the attachment spring.
- The heel wedge is designed such that no adhesive is required to affix it to the prosthetic foot.

4.2.3 Removing a heel wedge

INFORMATION

The heel wedge has a tendency to rotate during removal. Pulling the heel wedge in the direction of rotation helps make the removal easier.

- Push on the medial or lateral side of the heel wedge with your finger to free it from the mating stop (Fig. 3, Pos. K) on the attachment spring (Fig. 3, Pos. D).
- After the heel wedge is partially free grab it by the respective lateral or medial side and pull it free from the prosthetic foot.

4.3 Selecting and adjusting the functional ring

The 1C61 Triton Vertical Shock is delivered as a complete foot assembly. The stiffness category of the functional ring is synchronized with the weight category of the foot assembly. If during walking too much compression of the functional ring is observed, an exchange to a stiffer functional ring can be done.

NOTICE

Function loss as a result of improper use. The Functional Ring may only be selected/adjusted by prosthetists that have been authorized by Ottobock.

4.3.1 Selecting the right Functional Ring

Use the following selection table to select the right Functional Ring. The Functional Ring version is marked on the backside of the functional ring.

| | Body Weight | | | | | | | | | |
|-----------------------------------|-------------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| kg | 40-47 | 48-55 | 56-65 | 66-75 | 76-87 | 88-100 | 101-112 | 113-125 | 126-137 | 138-150 |
| lbs | 88-103 | 104-121 | 122-143 | 144-185 | 186-192 | 193-220 | 221-247 | 248-275 | 276-302 | 303-330 |
| Recommended Functional Ring 4X260 | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 |

Please note that the recommendations do not necessarily represent the optimal selection for the patient. When deviating from the selection data, the Functional Ring with the highest possible stiffness should be selected.

If the patient feels like they are “bottoming out”, a stiffer Functional Ring should be used.

4.3.2 Changing the functional ring

1. Remove the footshell (see Section 4.1.1)
2. Loosen the retaining screw on the distal end of the aluminum housing (Fig. 2, item E) by use of an allen wrench (**6 mm**). Remove the screw (Fig. 6, Step 1). This can be carried out from below or optionally by inserting an Allen key through the hole in the pyramid adapter.
3. Pull the shaft away from the Base (Fig. 6, Step 2).
4. Pull off the installed Functional Ring (Fig. 6, Step 3).
5. Select the desired Functional Ring (see Section 4.3.1).
6. Preparing the Functional Ring for initial use.

Before the first use, it is recommended by Ottobock that separately ordered functional rings are pre-compressed using, for example, the separately available 4X247 Pre-Compression Tool in order to guarantee trouble-free operation. Without pre-compression, new Functional Rings can be overly stiff. If using the preinstalled Functional Ring or an already pre-compressed unit please continue with step 7.

- a. Open the Pre-compression Tool.
- b. Place the Functional Ring in the Pre-compression Tool (Fig. 9)
- c. Completely compress the Functional Ring by tightening the screw completely or clamping in a vise. (Fig. 10)
- d. Compression time: Minimum **3 minutes** to a maximum of 15 minutes
- e. Open the Pre-compression Tool, remove the Functional Ring.
7. Slide the Functional Ring over the Shaft (Fig. 6, Step 4).
8. **Caution** - To prevent unintended loosening of the screwing and thus the risk of falling, clean the threads of the prosthetic foot and the retaining screw so that they are clean and free from grease and other contaminants.
9. Reassemble the prosthetic foot (Fig. 6, Step 5 & 6). Secure the Retaining Screw with Loctite 243 or 246 and tighten with a **6 mm** allen bit in a torque wrench to **15 Nm (133 in-lbs)**; This will need to be done by inserting the **6 mm** allen socket going through the top hole in pyramid.
10. Insert the foot into the footshell (see Section 4.1.1).

4.4 Foam cover

Both Footshell Connection Caps 2C20 (slim) and 2C19 (normal) serve as a detachable connection between the foam cover and the footshell. They snap into the proximal brim of the footshell.

Please observe the following recommendations:

1. For mounting to the footshell the connection caps of the sizes 21-25 have to be cut along the cutting mark (Fig. 8).
2. Cut the foam cover to length, allowing for distal compression about 0.4 in (10 mm) for transibial and about 1.6 in (40 mm) for transfemoral covers.
3. With the connection cap attached to the footshell, pull the foam cover over the prosthesis, reassemble the foot and mark the outer contours of the connection cap.
4. Remove the connection cap and clean the foam cover with Isopropyl Alcohol 634A58.
5. Apply 636N9 Ottobock Contact Adhesive or 636W17 Plastic Adhesive to the distal end of the foam cover and to the connection cap.
6. Allow the adhesive to dry (approx. 10 minutes) and cut out a space in the distal foam cover with a grinding and milling machine to allow for precise fit for the upper part of the adapter. Sufficient space should be available for the adapter to ensure proper function and eliminate noise during walking.
7. Assemble the foam cover and finish off the outer shape of the foam cover in usual way with carbon foot installed. Any compression of the foam cover by cosmetic stockings or SuperSkin needs to be taken into account.

4.5 Adjustment and final assembly

Due to the high forces in the adapter section, it is indispensable to tighten the adjustment screws of the tube adapter with a torque of **15 Nm (135 lb-in)**.

Use 710D1 Torque Wrench. Once you have attained the final setting for the prosthesis, use 636K13 Loctite to secure the set screws.

4.6 Maintenance instructions

- Foot assembly should be inspected after first 30 days (4 weeks) of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Footshell may require replacement if wear is excessive.

Ottobock recommends replacing the 4X260 Functional Ring after approximately 1 year of use because the stiffness and restoring force can decrease over time. Depending on activity and environmental conditions, this may become necessary at an earlier (or later) point in time. Reduced shock absorption is a typical indicator the patient may notice.

Note: New Functional Rings may feel stiffer for the patient on the initial days of use.

When replacing the Functional Ring at the recommended service interval, wipe all old grease from the shaft and bearings, apply new grease, and replace the two o-rings on the shaft. They are used to cushion undesirable but necessary movement in the bearing system.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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