6A53, 6A54

DE Gebrauchsanweisung	3	CS Návod k použití.....	44
EN Instructions for use	6	RO Instrucțiuni de utilizare.....	47
FR Instructions d'utilisation.....	9	HR Upute za uporabu	50
IT Istruzioni per l'uso	12	SL Navodila za uporabo	53
ES Instrucciones de uso	15	SK Návod na používanie	56
PT Manual de utilização.....	19	BG Инструкция за употреба	60
NL Gebruiksaanwijzing.....	22	TR Kullanma talimatı	63
SV Bruksanvisning	25	EL Οδηγίες χρήσης	66
DA Brugsanvisning	28	RU Руководство по применению	69
NO Bruksanvisning	31	JA 取扱説明書.....	73
FI Käyttöohje.....	34	ZH 使用说明书.....	76
PL Instrukcja użytkowania.....	38	KO 사용 설명서.....	78
HU Használati utasítás	41		

1

2

3

4

1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2021-06-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Design and Function

The 6A53 and 6A54 Sliding Adapters permit the distal components of the modular prosthesis to be shifted in parallel, regardless of the angle adjustment by the pyramid adapters. The adjustments can be made either in the frontal plane (medial or lateral) or in the sagittal

plane (anterior or posterior). The sliding adapters are equipped with different connections:

Product reference number	Connections
6A53	pyramid adapter/pyramid adapter
6A54	pyramid adapter/pyramid receiver

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

- The maximum approved body weight is specified in the technical data (see page 9).

2.3 Environmental conditions

Storage and transport

Temperature range $-20\text{ }^{\circ}\text{C}$ to $+60\text{ }^{\circ}\text{C}$ ($-4\text{ }^{\circ}\text{F}$ to $+140\text{ }^{\circ}\text{F}$), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature range: $-10\text{ }^{\circ}\text{C}$ to $+45\text{ }^{\circ}\text{C}$ ($14\text{ }^{\circ}\text{F}$ to $113\text{ }^{\circ}\text{F}$)

Moisture: relative humidity: 20% to 90%, non-condensing

Unacceptable environmental conditions

Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water

Solids: dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

3 Safety

3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

3.2 General safety instructions

Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 6).
- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.

Risk of product damage and limited functionality

- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Adapter	–
4	Set screw (only 6A54)	506G3=M8x12

5 Preparing the product for use

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

5.1 Adjusting adapter

The adapter includes a pyramid adapter adjustable in 2 mm increments. The pyramid adapter can be aligned proximally or distally.

> **Required materials:** 710D20 torque wrench, 636K13 Loctite 241

- 1) Align the sliding adapter in the desired displacement direction (**a – p** or **m – l**).

- 2) Lift the ring on the adjustable pyramid adapter slightly in order to displace the pyramid adapter (see fig. 3).
- 3) Set the desired displacement and press the ring back down until it engages.

TIP: Mark the final position with a pen.

- 4) Connect the pyramid adapter with the applicable connecting component (see fig. 4). This fastens the pyramid adapter in its position.

5.2 Connecting the pyramid adapter and pyramid receiver

The pyramid adapter is fixed with the set screws of the pyramid receiver.

> **Required materials:** torque wrench (e.g. 710D20), 636K13 Loctite 241

- 1) **Fitting:**
Screw in the set screws.
Use the torque wrench to tighten the set screws (**10 Nm**).
- 2) **Definitive assembly:**
Secure the set screws with Loctite.
Screw in the set screws.
Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten them (**15 Nm**).
- 3) Replace any set screws that are protruding or are recessed too much with appropriate set screws (see selection chart).

Selection table for set screws	
Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

6 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

7 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

8.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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