


## 4801N Epiflex

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use.....	7
<b>FR</b> Instructions d'utilisation .....	12
<b>IT</b> Istruzioni per l'uso.....	16
<b>ES</b> Instrucciones de uso .....	21
<b>NL</b> Gebruiksaanwijzing .....	26

	Größe / Size	Ellenbogenumfang / Elbow circumference	
		cm	inch
	S	24–25	9.51–9.8
	M	26–27	10.2–10.6
	L	28–30	11.0–11.8
	XL	31–33	12.2–13.0

<b>Material</b>	Kunststoffspange/plastic support: PC, PE, SBR (Neopren/neoprene)
	Gewebe/fabric: PC, EL

## INFORMATION

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### INFORMATION

Date of last update: 2020-04-03

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.

► Please keep this document for your records.

These instructions for use provide you with important information on the handling and intended use of the 4801N Epiflex epicondylitis brace.

## 2 Intended use

### 2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the upper limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 2.2 Indications

- Epicondylitis humeri radialis (tennis elbow)
- Epicondylitis humeri ulnaris (golfer elbow)

Indications must be determined by the physician.

### 2.3 Contraindications

#### 2.3.1 Absolute Contraindications

Not known.

#### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent, swollen scars, reddening and hyperthermia in the relevant area; pronounced varicose veins, especially with impaired blood return, lymph drainage disorders and unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the forearm and hand, e.g. associated with diabetic neuropathy.

### 2.4 Effects

The brace with a semi-rigid, moulded plastic support has an analgesic effect because it relieves strain on the muscle origins of the forearm musculature (depending on the positioning of the brace, the muscle origins of the hand and finger extensors or the hand and finger flexors). It allows the user to carry out everyday activities, despite elbow complaints.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION**

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

### 3.2 General safety instructions

**CAUTION**

#### **Use of a product with neoprene**

Allergic reactions due to neoprene or heat allergies

- ▶ Do not wear products with neoprene for more than 3 – 4 hours without interruption (unless explicitly instructed otherwise by the doctor).
- ▶ Do **not** use the product in case of a known neoprene or heat allergy.

**CAUTION**

#### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**CAUTION**

#### **Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**NOTICE****Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 4 Handling

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Size Selection

- 1) Measure the elbow circumference.
- 2) Determine the orthosis size (see size chart).

The brace is suitable for either side.

### 4.2 Application

**⚠ CAUTION****Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.
- 1) Undo the hook-and-loop closure. Do not pull the strap fully out of the loop as you do this.
  - 2) Pull the plastic support over the forearm. Position the large plate of the support on the outside (see fig. 1, item 1).

- 3) Position the support approx. **5 cm** below the elbow joint (see fig. 2).
- 4) Tighten the strap so that a finger still fits between the arm and the strap (see fig. 1, item 2).
- 5) Fasten the hook-and-loop closure.

## **5 Cleaning**

- 1) Fasten all hook-and-loop closures.
- 2) Recommendation: Use a laundry bag or net.
- 3) Wash the orthosis in warm water at **30 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## **6 Disposal**

Dispose of the product in accordance with national regulations.

## **7 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **7.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### **7.2 CE conformity**

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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