




Phoen 3R62=N, 3R62=1-N


DE Gebrauchsanweisung	7
EN Instructions for use	19
FR Instructions d'utilisation	31
IT Istruzioni per l'uso	43
ES Instrucciones de uso	56
PT Manual de utilização	68
NL Gebruiksaanwijzing	81
SV Bruksanvisning	93
DA Brugsanvisning	105
NO Bruksanvisning	117
FI Käyttöohje	128
PL Instrukcja użytkowania	140
HU Használati utasítás	152
CS Návod k použití	164
HR Upute za uporabu	176
TR Kullanma talimatı	188
RU Руководство по применению	200
JA 取扱説明書	214
ZH 使用说明书	225







7

INFORMATION

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INFORMATION

Last update: 2017-03-21

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 3R62* Prosthetic Knee Joint is a polycentric prosthetic knee joint. It is available in the following versions:

- 3R62=N, 3R62=1-N (with proximal pyramid adapter; in combination with an approved tube clamp adapter for hip disarticulation)

The maximum flexion angle of the prosthetic knee joint is **155°**.

In conjunction with the polycentric joint kinematics, the EBS (ergonomically balanced stride) feature secures the stance phase. At heel strike, prosthetic knee joint flexion is cushioned in the range from **0° to approx. 10°** by an adjustable spring damping element and stance phase extension is dampened. The EBS is adjusted to the user's body weight and activity level with the help of the adjusting wheel ② (see fig. 8).

In the swing phase, the spring action of the extension assist prevents the prosthetic foot from swinging too far posterior during flexion. The spring action of the extension assist helps the patient achieve full extension. The extension assist is adjusted with the help of the adjusting wheel ② (see fig. 10).

The swing phase is also influenced by an integrated friction brake, which prevents a hard extension stop in particular. The braking effect is adjusted by means of the adjusting screw ② (see fig. 12).

The upper joint section ① and lower joint section ② are connected by the two anterior links ③ and the posterior link ④ to form a four-bar linkage. In extension, the instantaneous centre of rotation ⑤ is located above the joint and behind the load line, thus providing knee stability in the stance phase (see fig. 2).

The distal portion of the prosthetic knee joint is fabricated as a tube clamp and accommodates the tube adapter.

1.2 Combination possibilities

INFORMATION

When fabricating a prosthesis for the lower limbs, all of the prosthetic components used have to meet the requirements in regard to the patient's body weight and activity level.

Prosthetic components that are especially well suited for combination with the product are listed below.

Designation	Reference number
Tube clamp adapter (with hip disarticulation or hemipelvectomy)	4R156* (up to 150 kg) 4R56* (up to 100 kg)
Tube adapter	2R49
Foam cover	3S107 3R27* (for hip disarticulation)

For additional combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

3R62=N

Recommended for mobility grades **1 and 2** (indoor walkers and restricted outdoor walkers). Approved for a body weight of **75 to 125 kg (165 to 275 lbs)**.

3R62=1-N

Recommended for mobility grades **1 and 2** (indoor walkers and restricted outdoor walkers). Approved for a body weight of **45 to 75 kg (100 to 165 lbs)**.

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification (see section "Area of application").

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

CAUTION

Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

CAUTION

Reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

⚠ CAUTION**Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Noticeable changes in functionality may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

4 Scope of delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

Fig.		Item no.	Quantity	Designation	Reference number
-	■	-	1	Instructions for use	647G876-1
1	-	①	1	Prosthetic knee joint	3R62*
1	■	②	1	Tube adapter	2R49 (299 mm)
1	■	③	1	Adjustment wrench	710H10=2X3
1	■	④	1	Pyramid Dome Protector Kit	2Z11=KIT
1	■	⑤	1	EBS adjusting aid	4G618

5 Preparation for use**5.1 Information on fabrication of a prosthesis****⚠ CAUTION****Incorrect alignment, assembly or adjustment**

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION**Initial use of the prosthesis by the patient**

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail) during initial standing and walking.
- ▶ Do not use devices with insufficient stability (e.g. crutches or walking sticks).

Fabricating a prosthesis that is optimally tailored to the patient demands technical knowledge of orthopaedic technology. The alignment and adjustments may only be carried out by a prosthetist.

The alignment and adjustments are mainly completed in the following steps:

1. Bench alignment
2. Static alignment
3. Dynamic trial fitting

These steps are first carried out with a test prosthesis in order to determine the best possible combination and positioning of the prosthetic components relative to each other. Once all settings have been tailored to the patient, the definitive prosthesis is produced. The same steps are repeated for this purpose.

In dynamic trial fitting, the patient must be given sufficient time to become familiar with the functions of the prosthesis, thereby learning safe use for everyday life.

5.2 Information on Using a Cosmetic Cover

CAUTION

Use of talcum

Risk of injury, damage to the product due to lack of lubrication

▶ Do not use talcum on the product or other prosthesis components.

- ▶ Only use the foam cover approved by the manufacturer (see section "Combination Possibilities" – see Page 20).
- ▶ To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.
- ▶ After fabricating the cosmetic cover, repeat the dynamic trial fitting since the cover affects the adjustments of the prosthesis (e.g. damping).

5.3 Bench Alignment

Bench alignment process	
+ = anterior placement / - = posterior placement (in relation to the alignment reference line)	
Item	see fig. 6
	Required equipment: Alignment apparatus (e.g., 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly), 743A80 50:50 gauge, 743S12 heel height measuring apparatus
	Position the prosthetic foot in the alignment apparatus in accordance with the following values:
①	Heel height: effective heel height (x) + 5 mm
②	a-p positioning, middle of the prosthetic foot to the sagittal alignment reference line ①: +30 mm
	External foot rotation: approx. 7°
	m-l positioning, middle of the prosthetic foot to the frontal alignment reference line ⑤: 0 mm medial
③	Check the prosthetic knee joint according to the section "Checking the factory settings" (see Page 27). Position the alignment reference point ② (proximal, anterior rotation axis) of the prosthetic knee joint in the alignment apparatus according to the following values:
	Height of the alignment reference point in the alignment apparatus: MTP-floor measurement + 20 mm
	a-p positioning of the alignment reference point to the alignment reference line: 0 mm
	m-l positioning, centre of the prosthetic knee joint to the frontal alignment reference line ⑤: 0 mm
	External rotation of the prosthetic knee joint: approx. 5°
④	INFORMATION: If the length-adjustable 2R45=34 tube adapter (max. body weight 125 kg) is used in the alignment of the test prosthesis, the required length for the definitive prosthesis is easier to determine.

Bench alignment process	
+ = anterior placement / - = posterior placement (in relation to the alignment reference line)	
Item	see fig. 6
	<p>When using the tube adapter, note the section "Adapting the length of the tube adapter" (see Page 24).</p> <p>Connect the prosthetic knee joint to the prosthetic foot and chosen adapters according to their instructions for use and the section "Installing the tube adapter" (see Page 25).</p>
5	<p>Position the 743A80 50:50 gauge on the lateral side of the prosthetic socket and mark a point respectively in the distal and proximal area.</p> <p>Using these points, mark the lateral centre line of the prosthetic socket.</p> <p>Mark the tuberosity height as the tuberosity reference point ④ on the centre line.</p> <p>Clearly mark the lateral socket reference point ③ 30 mm proximal to the tuberosity reference point.</p> <p>Position the tuberosity reference point ④ of the prosthetic socket at the height of the tuberosity-floor measurement.</p>
6	<p>CAUTION! The correct position and flexion of the prosthetic socket is crucial for the safety of the patient and functionality of the prosthetic knee joint, and avoids premature wear.</p> <p>INFORMATION: The 4R118 adapter plate (max. body weight 125 kg) can be used for further posterior placement.</p> <p>Position the prosthetic socket so the lateral socket reference point ③ lies on the sagittal alignment reference line ①.</p> <p>Adjust the socket flexion angle β depending on the hip flexion contracture (angle α) according to the following specifications:</p> <p>If the angle $\alpha = 0^\circ$, the optimum value of the angle $\beta = 3$ to 5° to the alignment reference line.</p> <p>If the angle $\alpha > 0^\circ$, the optimum value of the angle $\beta = \alpha + 5$ to 10°.</p>
7	<p>Position the socket adapter so that all prosthetic components between the socket and prosthetic knee joint can be connected to each other correctly.</p> <p>Join the prosthetic socket and the prosthetic knee joint using the adapters selected.</p> <p>Be sure to follow the adapter instructions for use during adjustment and installation.</p> <hr/> <p>Position the 743A80 50:50 gauge on the anterior side of the prosthetic socket and mark a point in the distal and proximal areas, respectively.</p> <p>Using these points, mark the anterior centre line of the prosthetic socket.</p> <p>Clearly mark the anterior socket reference point ⑥ (intersection of the socket edge and anterior centre line).</p> <p>Position the prosthetic socket so the anterior socket reference point ⑥ lies on the alignment reference line ⑤.</p> <p>Adjust the socket flexion angle γ according to the patient's adduction angle.</p>

Use 2Z11=KIT

INFORMATION: The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

5.3.1 Shortening the Tube Adapter

CAUTION

Incorrect processing of tube

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

> **Tools required:**

719R3 tube cutter, 718R1 tube deburrer

- 1) **CAUTION! The length of the tube must not be less than the required length for the patient, since length compensation in the insertion area of the prosthetic knee joint is prohibited.**

Use the tube cutter to shorten the tube to the required length (see fig. 3).

- 2) Carefully deburr the inside and outside of the tube with the tube deburrer (see fig. 4).

5.3.2 Installing the Tube Adapter

⚠ CAUTION

Incorrect mounting of the tube

Risk of injury due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthetic component when mounting.

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

Installation on the Prosthetic Knee Joint

> **Required tools and materials:**

710D4 Torque Wrench, degreasing cleaner (verify material compatibility)

- 1) Loosen the cap screw by **approx. 2 turns** (see fig. 5).
- 2) Clean the contact surfaces of the prosthetic knee joint and tube adapter using a degreasing cleaning agent.
- 3) Slightly rotate the tube as you insert it into the prosthetic knee joint to the stop (see fig. 5).
- 4) Using the torque wrench, tighten the cap screw (**10 Nm**) (see fig. 5).

Installation on the Prosthetic Foot

> **Required tools and materials:**

710D4 Torque Wrench, 636K13 Loctite®

- 1) Position the pyramid receiver of the tube adapter on the pyramid adapter of the prosthetic foot.
- 2) **For final assembly:**
Use Loctite® to secure the set screws of the pyramid receiver.
- 3) Screw in the set screws.
- 4) Use the torque wrench to tighten the set screws (**10 Nm**).
- 5) Replace any set screws that are protruding or recessed too much with suitable ones (see selection table).
- 6) **For final assembly:**
Use the torque wrench to retighten the set screws (**15 Nm**).

Selection table for set screws

Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

5.4 Static Alignment

Static alignment process	
+ = anterior placement / - = posterior placement (in relation to the load line)	
No.	see fig. 7
	Required equipment: 743L100 L.A.S.A.R. Posture
①	To determine the sagittal load line ①, position the patient on the L.A.S.A.R. Posture as follows: <ul style="list-style-type: none"> • Prosthetic foot (with shoe) on force measurement plate (sufficient load: > 35% body weight) • Other foot (with shoe) on height compensation plate • The tips of the shoes are in one line
②	Optimise the static alignment solely by changing the plantar flexion. Adjustment is performed only with the anterior and posterior set screws of the tube adapter on the prosthetic foot. <p>a-p positioning of the alignment reference point ② (distal, anterior rotation axis) to the sagittal load line ①: -10 mm</p> <p>m-l positioning, middle of the prosthetic foot to the frontal load line ④: 1 to 2 cm medial</p> <p>m-l positioning, centre of the prosthetic knee joint to the frontal load line ④: 1 to 2 cm medial</p> <p>m-l positioning of the prosthetic socket: Positioning of the anterior superior iliac spine (Spina iliaca anterior superior) ③ to the frontal load line ④: 1 to 2 cm lateral</p>

5.5 Dynamic Trial Fitting

⚠ CAUTION

Adjusting the Settings

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

INFORMATION

Established initial tension values are used for the settings of the prosthetic knee joint at the factory.

- ▶ If a reduction of the initial tension is required, ensure that the initial tension is always free of play.
- ▶ When initial tension for a setting is reduced to a minimum, ensure constant freedom of play for the setting (e.g. by establishing a maintenance interval to verify the adjustment values). This is because the adjustment values can change during use of the prosthesis.

During dynamic trial fitting, the alignment and settings of the prosthesis are checked and adapted for optimum walking according to the needs and abilities of the patient.

The patient has to learn the safe use of the prosthesis through intensive training.

The following subsections describe the adjustment possibilities for adapting the product to the patient.

The list below provides an overview of the sequence for working through the subsections:

- Checking the Factory Settings
- Adjusting the EBS
- Adjusting the extension assist
- Adjusting the swing phase friction

- Fine-tuning the settings

5.5.1 Checking the Factory Settings

Function	Means of adjustment	Factory setting	Meaning
EBS	Adjusting wheel (see fig. 9 - ②)	Spring with adjusting wheel, 3 mm initial tension	Low initial tension as the starting point for adjustment to the patient
Extension assist	Adjusting wheel (see fig. 11 - ②)	Spring with adjusting wheel, 1.5 mm initial tension	Low initial tension as the starting point for adjustment to the patient
Swing phase friction	Adjusting screw (see fig. 13 - ②)	Adjusting screw tightened to 2 Nm	Standard braking effect as the starting point for adaptation to the patient

1) **CAUTION! Do not adjust the factory settings prior to the dynamic trial fitting.**

- Note the position of the adjusting screw ② upon delivery as well as all changes so that the initial position can be restored if required. Use the markings on the adjusting screw and the plastic cover ① for this purpose (see fig. 13).
- 2) If the settings have been changed, restore the factory settings (see table) prior to the dynamic trial fitting.

5.5.2 Adjusting the EBS

INFORMATION

The EBS function must be explained to the patient in detail. Through intensive training (active loading and unloading of the prosthesis side while standing between parallel bars), the patient must become familiar with the functionality of the EBS and find the optimum setting by trial and error.

INFORMATION

The EBS adjusting aid is stuck on the side of the prosthetic knee joint.

The EBS adjusting aid is used to visually check EBS flexion during the training with the patient. With unloaded EBS, the EBS adjusting aid is roughly parallel to the longitudinal axis of the tube adapter. Under load the EBS adjusting aid turns away from the longitudinal axis of the tube adapter.

If the adjusting aid does not show any movement when there is load on the EBS and at heel strike, then the initial tension is too high and must be reduced.

If there is no movement visible despite minimum initial tension, then the heel strike is not correct and must be practised with the patient until he or she has learned how to load the EBS correctly.

Function	Means of adjustment	Setting	Meaning
EBS (controlled stance phase movement between 0° and approx. 10°)	Adjusting wheel (see fig. 8 - ②)	Rotation to the left ③	Decrease initial tension
		Rotation to the right ④	Increase initial tension

> **Required equipment and tools:**

Parallel Bars, 710H10=2x3 Adjustment Wrench

▶ **CAUTION! Explain the EBS function to the patient and familiarise the patient with the EBS between parallel bars.**

With the adjustment wrench, adjust the EBS so that the prosthetic knee joint goes into stance phase flexion at heel strike, not too fast and not too slow in accordance with the needs of the patient.

5.5.3 Adjusting the Extension Assist

Function	Means of adjustment	Setting	Meaning
Extension assist	Adjusting wheel (see fig. 10 - ②)	Rotation to the left ③	Decrease initial tension
		Rotation to the right ④	Increase initial tension

> **Tools required:**

710H10=2x3 Adjustment Wrench

▶ **CAUTION! The prosthetic foot must not swing through too far.**

Adjust the extension assist using the adjustment wrench so that the prosthetic foot swings through sufficiently according to the needs of the patient.

5.5.4 Adjusting the Swing Phase Friction

Function	Means of adjustment	Setting	Meaning
Swing phase friction	Adjusting screw (see fig. 12 - ②)	Rotate anticlockwise ③	Decrease braking effect Speed up extension Extension stop gets louder
		Rotate clockwise ④	Increase braking effect Slow down extension Extension stop gets quieter

> **Tools required:**

Allen key (4 mm)

▶ **CAUTION! The prosthetic knee joint has to reach full extension, even at slower walking speeds. Observe the information in the section "Fine-tuning the Settings" (see Page 28).**

Using the screwdriver, adjust extension damping so that the prosthetic knee joint does not swing too hard against the extension stop and so that the prosthetic foot swings through sufficiently according to the needs of the patient.

5.5.5 Fine-tuning the Settings

▶ **CAUTION! The various settings of the prosthetic knee joint cannot be adjusted entirely independently of each other. When the settings cannot be fully adjusted to the comfort needs of the patient, the settings should be established primarily based on safety aspects.**

Adapt the settings of the prosthetic knee joint to the patient by means of fine-tuning and exercises.

If the extension position is not reached with every step and at any walking speed, reduce the braking effect and/or increase the initial tension of the extension assist.

▶ Check the adjustment of the prosthesis settings during normal consultations and the annual safety inspection.

Advise the patient to check the prosthesis if changes in function occur.

5.6 Finishing the Prosthesis

> **Required tools:**

Torque wrench

- 1) Tighten all screws of the prosthetic components to the specified installation torque values.
- 2) Verify proper functioning of the prosthesis.

6 Use

6.1 Information for use

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

INFORMATION

At a flexion angle > **approx. 110°**, the extension assist automatically moves the prosthetic knee joint against the flexion stop. This function is activated in particular when kneeling.

- ▶ **Inform the patient.**

6.2 Cleaning

⚠ CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.

INFORMATION

- ▶ Clean the product if it gets dirty.
- ▶ For disinfection, only use disinfectants that do not affect the materials of the product. For further information, please consult the manufacturer.
- ▶ Observe the cleaning instructions for all prosthetic components.

Cleaning light soiling

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

Cleaning heavier soiling

- > **Required materials:** Cleaning cloth, 634A58 Isopropyl Alcohol, soft cloth
- 1) **NOTICE! Verify material compatibility! Do not remove lubricant from bearing points!**
Clean the product with a cleaning cloth and isopropyl alcohol.
- 2) Dry the product with a cloth.

Cleaning with Disinfectants

- > **Required materials:** colourless, alcohol-free disinfectant (verify material compatibility!), soft cloth
- 1) Disinfect the product with the disinfectant.
- 2) Dry the product with a cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

► Observe the following maintenance instructions.

- **NOTICE! Do not lubricate and grease the prosthetic joint.**
- **NOTICE! Repair work must be performed exclusively by manufacturer service.**
- Arrange regular maintenance intervals with the patient depending on the level of use.
- Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- The prosthetic components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.
- As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.
- Check the initial tension settings (EBS unit, extension assist) to make sure they are free of play. Increase the initial tension if needed.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

10 Technical data

Reference number	3R62=N, 3R62=1-N
Weight [g]	840
System height [mm]	142
Proximal system height up to the alignment reference point [mm]	-3
Distal system height up to the alignment reference point [mm]	145
Knee flexion angle	155°
Proximal connection	Pyramid adapter
Distal connection	Tube clamp (Ø 30 mm)

Reference number	3R62=N, 3R62=1-N
Max. body weight [kg]	3R62=1-N: 45 – 75 3R62=N: 75 – 125
Mobility grade	1 + 2

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Otto Bock HealthCare GmbH