


AxonRotation Adapter 9S501

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	8
FR Instructions d'utilisation (Personnel spécialisé)	12
IT Istruzioni per l'uso (Personale tecnico specializzato)	17
ES Instrucciones de uso (Personal técnico especializado)	22
PT Manual de utilização (Pessoal técnico)	27
NL Gebruiksaanwijzing (Vakmensen)	32
SV Bruksanvisning (Fackpersonal)	37
DA Brugsanvisning (Faguddannet personale)	42
NO Bruksanvisning (Fagpersonell)	47
FI Käyttöohje (Ammattihenkilöstö)	52
CS Návod k použití (Odborný personál)	56

INFORMATION

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INFORMATION

Last update: 2015-08-17

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

2 Product Description

2.1 Function

The AxonRotation Adapter 9S501 in combination with an Axon-Bus gripping component, e.g. the Michelangelo Hand 8E500, facilitates fitting with components in the modular Ottobock Axon-Bus prosthetic system.

The AxonRotation Adapter 9S501 allows passive pronation/supination. These functions support the patient during his or her daily activities and provide the highest rehabilitation value. The Axon-Rotation Adapter 9S501 promotes a physiologically correct body posture and reduces the need for unnatural, compensating body movements.

Passive, unlimited 360° pronation/supination of the Axon-Bus gripping component is possible from an activation moment > 2.0 Nm via the integrated ratchet system. There are 24 ratchet positions, one at every 15°, for this purpose.

3 Application

3.1 Intended Use

The AxonRotation Adapter 9S501 is to be used exclusively for upper limb exoprosthetic fitting and enables passive rotation of the Axon-Bus gripping component.

3.2 Use / Field of Application

In combination with myoelectrically controlled Axon-Bus gripping components in the modular Ottobock Axon-Bus prosthetic system, the AxonRotation Adapter 9S501 can be used from the transradial or transhumeral amputation level or, in case of dysmelia, for forearm or upper arm fittings.

3.3 Conditions of Use

The AxonRotation Adapter 9S501 was developed for everyday use and may not be used for unusual activities such as extreme sports (free climbing, paragliding etc.). The AxonRotation Adapter 9S501 is intended for use by one person only. Use of the product by another person is not approved by the manufacturer.

3.4 Qualification

Fitting a patient with the AxonRotation Adapter 9S501 may only be carried out by a prosthetist who has been authorised by Ottobock after completion of a corresponding training course.

4 Safety

4.1 Explanation of Warning Symbols

WARNING

Warnings regarding possible risks of severe accident or injury.

CAUTION

Warnings regarding possible risks of accident or injury.

NOTICE

Warnings regarding possible technical damage.

4.2 General Safety Instructions

CAUTION

Non-observance of safety instructions

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- ▶ Destruction of Axon-Bus components.
- Follow the safety instructions in this accompanying document.

CAUTION

Improper handling

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Instruct the patient in the proper use of the Axon-Bus prosthetic system.

4.3 Information to be Provided to the User

CAUTION

Mechanical overload caused by mechanical influences or external loads

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Do not subject the Axon-Bus components to mechanical vibrations or impacts.
- Check the Axon-Bus components for any visible damage before each use.

CAUTION

Independent manipulation of Axon-Bus components

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Manipulations on the Axon-Bus prosthetic system other than the tasks described in these instructions for use are not permitted.
- The Axon-Bus prosthetic system and any damaged Axon-Bus components may only be opened and repaired by certified Ottobock Myo-Service technicians.

CAUTION

Excessive strain or overloading due to unusual activities

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- The Axon-Bus prosthetic system and the Axon-Bus components have been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

- Careful handling of the Axon-Bus prosthetic system and the Axon-Bus components not only increases their service life, but is primarily in the interest of your own safety.
- If the Axon-Bus prosthetic system and the Axon-Bus components have been subjected to extreme loads (e.g., through a fall or similar occurrence), the Axon-Bus prosthetic system must be inspected by a prosthetist immediately for possible damage. The contact person for this is the responsible prosthetist, who will send the Axon-Bus prosthetic system to Ottobock Service if necessary.

⚠ CAUTION

Wearing in extreme ambient temperatures

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Avoid wearing the prosthesis in areas with temperatures outside of the permissible range. See the section "Technical Data".

⚠ CAUTION

Penetration of dirt and humidity

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Ensure that neither solid particles nor liquids can penetrate inside the Axon-Bus prosthetic system or Axon-Bus components.

5 Delivery/Storage

5.1 Scope of Delivery

- 1 AxonRotation Adapter 9S501
- 1 Lamination Ring
- 8 M3x8 Torx Screws
- 1 set Lamination Ring Dummy
- 1 647G934 Instructions for Use

5.2 Storage

Please see the specifications in the 'Technical Data' section.

6 Preparation for Use

6.1 Assembly on the Socket

Fabricating a socket with integrated lamination ring is a prerequisite for installation.

To fabricate the socket, please read through the 646T3=3.4 Technical Information 'Fabrication of a Transradial Prosthesis with a Michelangelo Hand'.

- 1) Connect the receptacle and plug of the Axon-Bus cable (see fig. 1).
- 2) Mount the oval plate of the AxonRotation Adapter on the integrated lamination ring (see fig. 2).
- 3) Screw together the oval plate and lamination ring with 4 M3x8 Torx screws (see fig. 3). Tighten the screws hand-tight and do not use a thread lock such as Loctite.

7 Handling

7.1 Disconnecting the AxonRotation Adapter from the Axon-Bus flexion component

When the unlock mechanism is activated the Axon-Bus flexion component, e.g. the AxonFlexion Adapter 9S500=*, can be disconnected from the AxonRotation Adapter 9S501. The AxonRotation Adapter 9S501 remains on the socket.

The Axon-Bus contacts are separated during detachment and automatically reconnected when the components are reattached.

CAUTION

Disconnecting or establishing electrical connections when switched on

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Power down the Axon-Bus prosthetic system by pressing the button in the charging receptacle before changing any Axon-Bus components (e.g. Axon-Bus gripping component).

8 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.
Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

9 Maintenance/Repair

9.1 Service Information

Since all moving mechanical parts are subject to wear and tear, regular service is required within the guarantee period. In the course of service, Ottobock Service (Ottobock Myo-Service) checks the entire product. Wear parts are replaced as required. Omitting a scheduled service inspection voids the guarantee.

10 Legal Information

10.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

11 Appendix

11.1 Symbols Used

Declaration of conformity according to the applicable European directives

SN YYY YW NNN

Serial number

Legal manufacturer

11.2 Technical data

Environmental conditions	
Storage and transport in original packaging	-20 °C/-4 °F to +40 °C/+104 °F
Storage and transport without packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 80 % relative humidity, non-condensing
Operation	-10 °C/+14 °F to +60 °C/+140 °F max. 80 % relative humidity, non-condensing

General	
Reference number	9S501
Product service life	5 years
Pronation	360° (24 ratchet positions)
Supination	360° (24 ratchet positions)
Weight	approx. 90 g/3.17 oz

11.3 Axon

The term "Axon" stands for **A**daptive **eX**change **O**f **N**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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