AxonArm Ergo Service 12K501=M-S

<input type="checkbox"/> DE	Gebrauchsanweisung (Fachpersonal)	3
<input type="checkbox"/> EN	Instructions for use (qualified personnel)	15

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DE

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Table of contents

EN

1	Foreword	16
2	Introduction	16
3	Purpose	16
3.1	Intended Use	16
3.2	Application	16
3.3	Conditions of Use	16
3.4	Qualification	16
4	Scope of Delivery, Tools and Instructions	17
4.1	Scope of Delivery	17
4.2	Accessories and Tools	17
4.3	Instructions	18
5	Preparation for Use	18
5.1	Adjustment, Fitting, Use	18
5.2	Disassembling the Patient's Prosthesis	18
5.3	Length Ranges	20
5.4	Adjusting the Length	21
5.5	Finishing	22
6	Service	24
7	Disposal	24
8	Legal Information	24
8.1	Liability	24
8.2	Trademarks	24
8.3	CE Conformity	24
9	Appendix	25
9.1	Technical data	25
9.2	Symbols Used	25

1 Foreword

INFORMATION

Last update: 2014-11-10

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The 12K501=M-S AxonArm Ergo (Service/Test) is hereinafter referred to as the AxonArm Ergo Service.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Introduction

This document explains all relevant steps from the disassembly of the patient prosthesis to the assembly and delivery of the finished AxonArm Ergo Service. It is intended for qualified professionals and assumes that this qualified personnel is trained in the handling of the various materials, machines and tools.

- ▶ **INFORMATION: This document does not replace the 647G790 Instructions for use of the 12K500=* AxonArm Hybrid/12K501=* AxonArm Ergo containing detailed information on the handling and safety of your product. Please read the 647G790 Instructions for Use carefully!**

3 Purpose

3.1 Intended Use

The AxonArm Ergo Service is to be used solely for the duration of the service and for trial fittings for exoprosthetic fittings of the upper limbs.

The AxonArm Ergo Service is compatible with the upper arm sockets of the 12K500=* AxonArm Hybrid and 12K501=* AxonArm Ergo.

3.2 Application

The AxonArm Ergo Service is suitable for unilateral and bilateral above-elbow amputees.

The AxonArm Ergo Service is to be used **exclusively** for the fitting of adults.

3.3 Conditions of Use


The AxonArm Ergo Service was designed for everyday activities and must not be used for unusual activities such as extreme sports (i.e. free climbing, paragliding, etc.).

3.4 Qualification

The fitting of patients with the AxonArm Ergo Service may only be carried out by a prosthetist.

4 Scope of Delivery, Tools and Instructions

4.1 Scope of Delivery

Item	Piece(s)	Product name	Article number
1	1	AxonArm Ergo Service	12K501=M-S
2	1	Axon adapter, short	13Z152
2	1	Counter nut, large	13Z111
3	1	Axon adapter, long	13Z153
3	1	Counter nut, small	13Z115
4	1	Adapter, short	13Z112
5	1	Adapter, long	13Z113

4.2 Accessories and Tools

The following accessories and tools are supplied for the assembly of the AxonArm Ergo:

Item	Tool	Article number
1	Disposable foam cover	13Z82
2	Strap wrench	13Z116
3	Hook wrenches	13Z109
4	Case key	

4.3 Instructions

The following accompanying documents necessary for assembly can be found on the inner side of the case under the foam cover:

Item	Description	Product number
not illustrated	Instructions for Use – 12K500=* AxonArm Hybrid/12K501=* AxonArm Ergo	647G790
not illustrated	Alignment and Assembly Instructions for 12K501=M-S AxonArm Ergo (Service/Test)	647G976
not illustrated	Information slip for reusable package	647H324

5 Preparation for Use

5.1 Adjustment, Fitting, Use

Read the 647G790 Instructions for Use before fitting the patient with the AxonArm Ergo. This contains important additional information about the use of the AxonArm Ergo.

5.2 Disassembling the Patient's Prosthesis

Before adapting the product to the patient, the patient's elbow component must be removed from the socket.


Example: patient socket, AxonArm and Michelangelo Hand

Turn off the AxonArm.

► **INFORMATION:** The AxonArm can only be turned off when the 13E500 AxonMaster is connected.

Pull off the prosthetic component.


► **INFORMATION:** For pulling off the prosthetic component, consult the corresponding instructions for use for AxonBus gripping components, e.g. 447D115 Instructions for Use for the 8E500 Michelangelo Hand.

Bend the AxonArm.

Release the lock using the pull cable if required.


Unscrew the AxonRotation adapter.


Set the dial to minimum compensation.


Extend the AxonArm.


Lock the elbow joint in this position using the pull cable to fix the removal and installation position.

Remove the screw on the strap clamp and remove the clamp.


Turn the screw. The humeral rotation feature is loosened and the patient socket can be unscrewed more easily.

Determine the length. Measure the distance between the front edge of the socket and the rotation centre of the joint. Document value 'X'.

5.3 Length Ranges


Length ranges from 190 mm to 302 mm can be covered by combining long and short Axon adapters.

Length range 190 mm – 210 mm


Required accessories:

1 pc. Axon adapter, short with counter nut (item 2)

**Length range 215 mm – 256 mm**

Required accessories:

- 1 pc. Axon adapter, long with counter nut (item 3)
- 1 pc. adapter, short (item 4)

**Length range 261 mm – 302 mm**

Required accessories:

- 1 pc. Axon adapter, long with counter nut (item 3)
- 1 pc. adapter, long (item 5)


5.4 Adjusting the Length

Based on the established length, the corresponding version can be configured and assembled.


The length range of 215 mm – 256 mm is illustrated here as an example.


The adapter, short and the Axon adapter, short with counter nut were already threaded into the AxonArm Ergo Service.

Screw the Axon adapter, short with counter nut in or out as needed and set the documented measurement "X". See the section "Disassembling the Patient's Prosthesis".


Mark the radial position on the Axon adapter, short or Axon adapter, long. The position between the 2 mounting screw holes must lie in an imaginary line to the charging receptacle.

► **INFORMATION:** The threaded tube must be screwed in at least to the point where the marker bore holes (arrow) are covered by the counter nut. If this is not possible, the next longer version has to be chosen.

Hand-tighten the adapter and counter nut of the threaded tube clockwise.


Thread the AxonArm Ergo Service into the strap wrench and fix it. First apply the hook wrench to the adapter (optional) and finally to the counter nut of the threaded tube. Tighten the nut clockwise until the prescribed torque is reached. Ensure that the radial position of the Axon adapter is not twisted.

The tightening torque of 16 Nm is reached as soon as the insides of the handles touch each other at the hook wrench.


5.5 Finishing

The process of adapting the AxonArm Ergo Service to the patient is completed by pulling on the foam cover and installing the AxonRotation adapter.


Pull the disposable foam cover over the AxonArm Ergo Service.


Mark the corresponding length according to the measurement 'X' determined previously (see the section "Disassembling the Patient's Prosthesis").


Use scissors to cut the disposable foam cover to length.

Use silicone grease to seal the plug connection.


Connect the AxonRotation adapter with the connecting cable.

Mount the AxonRotation adapter to the Axon adapter, short or the Axon adapter, long.

Observe the installation position when using the 9S503 AxonRotation adapter.

► **INFORMATION: Do not twist the Axon adapter, short or Axon adapter, long after installing the AxonRotation adapter. The cables could be damaged.**

Pull the disposable foam cover previously cut to length over the AxonArm Ergo Service.
 Mount the prosthetic component.
 Pull the cosmetic glove over the Michelangelo Hand if applicable.
 In case of a fitting with the 9S503 AxonRotation Adapter, make sure that the rotational movement is working smoothly.

- 1) Remove the lamination ring of the AxonArm Ergo Service (see the section "Disassembling the Patient Prosthesis").
- 2) Connect the elbow joint to the upper arm and configure the settings. See the 647G790 Instructions for Use, sections "Connecting the Elbow Joint to the Upper Arm", "Setting the Preflexion", "Adjusting the Humeral Rotation Feature", "Installing the Pull Cable" and "Checking the Symmetry".

6 Service

After each fitting, the AxonArm Ergo Service has to be sent to the authorised Ottobock Service Centre for inspection. There, the complete product is examined and readjusted if required. Wear parts are also replaced if required.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that does not comply with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

8 Legal Information

8.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

8.2 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

8.3 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

9 Appendix

9.1 Technical data

Description	Value
Reference number	12K501=M-S
Operating temperature	0 °C to 60 °C/32 °F to 140 °F
Storage and transport temperature	-20 °C to 40 °C/-4 °F to 104 °F
Relative humidity	max. 80 % relative humidity, non-condensing
Long-term storage temperature	-20 °C to 40 °C/-4 °F to 104 °F
Dimensions	approx. 345 mm x d55 mm to d80 mm (arm extended)
Weight without APS gripper	approx. 840 g/29.6 oz
Service life	5 years
Forearm length (as delivered)	185 mm
Maximum adjustable forearm length	302 mm
Minimum adjustable forearm length	190 mm
Flexion angle	15 ° - 145 °
Maximum holding force for forearm length of 305 mm	230 N
Protection class	IP20

9.2 Symbols Used

Symbol	Meaning
	This product must not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions provided by the responsible authorities in your country regarding return and collection processes.
	Declaration of conformity according to the applicable European directives
	Serial number of device
	Lot number
	Legal manufacturer

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.