

## 7E8

<b>DE</b> Gebrauchsanweisung .....	13	<b>NL</b> Gebruiksaanwijzing .....	50
<b>EN</b> Instructions for use .....	19	<b>SV</b> Bruksanvisning .....	56
<b>FR</b> Instructions d'utilisation .....	25	<b>PL</b> Instrukcja użytkowania .....	62
<b>IT</b> Istruzioni per l'uso .....	31	<b>HU</b> Használati utasítás .....	68
<b>ES</b> Instrucciones de uso .....	37	<b>TR</b> Kullanma talimatı .....	74
<b>PT</b> Manual de utilização .....	44	<b>EL</b> Οδηγίες χρήσης .....	80

1 **i** 5 Lieferumfang/  
Scope of delivery


501S71=M5X16


7Z68


647G98=all\_INT

709S10=2.5

709S10=1.5

7ES

DE Gebrauchsanweisung.....  
EN Instructions for use.....  
FR Instructions d'utilisation.....  
IT Istruzioni per l'uso.....  
ES Instrucciones de uso.....  
PT Manual de utilização.....

DE Gebrauchsanweisung.....  
EN Instructions for use.....  
FR Instructions d'utilisation.....  
IT Istruzioni per l'uso.....  
ES Instrucciones de uso.....  
PT Manual de utilização.....

2

**i** 6 Technische Daten/  
Technical data

①

≤45 kg  
≤99 lbs

②

≤145 cm


③


≤21 cm

215 g


Ø 22 mm


④

Systemhöhe (SH) **19 mm**  
System height (SH)



131 mm

125°

3

## i 7.2 Grundaufbau der Prothese/ Bench alignment of the prostheses


4


## 7.2 Grundaufbau der Prothese/ Bench alignment of the prosthesis

① Beckenkorb herstellen/Manufacturing the pelvic socket

②









③

## i 7.2 Grundaufbau der Prothese/ Bench alignment of the prosthesis


### ① Verstellbarer Extensionsanschlag/Adjustable extension stop


6 **i** 7.2 Grundaufbau der Prothese/  
Bench alignment of the prosthesis

① Rotation in der Achse des Rohrs/Rotation in the axis of the tube



1

2


3


7 **i** 7.2 Grundaufbau der Prothese/  
Bench alignment of the prosthesis

① Abduktion/Adduktion/Abduction/adduction


①

②

③

8

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### INFORMATION

Date of last update: 2021-10-21

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings.

## 1.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

## 2 Product description

The product (7E8) has the following key features:

- Monocentric prosthetic hip joint
- Product components for swing phase control:
  - Extensor (spring force adjustable)
- Material: Aluminium

## 3 Intended use

### 3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 3.2 Area of application

#### Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Only use the product according to its allowable field of application.

#### Allowable field of application

Children's prosthesis (Everyday prosthesis)

#### Allowable body weight

See item number ① in Figure ② "Technical data" on page 3

### Allowable body height

See item number ② in Figure 2 "Technical data" on page 3

### Allowable foot size

See item number ③ in Figure 2 "Technical data" on page 3

## 3.3 Combination possibilities

#### ⚠ CAUTION

#### Improper combination of prosthetic components

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's area of application.

#### INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

## 3.4 Environmental conditions

#### ⚠ CAUTION

#### Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

#### Allowable environmental conditions

Temperature range: -10 °C – +60 °C

#### Allowable environmental conditions

Use - Relative humidity: 0 % – 90 % (Non-condensing)

Storage/Transportation - Relative humidity: 100 %

Tyre pressure: 250 – 1100 mbar

Commercially available, solvent-free cleaning agents

#### Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

## 3.5 Reuse and lifetime

#### ⚠ CAUTION

#### Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

#### ⚠ CAUTION

#### Exceeding the lifetime

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the maximum lifetime defined in this section is not exceeded.

This product has been load-tested by the manufacturer. The maximum lifetime is 2 years.

## 4 General safety instructions

#### ⚠ CAUTION

#### Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

#### **⚠ CAUTION**

#### **Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

#### **Signs of changes in or loss of functionality during use**

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

#### **5 Scope of delivery**

The scope of delivery is shown in Figure 1 on page 2.

Only product components with reference numbers in the illustration can be reordered separately.

#### **6 Technical data**

The technical data are shown in Figure 2 on page 3.

- ④: System height (SH)

## **7 Preparing the product for use**

### **7.1 Information on fabrication of a prosthesis**

#### **⚠ CAUTION**

#### **Incorrect alignment, assembly or adjustment**

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

#### **⚠ CAUTION**

#### **Initial use of the prosthesis by the patient**

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

### **7.2 Bench alignment of the prosthesis**

#### **INFORMATION**

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

- ▶ **CAUTION!** To enable the patient to stand safely, align the prosthesis using the alignment images for the sagittal plane and frontal plane (Figure 3 on page 4) as well as the instructions for use of all prosthetic components that were used.

## Fabricating the pelvic socket (4 on page 5 – 1)

### CAUTION

#### Failure to observe instructions regarding the lamination plate and pelvic socket

Risk of falling due to breakage or deformation of the lamination plate and pelvic socket

- ▶ Follow all instructions in this section when fabricating the pelvic socket.
  - ▶ Do not modify the lamination plate (e.g. by shortening or bending).
  - ▶ Note that the section "Reuse and lifetime" also applies to the lamination plate.
- ▶ Manufacture the pelvic socket using the lamination plate in such a way as to ensure safe use with respect to the patient's body weight and level of activity.
- 1) Before laminating, close the threaded holes in the lamination plate with 636K6 plasticine.
  - 2) ②: For reinforcement, thread 699B1 fibreglass roving through the bores in the lamination plate and leave an overhang of about 4 cm.
  - 3) ③: Reinforce the pelvic socket in the area of the lamination plate using 616G12 carbon fibre cloth (up to three layers, depending on the weight and height of the patient).

## Installing the tube adapter (6 on page 7)

### CAUTION

#### Incorrect processing of tube

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

### CAUTION

#### Incorrect mounting of the tube

Risk of injury due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthetic component when mounting.

### CAUTION

#### Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

- ▶ Use the smallest value for the allowable insertion depth as the minimum insertion depth.
  - ②: Allowable insertion depth: 0 mm
  - ③: Unallowable insertion depth: > 0 mm

## Product adjustment options

- ▶ If necessary, use the adjustment options shown in illustrations (5 on page 6 – ①: adjustable extension stop; 6 on page 7 – ②: rotation in the axis of the tube; 7 on page 8 – ③: abduction/adduction).

## 7.3 Optimising the static alignment

- ▶ To optimise the static alignment, position the patient on the measuring device.
- ▶ In the sagittal plane, check the course of the load line relative to the marked alignment reference point and correct this if needed (see Figure 8 on page 9).

- ▶ In the frontal plane, check the course of the load line relative to the marked points (prosthetic foot, prosthetic knee joint, spina iliaca anterior superior) and correct if needed (see Figure 8 on page 9).

#### 7.4 Optimising during dynamic trial fitting

##### CAUTION

##### Adjusting the Settings

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

##### Setting the hip extensor

- ▶ Set the hip extensor according to Figure 9 on page 10 in such a way that hip extension is always achieved while walking, while ensuring it does not knock against the stop too hard (1 increasing spring pre-tensioning, 2 reducing spring pre-tensioning).

##### Fine-tuning the settings

- ▶ The various settings of the prosthesis cannot be adjusted entirely independently of each other. If the settings cannot be fully adjusted to the comfort needs of the patient, the settings should be established primarily based on safety aspects.
- ▶ Adapt the settings of the prosthetic components to the patient by means of fine-tuning and exercises.
- ▶ Check the adjustment of the prosthesis settings during normal consultations and as part of the annual safety inspections. Advise the patient to have the prosthesis inspected by qualified personnel if changes in function occur.

#### 7.5 Attaching the cosmetic cover

##### CAUTION

##### Use of highly hygroscopic particles (grease-absorbent substances such as talcum)

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not allow the product to come into contact with highly hygroscopic particles.
- ▶ To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.
- ▶ After affixing the cosmetic elements, verify proper functioning of the prosthesis.

#### 7.6 Finishing the prosthesis

- ▶ **CAUTION!** To avoid product damage and the risk of falling, finish the prosthesis by replacing set screws that are too short or too long, and by tightening all screw connections. In doing so, note the instructions for use for all prosthesis components regarding torque values and thread lock.

##### Securing the abduction/adduction position (10 on page 11)

- ▶ Secure the abduction/adduction position as shown in the illustration:
  - 1: Mark the position of the prosthetic hip joint
  - 2: Use the 7Z68 drilling template
  - 3: Bore through
  - 4: Tighten the set screw

### 8 Use

#### 8.1 Information for use

##### NOTICE

##### Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

## 8.2 Cleaning

### ⚠ CAUTION

#### Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
- ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.

#### Cleaning light soiling

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

#### Cleaning heavier soiling

- > **Required materials:** cleaning cloth, 634A58 isopropyl alcohol, soft cloth
- ▶ **NOTICE!** Verify material compatibility. Do not remove lubricant from bearing points.
- ▶ Clean the product with a cleaning cloth and isopropyl alcohol.
- ▶ Dry the product with a cloth.

#### Cleaning with Disinfectants

- > **Required materials:** colourless, alcohol-free disinfectant (verify material compatibility!), soft cloth
- ▶ Disinfect the product with the disinfectant.
- ▶ Dry the product with a cloth.
- ▶ Allow to air dry in order to remove residual moisture.

## 9 Maintenance

### ⚠ CAUTION

#### Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ **NOTICE!** Repair work must be performed exclusively by the manufacturer's technical service.
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Perform a safety inspection every 3 months.
- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.
- ▶ In case of changes in prosthesis function or relevant patient changes (such as the activity level, body weight, body height), perform unscheduled inspections and implement measures according to the situation.

#### Lubrication

Figure [11](#) on page 12 shows which lubricants to use.

- ▶ Use lubricants on the product sparingly.

## 10 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to

health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 11 Legal information

### 11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 11.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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