
4R68

DE	Gebrauchsanweisung	3
EN	Instructions for Use.....	5
FR	Instructions d'utilisation	8
IT	Istruzioni per l'uso	11
ES	Instrucciones de uso	13
PT	Manual de utilização.....	16
NL	Gebruiksaanwijzing.....	19
SV	Bruksanvisning	21
DA	Brugsanvisning.....	24
PL	Instrukcja użytkowania.....	26
HU	Használati utasítás	29
CS	Návod k použití.....	31
EL	Οδηγίες χρήσης	34
ZH	使用说明书.....	36

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English

Date of the last update: 2013-03-27

- Please read this document carefully.
- Follow the safety instructions.

1 Application and Description

The 4R68 Socket Adapter (aluminum) is to be **exclusively** used for the prosthetic fitting of amputations of the lower limb.

Field of application according to the **Ottobock Mobility System MOBIS:**

Approved for a patient with a **weight of up to 100 kg/220 lbs.**

The Socket Adapter is recommended for below-knee prostheses and is integrated as lamination anchor into the laminate of the prosthetic socket.

Note!

Do not bend the socket adapter. Risk of cracking in the aluminum.

**Attention!**

Bring the anchor arms in the correct **AP** and **ML** position.

**Attention!**

Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

2 Technical Data

	4R68
System height	-7 mm
Weight	70 g
Material	Aluminum
Max. patient weight	100 kg

3 Handling

3.1 Recommendations for the Layup

- 2 layers of 623T9 Nylglass Stockinette
- 1 layer of 616G12 Carbon Fiber Cloth
- 1 layer of 616G4 Fiberglass Matting
- 1 strip of 616G12 Carbon Fiber Cloth – lateral and medial (fig. 1)
- 1 strip of 616G12 Carbon Fiber Cloth – circularly applied below the patella
- 2 layers of 623T9 Nylglass Stockinette
- The hollow space of the socket adapter is filled with a mixture of 617H21 Orthocryl Sealing Resin and 617Z8 Microballoons. Following this the socket adapter is placed on the model. Then align the socket adapter perpendicularly according to the marking of the 4X3 Lamination Protection (fig. 2).
- 1 layer of 616G12 Carbon Fiber Cloth over the socket adapter (fig. 3)

- 1 ½ layers of 616G13 Fiberglass Stockinette (figs. 4, 5)
- 2 layers of 623T9 Nylglass Stockinette

The enclosed 4X3 Lamination Protection covers the adjustment pyramid during the lamination procedure with 617H19 Orthocryl Lamination Resin 80:20 or 617H55 C-Orthocryl and provides an accurately fitting brim of the laminate (fig. 6).

3.2 Maintenance Instructions

Note:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

4 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

5 Warranty

Ottobock warrants this product/these products from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local Ottobock distribution company for more information.

6 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH