4R103

DE Gebrauchsanweisung	3
EN Instructions for use	6
FR Instructions d'utilisation.....	10
IT Istruzioni per l'uso	14
ES Instrucciones de uso	18
PT Manual de utilização.....	22
NL Gebruiksaanwijzing.....	25
SV Bruksanvisning	29
DA Brugsanvisning	33
NO Bruksanvisning.....	36
FI Käyttöohje.....	40
PL Instrukcja użytkowania.....	44
HU Használati utasítás	47

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INFORMATION

Last update: 2016-09-30

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 4R103 sliding tube clamp adapter permits translation adjustments at the proximal end of a tube adapter. It therefore permits the distal components of the modular prosthesis to be shifted in parallel, regardless of the angle adjustment by the pyramid adapters. The

adjustments can be made either in the frontal plane (medial or lateral) or in the sagittal plane (anterior or posterior). This is also possible while wearing the prosthesis.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system and can be combined with other products of this modular system.

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Approved for a body weight of up to **85 kg**.

Ottobock recommends the use of structural components with a diameter of 34 mm when elevated loads are expected (e.g. for a transtibial prosthesis).

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the patient's level of activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 7).

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.

- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION

Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

4R103 sliding tube clamp adapter				
Fi-g.	Item	Quantity	Designation	Reference number
-	-	1 (■)	Instructions for use	647H66

4R103 sliding tube clamp adapter				
Fi-g.	Item	Quantity	Designation	Reference number
-	-	1 (▲)	Cap screw	501Z2=M4x-22
-	-	2 (▲)	Cap screw	501Z16
-	-	4 (▲)	Set screw	506G3=M8x-16
-	-	1 (▲)	Set screw	506G3=M5x-8

5 Preparation for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

⚠ CAUTION

Incorrect mounting of the tube

Risk of injury due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthetic component when mounting.

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

5.1 Mounting the adapter

- > **Required materials:** 710D4 torque wrench, degreasing cleaner (such as 634A58 isopropyl alcohol)
- 1) Loosen the tube clamp cap screw by **2 turns**.
- 2) **For final assembly:** Degrease the insertion range of the tube clamp adapter and tube adapter.
- 3) Fully insert the tube adapter into the tube clamp adapter up to the stop. Align the clamping slot:
Tube clamp adapter: **anterior**
Tube clamp adapter, sliding: **anterior** or **medial**
- 4) Tighten the cap screw on the tube clamp, choosing the torque according to the thread diameter (M4: **5 Nm**, M5: **10 Nm**, M6: **13 Nm**).

5.2 Installation in the modular prosthesis

- > **Required tools and materials:**
710D4 torque wrench, 636K13 Loctite®
- ▶ Install the other prosthetic components as described in the respective instructions for use.

5.3 Adjusting the adapter

The proximal section of the tube clamp adapter is movable in a dovetail guide. A set screw acts as a stop. It is fastened in the desired position.

Adjusting the stop

- > **Tools required:**
710D4 torque wrench
- ▶ **CAUTION! Only use the adapter with the set screw screwed in!**
Screw in the set screw at the desired position (installation torque: **5 Nm**).

Moving the adapter

- > **Tools required:**
710D4 torque wrench
- > No load on the prosthesis, or the prosthesis is not being worn.
- 1) Loosen the two cap screws of the dovetail guide (see fig. 1).
- 2) Move the proximal section of the tube clamp adapter in the dovetail guide. Note the maximum displacement (see fig. 2).
- 3) Tighten the two cap screws (installation torque: **15 Nm**).

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Reference number	4R103
Weight [g]	185
System height [mm]	51
Material	Titanium
Diameter [mm]	30
Max. displacement [mm]	±11
Max. body weight [kg]	85

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Otto Bock HealthCare GmbH