
28U9

DE Gebrauchsanweisung	3
EN Instructions for Use	7
FR Instructions d'utilisation	10
IT Istruzioni per l'uso	15
ES Instrucciones de uso	19
PT Manual de utilização	23
NL Gebruiksaanwijzing	27
SV Bruksanvisning	31
DA Brugsanvisning	35
PL Instrukcja użytkowania	38
HU Használati utasítás	43
CS Návod k použití	46
RO Instrucțiuni de utilizare	50
BG Инструкция за употреба	54
EL Οδηγίες χρήσης	59
RU Руководство по применению	63

	Artikelnummer / Article Number	Seite / Side	Schuhgröße / Shoe Size
	28U9=L35-37	links/left	L35-37
	28U9=L35-37	links/left	L35-37
	28U9=R35-37	rechts/right	R35-37
	28U9=R35-37	rechts/right	R35-37
	28U9=L37-39	links/left	L37-39
	28U9=L37-39	links/left	L37-39
	28U9=R37-39	rechts/right	R37-39
	28U9=R37-39	rechts/right	R37-39
	28U9=L39-41	links/left	L39-41
	28U9=L39-41	links/left	L39-41
	28U9=R39-41	rechts/right	R39-41
	28U9=R39-41	rechts/right	R39-41
	28U9=L41-44	links/left	L41-44
	28U9=L41-44	links/left	L41-44
	28U9=R41-44	rechts/right	R41-44
	28U9=R41-44	rechts/right	R41-44
	28U9=L44-47	links/left	L44-47
	28U9=L44-47	links/left	L44-47
28U9=R44-47	rechts/right	R44-47	
28U9=R44-47	rechts/right	R44-47	

Material	PA, PP
-----------------	--------

1 Vorwort

INFORMATION

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INFORMATION

Last update: 2013-05-21

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

These instructions for use provide important information on the fitting and application of your 28U9 ankle foot orthosis.

2 Application

2.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Dorsal flexor weakness (e.g. peroneal paralysis)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

- Moderate to severe spasticity in the lower leg
- Moderate to severe oedema
- Leg ulcers
- Moderate to severe foot deformities

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical aid will be applied; sensory and circulatory disorders in the body region being treated and in the foot and lower leg areas.

2.4 Further Usage Restrictions

The orthosis must not be used for the prevention of contractures (shortening) etc.

The orthosis must not be used for sports activities that involve jumping, sudden movements or rapid step sequences (e.g. basketball, badminton, riding sports). Sports activities should generally be discussed with the patient.

2.5 Effects

The ankle foot orthosis provides lifting of the foot and limits plantar flexion.

3 Safety

3.1 Explanation of Warning Symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

NOTICE

Warnings regarding possible technical damage.

3.2 General Safety Instructions

⚠ CAUTION**Using the product without following the instructions for use**

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

⚠ CAUTION**Reuse on other patients and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

⚠ CAUTION**Contact with heat, embers or fire**

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

⚠ CAUTION**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).
- ▶ **Inform the patient.**

NOTICE**Contact with salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

4 Handling

⚠ WARNING**Operation of motor vehicles**

Risk of accidents due to restricted body function

- ▶ Observe the applicable legal and insurance regulations for the operation of motor vehicles and have your driving ability examined and certified by an authorised agency.
- ▶ **Inform the patient.**

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

4.2 Adaptation

⚠ CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.
- ▶ **Inform the patient.**

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Use without suitable footwear

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.
- ▶ **Inform the patient.**

Adapting the orthosis size:

- 1) **If necessary:** Cut or sand down the orthosis to size according to shoe size.
- 2) Then, deburr the worked edges and polish out any notches.

Riveting of the 623Z3 hook and loop closure

- 1) **If necessary:** Thermoplastically re-form the orthosis by applying gentle heat to a limited extent.
NOTICE! Re-forming of the orthosis in cold conditions or at temperatures above 80°C is not permitted.
- 2) Gently heat the riveting points to avoid the formation of cracks.
- 3) Adapt the orthosis by cutting, milling and grinding.
- 4) Polish out any notches on the edges.

Final inspection:

- ▶ Conduct trial walking with the patient. Testing on inclines, ramps and stairs is mandatory. Fine tune as required.

4.3 Application

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).
- ▶ **Inform the patient.**

1. Open the hook-and-loop strap.
2. Insert the orthosis into the shoe.
3. Put the orthosis on together with the shoe.
4. Tighten and close the hook-and-look strap.

4.4 Cleaning

Clean the orthosis regularly:

Textile component:

- 1) Remove the textile component from the orthosis.
- 2) Close all hook-and-loop closures.
- 3) Hand wash the textile component in warm water at 30 °C with a standard mild detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

Composite component:

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH