
28L20

DE Gebrauchsanweisung	3
EN Instructions for use	5
FR Instructions d'utilisation	7
IT Istruzioni per l'uso	9
ES Instrucciones de uso	11
PT Manual de utilização	13
NL Gebruiksaanwijzing	15
SV Bruksanvisning	18
FI Käyttöohje	20
DA Brugsanvisning	22
HU Használati utasítás	24
TR Kullanma talimatı	26
EL Οδηγίες χρήσης	28
RU Руководство по применению	30
ZH 使用说明书	32

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English

Last update: 2015-11-26

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

The 28L20 Hip Abduction Orthosis is to be used **exclusively** for the treatment of congenital hip dysplasia in infants.

2 Indication and function

Treatment of hip dysplasia Type IIb (as categorized by Graf) as well as residual dysplasia after conservative reduction of hip dislocation

Specific indication should be determined by the prescribing physician.

This particular type of treatment (devised by Hoffman and Daimler) has been used effectively for decades. When the child begins to stand on its own, development of the acetabulum will occur until the child has reached the age of 18 months.

Note:

The 28L20 Hip Abduction Orthosis was designed for children being at the age of 6 to 18 months.

Note:

Please observe the indication period, since otherwise the stability of the components and thus the security of the patient will be endangered.

3 Components

The 28L20 Hip Abduction Orthosis is available in one size and consists of the following components (fig.1):

- (1) 29L48=L/RI Medial Condyle Band
- (2) 29L48=LA/RA Lateral Condyle Band
- (3) 29L51=40 Spreader Bar
- (4) 29Y27=M4×6 Hexagon Head Cap Screws M4 with Collar
- (5) 21Y136=38-7 Hook and Loop Strap with Roll Loop
- (6) 29L56 Rubber Collar
- (7) 29L103 Abduction Unit
- (8) 29L50=L/R Condyle Pads

- (9) 709S10=2.5 Allen Wrench, 2.5 mm
- (10) 501Z24=M4×4.5 Binding Head Hexagon Socket Screw 2.5 mm

4 Technical information

The anatomically formed thigh supports are divided into two and provided with parallel width adjustment. They guarantee exact adaptation to the thigh and good fit. The condyle wedge nearly excludes that the orthosis will slip down. The adjustment range of the condyle bands is 55 mm. To adjust the width of the bands, release the two binding head hexagon socket screws (2.5 mm). Later adjustment of the condyle bands is possible through thermoplastic forming. The condyle pads made of Pedilin® should be glued in place after initial fitting. A hook and loop closure serves to close the orthosis. Pads for the hollow of the knee may additionally be attached.

Two ball joints are mounted between the medial condyle bands and the abduction unit. The range of motion of these ball joints is limited to 45 degrees. Rubber collars serve to protect the ball joints. This medial bearing of the abduction unit represents an essential biomechanical feature of the orthosis and prevents torsion of the condyle band on the thigh.

The 28L20 Hip Abduction Orthosis is easy to use and waterproof.

5 Adjusting and applying the orthosis

The thigh condyle bands may be applied over bare skin, over cotton stockinette (623T1=6) or over thin trousers. The width adjustment of the condyle bands may be adjusted beforehand, according to taken measurements, or directly on the leg. For shifting the condyle bands (fig. 2) release the two binding head hexagon socket screws (Allen Wrench 2.5 mm).

To apply the 28L20 Hip Abduction Orthosis, place the supports on the front of the child's thighs. Take care to position the medial supracondylar profile correctly (fig. 3). Adjustment of the abduction unit is achieved by releasing the headless set screw (2.5 mm) (fig. 4). Initial fitting and later adjustment of the orthosis should only be carried out by a physician or by an orthotist.

The outer telescope has a hole so that the width setting is documented with a marker (Fig. 5) and markings, making it reproducible at any time after adjustments are made.

After adjustment of the orthosis, secure all screws with 636K13 Loctite!

6 Instructions for use

Parents should be instructed on the use of the orthosis at the time of fitting. Only the hook and loop closures are needed to apply or remove the orthosis (fig. 6). Highly stable shoes with good traction are suggested to help your child stand safely when wearing the orthosis. In addition, the outer edge of the shoes may be wedged by 3 mm to make it easier to stand. The ball joints in the orthosis minimize the torsion of the thigh supports, thus preventing a possible adduction at heel strike. Your child's doctor will determine a daily wearing schedule for the hip abduction orthosis. Since this orthosis is waterproof, it may remain in place during bathing and other activities. When the radiological findings show a satisfying result, the child is to be weaned from the orthosis. In general, the time to wean the child from the orthosis will extend over a period of 6 to 8 weeks. During this period, abduction will be reduced every two weeks by 1 cm.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document.

The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE conformity

This product meets the requirements of the European Directive 93 / 42 / EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive.

The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock Korea HealthCare Inc.
South East Asia Co., Ltd.
Otto Bock HealthCare GmbH