





9E369 / 9E370

DE 4 in1 Controller LS	3
EN 4 in1 Controller LS	7
FR Controler LS 4 en1	13
IT Controller LS 4 in1	17
ES 4 in1 Controlador LS.....	22
SV 4 in1 Controller LS, styrenhet.....	27
NL 4 in1 Controller LS	32

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Date of the last update: 2013-01-23

- Please read this document carefully.
- Follow the safety instructions.

1 Application

The 4in1 Controller LS 9E369 or 9E370 is to be used **solely** for control of the 8E51 Electrohand 2000.

1.1 Qualification

Patients may be fitted with the product only by qualified personnel authorised by Ottobock after completing the corresponding training.

2 Safety Instructions

Failure to follow the safety instructions mentioned below can lead to malfunction of the 4 in 1 Controllers LS and may cause a health risk for the patient.

Use the Patient Information section of this pamphlet (Section 3) to inform your patient about the function and use of the 4in1 Controller LS.

3 Patient Information

- Do not allow foreign particles or liquids to get into the 4in1 Controller LS.
- The 4in1 Controller LS should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Avoid staying near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores), as this can lead to malfunction of the 4in1 Controller and the Electric Hand.
- Modifications of the components without the written approval of the manufacturer are explicitly forbidden.

4 Scope of Delivery and Spare Parts

4.1 Scope of Delivery

Delivery of the 9E369 4in1 Controller LS for hand size 5 comprises:

- 1 pc. 9E371 Contact Plate (Fig. 2, item 1)
- 1 pc. 9E372 Pattern for Lamination (Fig. 2, item 2)
- 1 pc. 9E373 Controller for size 5 (Fig. 2, item 3)
- 1 pc. 647H209 Instructions for Use

The delivery does not include the 13E182 Coding Plug Set.

Delivery of the 9E370 4in1 Controller LS for hand sizes 5½, 6 and 6½ comprises:

- 1 pc. 9E371 Contact Plate (Fig. 2, item 1)
- 1 pc. 9E372 Pattern for Lamination (Fig. 2, item 2)
- 1 pc. 9E374 Controller for size 5½, 6 and 6½ (Fig. 2, item 3)
- 1 pc. 647H209 Instructions for Use

The delivery does not include the 13E182 Coding Plug Set.

4.2 Spare Parts

- 9E371 Contact Plate
- 9E372 Casting Pattern
- 9E373 Controller for size 5
- 9E374 Controller for size 5 ½, 6 and 6 ½

5 Description and Function

The 4in1 Controller LS offers further options in the fitting of children, particularly those with long below-elbow residual limbs or wrist disarticulation amputations.

Its flat, slightly vaulted structural shape allows the Controller to be mounted between inner and outer socket (Fig. 1). This minimizes the overall length of the prosthesis. The 9E371 Contact Plate connects to the lamination ring.

The Controller offers four different control options, indicated by coding plugs of different colors.

		9E369	9E370
1 white	L	EVO 1	digital
2 red	R	EVO 1	digital
3 green	L	DMC	DMC
4 blue	R	DMC	DMC
5 yellow	L	DMC „Low/ Input“	DMC „Low/ Input“
6 violet	R	DMC „Low/ Input“	DMC „Low/ Input“
7 orange	L	EVO 2	EVO 2
8 black	R	EVO 2	EVO 2

6 Fabrication of the Prosthetic Socket

Warning: Never apply carbon fiber layers in one piece from the wrist area to the electrode area, as this material is conductive and may interfere with the myoelectric signals.

7 Assembly

Lubricate all contact bushings of the cable and the 4in1 Controller LS using 633F11 Silicone Grease.

Insert the coding plug into the controller with the recess fitting into the outwardly curved side, as shown (Fig. 5). Insert the coding plug completely. Connect the battery.

To change control modes, the battery must be removed and reinserted (voltage cutoff) after a new coding plug has been inserted. Until this has been done, the 4in1 Controller LS will be operated by the previous control option.

Attach the "opening" electrode to 2 (■■■) and the "closing" electrode to 1 (■■■). Connect the connection cable of the contact plate and of the battery to the Controller as well (Fig. 3).

Place the 711M50-1 Adapter into the contact plate. The contact plate should audibly engage in the lamination ring (Fig. 4).

Attach the Controller inside the forearm shell using adhesive tape.

8 Control Options and Adjustment of the Electrodes

The best possible function of the prosthesis is achieved through optimal adjustment of the electrodes. This is easy with the MyoBoy.

8.1 Dynamic Mode Control

green coding plug = left Electrohand

blue coding plug = right Electrohand

Opening and closing of the Electrohand 2000 is controlled by two 13E200 Electrodes.

Grip speed and grip force are determined by the level of the muscle signal.

This provides a more natural gripping movement and gentle handling of small, fragile objects is enhanced.

Adjusting the Electrodes

Select the desired program version in the PAULA program. Adjust each electrode so that the patient can maintain the muscle signal for approximately two seconds at a point slightly above the value HIGH (Fig. 6).

8.2 Dynamic Mode Control – „Low/ Input“

yellow coding plug = left Electrohand

violet coding plug = right Electrohand

In this mode, maximum grip speed and grip force are reached more quickly. This control should be used **only** for children with a low muscle potential (below 20 μ V).

Opening and closing of the Electrohand 2000 is controlled by two 13E200 Electrodes.

The grip speed and grip force are controlled by the level of the muscle signal.

This provides a more natural gripping movement and gentle handling of small, fragile objects is enhanced.

Adjusting the Electrodes

Select the desired program version in the PAULA program. Adjust each electrode so that the patient can maintain the muscle signal for approximately two seconds at a point slightly above the value LOW (Fig. 7).

8.3 Digital/EVO 1 (Electronic Voluntary Opening)

white coding plug = left Electrohand

red coding plug = right Electrohand

Two-electrode control –Digital– (only in combination with 9E370 4in1 Controller LS)

Opening and closing of the Electrohand 2000 is controlled by two 13E200 Electrodes.

When sufficient muscle signal is generated, the Electrohand opens and closes at a constant speed.

One-electrode control –EVO 1– (only in combination with 9E369 4in1 Controller LS)

The EVO mode has been developed to facilitate control of the Electrohand for children who cannot differentiate between two muscle signals.

Opening and closing of the Electrohand 2000 is controlled by only one 13E200 Electrode.

When sufficient muscle signal is generated, the Electrohand opens at a constant speed, and when the myo-signal stops, the Electrohand closes at a constant speed. The electrode is usually placed over remnant extension muscles.

If this is not practical, the Electrohand can be controlled by any other available myoelectric signals.

The electrode cable must be plugged into the electrode connector plug 1 or 2 of the 4in1 Controller LS.

Adjusting the Electrodes

Select the desired program version in the PAULA program. Adjust each electrode so that the patient can maintain the muscle signal for approximately two seconds at a point slightly above the value ON (Fig. 8).

8.4 EVO 2 (Electronic Voluntary Opening)

range coding plug = left Electrohand

black coding plug = right Electrohand

The EVO mode has been developed to facilitate control of the Electrohand for children who cannot differentiate between two muscle signals.

Opening of the Electrohand 2000 is controlled by one 13E200 Electrode.

The stronger the muscle signal, the faster the Electrohand opens.

When the myo-signal stops, the motor automatically closes the Electrohand.

The electrode is usually placed over remnant extension muscles.

If this is not practical, the Electrohand can be controlled by any other available myoelectric signals.

Grip force in the EVO mode can be adjusted.

Connecting the Electrode

The electrode cable must be plugged into the electrode connector plug 1 or 2 of the 4in1 Controller LS.

Adjusting the Electrode

Select the desired program version in the PAULA program. Adjust each electrode so that the patient can maintain the muscle signal for approximately two seconds at a point slightly above the value LOW (Fig. 7).

Attention:

Allow the patient to rest during the adjustment of the electrodes, otherwise muscle fatigue will lead to inconsistent results and the therapist will tend to adjust the electrodes with too high a sensitivity.

The electrodes are to be placed only on intact skin and with as much electrode-skin interface as possible. If strong interference by electrical appliances still occurs, please check the positioning of the electrodes and – if need be – reposition them. If the interference cannot be eliminated, please contact an Ottobock Service Facility.

9 Handling the Batteries

For operation of the Electrohand 2000, use only fully charged 757B13 batteries. To increase the useful life of the batteries, the 4in1 Controller LS automatically turns itself off whenever the battery is discharged.

An intelligent battery management alerts the patient to the decreasing charging condition of the battery as follows:

Stage 1: Operation of the hand is very slow, and even with maximum electrode voltage only minimum grip force is achieved (maximum myoelectric signal).

Stage 2: The hand can be opened, but once it stops no further movement is possible.

10 Disposal

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

11 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

12 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

13 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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HealthCare Inc.
South East Asia Co., Ltd.