

17B105, 17B106

DE	Gebrauchsanweisung	3
EN	Instructions for use	7
FR	Instructions d'utilisation	11
IT	Istruzioni per l'uso	15
ES	Instrucciones de uso	19
PT	Manual de utilização	24
NL	Gebruiksaanwijzing	28
SV	Bruksanvisning	32
DA	Brugsanvisning	36
NO	Bruksanvisning	40
PL	Instrukcja użytkowania	44
HU	Használati utasítás	49
CS	Návod k použití	53
HR	Upute za uporabu	57
TR	Kullanma talimatı	61
EL	Οδηγίες χρήσης	65
ZH	使用说明书	70

INFORMATION

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INFORMATION

Last update: 2015-10-06

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

Explanation of symbols

⚠ WARNING Warnings regarding possible risks of severe accident or injury.

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

1 Application

The System Knee Joints with Wedge Lock and Remote Cable-Release must be used in pairs, **solely** for the orthotic fitting of the lower limb.

INFORMATION

These Titanium System Knee Joints are designed to be used together with free-motion or dorsiflexion assist ankle joints. Please note: combining this product with ankle joints having a dorsiflexion stop increases stresses on the entire orthosis and results in increased wear and tear on the components.

The remote cable-release should only be used in combination with the integral Cable Housing (as delivered.) Refer to 5. Maintenance.

2 Indications for use

Partial paralysis or complete paralysis of the lower limb musculature. Specific indications must be determined by the physician.

3 Construction / Components

3.1 Components (fig. 1)

(1)	30X75=*	Upper Joint Section
(2)	30U113=*	Lower Joint Section
(3)	513D19=0.63x4x20	Compression Spring
(4)	17Y122	Wedge Lock
(5)	4Z80=*	Bearing Washer
(6)	17Y17=*	Brass Bushing
(7)	17Y93=*	Cotter Bolt
	and 501S32=*	Joint Screw
(8)	17Y121=*	Cover
(9)	17Y123=M4x10	Countersunk Head Screw
(10)	21A17=*	Cable Housing
(11)	21A12	Threaded Sleeve
(12)	21A7	Coupling Piece
(13)	21A25	Spring

The upper joint section (1), lower joint section (2), and cover (8) are made of titanium. The 17B105/17B106 System Knee Joints have both axial bearings (bearing washers) and radial bearings (brass bushings) and have hardened cotter bolts. The 17Y122 Wedge Lock is mounted in the 30U113 Lower Joint Section along with the 513D19=0.63x4x20 Compression Spring. The 21A18=2 Perlon Cable is led through the 21A17 Cable Housing and allows unlocking of the joint.

4 Processing

INFORMATION

Please note that the joints are exposed to higher structural loading as a result of individual clinical parameters such as forefoot lever length or angular limb deformities. These must be accommodated in the alignment and in the structural design of the orthosis.

Insert the appropriate side bars into the connection channels. Options include the 17B38=T (titanium), 17B38, 17B39, 17B40, 17B41, 17B84, 17B85 System Side Bar Sets or the 17Y105, 17Y107 System Lamination Bars.

The 636W28 Ottobock Special Adhesive is required to stabilize the bars in the connection channels. To do this, degrease the attachment surfaces using 634A1 Ottobock Thinner and Solvent or equal. Mix the contents of the two tubes containing special adhesive and hardener (included), and apply to the joint's connection channels. Insert the System Bars and screw in the 501T7=7.5x9xM5 Phillips Head Screws firmly. Allow three to four hours for the adhesive to dry prior to use. Final bonding strength will be achieved after about 16 hours.

NOTICE

Do not apply heat in excess of H 300 °C (H 572 °F) because the material's temper will be destroyed.

⚠ CAUTION

Before delivering the orthosis to the patient, make sure the joint systems are precisely parallel and that both locks are parallel and engaging simultaneously. Instruct the patient on use and maintenance of the orthosis.

Note that the wedge locks will not disengage when there is a flexion moment on the joints. The patient must keep the leg fully extended to unlock the orthosis safely.

5 Maintenance

For cleaning and other maintenance of the upper joint section (1), the cover (8) can be removed by removing the countersunk head screws (9).

When the brass bushings (6) show traces of wear, they can be replaced. To assemble, place the lubricated bearing washers (5) on the protruding bushings.

The 21A17 Cable Housing can be attached to the orthosis using the 21A5 Holding Device, which helps guide the 21A18=2 Perlon Cable. The Cable Housing can be shortened to within 10 mm [3/8 inch] of the opening of the upper joint section and must be screwed 3 mm [1/8 inch] into the upper joint section for secure attachment.

The Cable Housing must be secured with 636K13 Loctite.

INFORMATION

Screwing the cable housing more than 3 mm [1/8 inch] into the upper joint section compresses the compression spring and thus makes it more difficult to release the lock.

The coupling piece (12) with inner spring (13) can be used to connect the two perlon cables (fig. 2).

Use only 633F7 Ottobock Special Lubricant. To ensure proper function, make sure no dirt particles enter the joint mechanism. Use appropriate protection devices (e.g. joint covers), as necessary.

CAUTION

Visible wear on the lock mechanism may affect the patient's safety. Replace the wedge lock and/or the lower joint section, according to each individual's use pattern.

INFORMATION

It is recommended that the function and condition of this orthosis be checked every six months. The effectiveness of the 513D19=0.63x4x20 Compression Springs should be examined at this time, and they should be replaced as necessary.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH