

**SensorHand Speed 8E38=8\*, 8E39=8\*, 8E41=8\***  
**MyoHand VariPlus Speed 8E38=9\*, 8E39=9\*, 8E41=9\***

<b>DE</b> Patienteninformation .....	3
<b>EN</b> Patient information .....	24
<b>FR</b> Informations destinées au patient .....	43
<b>IT</b> Informazioni per il paziente .....	63
<b>ES</b> Información para el paciente .....	83
<b>PT</b> Informação para o paciente .....	103
<b>NL</b> Patiënteninformatie .....	123
<b>SV</b> Brukarinformation .....	143
<b>EL</b> Πληροφορίες για τον ασθενή .....	163

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---

**Date of the last update: 2021-04-14**

- Please read this document carefully before using the product and observe the safety notices.
- Obtain instruction from the qualified personnel in the safe use of the product.
- Please contact the qualified personnel if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

## Introduction

Dear Patient. Please note that this information brochure has been prepared for the Ottobock products SensorHand Speed and MyoHand VariPlus Speed. For clear distinction, the chapter headings therefore contain the respective product names. Only start up the product in accordance with the information contained in the accompanying documents supplied.

## SensorHand Speed with SUVA\*-Sensor Technology

During the fitting, your prosthetist considered six control modes to select one most suitable for your specific requirements. All six control modes are described in the appendix with examples.

To help you learn how the SensorHand Speed works, this booklet summarises the basic features, based on the most common type of control, DMC or Dynamic Mode Control:

- You can control your hand's opening or closing speed based on the strength of your muscle contraction. This makes gripping easy and very natural.
- When you hold an object in the SensorHand Speed, the strength of your muscle contraction controls the degree of force used to grip the object. This control makes it easier to handle small or fragile objects.

### Auto-Grasp SUVA-Sensor Technology (Fig. 1)

If an object begins to slip out of your prosthetic hand, the SensorHand Speed automatically increases the gripping force until the object is held securely.

To function in this way, sensors integrated in the thumb (See arrow in Fig. 1) and finger area transmit signals which are evaluated by special electronics. The electronics then tell the motor the required gripping force. You can override the Auto-Grasp feature at any time by making a brief opening signal.

With the SensorHand Speed, you don't need to constantly control your grip: you can pay attention to other things instead. Suddenly everyday situations – which might require concentration when using other prosthetic systems – are easier to tackle.

- Plastic grocery bags, for example, will no longer tend to slip out of your hand.
- Tying your shoelaces becomes easier, too.
- Delicate drinkware no longer slips when you pour liquids into them.
- Etc.

A brief opening signal will immediately stop the SensorHand Speed's Auto-Grasp response at any time. Also, the maximum value of the automatic gripping force is limited, and the Auto-Grasp feature is available only for the control modes 1 to 4.

For grasping flexible foam materials or tools, control mode 6 offers you the possibility to manually deactivate the SUVA-Sensor.

\* Developed by Otto Bock HealthCare Products GmbH in cooperation with the Schweizerische Unfall Versicherungs Anstalt, SUVA (Swiss Insurance Agency)

**Attention:**

The Auto-Grasp feature of the SUVA-Sensor Technology, which is primarily initiated by contact with the sensor, should be monitored in certain situations, such as when shaking hands, to avoid applying too much grip force.

Please always leave your prosthetic hand in slightly open position when not in use. This will protect the hand's sensors and mechanism from possible damage caused by high permanent load.

**FlexiGrip**

Thanks to the FlexiGrip feature, you can change the position of an object within the Electric Hand without bothering to open and close the hand again. Just shift the object in the prosthetic hand, and the fingers will automatically re-grasp the object – much like a natural hand.

**MyoHand VariPlus Speed**

The MyoHand VariPlus Speed is a myo-electrically controlled prosthetic hand that is characterised by a high grip speed in connection with an innovative and sensitive control concept.

Different proportional and digital control modes with 1 or 2 electrodes, linear transducer or switch allow individual adaptation to the requirements of the respective patient.

The MyoHand VariPlus Speed is based on the Ottobock DMC System (DMC = Dynamic Mode Control). This system, which has been developed by Ottobock, uses two independent measurement and control systems in order to optimally adapt the grip speed and grip force to the patient's muscle signal.

The proportional DMC control allows the amputee to control grip speed and grip force proportionally to his muscle signal. As soon as the strength of the muscle signal changes, grip speed and grip force immediately adapt to the new muscle signal.

During the fitting, your prosthetist considered six control modes to select one most suitable for your specific requirements. All six control modes are described in the appendix with examples.

**General Safety Instructions**

**Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the stated precautions in this accompanying document.**

**Driving Motor Vehicles**

An upper extremity amputee's ability to drive a motor vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre. For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork). It is imperative that the driver is able to operate the vehicle without any risk with the System Electric Hand turned off. Driving with the System Electric Hand turned ON may present a risk if the hand inadvertently moves due to unintentional muscle contraction or other causes.

**Caution:**

Before an arm prosthesis with a quick-disconnect wrist can be used to drive a vehicle, the terminal device must be positioned in such a way that a slight turn of it cannot disconnect the hand from the prosthesis.

### Distance to communication devices

If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WLAN devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour. It is recommended to maintain the following minimum distances:

- Mobile phone GSM 850/GSM 900: 0.99m
- Mobile phone GSM 1800/GSM 1900/UMTS: 0.7m
- DECT cordless phones incl. base station: 0.35m
- WLAN (routers, access points,...): 0.22m
- Bluetooth devices (third-party products not approved by Ottobock): 0.22m

- The System Electric Hand and the connected components might malfunction when near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). Setting the electrode gains as low as possible will minimise this risk. If such malfunctions occur repeatedly, please have the electrode settings checked by your prosthetist.
- Do not let foreign particles or liquids enter the System Electric Hand. The System Electric Hand should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- The System Electric Hand has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Myo-Service.
- The System Electric Hand as well as damaged components of the hand may be opened or repaired only by the certified Ottobock Myo-Service.
- Before interrupting or establishing electrical connections (such as when removing the hand from the prosthesis), the system must be disconnected from the power supply. For this purpose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.
- Do not use silicone spray when putting on the prosthetic glove; silicone spray may compromise secure fit of the glove. As a mounting aid the 633S2 Ottobock Procomfort Gel is recommended.
- Use the product only in combination with components approved by Ottobock (combination options). The use of unsuitable prosthetic components can cause the user to fall due to malfunctions or breakage of load-bearing components. Ottobock shall not accept liability in the event that the product is used with components other than those Ottobock specifies.
- Any manipulations to the prosthesis not described in the tasks given in these Instructions for Use are not permitted.
- The battery may only be handled by Ottobock Service Centres (batteries are not user-serviceable).
- When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

### ***Entry of Water and Humidity***

The electrical and mechanical systems of your System Electric Hand are not water-resistant. You must prevent water from entering the System Electric Hand. You may wash the System Electric Hand, do the dishes and similar activities only when wearing a prosthetic glove which has no tears or punctures. Be careful not to let water run over the top of the glove and enter the hand mechanism. If water enters the inside of the hand for any reason, please contact your prosthetist immediately to avoid further damage.

## Hand Switch (Fig. 2 and Fig. 3)

The System Electric Hand has an inconspicuous on/off switch inside the inner hand.

Figure 2: Hand switch position "on", the System Electric Hand is switched on.

Figure 3: Hand switch position "off", the System Electric Hand is switched off.

### ***When should I use the hand switch?***

When you'd like to maintain a certain hand position for a period of time, such as while holding silverware or while writing.

Turning the hand off prevents unintentional opening of the hand due to accidental muscle tension or electrical interference. It also extends the operating time of the battery in your prosthesis.

The hand switch can be pressed easily with the other hand, or more subtly on an armrest or against your leg.

## Quick-Disconnect Wrist (Fig. 4)

If your System Electric Hand is supplied with a quick-disconnect wrist (allowing passive pronation and supination), you can rotate it to a better position by turning the hand. This makes many objects easier to grasp. You can also easily separate the System Electric Hand from the socket (to replace it with a System Electric Greifer) by rotating the hand in either direction (approx. 360°) until you feel slight resistance (Fig. 4). Once you overcome the resistance, the System Electric Hand will disconnect from the wrist.

## Replacing the Terminal Device

Insert the quick-disconnect wrist into the lamination ring, press it in firmly and turn the System Electric Hand or Electric Greifer slightly to the right or left. The pressure activates the locking mechanism and anchors the Electric Hand or Electric Greifer into the lamination ring.

**Make sure the Electric Hand or Electric Greifer is properly attached before picking anything up.**

## Rechargeable Batteries

Only a fully charged Ottobock battery should be used to operate the System Electric Hand. Always keep a charged, back-up Ottobock battery available as a replacement.

An intelligent „battery manager“ alerts you to low voltage conditions: The hand will become slower and the grip force will decrease. Finally, when there is very little battery capacity remaining, the hand will switch off automatically to protect the battery against harmful total discharge.

Ottobock recommends that you charge or replace the battery as soon as you notice that the battery capacity is decreasing. Simply push the locking lever away from the battery, remove the depleted battery and replace it with a new one. When re-inserting a battery, the locking lever automatically snaps in to secure it.

 We recommend operating the System Electric Hand with the 757B20/757B21 EnergyPack, in order to be able to take advantage of the Hand's full capacity over an extended period of time. For detailed instructions on how to handle the batteries, please refer to the instructions for use enclosed with the batteries.

To recharge the batteries, use only the Ottobock Battery Chargers mentioned in the instructions for use.

### **Attention!**

As soon as the System Electric Hand recognises a fully charged 757B20/757B21 Energy-Pack or 757B35=\* MyoEnergy Integral, it automatically switches to Li-Ion battery technology. Once this has happened, only this technology should be used to operate the system. If an 757B15 Ottobock X-ChangePack is employed after the control has switched to Li-Ion battery technology, the battery can not be used to its full capacity. The 757T13 Ottobock MyoSelect can be used, however, to reset the System Electric Hand control for use with X-ChangePack. For operational security and reliability, only 757B20, 757B21, 757B35=\* or 757B15 Ottobock Batteries are to be used.

As delivered, the System Electric Hand's default setting is optimised for the X-ChangePack.

## Repair/Maintenance

Instructions for care of the prosthetic glove are enclosed with the glove. For further information, please contact your prosthetist.

This product may not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

## Description of Control Modes – SensorHand Speed

**Notice: A brief "OPENING" signal will immediately stop both the SensorHand Speed's Auto-Grasp response and the FlexiGrip function at any time.**

### Control Mode 1: *DMC plus Auto-Grasp*

**Requires: 2 electrodes**

#### ► Control with two electrodes

This control mode corresponds to the DMC plus control with integrated "virtual hand switch" (Dynamic Mode Control), and additionally features the Auto-Grasp SUVA Sensor Technology. The grip speed and grip force are determined by the strength of the EMG muscle signal. After gripping once with maximum force, the EMG signal required to open the hand will be set to a higher value ("virtual hand switch"). This reduces the risk of the hand opening due to unintentional muscle contraction. Gripping security for situations such as when holding silverware will be increased.

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** Proportionally controlled by the CLOSING electrode.

**Example 1:** A gentle muscle signal applies the smallest grip force (10 N) to the object. If the sensor recognises that the object is slipping, it automatically increases the grip force up to 50% higher than the initial grip force (15 N), as necessary.

FlexiGrip becomes active when a 20 N force is applied when moving an object within the hand. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

**Example 2:** A stronger muscle signal creates a more powerful grip force. If the sensor recognises that the object is slipping, it automatically increases the grip force until the object is held firmly (up to max. 130 N).

FlexiGrip becomes active when a 130 N force is applied when moving an object within the hand. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

### Control Mode 1: *DMC plus Auto-Grasp: Two electrodes*

Opening	Closing
Myo-signal to the electrode Opening speed Proportional 15 mm/s - 300 mm/s	Myo-signal to the electrode Closing speed Proportional 15 mm/s - 300 mm/s

Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
--------------------	--------------------------------	--------------------

Proportional: 0 N - 100 N	Proportional: up to 50% higher than initial grip force e.g. initial grip force 10 N grip force increases as necessary, up to 15 N	Occurs when repositioning force is higher than the present grip force Effective from 20 - 130 N
------------------------------	--	---

The patient can always overcome the Auto-Grasp feature with a stronger EMG muscle signal and increase the grip force up to the maximum of 100 N.

## **Control Mode 2: *AutoControl - LowInput***

**Requires: 2 electrodes,  
1 electrode and 1 switch or  
1 switch**

The hand closes at maximum speed and grasps an object with the minimal amount of grip force (10 N). If the sensor recognises that the object is slipping, it automatically increases the grip force up to the maximum (100 N), as necessary.

FlexiGrip becomes active when the maximum grip force is applied when moving an object within the hand. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

### **► Control with two electrodes**

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** At highest speed by short muscle signal of any level above the ON-threshold to the CLOSING electrode.

### **Control Mode 2: *AutoControl - LowInput: Two electrodes***

<b>Opening</b>	<b>Closing</b>
Myo-signal to the electrode	Myo-signal to the electrode digitally (short signal at any level)
Opening speed Proportional 15 mm/s - 300 mm/s	Closing speed Constant 300 mm/s

<b>Initial grip force</b>	<b>Auto-Grasp grip force increase</b>	<b>FlexiGrip function</b>
10 N	Up to 130 N	Effective at maximum grip force

### **► Control with one electrode and one switch**

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** At highest speed by short operation of switch.

## Control Mode 2: *AutoControl - LowInput: One electrode and one switch*

Opening		Closing
Myo-signal to the electrode		Signal from the CLOSE-side of the switch. SensorHand Speed closes.
Opening speed Proportional 15 mm/s - 300 mm/s		Closing speed Constant 300 mm/s
Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	Up to 130 N	Effective at maximum grip force

### ► Control with one switch

This control mode can be used with any Ottobock MyoBock switch.

**Opening:** At highest speed, as long as the OPEN-side of the switch is operated. The hand remains open.

**Closing:** At highest speed by operating the CLOSE-side of the switch.

## Control Mode 2: *AutoControl - LowInput: One switch*

Opening		Closing
Hand opens as long as the OPEN-side of the switch is operated.		Signal from the CLOSE-side of the switch: Hand closes.
Opening speed Constant 300 mm/s		Closing speed Constant 300 mm/s
Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	Up to 130 N	Effective at maximum grip force

## Control Mode 3: *AutoControl*

---

**Requires:** 1 electrode or  
1 linear transducer or  
1 switch

The hand closes at maximum speed and grasps an object with the minimal amount of grip force (10 N). If the sensor recognises that the object is slipping, it automatically increases the grip force up to the maximum (130 N), as necessary.

FlexiGrip becomes active when a 130 N force is applied when moving an object within the hand. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

### ► Control with one electrode

**Opening:** At maximum speed with a quick, constant muscle signal to the electrode.

**Closing:** At maximum speed by quick **relaxation** of muscle.

**Stopping:** By very slow muscle relaxation to the electrode, the hand stops in an opened position.

Example 1: After opening the hand, the patient very slowly relaxes the muscle.

The opening position remains unchanged.

Example 2: After opening the hand, the patient relaxes the muscle at maximum speed. The hand automatically closes at maximum speed and begins to grasp the object with 10 N grip force.

If the sensor recognises that the object is slipping, it automatically increases the grip force up to the maximum (130 N), as necessary.

### Control Mode 3: *AutoControl*: One electrode

Opening	Closing
Quick, sustained myo-signal to the electrode.	<u>Very slow</u> muscle relaxation to the electrode: Hand stops in open position.
	Quick muscle relaxation to the electrode: Hand closes.
Opening speed Constant 300 mm/s	Closing speed Constant 300 mm/s

Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	Up to 130 N	Effective at maximum grip force

### ► Control with one linear transducer

**Opening:** At maximum speed by a quick pull on the linear transducer.

**Closing:** At maximum speed by quick release of the linear transducer.

**Stopping:** By very slow release of the linear transducer, the hand stops in an opened position.

Example 1: After opening the hand, the patient very slowly releases the linear transducer.

The opening position remains unchanged.

Example 2: After opening the hand, the patient releases the linear transducer at maximum speed.

The hand automatically closes at maximum speed and begins to grasp the object with 10 N grip force.

If the sensor recognises that the object is slipping, it automatically increases the grip force up to the maximum (130 N), as necessary.

Opening	Closing
Pull on the linear transducer with high speed.	<u>Very slow</u> release of the linear transducer. Hand stops in open position.
	Quick release of the linear transducer. Hand closes.
Opening speed Constant 300 mm/s	Closing speed Constant 300 mm/s

Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	Up to 130 N	Effective at maximum grip force

## ► Control with one switch

**Opening:** At highest speed, as long as the switch is operated.

**Closing:** After the switch is released, the hand automatically closes at maximum speed and begins to grasp the object with 10 N grip force.

### Control Mode 3: *AutoControl*: One switch

Opening	Closing
Opens as long as the switch is operated.	Automatically closes as soon as the switch is released.
Opening speed Constant 300 mm/s	Closing speed Constant 300 mm/s

Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	Up to 130 N	Effective at maximum grip force

### Control Mode 4: *VarioControl*

---

**Requires:** 1 electrode or 1 linear transducer

## ► Control with one electrode

With this control mode, the opening speed is determined by the level and speed of the muscle signal. The closing speed depends on the reduction of muscle tension.

FlexiGrip becomes active when the maximum grip force is applied when moving an object within the hand. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

**Opening:** Proportional. The speed and strength of the **muscle tension** determine the opening speed of the hand.

**Closing:** Proportional. The speed and amount of the **muscle relaxation** determine the closing speed of the hand. This also affects the level of the maximal grip force.

**Stopping:** By very slow muscle relaxation to the electrode, the hand stops in an opened position.

Example 1: After opening the hand, the patient slowly relaxes his muscle. The hand closes slowly, analogous to the duration of the muscle relaxation. The hand will grasp gently at 10 N and there is no other Auto-Grasp function.

Example 2: After opening the hand, the patient relaxes the muscle at maximum speed. The hand closes at maximum speed and begins to grasp an object with 10 N grip force. If the sensor recognises that the object is slipping, it automatically increases the grip force up to the maximum (130 N), as necessary.

#### Control Mode 4: *VarioControl*: One electrode

Opening	Closing
By speed and strength of the muscle tension to the electrode	By speed and amount of the muscle relaxation to the electrode
Opening speed Proportional 15 mm/s - 300 mm/s	Closing speed Proportional 15 mm/s - 300 mm/s

Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	At low to moderate closing speed  None	Effective at 20 N
10 N	At moderate to high closing speed  Up to 130 N	Effective at maximum grip force

#### ► Control with one linear transducer

**Opening:** Proportional. The speed and strength of the pull on the linear transducer determine the opening speed of the hand.

**Closing:** Proportional. The speed of the release of the linear transducer determines the closing speed. This also affects the level of the maximal grip force.

**Stopping:** By very slow release of the linear transducer, the hand stops in an opened position.

Example 1: After opening the hand, the patient slowly releases the linear transducer.

The hand closes slowly, analogous to the duration of the muscle relaxation. The hand will grasp gently at 10 N and there is no other Auto-Grasp function.

Example 2: After opening the hand, the patient releases the linear transducer at maximum speed.

The hand closes at maximum speed and begins to grasp an object with 10 N grip force. If the sensor recognises that the object is slipping, it automatically increases the grip force up to the maximum (130 N), as necessary.

#### Control Mode 4: *VarioControl*: One linear transducer

Opening	Closing
By speed and strength of the pull on the linear transducer	By speed of release of the linear transducer
Opening speed Proportional 15 mm/s - 300 mm/s	Closing speed Proportional 15 mm/s - 300 mm/s

Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	At low to moderate closing speed None	At low to moderate closing speed Effective at 20 N
10 N	At moderate to high closing speed Up to 130 N	At moderate to high closing speed Effective at maximum grip force

## Control Mode 5: *VarioDual*

Requires: 2 electrodes

### ► Control with two electrodes

With this control mode, the opening speed is determined by the strength and speed of the **muscle contraction**. The closing speed until the minimal grip force of approx. 10 N has been reached depends on the speed of the **muscle relaxation**. The following or simultaneous muscle signal to the second electrode determines the grip force.

FlexiGrip occurs when the repositioning force is higher than the present grip force. Once the object has been repositioned, the SensorHand Speed reverts to the previous grip force.

#### **Electrode 1:**

**Opening:** Proportional. The speed and strength of the **muscle tension** determine the opening speed of the hand.

**Closing:** Proportional. The speed and amount of the **muscle relaxation** determine the closing speed of the hand. The grip force is approx. 10 N.

**Stopping:** By **very slow muscle relaxation** to the electrode, the hand stops in an opened position.

#### **Electrode 2:**

**Grasping:** **The level of the muscle signal to the second electrode determines the build-up of the grip force.** The maximal grip force is approx. 100 N.

Example 1: After opening the hand, the patient relaxes the muscle at any speed. The hand closes proportionally to the speed of the muscle relaxation. The hand will grasp an object with minimal grip force (10 N). FlexiGrip becomes active when a 20 N force is applied when moving an object within the hand. Once the object has been repositioned, the SensorHand Speed reverts to minimal grip force.

Example 2: After grasping as described in example 1, the hand shall grasp the object with a higher grip force. Create a muscle signal to the second electrode. Proportionally, a grip force between 10 N and 100 N can be built up. If the object is slipping, the grip force increases up to 50% higher than the preset grip force. FlexiGrip becomes active when approx. twice the value of the preset grip force is applied when moving an object within the hand (max. 130 N). Once the object has been repositioned, the SensorHand Speed reverts to the original grip force.

## Control Mode 5: *VarioDual*: Two electrodes

Opening	Closing
By speed and strength of the muscle tension to the 1st electrode	By speed and amount of the muscle relaxation to the 1st electrode
Opening speed Proportional 15 mm/s - 300 mm/s	Closing speed Proportional 15 mm/s - 300 mm/s
	Build-up of grip force: The grip force depends on the level of the muscle signal to the 2nd electrode. Grip force: Proportional 10 N - 100 N

Initial grip force	Auto-Grasp grip force increase	FlexiGrip function
10 N	While closing: No grip force increase	At 20 N
Proportional: 10 N - 100 N	While building up grip force: Proportional, up to 50% higher than the initial grip force	Occurs when repositioning force is higher than the present grip force Effective from 20 - 130 N

## Control Mode 6: *DMC plus Auto-Grasp, can be turned off*

**Requires: 2 electrodes**

### ► Control with two electrodes

This control mode corresponds to control mode 1; the SUVA-Sensor Technology and FlexiGrip function can be deactivated, however.

### ***Turning the SUVA-Sensor and FlexiGrip function on and off***

To grasp very soft and flexible objects such as very flexible foam materials or tweezers, the SUVA-Sensor can be deactivated. Open the SensorHand Speed as far as possible and keep open with a muscle signal of any strength. At the same time, apply slight pressure on the SUVA-Sensor (Fig. 1), e.g. pressure against a table edge. A short vibration signal will confirm that the SUVA-Sensor has been switched off.

To turn on the SUVA-Sensor, repeat the procedure. Two short vibration signals will confirm reactivation of the SUVA-Sensor.

**Notice:** Please note, that with the SUVA-Sensor turned off, the grip force will not increase automatically if an object slips off your hand.

After inserting the battery, vibration signals will inform you on the current state.

One vibration signal: Sensor is switched off

Two vibration signals: Sensor is switched on

# Description of Control Modes – MyoHand VariPlus Speed

## Control Mode 1: *DMC plus*

---

**Requires: 2 electrodes**

### ► Control with two electrodes

The grip speed and grip force are determined by the strength of the EMG muscle signal. With another, stronger electrode signal, the amputee can always increase the grip force up to the maximum of approx. 100 N.

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** Proportionally controlled by the CLOSING electrode.

Example 1: A gentle electrode signal applies the smallest grip force to the object.

Example 2: A stronger electrode signal creates a more powerful grip force. Another, stronger electrode signal increases the grip force up to the maximum of approx. 100 N.

### Control Mode 1: *DMC plus*: Two electrodes

Opening	Closing
Sustained electrode signal	Sustained electrode signal
	Gripping once with maximum force prevents the hand from opening due to unintentional electrode signals.
<b>Speed:</b> Proportional	<b>Speed:</b> Proportional

## Control Mode 2: *AutoControl - LowInput*

---

**Requires: 2 electrodes,  
1 electrode and 1 switch or  
1 switch**

The hand closes at constant speed and the grip force increases proportionally to the duration of gripping.

### ► Control with two electrodes

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** At constant speed by a muscle signal of any level above the ON-threshold to the CLOSING electrode.

### Control Mode 2: *AutoControl - LowInput*: Two electrodes

Opening	Closing
Sustained electrode signal	Sustained electrode signal
	Gripping once with maximum force prevents the hand from opening due to unintentional electrode signals.
<b>Speed:</b> Proportional	<b>Speed:</b> Constant

## ► Control with one electrode and one switch

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** At constant speed by operation of the MyoBock switch.

### Control Mode 2: *AutoControl – LowInput: One electrode and one switch*

Opening	Closing
Sustained electrode signal	Signal from the switch
	Gripping once with maximum force prevents the hand from opening due to unintentional electrode signals.
<b>Speed:</b> Proportional	<b>Speed:</b> Constant

## ► Control with one switch

This control mode can be used with any MyoBock switch.

**Opening:** At constant speed, as long as the OPEN-side of the switch is operated. The hand remains open.

**Closing:** At constant speed by operating the CLOSE-side of the switch.

### Control Mode 2: *AutoControl - LowInput: One switch*

Opening	Closing
Hand opens as long as the OPEN-side of the switch is operated.	Hand closes as long as the CLOSE-side of the switch is operated.
<b>Speed:</b> Constant	<b>Speed:</b> Constant

## Control Mode 3: *VarioControl*

---

**Requires: 1 electrode or  
1 linear transducer**

With this control mode, the opening speed is determined by the strength and speed of the **muscle contraction**. The closing speed and grip force depend on the reduction of **muscle tension**.

## ► Control with one electrode

**Opening:** Proportional. The opening speed is determined by the speed and strength of the **muscle contraction**.

**Closing:** Proportional. The closing speed is determined by the speed and amount of the **muscle relaxation**.

**Stopping:** By very slow **muscle relaxation** to the electrode, the hand stops in an opened position.

Example 1: After opening the amputee slowly relaxes the muscle. Closing takes place proportional to the decline of the electrode signal. The object is grasped with low force.

Example 2: After opening the amputee relaxes the muscle at maximum speed. The hand closes at maximum speed and begins to grasp an object with maximum grip force (approx. 100 N).

### Control Mode 3: *VarioControl*: One electrode

Opening	Closing
Increasing electrode signal	Declining electrode signal
	<b>Grip force:</b> Proportional to the decline of the electrode signal
	Gripping once with maximum force prevents the hand from opening due to unintentional electrode signals.
<b>Speed:</b> Proportional	<b>Speed:</b> Proportional

#### ► Control with one linear transducer

**Opening:** Proportional. The speed and strength of the pull on the linear transducer determine the opening speed of the MyoHand VariPlus Speed.

**Closing:** Proportional. The speed and amount of the release of the linear transducer determine the closing speed.

**Stopping:** By very slow release of the linear transducer, the hand stops in an opened position.

Example 1: After opening the amputee slowly releases the pull on the linear transducer. The object is grasped with low force.

Example 2: After opening the amputee releases the pull on the linear transducer at maximum speed. The hand closes at maximum speed and grasps an object with of grip force of approx. 100 N.

### Control Mode 3: *VarioControl*: One linear transducer

Opening	Closing
Increasing electrode signal	Declining electrode signal
	<b>Grip force:</b> Proportional to the decline of the electrode signal
	Gripping once with maximum force prevents the hand from opening due to unintentional electrode signals.
<b>Speed:</b> Proportional	<b>Speed:</b> Proportional

### Control Mode 4: *VarioDual*

---

**Requires: 2 electrodes**

#### ► Control with two electrodes

With this control mode, the opening speed is determined by the strength and speed of the **muscle contraction**. The closing speed depends on the speed of the **muscle relaxation**. The following or simultaneous muscle signal to the second electrode determines the grip force.

**Electrode 1:**

**Opening:** Proportional. The opening speed is determined by the speed and strength of the **muscle contraction**.

**Closing:** Proportional. The closing speed is determined by the speed and amount of the **muscle relaxation**.

**Stopping:** By very slow muscle relaxation to the electrode, the hand stops in an opened position.

**Electrode 2:**

**Grasping:** The level of the muscle signal to the second electrode determines the build-up of the grip force.

Example 1: After opening the hand, the patient relaxes the muscle at any speed. The hand closes proportionally to the speed of the muscle relaxation. The object is grasped with minimal force.

Example 2: After grasping as described in example 1, the hand shall grasp the object with a higher grip force. For this purpose, the amputee must create an electrode signal to the second electrode. Proportionally, a grip force of up to approx. 100 N can be built up.

**Control Mode 4: VarioDual: Two electrodes**

Opening	Closing
Increasing electrode signal through muscle contraction to the first electrode	Declining electrode signal through muscle relaxation to the first electrode
	<b>Grip force:</b> Proportional to the signal strength to the second electrode
	Gripping once more with maximum force to the 2nd electrode prevents the hand from opening due to unintentional electrode signals.
<b>Speed:</b> Proportional	<b>Speed:</b> Proportional

**Control Mode 5: DigitalControl**

**Requires: 2 electrodes,  
1 electrode and 1 switch or  
1 switch**

The hand closes at constant speed, and the grip force increases proportionally to the duration of gripping.

► **Control with two electrodes**

**Opening:** Digital. Controlled by the OPENING electrode. The hand opens at constant speed.

**Closing:** Digital. Controlled by the CLOSING electrode. The hand closes at constant speed.

**Control Mode 5: DigitalControl: Two electrodes**

Opening	Closing
Sustained electrode signal	Sustained electrode signal
	<b>Grip force:</b> Depends on the duration of the signal
<b>Speed:</b> Constant	<b>Speed:</b> Constant

► **Control with one electrode and one switch**

**Opening:** Electrode signal to the OPENING electrode.

**Closing:** Signal from the switch The hand closes at constant speed.

### Control Mode 5: *DigitalControl*: One electrode and one switch

Opening	Closing
Sustained electrode signal	Signal from the switch
	<b>Grip force:</b> Depends on the duration of the signal
<b>Speed:</b> Constant	<b>Speed:</b> Constant

#### ► Control with one switch

This control mode can be used with any MyoBock switch.

**Opening:** At constant speed, as long as the OPEN-side of the switch is operated. The hand remains open.

**Closing:** At constant speed, as long as the CLOSE-side of the switch is operated.

### Control Mode 5: *DigitalControl*: One switch

Opening	Closing
Signal from the switch	Signal from the switch
	<b>Grip force:</b> Depends on the duration of the signal
<b>Speed:</b> Constant	<b>Speed:</b> Constant

### Control Mode 6: *Double Channel Control*

---

**Requires:** 1 electrode

With this control mode, a quick and strong signal opens the hand, and a slow, gently signal closes it.

#### Control Mode 6: *Double Channel Control*: One electrode

Opening	Closing
Quick, strong electrode signal	Slow, gentle electrode signal
	<b>Grip force:</b> Depends on the duration of the signal
<b>Speed:</b> Constant	<b>Speed:</b> Constant

## Technical Data

Static current	2 mA
Operating temperature	0-70°C
Opening width	100 mm
Proportional speed	15-300 mm/s
Proportional gripping force	0 to approx. 100 N
Service life of the terminal device	5 years
Battery service life	2 years
Battery power supply:	
757B20/757B21 EnergyPack	7.2 V
757B35=* MyoEnergy Integral	7.4 V
757B15 X-ChangePack	6 V

## Ambient conditions

Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing

## Symbols used

**MD** Medical device

## Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product. The grace period is no more than one month before or three months after maintenance is due. In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover. Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty. The following components must always be sent in for maintenance and repairs: The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

## Cleaning and Care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Ensure that no liquid penetrates into the system components.

Then, dry the prosthesis component with a soft cloth.

**Liability**

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in component combinations that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

**Trademarks**

All product names mentioned in this accompanying document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners. Should trademarks used in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

**CE Conformity**

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices. The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices. The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>.

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