8E38=8\*, 8E39=8\*, 8E43=8\*

8E38=9\*, 8E39=9\*, 8E41=9\*

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	25
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	47
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	70
<b>ES</b> Instrucciones de uso (Usuario) .....	92
<b>PT</b> Manual de utilização (Usuário) .....	115
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	138
<b>SV</b> Bruksanvisning (Användare) .....	160
<b>HU</b> Használati utasítás (felhasználó) .....	182
<b>EL</b> Οδηγίες χρήσης (Χρήστης) .....	203

## INFORMATION

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### INFORMATION

Date of last update: 2022-03-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 8E38=8\*, 8E38=9\*, 8E39=8\*, 8E39=9\*, 8E43=8\* and 8E41=9\* System Electric Hands are referred to as the product/terminal device in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 SensorHand Speed with SUVA sensors

#### 2.1.1 Grip stabilisation system

The grip stabilisation system of the SensorHand Speed prevents the object being held from slipping by automatically increasing the gripping force. This process is interrupted as soon as the object is secured. Sensors integrated into the SensorHand Speed make this possible.

Grip stabilisation is only possible up to the maximum value for automatic gripping force readjustment and is intended for programs 1–4. Program 6 is intended for gripping soft foams or tools.

Automatic gripping force readjustment can be stopped at any time with a brief “OPEN” impulse.


#### 2.1.2 FlexiGrip

The FlexiGrip function allows a grasped object to be turned or shifted in the hand without having to loosen the grip via electrode signals and subsequently re-grip the object.

The SensorHand Speed follows the changes in the position of the gripped objects like a natural hand would. This makes the grip flexible.

The FlexiGrip function can be stopped at any time with a brief “OPEN” impulse.

#### 2.1.3 Unlocking/locking the wrist joint (8E43=\*)

The individual flexion and extension of the wrist joint can be locked in five different positions (in 20° increments each).

- 1) Press the release button in the direction indicated by the arrow.
- 2) Move the terminal device to the desired position while holding the release button. From the relaxed hand position, it engages at 20° and 40° in each direction.
- 3) Let go of the release button to lock the terminal device in the respective position.

#### 2.1.4 Program description

Six different control variants can be configured by the O&P professional with this terminal device.

##### Overview of control variants

- Program 1: DMC plus Sensorik (control with two electrodes)
- Program 2: AutoControl LowInput (control with two electrodes, with one electrode and one switch, or with one switch only)
- Program 3: AutoControl (control with one electrode, with one linear control element or with one switch)
- Program 4: VarioControl (control with one electrode or with one linear control element)
- Program 5: VarioDual (control with two electrodes)
- Program 6: DMC plus Sensorik, can be deactivated (control with two electrodes)

##### 2.1.4.1 Program 1: DMC plus Sensorik

Program 1 is most frequently used. The key functional characteristics are described below.

- The product’s opening and closing speed is controlled by the strength of the muscle contractions. Gripping is straightforward and appears natural as a result.
- The strength of the muscle contractions determines the force used to hold an object in the product. Small and fragile objects can therefore be grasped with ease.

**Control with two electrodes**

This control variant corresponds to DMC plus control with the integrated “virtual hand switch” (dynamic mode control) but also features the “SUVA sensors” grip stabilisation system. The level of the gripping speed or gripping force is determined by the strength of the electrode signal (resulting from the muscle contraction). After gripping with maximum force, the activation threshold in the OPEN direction is raised to a higher value (“virtual hand switch”). This increase reduces the risk of opening the hand with unintended muscle signals. The grip is more reliable, for example, when holding cutlery.

Opening: Proportional via the OPEN electrode  
 Closing: Proportional via the CLOSE electrode

Example 1: With a low muscle signal, the lowest gripping force (10 N) is built up to grasp an object. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed at up to 1.5 times the initial gripping force (15 N) as needed. FlexiGrip takes effect starting at 20 N. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Example 2: A higher muscle signal produces a higher gripping force and, when the position of the gripped object changes, a readjustment is performed up to the maximum gripping force (130 N) as needed. FlexiGrip takes effect when the load on the hand exceeds 130 N (gripping force of the hand and externally applied force). When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Open	Close
Muscle signal via the electrode	Muscle signal via the electrode
Speed: Proportional, 15 mm/s to 300 mm/s	Speed: Proportional, 15 mm/s to 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
Proportional: 0 N bis 100 N	Proportional: Up to max. 1.5 times the initial gripping force, e.g. if the initial gripping force is 10 N, gripping force readjustment up to max. 15 N	Depending on the initial gripping force, takes effect slightly above the maximum gripping force readjustment in each case Min. from 20 N Max. from 130 N

The gripping force can be increased up to the maximum gripping force (100 N) at any time as needed via a stronger muscle signal, independently of the automatic gripping force readjustment.

**2.1.4.2 Program 2: AutoControl - LowInput**

**Control with: Two electrodes, one electrode and one switch or one switch**

The hand closes at the fastest speed and grips an object with the lowest gripping force (10 N). When the sensors recognise that the position of the object is changing, a readjustment is auto-

matically performed up to the maximum gripping force (130 N) as needed. FlexiGrip takes effect starting at the maximum gripping force. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

### Control with two electrodes

Opening: Proportional via the OPEN electrode.

Closing: At maximum speed with short muscle signal at any level over the ON threshold on the CLOSE electrode.

Open	Close
Muscle signal via the electrode	Muscle signal via the electrode
Speed: Proportional, 15 mm/s to 300 mm/s	Speed: Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

### Control with one electrode and one switch

Opening: Proportional via the OPEN electrode.

Closing: At maximum speed by briefly activating the switch.

Open	Close
Muscle signal via the electrode	Signal via the CLOSE side of the switch
Speed: Proportional, 15 mm/s to 300 mm/s	Speed: Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

### Control with one switch

This program can be used in combination with any MyoBock switch.

Opening: At maximum speed as long as the OPEN side of the switch is activated. The hand then remains open.

Closing: At maximum speed by activating the CLOSE contact of the switch.

Open	Close
The hand opens as long as the OPEN side of the switch is activated	Signal via the CLOSE side of the switch
Speed: Proportional, 15 mm/s to 300 mm/s	Speed: Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function

10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force
------	------------------	-----------------------------------------------------

### 2.1.4.3 Program 3: AutoControl

#### Control with: One electrode or one linear control element or one switch

The hand closes at the fastest speed and grips an object with the lowest gripping force (10 N). When the sensors recognise that the position of the object is changing, a readjustment is automatically and continuously performed up to the respective gripping force required (max. 130 N). FlexiGrip takes effect when the load on the hand exceeds 130 N (gripping force of the hand and externally applied force). When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

#### Control with one electrode

Opening: At maximum speed with fast, sustained muscle signal via the electrode.

Closing: At maximum speed with quick relaxation of the muscle.

Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.

Example 1: After opening, relax the muscle at very slow speed. The opening position remains unchanged.

Example 2: After opening, relax the muscle at the fastest speed. The hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

Open	Close
Fast, sustained muscle signal via the electrode.	Relaxing the muscles very slowly via the electrode:
Speed:	Hand stops in the open position.
Constant, 300 mm/s	Relaxing the muscles quickly via the electrode:
	Hand closes
	Speed:
	Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

#### Control with one linear control element

Opening: At maximum speed with a quick pull on the linear control element.

Closing: At maximum speed by quickly relaxing the pull on the linear control element.

Stopping: Relaxing the pull on the linear control element very slowly causes the hand to stop in the open position.

Example 1: After opening, relax the pull on the linear control element very slowly. The opening position remains unchanged.

Example 2: After opening, relax the pull on the linear control element at the fastest speed. The hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the

object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

Open	Close
Pull on the linear control element at high speed.  Speed: Constant, 300 mm/s	Relaxing the pull on the linear control element very slowly: Hand stops in the open position Relaxing the pull on the linear control element quickly: Hand closes Speed: Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

### Control with one switch

- Opening: At maximum speed as long as the switch is activated.  
 Closing: After releasing the switch, the hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N.

Open	Close
Opens as long as the switch is activated. Speed: Constant, 300 mm/s	Closes automatically as soon as the switch is released. Speed: Constant, 300 mm/s

### 2.1.4.4 Program 4: VarioControl

#### Control with: One electrode or one linear control element

With this program, the opening speed is determined by the level and speed of the muscle contraction. The closing speed depends on the reduction of muscle tension. FlexiGrip takes effect starting at the maximum gripping force. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

#### Control with one electrode

- Opening: Proportional. The opening speed is determined by the speed and strength of the muscle contraction.  
 Closing: Proportional. The closing speed is determined by the speed and strength with which the muscles are relaxed. This also specifies the level of the maximum gripping force readjustment.  
 Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.  
 Example 1: After opening, relax the muscle at slow speed. The hand is closed analogue to the duration of muscle relaxation, at slow speed. The object is gripped with low force (10 N). There is no automatic gripping force readjustment.  
 Example 2: After opening, relax the muscle at the fastest speed. The hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

<b>Open</b>	<b>Close</b>
Via speed and strength of the muscle contraction on the electrode. Speed: Proportional, 15 mm/s to 300 mm/s	Via speed and strength of muscle relaxation on the electrode. Speed: Proportional, 15 mm/s to 300 mm/s

<b>Initial gripping force</b>	<b>Automatic gripping force readjustment</b>	<b>FlexiGrip function</b>
10 N	At low to moderate closing speed: None	Effective from 20 N
10 N	At moderate to high closing speeds: Up to max. 130 N	Takes effect starting at the maximum gripping force

### **Control with one linear control element**

- Opening:** Proportional. The opening speed is determined by the speed and strength of the pull on the linear control element.
- Closing:** Proportional. The closing speed is determined by the speed at which the pull on the linear control element is relaxed. This also specifies the level of the maximum gripping force readjustment.
- Stopping:** Relaxing the pull on the linear control element very slowly causes the hand to stop in the open position.
- Example 1:** After opening, relax the pull on the linear control element slowly. The hand is closed analogue to the duration of muscle relaxation, at slow speed. The object is gripped with low force (10 N). There is no automatic gripping force readjustment.
- Example 2:** After opening, relax the pull on the linear control element at the fastest speed. The hand closes at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

<b>Open</b>	<b>Close</b>
Via speed and strength of the pull on the linear control element Speed: Proportional, 15 mm/s to 300 mm/s	Via speed at which the pull on the linear control element is relaxed Speed: Proportional, 15 mm/s to 300 mm/s

<b>Initial gripping force</b>	<b>Automatic gripping force readjustment</b>	<b>FlexiGrip function</b>
10 N	At low to moderate closing speed: None	At low to moderate closing speed: Effective from 15 N
10 N	At moderate to high closing speeds: Up to max. 130 N	At moderate to high closing speeds: Effective starting at the maximum gripping force

### **2.1.4.5 Program 5: VarioDual**

#### **Control with two electrodes**

With this program, the opening speed is determined by the level and speed of the muscle contraction. The closing speed up to reaching the minimum gripping force of approx. 10 N depends on the speed of muscle relaxation. The gripping force is determined by the subsequent or simultaneous muscle signal on the second electrode. Depending on the initial gripping force, FlexiGrip

takes effect slightly above the maximum gripping force readjustment in each case. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

### Electrode 1

- Opening: Proportional. The opening speed is determined by the speed and strength of the muscle contraction.
- Closing: Proportional. The closing speed is determined by the speed and strength with which the muscles are relaxed. The gripping force is approx. 10 N.
- Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.

### Electrode 2

- Gripping: The build-up of the gripping force is determined by the level of the muscle signal on the second electrode. The maximum gripping force is approx. 100 N.
- Example 1: After opening, relax the muscle at any speed. Closing is proportional to the speed of muscle relaxation. The object is grasped with the minimum gripping force (10 N). FlexiGrip takes effect starting at 20 N. When there is no longer a load, the SensorHand Speed grips again with the minimum gripping force.
- Example 2: After grasping as in example 1, the object is to be gripped with a higher gripping force. To do so, produce a muscle signal on the second electrode. A gripping force between 10 N and 100 N can be built up proportionally. When the position of the gripped object changes, the gripping force is increased by up to approx. 1.5 times the specified gripping force. FlexiGrip takes effect at about twice the value of the specified gripping force up to max. 130 N; when there is no longer a load, the SensorHand Speed grips again with the original gripping force.

Open	Close
Via speed and strength of the muscle contraction on electrode 1 Speed: Proportional, 15 mm/s to 300 mm/s	Via speed and strength of muscle relaxation on electrode 1 Speed: Proportional, 15 mm/s to 300 mm/s Build-up of gripping force: Gripping force depending on the level of the muscle signal on electrode 2. Gripping force: Proportional, 10 N to 100 N

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N Proportional: 10 N to 100 N	During closing: No gripping force readjustment During build-up of gripping force: Proportional, up to max. 1.5 times the initial gripping force	From 20 N depending on the initial gripping force, takes effect slightly above the maximum gripping force readjustment in each case Min. from 20 N Max. from 130 N

#### 2.1.4.6 Program 6: DMC plus Sensorik, can be deactivated

##### Control with two electrodes

This control variant corresponds to program 1, but the “SUVA sensors” and FlexiGrip function can be temporarily deactivated.

##### Activating and deactivating the “SUVA sensors” and FlexiGrip function

The “SUVA sensors” can be deactivated for gripping very soft and pliable objects, such as very soft foams or tweezers. To do so, open the SensorHand Speed to the stop and keep it open with

a muscle signal at any level. At the same time, apply some pressure to the “SUVA sensors” (Fig. 1), e.g. by pushing on the edge of a table. A short vibration signal confirms deactivation. Repeat the same process to activate the “SUVA sensors”. Two short vibration signals confirm activation of the “SUVA sensors”.

**INFORMATION:** Note that when the “SUVA sensors” are deactivated, the gripping force is not readjusted automatically and gripped objects may slip from the grasp. After the battery is inserted, vibration signals inform you of the current mode.

One vibration signal: Sensors are deactivated

Two vibration signals: Sensors are activated

## 2.2 MyoHand VariPlus Speed

The MyoHand VariPlus Speed is a prosthetic hand with myoelectric control, featuring a particularly high gripping speed combined with an innovative, sensitive control concept. Various proportional and digital programs for control are available, e.g. with one or two electrodes, a linear control element or a switch. These permit individual adaptation to the needs of the respective patient. The terminal device is based on the Ottobock DMC system (DMC = dynamic mode control). This system uses two independent measurement and control systems to optimally adapt the gripping speed and gripping force to the muscle signal. Proportional control with DMC allows the gripping speed and gripping force to be regulated in proportion to the level of the muscle signal. When the strength of the muscle signal changes, the gripping speed and gripping force immediately adjust to the changed muscle signal.

### 2.2.1 Program description

Six different control variants can be configured by the O&P professional with this terminal device.

#### Overview of control variants

- Program 1: DMC plus (control with two electrodes)
- Program 2: AutoControl LowInput (control with two electrodes, with one electrode and one switch, or with one switch only)
- Program 3: VarioControl (control with one electrode or with one linear control element)
- Program 4: VarioDual (control with two electrodes)
- Program 5: DigitalControl (control with two electrodes, with one electrode and one switch, or only with one switch)
- Program 6: Double Channel Control (control with one electrode)

#### 2.2.1.1 Program 1: DMC plus

##### Control with two electrodes

The level of the gripping speed or gripping force is determined by the strength of the electrode signal (resulting from the muscle contraction). The gripping force can be increased up to the maximum gripping force (approx. 100 N) at any time as needed via a stronger electrode signal.

Opening: Proportional via the OPEN electrode.

Closing: Proportional via the CLOSE electrode.

Example 1: With a low electrode signal, the lowest gripping force is built up to grasp an object.

Example 2: A higher gripping force is generated with a higher electrode signal. Re-gripping with a stronger electrode signal increases the gripping force up to the maximum of approx. 100 N.

Open	Close
Sustained electrode signal	Sustained electrode signal
	A maximum gripping force prevents the hand from opening due to unintentional muscle signals.

Speed: Proportional	Speed: Proportional
---------------------	---------------------

### 2.2.1.2 Program 2: AutoControl - LowInput

#### Control with: Two electrodes, one electrode and one switch or one switch

The hand closes at constant speed, and the gripping force increases proportionally with the duration of gripping.

#### Control with two electrodes

Opening: Proportional via the OPEN electrode.

Closing: At constant speed with a muscle signal at any level over the ON threshold on the CLOSE electrode.

Open	Close
Sustained electrode signal	Sustained electrode signal
	A maximum gripping force prevents the hand from opening due to unintentional muscle signals.
Speed: Proportional	Speed: Constant

#### Control with one electrode and one switch

Opening: Proportional via the OPEN electrode.

Closing: At constant speed by activating the MyoBock switch.

Open	Close
Sustained electrode signal	Signal via the switch
	A maximum gripping force prevents the hand from opening due to unintentional muscle signals.
Speed: Proportional	Speed: Constant

#### Control with one switch

This program can be used in combination with any MyoBock switch.

Opening: At constant speed as long as the OPEN side of the switch is activated. The hand then remains open.

Closing: At constant speed by activating the CLOSE side of the switch.

Open	Close
The hand opens as long as the OPEN side of the switch is activated.	The hand closes as long as the CLOSE side of the switch is activated.
Speed: Constant	Speed: Constant

### 2.2.1.3 Program 3: VarioControl

#### Control with: One electrode or one linear control element

With this program, the opening speed is determined by the speed and strength of the muscle contraction. The closing speed and gripping force depend on the reduction of **muscle tension**.

#### Control with one electrode

Opening: Proportional. The opening speed is determined by the speed and strength of the muscle contraction.

Closing: Proportional. The closing speed is determined by the speed and strength with which the muscles are relaxed.

- Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.
- Example 1: After opening, relax the muscle at slow speed. Closing is proportional to the decline of the electrode signal. The object is gripped with low force.
- Example 2: After opening, relax the muscle at the fastest speed. The hand closes at the fastest speed and starts grasping the object with the maximum gripping force (approx. 100 N).

Open	Close
Increasing electrode signal	Decline of the electrode signal
	Gripping force: Proportional to the decline of the electrode signal.
	A maximum gripping force prevents the hand from opening due to unintentional muscle signals.
Speed: Proportional	Speed: Proportional

### Control with one linear control element

- Opening: Proportional. The opening speed is determined by the speed and strength of the pull on the linear control element.
- Closing: Proportional. The closing speed is determined by the speed and strength with which the pull on the linear control element is relaxed.
- Stopping: Relaxing the pull on the linear control element very slowly causes the hand to stop in the open position.
- Example 1: After opening, relax the pull on the linear control element slowly. The object is gripped with low force.
- Example 2: After opening, relax the pull on the linear control element at the fastest speed. The hand closes at the fastest speed and starts grasping the object with the maximum gripping force (approx. 100 N).

Open	Close
Increasing electrode signal	Decline of the electrode signal
	Gripping force: Proportional to the decline of the electrode signal.
	A maximum gripping force prevents the hand from opening due to unintentional muscle signals.
Speed: Proportional	Speed: Proportional

#### 2.2.1.4 Program 4: VarioDual

##### Control with: Two electrodes

With this program, the opening speed is determined by the level and speed of the muscle contraction. The closing speed depends on the speed of muscle relaxation. The gripping force is determined by the subsequent or simultaneous muscle signal on the second electrode.

##### Electrode 1:

- Opening: Proportional. The opening speed is determined by the speed and strength of the muscle contraction.
- Closing: Proportional. The closing speed is determined by the speed and strength with which the muscles are relaxed.

Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.

**Electrode 2:**

Gripping: The build-up of the gripping force is determined by the level of the muscle signal on the second electrode.

Example 1: After opening, relax the muscle at any speed. Closing is proportional to the speed of muscle relaxation. The object is grasped with the minimum gripping force.

Example 2: After grasping as in example 1, the object is to be gripped with a higher gripping force. To do so, the patient generates an electrode signal on the second electrode. A gripping force of up to approx. 100 N can be built up proportionally.

Open	Close
Increasing electrode signal via muscle contraction on electrode 1	Declining electrode signal via muscle relaxation on electrode 1
	Gripping force: Proportional to the signal strength on electrode 2
	Re-gripping with the maximum gripping force on electrode 2 prevents the hand from opening due to unintentional electrode signals.
Speed: Proportional	Speed: Proportional

**2.2.1.5 Program 5: DigitalControl**

**Control with: Two electrodes, one electrode and one switch or one switch**

The hand closes at constant speed, and the gripping force increases proportionally with the duration of gripping.

**Control with two electrodes**

Opening: Digital. Via the OPEN electrode. The hand opens at constant speed.

Closing: Digital. Via the CLOSE electrode. The hand closes at constant speed.

Open	Close
Sustained electrode signal	Sustained electrode signal
	Gripping force: Depending on the duration of the signal
Speed: Constant	Speed: Constant

**Control with one electrode and one switch**

Opening: Electrode signal via the OPEN electrode.

Closing: Signal via the switch. The hand closes at constant speed.

Open	Close
Sustained electrode signal	Signal via the switch
	Gripping force: Depending on the duration of the signal
Speed: Constant	Speed: Constant

**Control with one switch**

This program can be used in combination with any MyoBock switch.

Opening: At constant speed as long as the OPEN side of the switch is activated. The hand then remains open.

Closing: At constant speed as long as the CLOSE side of the switch is activated.

Open	Close
Signal via the switch	Signal via the switch
	Gripping force: Depending on the duration of the signal
Speed: Constant	Speed: Constant

### 2.2.1.6 Program 6: Double Channel Control

#### Control with: One electrode

In this program, the hand is opened with a fast and strong signal and closed with a slow, gentle signal.

Open	Close
Fast, strong electrode signal	Slow, gentle electrode signal
	Gripping force: Depending on the duration of the signal
Speed: Constant	Speed: Constant

## 3 Intended use

### 3.1 Indications for use

The product is to be used exclusively for exoprosthetic fittings of the upper limbs.

### 3.2 Conditions of use

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

Permissible ambient conditions are described in the technical data (see page 46).

### 3.3 Indications

- Amputation level below-elbow, above-elbow and shoulder disarticulation
- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm
- The user must be able to understand usage instructions and safety notices and put them into practice.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

### 3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

## 4 Safety

### 4.1 Explanation of warning symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 **WARNING**

#### The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 4.3 General safety instructions

 **WARNING**

#### Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

 **WARNING**

#### Operating a vehicle with the prosthesis

Accidents due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

 **WARNING**

#### Using the prosthesis while operating machines

Injury due to unexpected prosthesis behaviour.

- ▶ The prosthesis should not be used to operate industrial machines or motor-driven equipment.

 **WARNING**

#### Operating the prosthesis near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**⚠ WARNING**

**Skin contact with leaking lubricants due to mechanical defects**

Injury due to skin irritation.

- ▶ Avoid contact between leaking lubricants and the mouth, nose and eyes.
- ▶ The product must be inspected by an authorised Ottobock Service Center.

**⚠ CAUTION**

**Penetration of product with dirt and humidity**

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

**⚠ CAUTION**

**Independent manipulation of the product**

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

#### 4.4 Information on the Power Supply/Battery Charging

**⚠ WARNING**

**Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

**⚠ WARNING**

**Charging the prosthesis without taking it off**

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

**⚠ CAUTION**

**Prosthesis use when battery charge level is too low**

Injury due to unexpected prosthesis behaviour

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the prosthesis may be reduced at low ambient temperatures or due to ageing of the battery.
- ▶ Note that the actions/reactions of the terminal device become slower when the battery voltage gets very low.
- ▶ Note that only a few grips or actions will still be possible with the terminal device when the battery voltage gets very low.
- ▶ A small opening width may indicate a low battery voltage.

## 4.5 Information on Proximity to Certain Areas

### CAUTION

#### **Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

### CAUTION

#### **Operating the product in very close proximity to other electronic devices**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

### CAUTION

#### **Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

### CAUTION

#### **Wearing in extreme ambient temperatures**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see page 46).

## 4.6 Information on Use

### CAUTION

#### **Mechanical stress on the product**

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

### CAUTION

#### **Improper handling**

Injury due to faulty operation or malfunction of the product.

- ▶ Seek expert guidance in the proper use of the product.

**⚠ CAUTION**

**Improper product care**

- > Injuries due to faulty control/malfunction of the product or damage to the mechanical components
- > Damage or breakage due to brittleness of plastics caused by the use of acetone, petrol or similar solvents.
- ▶ Clean the product only as described in the section "Cleaning and care" (see page 45).
- ▶ Do not clean the product under running water.
- ▶ When using a prosthetic glove, also note the instructions for use of the glove.

**⚠ CAUTION**

**Grasping objects with incorrect gripping forces**

Injury due to unexpected product behaviour.

- ▶ Note that the gripping has to be controlled manually depending on the consistency (soft/hard) of the object being grasped.

**⚠ CAUTION**

**Overloading due to unusual activities**

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Risk of pinching between the fingertips**

Injury due to pinching parts of the body.

- ▶ Ensure that no body parts are between the fingertips when using the product.
- ▶ Ensure that no body parts are between the fingertips when closing the hand.
- ▶ Make sure the product is switched off for cleaning.

**⚠ CAUTION**

**Insufficient distance from sources of high heat**

Inflammation of the product.

- ▶ Do not expose the product to sources of high heat (fire, stovetop, space heater, radiator etc.).
- ▶ Do not grasp or hold any red-hot objects with the product.

**⚠ CAUTION**

**Unintentional unlocking of the terminal device**

Injury due to releasing the terminal device from the forearm (e.g. while carrying objects).

- ▶ When connecting the hand to the prosthetic socket or components, ensure that the connection is carried out correctly.

## CAUTION

### Insufficient skin contact of the electrodes

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.
- ▶ Make sure that skin contact of the electrodes continues to be given when carrying heavy loads.
- ▶ If the product cannot be controlled correctly using muscle signals, turn off the entire prosthesis and contact your O&P professional.

## CAUTION

### Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

## 5 Scope of delivery

- 1 pc. 8E38=8\*, 8E39=8\*, 8E43=8\* SensorHand Speed
- or
- 1 pc. 8E38=9\*, 8E39=9\*, 8E41=9\* MyoHand VariPlus Speed
- 1 pc. Instructions for use (user)

## 6 Use

### 6.1 Terminal device on/off

The ON/OFF button is located directly under the palm of the hand.

The ON/OFF button can be operated by pressing on the prosthetic glove.

Area	Function
Back of the hand	ON (see fig. 1)
Thumb	OFF (see fig. 2)

### 6.2 Quick-disconnect wrist

The quick-disconnect wrist makes switching between the product and the System Electric Greifer quick and easy.

If the product is equipped with a quick-disconnect wrist, the following options are available:

- Bring the terminal device into the most favourable gripping position by rotating it (passive pronation and supination).
- Easily disconnect the terminal device from the prosthetic socket to replace it with a System Electric Greifer.

To do so, turn the terminal device once around its own axis (360°), clockwise or anticlockwise, until a slight resistance is felt. The terminal device can be removed after overcoming this resistance. (see fig. 4)

### Engaging the quick-disconnect wrist

- 1) Insert the quick-disconnect wrist of the prosthesis into the lamination ring and press it in firmly.
  - 2) Turn the prosthesis slightly to the left or right.
  - 3) Pressure triggers the locking mechanism so the hand is firmly anchored in the lamination ring.
- INFORMATION: After it engages, check that the prosthesis or System Electric Greifer is firmly attached.**

### 6.3 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

#### 6.3.1 757L20 battery charger and 757B2\* EnergyPack

##### NOTICE

##### Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.  
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ●.
- 5) Remove the rechargeable battery after the charging process is complete.

##### 6.3.1.1 Display of the current charge level

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

LED indicator	Event
●	Charge level above 50% (illuminated in green)
● and ●	Charge level below 50% (flashing alternating green and orange)
●	Charge level below 5% (illuminated in orange)

##### 6.3.1.2 Troubleshooting

Charging socket	LED for the charging socket	Event
Empty	● or ● or ●	<b>The battery charger is defective</b> The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
Empty or rechargeable battery inserted	●	
Rechargeable battery inserted	●	<b>The rechargeable battery is charging</b>
	●	<b>The rechargeable battery is fully charged</b>

Charging socket	LED for the charging socket	Event
Rechargeable battery inserted		<b>The temperature of the rechargeable battery is too high</b> Remove the rechargeable battery and let it cool down
		<b>The inserted rechargeable battery is defective</b> The rechargeable battery must be inspected by an authorised Ottobock Service Centre.

### 6.3.2 757L35 battery charger and 757B35=\* rechargeable battery

#### 6.3.2.1 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.  
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.  
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

#### 6.3.2.2 Safety Shutoff

The purpose of the prosthesis safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

#### 6.3.2.3 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level.

Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50% charged (illuminated in yellow)
	Battery drained (illuminated in orange)

#### 6.3.2.4 Beep signals

Beep signal	Additional display	Event
1 x long	-	<ul style="list-style-type: none"> <li>• Turn off prosthesis on charging receptacle</li> <li>• Charging begins (charging plug connected to charging receptacle)</li> <li>• Charging ends (charging plug disconnected from charging receptacle)</li> </ul>

Beep signal	Additional display	Event
2 x short	LED on charging receptacle lights up briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis shuts off automatically

## 7 Storage

If the System Electric Hand will not be used, store the System Electric Hand in the open position to protect the sensors and mechanism.

## 8 Cleaning and Care

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1-N Derma Clean) when needed.  
Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

### INFORMATION

Avoid exposing the product without a prosthetic glove to direct sunlight or UV light (solarium) for extended periods.

Please see the information included with the glove for prosthetic glove care instructions. Further information is available from your O&P professional.

## 9 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

## 10 Legal information

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 11 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

Reference number	8E38=8*, 8E39=8*, 8E43=8*	8E38=9*, 8E39=9*, 8E41=9*
Opening width	100 mm	
Proportional speed	15–300 mm/s	
Proportional gripping force	0–100 N	
Lifetime	5 years	

Rechargeable battery reference number	757B20	757B21	757B35 =1	757B35 =3	757B35 =5
Temperature range for charging [°C]	+5 to +40				
Capacity [mAh]	900 <sup>1</sup> / 950 <sup>1</sup>	680 <sup>1</sup> / 800 <sup>1</sup>	600	1150	3450
Rechargeable battery lifetime [years]	2				
Product behaviour during the charging process	The product is non-functional				
Operating time of the product with fully charged battery [grip cycles]	Approx. 2500– 3000	Approx. 2000– 2500 (depend- ing on the bat- tery capacity)	Approx. 1000– 2000	Approx. 4000	Approx. 10000
Charging times (rechargeable battery fully dis- charged) [hours]	Approx. 3.5	Approx. 3	Approx. 2.5		Approx. 3
Nominal voltage [V]	Approx. 7.4				
Battery technology	Li-Ion		Li-Po	Li-Ion	
Approved battery chargers	757L20		757L35		

<sup>1</sup> See labeling on battery

## 12 Appendix

### 12.1 Symbols Used

Manufacturer

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

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