







757Z184=2, 757Z190=2

DE	Gebrauchsanweisung	4
EN	Instructions for use	8
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	16
ES	Instrucciones de uso	20
PT	Manual de utilização.....	24
NL	Gebruiksaanwijzing	28
SV	Bruksanvisning	32
HU	Használati utasítás	36

Information!

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**Information!**

Date of the last update: 2015-06-26

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

1 Application

The 757Z184=2/757Z190=2 Battery Mounting Set is to be used **solely** for exoprosthetic upper extremity fittings and installation in arm prostheses.

2 Safety instructions

**Attention!**

Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

- Use the appropriate safety gear during lamination and installation.
- Use the Patient Information section of this pamphlet (Section 3) to inform your patient about the function and use of the Battery Receptacle.

3 Patient information

- Do not let foreign particles or liquids get into the Battery Receptacle. The Battery Receptacle should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- The electrical contacts of the Battery Receptacle must be clean and free of grease to achieve optimal function, and can be cleaned with a dry cotton swab.

4 Scope of delivery (Fig. 1)

	757Z184=2 for EnergyPack 757B20	757Z190=2 for EnergyPack 757B21
Item 1, 1 pc.	757Z185=2 Battery Receptacle	757Z191=2 Battery Receptacle
Item 2, 1 pc.	757Z186 Engaging Lever	757Z186 Engaging Lever
Item 3, 1 pc.	757Z187=1 Lamination Pattern	757Z192=1 Lamination Pattern
Item 4, 1 pc.	757Z189=1 Pattern	757Z194=1 Pattern

5 Description (Fig. 2)

The 757Z184=2 und 757Z190=2 Battery Mounting Sets are used to mount battery receptacles in either a laminated socket or an Ottobock elbow unit, such as an ErgoArm elbow.

The Battery Receptacles are designed to accommodate an EnergyPack 757B20 or 757B21 to provide electric power for the prosthesis.

The integrated snap connection provides a quick way to remove the Battery Receptacle from the socket in order to, for example, do repairs or maintenance work inside the socket.

By connecting the Battery Receptacle to the 757T11 MyoCom it is possible to make customer specific adjustments on Ottobock System Electric Hands and System Electric Greifers with lamination ring, lamination plate or threaded connection while the patient is wearing the prosthesis. The System Electric Hand or Greifer can remain attached to the socket (Fig. 2). By connecting the Battery Receptacle to the 757M11 MyoBoy the muscle signals can be evaluated as well while the patient is wearing the prosthesis.

6 Installation

6.1 Installation into a laminated socket (Fig. 3 – Fig. 9)

Prepare the inner socket following the appropriate procedures, foam up and rough-grind. Set the lamination pattern into the rigid foam so that the edge rises **at least 2 mm** above the foam.

**Attention!**

The surface beneath the lamination pattern must be plane.

Create a contoured edge on the socket by applying Plastilin around the lamination pattern. Be careful not to damage the edges of the lamination pattern (Fig. 3).

Pull on a PVA bag and laminate the forearm shell in the usual manner. Once the lamination has hardened, sand down the visible edge of the lamination pattern so that it is flush (Fig. 4) with the socket.

Press in the two cut-outs in the middle of the lamination pattern with a screwdriver or similar tool (Fig. 5) and remove the lamination pattern with a pair of pliers (Fig. 6).

Connect the two electrode cables and the 4-pole cable from the terminal device to the battery receptacle (Fig. 7). Set the battery receptacle into the socket in such a way that the engaging lever is positioned proximally. Insert the battery receptacle into the socket (Fig.8) and then press into the socket in place so that it is flush.

Attach the battery receptacle to the socket by screwing in both Allen Head bolts with a SW 1.5 mm Allen Wrench (Fig. 9).

**Attention!**

The battery receptacle must not be glued into the socket!

6.2 Installation into the forearm of an Ottobock elbow component

(Fig. 10 – Fig. 11)

Use the included pattern (757Z189=1/757Z194=1) to determine the position of the battery receptacle.

**Attention!**

When using an elbow component with AFB (Automatic Forearm Balance) in conjunction with an electric wrist rotator, allow a minimum distance of 103 mm between the elbow joint and battery receptacle to avoid possible interference (Fig. 10).

Drill four holes inside the frame of the pattern and saw out the space for the battery receptacle (Fig. 11).

Connect the two electrode cables and the 4-pole connection cable from the terminal device to the battery receptacle (Fig. 7). Set the battery receptacle into the socket in such a way that the engaging lever is positioned proximally. Insert the battery receptacle into the socket (Fig. 8) and then press into the socket in place so that it is flush.

Attach the battery receptacle to the socket by screwing in both Allen Head bolts with a SW 1.5 mm Allen Wrench (Fig. 9).

Attention!

The battery receptacle must not be glued into the socket!

7 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise.

The product may only be opened and repaired by authorised Ottobock technicians.

8 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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