C-Brace joint unit 17K01=*

DE Gebrauchsanweisung (Benutzer)	3
EN Instructions for use (user)	51
FR Instructions d'utilisation (Utilisateur)	97
IT Istruzioni per l'uso (Utilizzatore)	147
ES Instrucciones de uso (Usuario)	195
PT Manual de utilização (Usuário)	243
NL Gebruiksaanwijzing (Gebruiker)	291
TR Kullanma talimatı (Kullanıcı)	339
JA 取扱説明書(ユーザー用)	385

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1	Foreword	54
2	Product description	54
2.1	Design	54
2.1.1	Control panel on the joint unit	55
2.2	Function	55
3	Application	56
3.1	Indications for use	56
3.2	Conditions of use	56
3.3	Indications	56
3.4	Contraindications	56
3.4.1	Absolute Contraindications	56
3.4.2	Relative Contraindications	56
3.5	Qualifications of the therapist or nursing staff	56
4	Safety	56
4.1	Explanation of warning symbols	56
4.2	Structure of the safety instructions	56
4.3	General safety instructions	57
4.4	Information on the Power Supply/Battery Charging	59
4.5	Battery charger information	59
4.6	Notes on applying the product	60
4.7	Information on Proximity to Certain Areas	60
4.8	Information on Use	61
4.9	Notes on the movement patterns	62
4.10	Notes on the safety modes	64
4.11	Information on the use of a mobile device with the cockpit app	65
5	Scope of Delivery and Accessories	66
5.1	Scope of delivery	66
5.2	Accessories	66
6	Charging the battery	66
6.1	Connecting the power supply and battery charger	67
6.2	Connect battery charger to the product	67
6.3	Display of the current charge level	68
6.3.1	Display of battery charge level without additional devices	68
6.3.2	Display of the current charge level using the Cockpit app	68
6.3.3	Display of the current charge level during the charging process	68
7	Cockpit app	69
7.1	System Requirements	69
7.2	Initial connection between cockpit app and prosthesis	69
7.2.1	Starting the cockpit app for the first time	70
7.3	Control elements for cockpit app	71
7.3.1	Cockpit app navigation menu	71

7.4	Managing components	71
7.4.1	Adding component	72
7.4.2	Deleting a component	72
8	Use	72
8.1	Application	73
8.2	Removal	73
8.3	Movement patterns in basic mode (mode 1)	74
8.3.1	Standing	74
8.3.1.1	Stance function	74
8.3.2	Walking.....	74
8.3.3	Sitting down	75
8.3.4	Sitting	75
8.3.5	Standing up.....	75
8.3.6	Walking down stairs	76
8.3.6.1	Orthosis construction with rigid ankle joint or dorsal spring element	76
8.3.6.2	Orthosis construction with movement in the ankle joint.....	76
8.3.7	Walking up stairs.....	76
8.3.8	Walking up a ramp	77
8.3.9	Walking down a ramp	77
8.3.10	Walking down flat steps	77
8.3.11	Kneeling.....	77
8.4	Changing orthosis settings	77
8.4.1	Changing the orthosis settings using the Cockpit app.....	78
8.4.2	Overview of adjustment parameters in basic mode.....	78
8.4.3	Overview of adjustment parameters in MyModes	79
8.5	Switching the product on/off.....	80
8.6	Switching Bluetooth of the component on/off	80
8.6.1	Switching off Bluetooth	81
8.6.2	Switching on Bluetooth	81
8.7	Querying the status of the component	81
8.7.1	Query status through cockpit app.....	81
8.7.2	Status display in the cockpit app	81
8.8	Recommendations for air travel	81
9	MyModes	81
9.1	Switching MyModes with the cockpit app	81
9.2	Basic mode	82
9.3	MyMode "Training mode"	82
9.4	MyMode "Freeze position".....	83
9.5	MyMode "User defined"	83
10	Additional operating states (modes)	83
10.1	Empty battery mode	83
10.2	Mode for charging the product.....	83
10.3	Safety mode.....	83
10.4	Overheating mode	84
11	Cleaning	84

- 12 Maintenance84**
- 13 Legal information85**
 - 13.1 Liability85
 - 13.2 Local Legal Information85
 - 13.3 CE conformity86
 - 13.4 Trademarks86
- 14 Technical data86**
- 15 Appendices88**
 - 15.1 Symbols Used.....88
 - 15.2 Operating states/error signals89
 - 15.2.1 Status display on the control panel89
 - 15.2.2 Error messages while establishing a connection with the cockpit app93
 - 15.2.3 Error while charging the product94
 - 15.3 Directives and manufacturer’s declaration94
 - 15.3.1 Electromagnetic environment94

1 Foreword

INFORMATION

Date of last update: 2020-07-13

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The product "C-Brace joint unit 17KO1=* " is referred to as the product/component/orthosis/joint unit below.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Design


The fabricated orthosis and the course of the orthosis shell edges can be individually designed. For this reason, the following illustration and the illustration on the cover page of these instructions for use only show one possible version:

1. Thigh closure straps
2. Thigh shell
3. 17KO1=* C-Brace joint unit with control panel
4. Medial joint support
5. Lower leg closure straps
6. Lower leg shell
7. Connecting element to ankle joint.
8. Ankle joint
9. Foot component

The connection can also be established with a dorsal spring element.

2.1.1 Control panel on the joint unit

1. Long button press: turns the component on/off (see Page 80)
Short button press: queries the status (see Page 89)
2. Displays the charge level of the integrated battery (see Page 68)
3. Maintenance required (see Page 89)
4. Long button press: turns the Bluetooth function on/off (see Page 80).
Short button press: if the Bluetooth function is on, the component becomes "visible" for about 2 minutes so it can be detected by another device such as a smartphone.

Under the cover at the end of the control panel: Charging receptacle to connect the battery charger (see Page 66)

2.2 Function

The product is an orthopaedic technology device that makes everyday activities such as standing and walking possible or easier for the user. If the user's activity level increases in the course of use, adjustments are required. Changes will also have to be made if the user's physical condition worsens and more support from the product is needed.

This product features microprocessor control of the stance and swing phase (SSCO).

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the extension and flexion resistance of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

Thanks to the microprocessor-controlled stance and swing phase, the product can be individually adapted to your needs.

In order to do so, the product is adjusted by the O&P professional using an adjustment app.

The product features MyModes for special motion types (e.g. cycling ...). These are pre-configured by the O&P professional using the adjustment software and can be activated using the Cockpit app (see Page 81).

In case of an error in the sensor system, hydraulic control or when the battery is empty, safety mode makes restricted operation possible. Resistance parameters that are predefined by the product are configured for this purpose (see Page 83).

The microprocessor-controlled hydraulic unit offers the following advantages

- Approximation of the physiological gait pattern
- Stability while standing and walking
- Adaptation of product characteristics to various surfaces, inclines, gait situations and walking speeds

Essential performance of the product

- Stability in the stance phase

3 Application

3.1 Indications for use

The product is intended **solely** for orthotic fittings of the lower limbs.

3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sports activities that involve jumping, sudden movements or fast sequences of steps (e.g. basketball, badminton, sporty riding).

Possible sports activities must be discussed with the O&P professional.

Permissible ambient conditions are described in the technical data (see Page 86).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

3.3 Indications

- Unilateral or bilateral lower limb paresis or flaccid paralysis, e.g. due to post-polio syndrome, traumatic paresis, incomplete paraplegia.
- Physical prerequisites such as muscle status, joint mobility and possible axis deviations are crucial and proper control of the orthosis must be guaranteed.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations
- The existing muscle strength of the hip extensors and flexors must permit the controlled swing-through of the limb (compensation using the hip is possible).

3.4 Contraindications

3.4.1 Absolute Contraindications

- Flexion contracture in the knee and/or hip joint in excess of 10°
- Knee varus/valgus malposition in excess of 10°
- Severe spasticity
- Body weight over 125 kg / 275 lbs

3.4.2 Relative Contraindications

- Moderate spasticity

3.5 Qualifications of the therapist or nursing staff

The therapists or nursing staff must be trained in handling the product. Training must be carried out by the authorised orthopaedic technician.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

WARNING

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

WARNING

Using the orthosis while operating a vehicle

Accident due to unexpected orthosis behaviour.

- ▶ All users are required to observe their country's national and state driving laws when operating vehicles with the orthosis. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.
- ▶ The leg on which the orthosis is worn may not be used to control the vehicle or its peripheral components (e.g. clutch pedal, brake pedal, gas pedal...).

WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

CAUTION

Failure to observe warning/error signals

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (Warnings/error signals) and corresponding change in damping settings must be observed.

CAUTION

Penetration of dirt and moisture into components of the orthosis

- > Falling due to unexpected behaviour of the orthosis as a result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that neither solid particles, foreign bodies nor liquids penetrate the orthosis components.
- ▶ The orthosis components are protected against splashed water from all directions.
- ▶ However, the orthosis components are not protected against submersion, jets of water and steam.
- ▶ If water has penetrated the components of the orthosis, remove the straps and pads if possible and allow the components to dry.
- ▶ In case of a malfunction after drying, the orthosis must be inspected by an authorised Ottobock Service Centre. The O&P professional is your contact person.
- ▶ Always close the charging receptacle with the protective cap after disconnecting the charging cable.

CAUTION

Independent manipulations on the joint unit and orthosis components

Falling due to the breakage of load-bearing components or malfunction of the orthosis.

- ▶ Manipulations on the joint unit and orthosis components other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The joint unit and orthosis components may only be opened and repaired by authorised, qualified Ottobock personnel.

CAUTION

Use of the product when battery charge level is too low

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the product if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

CAUTION

Mechanical stress on the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

CAUTION

Signs of wear and tear on the product components

Falling due to damage or malfunction of the product.

- ▶ Regular service inspections (maintenance) are mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty.

CAUTION

Use of unapproved accessories

- > Falling due to product malfunction as a result of reduced interference resistance.
- > Interference of other electronic devices due to increased emissions.
- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the sections "Scope of delivery" (see Page 66) and "Accessories" (see Page 66).

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product with a damp cloth only (fresh water).

INFORMATION

Movement noises of the joint unit and other joints (medial support, ankle joint)

When using orthotic joints, servo motor, hydraulic, pneumatic or brake load dependent control functions can cause movement noises. This kind of noise is normal and unavoidable. It generally does not indicate any problems. If movement noise increases noticeably during the lifecycle of

the orthosis, the orthosis should be inspected by the orthotics and prosthetics professional (O&P professional) immediately.

4.4 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the product with damaged power supply unit/charger/charger cable

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

⚠ CAUTION

Charging while wearing the product

- > Falling due to walking and getting caught on the connected battery charger.
- > Falling due to unexpected product behaviour because of changed damping behaviour.
- ▶ Only charge the product in a sitting position while wearing it.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

4.5 Battery charger information

NOTICE

Penetration of dirt and humidity into the product

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Mechanical stress on the power supply/battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

NOTICE

Operating the power supply unit/charger outside of the permissible temperature range

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see Page 86).

NOTICE

Independent changes or modifications carried out to the battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

INFORMATION

Depending on the battery charge level, the battery charger can warm up considerably during the charging process. This is not a malfunction.

4.6 Notes on applying the product

⚠ CAUTION

Foreign objects between the leg and orthosis shells

Pressure points on the leg due to foreign objects at the contact points between the leg and orthosis shells.

- ▶ Smooth out wrinkles in the padding material and clothing.
- ▶ Check the leg for pressure points.

⚠ CAUTION

Improper fit of the orthosis

Falling/skin irritation due to insufficient support/hold by the orthosis shells.

- ▶ Take the orthosis off immediately and put it back on.
- ▶ Observe the instructions for applying and removing.

⚠ CAUTION

Pinching skin in the area of the closures

Injuries and swelling due to disruption of blood circulation to the skin in the area of the closures.

- ▶ Do not over-tighten the closures when applying the product.

⚠ CAUTION

Volume fluctuations of the leg or problems with the fit of the orthosis

Injuries, friction and local pressure caused by improper fit of the orthosis shells (too firm/too loose).

- ▶ Sensitivity problems and skin damage require special attention with regard to the fit. Check skin areas daily.
- ▶ Consult a physician or O&P professional in case of any signs of skin damage, even if they are only minor.
- ▶ If the improper fit is the result of weight gain or loss, have new thigh shells and lower leg shells fabricated based on a new plaster cast.
- ▶ Check the leg for pressure points.

4.7 Information on Proximity to Certain Areas

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

⚠ CAUTION

Operating the product in very close proximity to other electronic devices

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for people (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations...).
- ▶ If this cannot be avoided, make sure to at least have a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected changes in the damping behaviour of the product.
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.

⚠ CAUTION

Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of a strong magnetic field.
- ▶ Take off the product before entering a room or area with strong magnetic fields and store the product outside this room or area.
- ▶ Damage to the product caused by exposure to a strong magnetic field cannot be repaired.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see Page 86).

4.8 Information on Use

⚠ CAUTION

Improper mode switching

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.

- ▶ Take the weight off the product and correct the switching, if required.

⚠ CAUTION

Risk of pinching in the joint flexion area

- > Injury due to pinching of body parts.
- > Damage to garments due to pinching in the mechanism of the joint unit or medial support.
- ▶ Ensure that no body parts or garments are in this area when flexing the product.

⚠ CAUTION

Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- ▶ Be sure to pay attention when pulsating vibration signals start (Warnings/error signals). They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you have to reduce the activity level so the hydraulic unit can cool down.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.
- ▶ If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. In this case, the product has to be inspected by an authorised Ottobock Service Centre. The O&P professional is your contact person.

⚠ CAUTION

Overloading due to unusual activities

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sports activities that involve jumping, sudden movements or fast sequences of steps (e.g. basketball, badminton, sporty riding).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional.

4.9 Notes on the movement patterns

⚠ CAUTION

Improper use of the stance function

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using the stance function and check the lock of the knee joint before placing your full weight on the orthosis.
- ▶ Make sure you have been instructed in the proper use of the stance function by an O&P professional or therapist. Information on the stance function see Page 74.

⚠ CAUTION

Walking up stairs

- > Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.
- > Falling due to breakage of load-bearing components.
- ▶ Always use the handrail when walking up stairs and place most of the sole of your foot on the stair surface.
- ▶ Always step up with the less affected leg, followed by the leg with the product.
- ▶ If there is a connection to a foot component, support by bouncing with the product extended while walking up stairs is not allowable and therefore has to be avoided.
- ▶ Particular caution is required when carrying children up the stairs.

⚠ CAUTION

Walking down stairs with moveable orthosis foot component

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and place most of the sole of your foot on the stair surface.
- ▶ Observe the warnings/error signals (see Page 89).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warning and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

⚠ CAUTION

Walking down stairs with rigid orthosis foot component

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and roll over the edge of the step with the middle of the shoe.
- ▶ Observe the warnings/error signals (see Page 89).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warnings and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

⚠ CAUTION

Walking down stairs and ramps

Falling due to unexpected, high stance phase damping during the transition from a level surface to stairs or ramps.

- ▶ Pay attention to the changed product behaviour.
- ▶ Verify changed stance phase damping before walking on stairs or ramps.

⚠ CAUTION

Lack of support from the product when walking down stairs

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Check whether corresponding flexion of the joint is possible before walking down stairs. If this is not possible, either the Cockpit app or turning the product off/on must be used to switch back to basic mode.

⚠ CAUTION

Improper use of the MyMode "Training mode"

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using this MyMode and check the lock of the knee joint before placing your full weight on the orthosis.
- ▶ Note that the knee joint is locked in the flexion direction in this MyMode.
- ▶ Make sure you have been instructed in the correct use of this MyMode by the O&P professional and/or therapist. For information about this mode, see Page 82.
- ▶ Switch back to basic mode once the activities in this MyMode have been completed.

⚠ CAUTION

Improper use of the MyMode "Freeze position"

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using this MyMode and check the lock of the knee joint before placing your full weight on the orthosis.
- ▶ Note that the knee joint is locked in the both the flexion and extension direction in this MyMode.
- ▶ Make sure you have been instructed in the correct use of this MyMode by the O&P professional and/or therapist. For information about this mode, see Page 83.
- ▶ Switch back to basic mode once the activities in this MyMode have been completed.

⚠ CAUTION

Lack of support from the product while sitting down

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Check whether corresponding flexion of the joint is possible before sitting down. If this is not possible, either the Cockpit app or turning the product off/on must be used to switch back to basic mode.

⚠ CAUTION

Quickly pushing the hip forward with the orthosis extended

- > Falling due to unexpected activation of the swing phase.
- ▶ Note that the joint may flex unexpectedly when the hip is pushed forward quickly while the orthosis is extended.
- ▶ Therefore, familiarise yourself with swing phase activation in such situations under secure conditions (e.g. while holding on to parallel bars) and with the instruction of trained, qualified personnel.

⚠ CAUTION

Changing the gait pattern influences activation of the swing phase

- > Falling due to unexpected activation of the swing phase.
- ▶ Changing the gait pattern can influence activation of the swing phase. Adjustment by the O&P professional has to be repeated in this case.

4.10 Notes on the safety modes

⚠ CAUTION

Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Do not continue using the defective product.
- ▶ Consult the O&P professional promptly.

CAUTION

Safety mode cannot be deactivated

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ If you cannot deactivate safety mode by recharging the battery, a permanent malfunction has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected immediately by an O&P professional.

CAUTION

Safety signal occurs (ongoing vibration)

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warnings/error signals (see Page 89).
- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. In this case, the contact person is your prosthetist.

CAUTION

Using the product in safety mode

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see Page 89) have to be observed.
- ▶ Particular caution is necessary when using a bicycle without a freewheel (with a fixed gear).

4.11 Information on the use of a mobile device with the cockpit app

CAUTION

Improper use of the mobile device

Falling due to changed damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Make sure you have been instructed in the proper use of the mobile device with the Cockpit app.

CAUTION

Independently applied changes or modifications made to the mobile device

Falling due to altered damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Do not make any independent changes to the hardware of the mobile device on which the app is installed.
- ▶ Do not make any independent changes to the software/firmware of the mobile device that are not included in the update function of the software/firmware.

CAUTION

Improper mode switching with the device

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter (beeper) and the display on the device.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.

NOTICE

Destruction of the mobile device due to falling or penetration of water

Malfunction of the mobile device.

- ▶ If necessary, let the mobile device dry at room temperature (at least 1 day).
- ▶ Should it no longer be possible to switch back from a MyMode to basic mode, the component can only be switched back to basic mode by turning it off/on (see Page 80).

NOTICE

Failure to observe the system requirements for the installation of the Cockpit app

Malfunction of the mobile device.

- ▶ Install the Cockpit app only on the operating systems listed in the section "System requirements" (see Page 69).

INFORMATION

The illustrations in these instructions for use are only examples and may deviate from the respective mobile device being used and the version.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. 757L16-4 power supply
- 1 pc. 4E50-2 C-Leg battery charger
- 1 pc. cosmetic case for battery charger and power supply
- 1 pc. 646C107 Bluetooth PIN card
- 1 pc. Orthotic passport
- 1 pc. 17KO1=L C-Brace joint unit left, or
17KO1=R C-Brace joint unit right
- 1 pc. Instructions for use (user)

5.2 Accessories

Cockpit app for download from the website: <https://www.ottobock.com/cockpitapp>

- "Cockpit 4X441-IO5=V*" app for iOS
- "Cockpit 4X441-ANDR=V*" app for Android

6 Charging the battery

The following points must be observed when charging the battery:

- Use the 757L16-4 power supply and 4E50-2 battery charger to charge the battery.
- With uninterrupted walking, the capacity of the fully charged battery is sufficient for at least 18 hours. It lasts about 2 days with average use.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.
- For the maximum operating time with one battery charge, turning the product off when it is not being used is recommended.
- Switching the MyModes and changing setting parameters using the Cockpit app is not possible during the charging process.
- After disconnecting the battery charger, the orthosis is in the same state it was in before connecting the battery charger. For example, if the orthosis was switched off before connecting the battery charger, it will also be switched off after disconnecting the battery charger.

6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
 - 2) Connect the round, **four-pin** plug of the charging cable to the **OUT** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
 - 3) Connect the round, **three-pin** plug of the power supply to the **12 V** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
 - 4) Plug the power supply unit into the wall socket.
 - The green LED on the back of the power supply (see fig. 3) and the green LED on the battery charger light up.
- If the green LED on the power supply and the green LED on the battery charger do not light up, there is an error (see Page 94).

6.2 Connect battery charger to the product

- 1) Open the charging receptacle cover.
- 2) Connect the charging plug to the charging receptacle of the product.
INFORMATION: Take note of the plug direction!
 - After the displays are tested, a short beep is emitted followed by three short vibration signals.
 - If the **!** symbol lights up, this means an error was identified during the self-test (see Page 89).
- 3) The charging process starts.
 - To check the charge level while the battery charger is connected, briefly press the **!** button on the control panel.
- 4) Disconnect the product after the charging process is complete.
 - After disconnecting the battery charger, a vibration signal is emitted followed by a short beep, and the current status is displayed for approx. 5 seconds (see Page 89).
 - When the symbol in the **!** button is lit up green **!**, the product is switched on and ready for operation.
- 5) Close the charging receptacle cover.

INFORMATION

After disconnecting the battery charger, the orthosis is in the same state it was in before connecting the battery charger. For example, if the orthosis was switched off before connecting the battery charger, it will also be switched off after disconnecting the battery charger. A descend-

ing sequence of beeps is emitted when the battery charger is disconnected while the orthosis is switched off .

INFORMATION

No display after connecting the battery charger

If no symbols light up on the control panel after connecting the battery charger, the battery may be deep discharged. Leave the battery charger connected for at least 15 minutes and check the charge level while charging by disconnecting/connecting the battery charger.

6.3 Display of the current charge level

6.3.1 Display of battery charge level without additional devices

The current charge level of the installed battery can be displayed by briefly pressing the  button:

Symbol	Battery charge level
	Charge level 67% to 100%
	Charge level 34% to 67%
	Charge level 10% to 34%
	Charge level 5% to 10%
	Charge level 1% to 5%
	Battery empty


6.3.2 Display of the current charge level using the Cockpit app

Once the Cockpit app has been started, the current charge level is displayed in the bottom line of the screen:

-  38% – Charge level of battery for currently connected component

6.3.3 Display of the current charge level during the charging process

After connecting the battery charger or pressing the  button while the battery charger is connected, the current charge level is shown by an animated symbol on the control panel in addition to the display on the battery charger (, , ).

Control panel	Battery charger	
	   	Battery is being charged, battery charge level is less than 34%
		Battery is being charged, battery charge level is 34% to 50%

Control panel	Battery charger	
		Battery is being charged, battery charge level is 50% to 67%
	The on time of the yellow LED gets shorter as the charge level increases. It only flashes briefly at the end of the charging process.	Battery is being charged, battery charge level is 67% to 99%
		Battery is fully charged

7 Cockpit app

The Cockpit app enables switching from basic mode to the pre-configured MyModes. In addition, information about the product (step counter, battery charge level, etc.) can be retrieved.

The everyday behaviour of the product can be changed to a certain extent using the app (e.g. while becoming accustomed to the product). The O&P professional can use the adjustment software to track these changes at the next appointment.

Information on the Cockpit app

- The Cockpit app can be downloaded free of charge from the respective online store. For more information, please visit the following website: <https://www.ottobock.com/cockpitapp>. To download the Cockpit app, the QR code on the supplied Bluetooth PIN card can also be read with the mobile device (requirement: QR code reader and camera).
- The language of the Cockpit app user interface can only be changed by the O&P professional using the adjustment app.
- The serial number of the component to be connected has to be registered with Ottobock the first time it is connected. If the registration is not accepted, use of the Cockpit app for this component will be limited.
- Bluetooth on the prosthesis must be turned on in order to use the Cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards ().
- Keep the mobile app up to date at all times.
- Please contact the manufacturer if you suspect cybersecurity problems.

7.1 System Requirements

The functioning of the Cockpit app is assured on mobile devices that support the following operating systems:

- **iOS (for iPhone, iPad, iPod):** version 10.0 or higher
- **Android:** version 5.0 or higher

7.2 Initial connection between cockpit app and prosthesis

The following points need to be observed before establishing the connection:

- Bluetooth of the component must be switched on (see Page 80).
- Bluetooth on the mobile device must be switched on.
- The mobile device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The mobile device must be connected to the Internet.**

- The serial number and Bluetooth PIN of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".

INFORMATION

If the Bluetooth PIN card with the Bluetooth PIN and serial number of the component is lost, contact your O&P professional.

7.2.1 Starting the cockpit app for the first time

- 1) Tap the symbol of the Cockpit app ().
 - The end user license agreement (EULA) is displayed.
 - 2) Accept the end user license agreement (EULA) by tapping the **Accept** button. If the end user license agreement (EULA) is not accepted, the Cockpit app cannot be used.
 - The welcome screen appears.
 - 3) Briefly press the  button on the control panel of the component to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
 - 4) Tap the **Add component** button.
 - The Connection Wizard opens and guides you through the process of establishing a connection.
 - 5) Follow the subsequent instructions on the screen.
 - 6) After the Bluetooth PIN is entered, a connection to the component is established.
 - While the connection is being established, 3 beep signals sound and the  symbol appears in the bottom right corner of the screen.
 - The  symbol is displayed when the connection has been established.
- Once the connection has been established, the data are read from the component. This process may take up to a minute.
Then the main menu appears with the name of the connected component.

INFORMATION

After the initial connection to the component has been established successfully and when the Bluetooth function is activated, the app will connect automatically each time it is started. No further steps are required.


INFORMATION

After activating the "visibility" of the component (briefly pressing the  button), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case the  button on the control panel must be briefly pressed again.

INFORMATION


If the connection to the component could not be established or the connection process is terminated, reset the component by connecting and disconnecting the battery charger.

7.3 Control elements for cockpit app

1. ☰ Access the navigation menu (see Page 71)
2. Name of the component. This name can only be changed using the adjustment app.
3. If connections to more than one component have been saved, you can switch between the saved components by tapping the **change** item (see Page 71).
4. MyModes configured with the adjustment app. Switching the mode by tapping the corresponding symbol and confirming by tapping "OK".
5. Currently selected mode
6. Charge level of the component.
 - 🔋 Component battery fully charged
 - 🔌 Component battery empty
 - 🔌 Component battery is being charged
 The current charge level is also displayed in %.
7. Display and name of the currently selected mode (e.g. **1. Basic Mode**)
8.
 - 📶 Connection to component has been established
 - 📶 Connection to component has been interrupted. The app is attempting to re-establish the connection automatically.
 - 📶 No existing connection to the component.

7.3.1 Cockpit app navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

Product

Name of the connected component

MyModes

Return to the main menu to switch MyModes

Settings

Change settings of the currently selected mode (see Page 77)

Status

Query status of the connected component (see Page 81)

Manage components

Add or delete components (see Page 71)

Imprint/Info

Display information/legal notice for the Cockpit app and the connected component

7.4 Managing components

Connections with up to four different components can be stored in the app. However, a component can only be connected to one mobile device at a time.

INFORMATION

Before establishing the connection, observe the points in the section "Initial connection between Cockpit app and component" (see Page 69).

7.4.1 Adding component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Briefly press the ⌘ button on the control panel of the joint unit to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the "+" button.
→ The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.
- 6) After the Bluetooth PIN is entered, a connection to the component is established.
→ While the connection is being established, 3 beep signals sound and the (i) symbol appears.
The (i) symbol is displayed when the connection has been established.
→ Once the connection has been established, the data are read from the component. This process may take up to a minute.
Then the main menu appears with the name of the connected component.

INFORMATION

If establishing a connection to a component is not possible, perform the following steps:

- ▶ Delete the component from the Cockpit app if applicable (see the section "Deleting a component").
- ▶ If present, delete the component from the Bluetooth settings on the device (unpair).
- ▶ Rest the component by connecting and disconnecting the battery charger.
- ▶ Add the component again in the Cockpit app (see the section "Adding a component").

INFORMATION

After activating the "visibility" of the component (briefly pressing the ⌘ button), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case the ⌘ button on the control panel must be briefly pressed again.

7.4.2 Deleting a component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Tap the "Edit" button.
- 4) Tap the ☒ symbol under the component you want to delete.
→ The component is deleted.

8 Use

INFORMATION

Before each use, check the product for functional reliability and for possible wear or damage.

Handling and use of the product must be learned before using it for the first time. Applying and removing, sitting down and standing up as well as walking must be practised.

In general, major discomfort should not be experienced when wearing the product. The affected limb should be checked periodically for areas of excessive pressure. If pressure points are observed, stop using the product and consult the O&P professional as soon as possible.

8.1 Application

INFORMATION

The individual approach to put on and take off the product should be developed in cooperation with the orthopaedic technician and/or therapist.

- 1) Open all orthosis shell closures.
- 2) Remove the shoe.
- 3) Sit on the front edge of a chair.
- 4) Flex the orthotic joint.
- 5) Insert the foot into the foot component. While doing so, position the heel and lower leg in the shell.
- 6) Slightly extend the leg and apply the orthosis to the thigh and lower leg.
- 7) Fasten the ankle closure if applicable.
- 8) Fasten the lower leg shell closure.
- 9) Fasten the thigh shell closure.
- 10) Retighten the top closure.
- 11) Put on the shoe.
- 12) Stand up from the chair and retighten all closures.
- 13) Check the correct fit of the orthosis.

⚠ CAUTION

Pinching skin in the area of the closures

Injuries and swelling due to disruption of blood circulation to the skin in the area of the closures.

- ▶ Do not over-tighten the closures when applying the product.

8.2 Removal


INFORMATION

Before sitting down, check whether adequate resistance is provided for support while sitting down. If a MyMode with very high or very low stance phase flexion resistance is activated, the Cockpit app must be used to switch back to basic mode before sitting down. Basic mode can also be activated by switching the product off and back on again.

- 1) Sit on a chair.
- 2) Open the orthosis shell closures.
- 3) Take off the orthosis.
- 4) Fasten the orthosis shell closures.
- 5) Store the orthosis and charge the battery if possible.

8.3 Movement patterns in basic mode (mode 1)

8.3.1 Standing

Knee control through high flexion resistance and static alignment.

The knee joint of the product has no locking function, so that slow bending is possible under flexion load. To restore the stable standing position, place the leg back under the body and put weight on the heel.

A stance function can be enabled by the O&P professional. Please see the following section for further information on the stance function.

8.3.1.1 Stance function

INFORMATION

To use this function, it must be enabled by the O&P professional. It also has to be activated using the Cockpit app (see Page 78).

The stance function is a functional supplement to the basic mode. This function makes it easier for the user to stand on an inclined surface for a longer time. The joint is fixed in the flexion direction at a flexion angle between 5° and 65°.

The standing function must be enabled by the O&P professional. In addition, the type of joint locking (intuitive/manual) has to be established by the O&P professional. Changing the locking type using the Cockpit app is not possible.

Intuitive locking of the joint

The intuitive stance function recognises any situation that puts strain on the orthosis in the flexion direction but where flexion is not permitted. Examples of this include standing on uneven or sloping surfaces. The joint is always locked in the flexion direction when the leg with the orthosis is not fully extended and is kept still for a brief moment.

The joint is not locked when the above conditions are met and a sitting position is assumed.

Intuitive unlocking of the joint

► Upon forward or backward rollover or extension, the flexion resistance is immediately switched from high to stance phase resistance again.

Manual joint lock


► Flex the joint between 5° and 60°.

→ Now the blocked joint can be loaded in the flexion direction.

Releasing the manual joint lock

► The Manual Stance function is automatically deactivated again by extending the knee or by repositioning the leg (e.g. taking a step).

8.3.2 Walking

Initial attempts at walking with the product always require instruction from trained, qualified personnel.

The hydraulics stabilise the knee joint in the stance phase and release the knee joint in the swing phase so that the leg can swing forward freely.

Switching to the swing phase requires rolling the product over to the front out of the stride position.

An acoustic feedback signal when initiating the swing phase can be activated by the O&P professional.

INFORMATION

Changing the gait pattern influences activation of the swing phase. Adjustment by the O&P professional has to be repeated in this case.

8.3.3 Sitting down

INFORMATION

Before sitting down, check whether adequate resistance is provided for support while sitting down. If a MyMode with very high or very low stance phase flexion resistance is activated, the Cockpit app must be used to switch back to basic mode before sitting down. Basic mode can also be activated by switching the product off and back on again.

The resistance in the product's knee joint while sitting down ensures even lowering into the sitting position.


- 1) Position both feet side by side at the same level.
- 2) While sitting down, distribute your weight evenly between both legs and use the arm supports where applicable.
- 3) Move the buttocks in the direction of the back support and lean the upper body forward.

INFORMATION: Resistance while sitting down can be changed with the Cockpit app via the parameter "Stance flexion resistance" (see Page 78).

8.3.4 Sitting

INFORMATION

While sitting, the knee joint also switches to energy saving mode. This energy saving mode is activated regardless of whether the sitting function is activated or not.

When a sitting position is maintained for more than two seconds (i.e. the thigh is close to horizontal and there is no load on the leg), the product switches the resistance in the extension direction to a minimum.

If the sitting function has been activated by the O&P professional and switched on via the Cockpit app (see Page 78), the resistance in the flexion direction is also reduced.

8.3.5 Standing up

INFORMATION

If a MyMode with a very high extension resistance (component remains flexed while standing up) or very low flexion resistance (no support as expected) is activated, it is necessary to switch back to basic mode. Basic mode can also be activated by switching the product off and back on again.

- 1) Place the feet at the same level.
- 2) Lean the upper body forward.
- 3) Put the hands on armrests, if available.
- 4) Stand up with support from the hands, while keeping weight evenly distributed on the feet.


8.3.6 Walking down stairs

INFORMATION

Before walking down stairs, check whether adequate resistance is provided for walking down stairs. If a MyMode with very high or very low stance phase flexion resistance is activated, the Cockpit app must be used to switch back to basic mode before walking down stairs. Basic mode can also be activated by switching the product off and back on again.

Depending on the construction of the orthosis (rigid or with movement in the ankle joint possible), the movement pattern is carried out as follows.

8.3.6.1 Orthosis construction with rigid ankle joint or dorsal spring element

This function must be practised and executed consciously. Only when the sole is properly positioned can the product react correctly and permit controlled flexion.


- 1) Hold the handrail with one hand.
- 2) Position the leg with the product on the step so that the foot projects halfway over the edge of the step.
→ This is the only way to ensure a secure rollover.
- 3) Roll the foot over the edge of the step.
→ This flexes the product slowly and evenly at the knee joint.
- 4) Place the foot of the less affected leg onto the next step.

8.3.6.2 Orthosis construction with movement in the ankle joint

This function must be practised and executed consciously. Only by properly stepping down with the sole can the product respond correctly and permit controlled flexion.

- 1) Hold the handrail with one hand.
- 2) Position the leg with the product on the step so that as much of the sole of the foot as possible is on the step.
- 3) Place the foot of the less affected leg onto the next step.

8.3.7 Walking up stairs


Walking up stairs one step at a time

- 1) Hold the handrail with one hand.
- 2) Place the foot of the less affected leg onto the first step.
- 3) Pull up the leg with the product.

Walking up stairs step-over-step


The product does not have an active drive system to support walking up stairs step-over-step. This is possible given certain physical prerequisites (residual muscular function present in the affected leg) and with corresponding training.

8.3.8 Walking up a ramp

- 1) Hold the handrail with one hand.
- 2) Place the foot of the less affected leg onto the ramp.
- 3) Move the leg with the product forward.

8.3.9 Walking down a ramp

Use the handrail if possible.


Slight slope (< 5-10%)

The sequence of movements is equivalent to walking on the level. Swing phase initiation is possible.

Moderate/steep slope (> 5-10%)

The sequence of movements is similar to walking down stairs. Place a load on the orthosis, permit knee flexion against the flexion resistance and roll over the forefoot.

8.3.10 Walking down flat steps

A (possibly surprising) swing phase may be triggered when walking down a flat step (such as a curb). The user can either use the stance phase flexion resistance or initiate a swing phase.

8.3.11 Kneeling

The individual approach to kneeling down and getting up again should be developed in cooperation with the O&P professional and/or therapist. Kneeling down is supported by increased flexion resistance, which permits controlled flexion of the knee joint.

8.4 Changing orthosis settings

Once an active connection to a component has been established, the settings **of the respective active mode** can be changed using the Cockpit app.

INFORMATION


Bluetooth of the component must be switched on to change the orthosis settings. Briefly press the  button on the control panel to check.

If the  symbol is not lit up, Bluetooth has to be activated by a long press of the  button on the control panel.

Information for changing the orthosis settings

- Before changing settings, always check the main menu of the Cockpit app to make sure the correct component has been selected. Otherwise parameters could be changed for the wrong component.
- It is not possible to change orthosis settings nor to switch to a different MyMode while the orthosis battery is being charged. Only the status of the orthosis can be called up. Instead of the  symbol, the  symbol appears in the bottom row of the screen in the Cockpit app.
- If the settings of a MyMode are to be modified, one must first switch to this MyMode.

8.4.1 Changing the orthosis settings using the Cockpit app

- 1) Once the component is connected and in the desired mode, tap the  symbol in the main menu.
→ The navigation menu opens.
- 2) Tap the menu entry "**Settings**".
→ A list appears with the parameters for the currently selected mode.
- 3) Change the setting of the desired parameter by tapping the "<", ">" symbols.

INFORMATION: The prosthetist's setting is marked and, after the setting has been changed, can be restored by tapping the "Standard" button.

8.4.2 Overview of adjustment parameters in basic mode

The parameters in basic mode describe the dynamic behaviour of the orthosis in the gait cycle. These parameters act as basic settings for automatically adjusting the damping behaviour to the current motion situation (e.g. ramps, slow walking speed, etc.).

The stance function and/or the sitting function can also be activated/deactivated. Further information on the stance function (see Page 74). Further information on the sitting function (see Page 75).

The following parameters can be modified:

Parameter	Adjustment range of the O&P professional via the adjustment app	Cockpit app adjustment range	Meaning
Stance flexion resistance	120 to 190	+/- 10 of the configured value	Resistance against flexion motion e.g. when walking down stairs or while sitting down

Parameter	Adjustment range of the O&P professional via the adjustment app	Cockpit app adjustment range	Meaning
Stance function		0/Off - deactivated 1/On - activated	Activation/deactivation of the stance function. This function needs to be enabled in the adjustment app.
Sitting function		0/Off - deactivated 1/On - activated	Activation/deactivation of the sitting function. This function needs to be enabled in the adjustment app.
Volume of acoustic signal	0 to 100	0 to 100	Volume of beep signal for confirmation tones (e.g. when switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals are still generated if errors occur.
Pitch of acoustic signal	1500 Hz — 3000 Hz	1500 Hz — 3000 Hz	Pitch of beep signal for confirmation tones

INFORMATION

Confirmation of the successful transfer of the parameters

A beep and vibration signal is emitted by the orthosis while changing the parameters using the Cockpit app. If the parameter "**Volume of acoustic signal**" is set to "0", only a vibration signal is emitted.

8.4.3 Overview of adjustment parameters in MyModes

The parameters in the MyMode "**User defined**" describe the static behaviour of the orthosis for a specific movement pattern such as cycling. There is no automatically controlled adjustment of the extension and flexion resistance.

The parameters of the MyModes "**Training mode**" and "**Freeze position**" are preconfigured and cannot be changed.

The following parameters can be modified in MyModes:

Parameter	Adjustment app setting range	Cockpit app adjustment range	Meaning
Flexion resistance increase	0 to 100	+/- 10 of the configured value	Value for speed at which flexion resistance increases as the knee angle increases This parameter can only be changed in the MyMode " User defined ".
Volume of acoustic signal	0 to 100	0 to 100	Volume of beep signal for confirmation tones (e.g. when switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals are still generated if errors occur.

Parameter	Adjustment app setting range	Cockpit app adjustment range	Meaning
Pitch of acoustic signal	1500 Hz — 3000 Hz	1500 Hz — 3000 Hz	Pitch of beep signal for confirmation tones

INFORMATION

Confirmation of the successful transfer of the parameters

A beep and vibration signal is emitted by the orthosis while changing the parameters using the Cockpit app. If the parameter "**Volume of acoustic signal**" is set to "0", only a vibration signal is emitted.

8.5 Switching the product on/off

In certain cases, e.g. for storage or transportation, the product can be purposely switched off.

CAUTION

Using the product while switched off

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Prior to using the product, briefly press the  button on the control panel to check whether it is switched on. The button symbol lights up green  when the product is switched on.

Switching off

- 1) Briefly press the  button on the control panel to check whether the product is switched on. The symbol in the button lights up green  and the current charge level is displayed (see Page 89).
- 2) With the product switched on, press and hold the  button on the control panel until the symbol in the button turns off and a descending series of beeps (shut-down melody)  is emitted.

Switching on

- ▶ Press the  button on the control panel to switch the product on.
- A long vibration signal followed by a short beep is emitted, and the current status is displayed for about 5 seconds (see Page 89).
- When the symbol in the  button is lit up green , the product is switched on and operational.
- Basic mode is activated after switching on.

INFORMATION

No display after switching on

If the  symbol does not light up after a long press of the  button, the battery may be deep discharged. Charging for at least 15 minutes is required in this case.

8.6 Switching Bluetooth of the component on/off

INFORMATION

Bluetooth on the component must be turned on in order to use the Cockpit app. Briefly press the  button on the control panel to check. In the  button, the symbol must be lit up blue .

If the symbol in the button is not lit up, Bluetooth has to be switched by a long press of the  button.

8.6.1 Switching off Bluetooth

- ▶ With the Bluetooth function activated, press and hold the  button on the control panel until a vibration signal is emitted and the symbol in the button turns off.
 - Bluetooth is switched off.
- To check whether the Bluetooth function was deactivated correctly, query the status by pressing the  button (see Page 89).

8.6.2 Switching on Bluetooth

- ▶ With the Bluetooth function deactivated, press and hold the  button on the control panel until a short vibration signal is emitted and the symbol in the button lights up blue .
- Bluetooth is switched on.
- To check whether the Bluetooth function was activated correctly, query the status by pressing the  button (see Page 89).

8.7 Querying the status of the component

8.7.1 Query status through cockpit app

- 1) When the component is connected, tap the  symbol in the main menu.
- 2) In the navigation menu, tap the entry "**Status**".

8.7.2 Status display in the cockpit app

Menu option	Description	Possible actions
Steps per day: 25	Daily step counter	Reset the counter by tapping the " Reset " button.
Overall steps: 1745	Total steps since most recent maintenance	For informational purposes only
Batt.: 68	Current charge level of the component in percent	For informational purposes only

8.8 Recommendations for air travel

Observing the following points is recommended prior to air travel or on board the airplane:

- Take along the 647F558 orthotic passport so you can show it on request or in case of questions.
- Turn off the component's Bluetooth function on board the aircraft if required (see Page 80).
- Take along the corresponding power supply adapter according to your destination. The power supply is suitable for connecting to alternating voltages from 100 V to 240 V at a network frequency of 50 Hz to 60 Hz.

9 MyModes

In addition to the MyModes "Training mode" and "Freeze position", the O&P professional can configure a MyMode "User defined" using the adjustment app. These MyModes can be selected using the Cockpit app.

9.1 Switching MyModes with the cockpit app


Once a connection to an orthosis has been established, the Cockpit app can be used to switch between the MyModes.

Information on switching

- It is always possible to switch back to basic mode by turning the product off and back on (see Page 80) and by connecting/disconnecting the battery charger.
- Before the first step, always check whether the selected mode corresponds to the required motion type.
- Check whether the battery charger is connected. Changing the mode is not possible when the battery charger is connected, even though the  symbol is lit up.

- On the device, check whether a connection to the component has been established. The (••) symbol must be lit up.
- On the device, check whether the correct component has been selected.

Switching process

- 1) Tap the symbol of the desired MyMode (1) in the main menu of the Cockpit app.
→ A security question for changing the MyMode appears.
- 2) If you want to change the mode, tap the "OK" button.
→ A beep signal sounds to confirm the switch.
- 3) After switching, a symbol (2) is displayed to identify the active mode.
→ The current mode is also indicated by the name on the lower edge of the screen (3).

9.2 Basic mode


This mode is intended for daily use.

9.3 MyMode "Training mode"

(if the knee joint is locked during the stance phase and swing phase release is possible)

The knee joint is locked in the stance phase. Swing phase initiation is possible as in basic mode. In this mode, walking down stairs or a ramp is only possible one step at a time.

Switching process


- 1) On the screen of the device, tap the corresponding symbol to activate training mode.
→ A confirmation signal sounds to indicate the switch to training mode.
- 2) The hydraulics stabilise the joint with high flexion resistance in the stance phase and release the joint in the swing phase so that the leg can swing forward freely.

Exit training mode by selecting a different MyMode or by turning the component off and then on again.

9.4 MyMode "Freeze position"

(if the knee joint is permanently locked)

Switching process

In this MyMode, the joint is locked in its current position and cannot be moved in the flexion nor the extension direction.


- 1) Flex or extend the knee joint to the angle where it should be locked.
- 2) On the screen of the device, tap the corresponding symbol to activate the lock.

→ A confirmation signal sounds to indicate activation of the lock.

- 3) The orthosis is locked in the flexion and extension directions.

Deactivate the lock by selecting a different MyMode/basic mode or by turning the component off and then on again.

9.5 MyMode "User defined"

The MyMode "**User defined**" is intended for specific motion patterns or postures (e.g. cycling ...). This MyMode can be configured individually by the O&P professional using the adjustment app.

Settings can be adjusted using the Cockpit app (see Page 79).

10 Additional operating states (modes)

10.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level drops to 5% (see Page 89). The settings are changed to the safety mode parameters during this time. Then the product is switched off.

INFORMATION

After disconnecting the battery charger, the orthosis is in the same state it was in before connecting the battery charger. For example, if the orthosis was switched off before connecting the battery charger, it will also be switched off after disconnecting the battery charger. A descending sequence of beeps is emitted when the battery charger is disconnected while the orthosis is switched off .

10.2 Mode for charging the product

The product is non-functional during charging.

To switch to basic mode, the battery charger for the product must be disconnected after the battery is charged.

10.3 Safety mode

The product automatically switches to safety mode if a critical fault occurs (e.g. failure of a sensor signal). Safety mode remains in effect until the error has been rectified.

Flexion resistance preconfigured by the O&P professional is activated in safety mode, corresponding at least to the stance phase resistance. The extension resistance is low and cannot be changed. This makes limited walking possible for the user and allows the user to sit down, even though the sensor system is not active.

The switch to safety mode is indicated by the red **!** symbol on the control panel and by beeps and vibration signals immediately prior to switching (see Page 89).

The current safety mode can be reset by switching the product off/on (see Page 80) and by connecting then disconnecting the battery charger. If the product switches into safety mode again or at a later time, this means a permanent error exists. The product must be inspected immediately by an O&P professional.

10.4 Overheating mode

When the hydraulic unit overheats due to uninterrupted, increased activity (e.g. extended walking downhill), the flexion resistance is increased along with the rising temperature in order to counteract the overheating. When the hydraulic unit cools down, the product switches back to the settings that existed before the overheating mode.

Overheating mode is not activated in the MyMode "User defined".

Overheating mode is indicated by a long vibration every 5 seconds and additionally by the orange **!** symbol lighting up.

The following functions are deactivated in overheating mode:

- Sitting function
- Switching to a MyMode
- Changing the orthosis settings

11 Cleaning

1) Clean the product with a damp cloth and mild soap when needed.

2) Dry the product with a lint-free cloth and allow it to air dry fully.

Orthosis shell padding

1) Remove the padding materials from the orthosis shells.

2) Fasten all hook-and-loop closures, if any.

3) If the padding materials "623P3 terry cloth padding fabric" or "623F62 SpaceTex padding fabric" are used, hand wash these in warm water at 30 °C/86 °F with a standard mild detergent.

If other padding materials are used, observe the care symbols and cleaning instructions for those materials.

4) Rinse out detergent residues thoroughly.

5) Allow to air dry. Do not expose to sources of direct heat (e.g. sunlight, stove or radiator).

6) After the padding material is completely dry, reattach it to the orthosis shells on the correct sides.

12 Maintenance

Regular maintenance (service inspections) is mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty, maintain basic safety and the essential performance characteristics, and ensure safety in regards to EMC.

When maintenance is due, this is indicated by feedback after connecting the battery charger (see the section "Operating states/error signals", see Page 89). The manufacturer grants a grace period of no more than one month before, or two months after, the due date.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components always have to be submitted to the O&P professional for maintenance and repairs:

The orthosis, battery charger and power supply unit.

17KO1=* C-Brace joint unit

The manufacturer has prescribed a service inspection for the product no later than every 24 months or upon reaching one million steps at the latest, whichever comes first.

Additional orthosis components used, e.g. joints

Observe the service intervals and maintenance instructions for all orthosis components used.

13 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

13.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

13.2 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA
Phone + 1-801-956-2400
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.

5470 Harvester Road

L7L 5N5 Burlington, Ontario

Canada

Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

13.3 CE conformity

Only applicable for "17KO1=L C-Brace joint unit"/"17KO1=R C-Brace joint unit"

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

13.4 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

14 Technical data

Ambient conditions	
Transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage in the original packaging (≤3 months)	-20 °C/-4 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Long-term storage in the original packaging (>3 months)	-20 °C/-4 °F to +25 °C/+77 °F Max. 93% relative humidity, non-condensing

Ambient conditions	
Transport and storage between applications (without packaging)	-25 °C/-13 °F to +35 °C/95 °F Max. 93% relative humidity, non-condensing +35°C/95°F to +70°C/158°F Water vapour pressure up to 50 hPa
Operation	-10 °C/+14 °F to +40 °C/+104 °F Relative humidity 15% to 93%, non-condensing, at a water vapour pressure up to 50 hPa Air pressure: 606.3 hPa (up to 4000 m without pressure equalisation)
Maximum temperature that can be reached on the surfaces of the orthosis during operation	+44 °C/+111 °F
Time for warming to the operating temperature after storage between applications, from -25 °C/-13 °F at an ambient temperature of +20 °C/+68 °F	30 minutes
Time for cooling to the operating temperature after storage between applications, from +70 °C/+158 °F at an ambient temperature of +20 °C/+68 °F	30 minutes
Charging the battery	+10 °C/+50 °F to +40 °C/+104 °F

General information	
Reference number	17KO1=L C-Brace joint unit left/17KO1=R C-Brace joint unit right
Joint unit weight [g/oz]	Approx. 1000 / 35
Maximum user body weight [kg/lbs]	125 / 276
Product service life [years]	6
Information on the product's ruleset and firmware version	Accessible via the Cockpit app navigation menu and the menu item "Imprint/Info"

Data transfer	
Wireless technology	Bluetooth Smart Ready
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK, $\pi/4$ DQPSK, 8DPSK
Data rate (over the air)	2178 kbps (asymmetrical)
Maximum output power (EIRP):	+8.5 dBm

Battery charger	
Reference number	4E50-2
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F Max. 93% relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Input voltage	12 V \equiv
Lifetime	8 years

Power supply unit	
Reference number	757L16-4
Type	FW8001M/12
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Operation	0 °C/+32 °F to +50 °C/+122 °F Max. 95% relative humidity Air pressure: 70–106 kPa (up to 3,000 m without pressure equalisation)
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ==

Orthosis battery	
Battery type	Li-Ion
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500
Charge level after 1 hour charging time	30 %
Charge level after 2 hours charging time	50 %
Charge level after 4 hours charging time	80 %
Charge level after 8 hours charging time	Fully charged
Behaviour of the orthosis while being charged	The orthosis has no function
Operating time of the orthosis with new, fully charged battery at room temperature	At least 18 hours of uninterrupted walking Approx. 2 days with average use

Cockpit app	
Reference number	4X441-IOS=*/4X441-Andr=V* Cockpit
Supported operating system	iOS 10.0/Android 5.0 or higher
Website for download	https://www.ottobock.com/cockpitapp

15 Appendices

15.1 Symbols Used

Please note the instructions for use

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Type BF applied part

Compliance with the requirements according to "FCC Part 15" (USA)

Non-ionising radiation

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Declaration of conformity according to the applicable European directives

Manufacturer

IP54

Protected against dust, protected against splashing water

Serial number (YYYY WW NNN)

YYYY – year of manufacture

WW – week of manufacture

NNN – sequential number

Article number


Medical device

15.2 Operating states/error signals

The orthosis indicates operating states and error messages with symbols on the control panel and through beeps and vibration signals.

15.2.1 Status display on the control panel

The current status of the component is displayed on the control panel for 5 seconds in the following cases:

- The  button on the control panel was pressed briefly.
- The component was switched on by pressing the  button.
- The battery charger was disconnected from the component.
- Battery charger was connected to the component.
- An error was identified during use.

INFORMATION

No status display due to battery deep discharge

If the status of the component is not displayed on the control panel, the battery may be deep discharged. In this case it must be charged for at least 15 minutes before the status can be queried again.

Symbols on the control panel	Beep signal	Vibration signal	Event	Required action
All symbols on the control panel light up alternately	—	—	Test of the indicators (LEDs) after connecting the battery charger	Check whether all symbols (LEDs) light up alternately and in the corresponding colours. If a symbol (LED) does not light up in a colour, the product should be inspected by an O&P professional. If no symbols (LEDs) light up, the battery may be deep discharged. Leave the battery charger connected for at least 15 minutes and then repeat this test by disconnecting/connecting the battery charger.
	—	—	The product is switched on and ready for operation	
	1 x short	1 x long and 1 x short	The product was switched on by pressing the  button on the control panel.	
	—	1 x long at interval of approx. 5 seconds	Overheated hydraulics	Reduce activity.
	—	—	A self-test error was detected when connecting the battery charger.	<ul style="list-style-type: none"> Conduct the self-test again by connecting/disconnecting the battery charger. If the  symbol lights up again, the product must be inspected by an O&P professional.

Symbols on the control panel	Beep signal	Vibration signal	Event	Required action
	30 x long	30 x long	Severe error/indication of safety mode activation (see Page 83) Switch to safety mode if possible.	Walking possible with restrictions. Please note the possible change in flexion/extension resistance. 1. Attempt to reset this error by switching the product off/on (see Page 80). 2. If the beep/vibration signal recurs, attempt to reset this error by connecting/disconnecting the battery charger. 3. If the beep/vibration signal recurs, use of the product is prohibited. The product must be inspected by an O&P professional.
	—	—	Charge level 10% to 34%	
	—	—	Charge level 34% to 67%	
	—	—	Charge level 67% to 100% Indicates that the battery is fully charged during the charging process.	
	—	—	Battery is being charged, charge level is less than 34%	
	—	—	Battery is being charged, charge level is 34% to 67%	
	—	—	Battery is being charged, charge level is 67% to 99%	
	3 x long	3 x long	Charge level between 5% and 10%	Charge battery soon. Remaining operating time approx. 1 hour.
	5 x long	5 x long	Charge level between 0% and 5%	Charge battery immediately; the product will be switched off after the next warning signal.
	10 x long	10 x long	Charge level 0% After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.

Symbols on the control panel	Beep signal	Vibration signal	Event	Required action
	4 x short at intervals of approx. 65 sec. (continuously)	—	Charging the battery outside the allowable temperature range	Check whether the specified ambient conditions for charging the battery were met (see Page 86).
	—	—	Maintenance required because the maintenance date has already been reached or exceeded.	You should visit the O&P professional soon. They will forward the component to an authorised Ottobock Service Center.
	—	—	Bluetooth activated	
—	1 x long	—	Battery charger is connected	
—	—	3 x short	Charging mode started (3 sec. after connecting battery charger)	
—	1 x short	1 x short	Mode switching or changing adjustment parameters is performed using the Cockpit app. If the parameter " Volume of acoustic signal " is set to "0" in the Cockpit app, only a vibration signal is emitted.	
—		—	The product switches itself off. This occurs automatically in the following cases: <ul style="list-style-type: none"> The  button on the control panel was pressed and held for longer than approx. 3 seconds while product was switched on. After the  symbol lights up. After disconnecting the battery charger, if the product was already switched off before connecting the battery charger. 	<ul style="list-style-type: none"> Charge the battery. If desired, switch on the product with the  button.

Symbols on the control panel	Beep signal	Vibration signal	Event	Required action
—	—	Ongoing	Total failure Electronic control no longer possible. Safety mode active or undetermined valve state. Undetermined product behaviour.	<ol style="list-style-type: none"> 1. Press and hold the  button on the control panel until the vibration signal stops (approx. 10 seconds), completely switching off the product. 2. If the vibration signal recurs after switching on, attempt to reset this error by connecting/disconnecting the battery charger. 3. If the vibration signal recurs, use of the product is prohibited. The product must be inspected by an O&P professional.

15.2.2 Error messages while establishing a connection with the cockpit app

Error message	Cause	Correction
Component was connected to another device. Establish connection?	The component was connected to another device	To disconnect the original connection, tap the " OK " button. If the original connection is not to be disconnected, tap the " Cancel " button.
Mode change failed	An attempt was made to switch to a different MyMode while the component was in motion (e. g. while walking)	For safety reasons, switching MyModes is only permitted when components are at rest, e.g. while standing or sitting.
	A current connection to the component was interrupted	<p>Check the following points:</p> <ul style="list-style-type: none"> • Distance from the component to the device • Charge level of the component's battery • Bluetooth of the component switched on? (see Page 80) • Briefly press the  button on the control panel to make the component "visible" for 2 minutes. • Component switched on? (see Page 80) • If multiple components were stored, was the correct component selected? • Is the component still connected to another device, and is this connection still active?

15.2.3 Error while charging the product

LED on power supply	LED on battery charger	Error	Resolution
○	 ○ ○ 	Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.
		Non-functional socket	Check socket with another electric device.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	 ○ ○ 	No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	 ○ ● 	Battery is fully charged (or connection with product is interrupted).	<p>Take note of the confirmation signal for differentiation.</p> <p>When the battery charger is connected, a self-test is conducted and confirmed by beeps/vibration signals.</p> <p>When these signals are emitted, the battery is fully charged.</p> <p>If no signals are emitted, the connection to the product is interrupted.</p> <p>If the connection to the product is interrupted, an authorised Ottobock Service Centre must inspect the product, battery charger and power supply.</p>

15.3 Directives and manufacturer's declaration

15.3.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (see Page 60).

Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage drops	IEC 61000-4-11	0% U _T ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% U _T ; 1 period and 70% U _T ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% U _T ; 250/300 periods

Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/90-0, TETRA 800, iDEN 820, CDMA 850, GSM 800/90-0, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 to 2,570	Bluetooth WLAN 802.1- 1 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.1- 1 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

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Legal Manufacturer of Custom Made Device

The product "C-Brace joint unit 17KO1=" is covered by the following patents:

USA: Patent US 9,022,965
European Patent: EP 2276433 in CH, DE, FR, IT, NL, SE
EP 2380529 in CH, DE, FR, IT, NL, SE
Germany DE 10 2018 111 234, DE 10 2017 13139

Patents pending in Germany, Brazil, USA, China, South Korea, Israel, Japan, Russia, Europe and as International Application.

The product "C-Brace joint unit 17KO1=" is covered by the following registered designs and design patents:

Australia: 201717600; 201810549
China: ZL 201730629343.9
European Design: No.004043412
Russia: 111643
Switzerland DM/098883
Turkey: DM/098883
India: 62940
Hongkong: 1702435.3
USA: Patent US D871,595

Legal Manufacturer of C-Brace joint unit 17KO1=
Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com